


# BD Integra™ Syringe

with retracting BD PrecisionGlide™ Needle

## Retract and protect with speed



## The latest innovation in retracting technology from BD

**The BD Integra™ Syringe is designed to have low waste space.  
Low waste space may result in reduced medication waste.\***

- The BD Integra™ Syringe provides one-handed activation.
- Can be activated before or after the needle is withdrawn from the patient.

Scan this code to watch a video on  
how to use the **BD Integra™ Syringe**:



\* Compared to a traditional syringe and needle



# BD Integra™ Syringe

with retracting BD PrecisionGlide™ Needle

BD Integra™ 3 mL Syringes with detachable needle		
Product number	Description	Packaging
305269	25 G x 1/2 in.	100/box, 400/case
305270	25 G x 1 in.	100/box, 400/case
305271	23 G x 1 in.	100/box, 400/case
305273	21 G x 1 in.	100/box, 400/case
305274	21 G x 1 1/2 in.	100/box, 400/case

Note: BD Integra™ Needle and Syringe do not have a standard luer-lock connection and thus conventional luer-lock needles and syringes cannot be used with BD Integra™ Needles and Syringes.

For more information on our products, please contact Medical Support Services at **800.555.7422** or **medical.services@bd.com**

To speak with your BD sales representative, call **888.237.2762**

**bd.com**

