

 HENRY SCHEIN®
MEDICAL



VALUE

SELECTION



QUALITY

Rely On
HENRY SCHEIN
BRAND PRODUCTS
FOR YOUR AESTHETIC PRACTICE

Rely on Us



{ *Quality, Selection and Value You Can Rely On* }

Now, more than ever, you require value-priced products for your practice. At Henry Schein Medical, you can rely on us to provide diverse product options that will help you achieve maximum profitability while maintaining patient satisfaction throughout your practice. We are dedicated to constantly improving and expanding our selection of items to meet the growing needs of our customers.

In this special publication of our Henry Schein Brand products, you will find a wide variety of choices focused on everyday essentials, including infection control supplies for instrument reprocessing and overall protection; exam room disposables from table paper to drapes, capes, and gowns; surgical and exam gloves featuring the popular Criterion line; and wound care items such as gauze, dressings, and bandages.

Henry Schein is a name you can trust. The Henry Schein Brand product line has met our high standards of quality and effectiveness—earning the Henry Schein Seal of Excellence—your guarantee of satisfaction. We invite you to try our line and experience the value, selection, and quality that we stand behind with an 80-year history of serving health care professionals around the world.



Celebrating 80 Years of Excellence

Henry Schein is proud of its history. From a New York pharmacy to a Fortune 300 company—we are dedicated to serving our customers, supplier partners, investors, society, and Team Schein.





Henry Schein Medical is your reliable network of trusted advisors.

Rely on Us for...

Expert Guidance

Our knowledgeable team of trusted advisors is dedicated to collaborating with you to meet specific goals related to quality patient care; staff and team growth; and profitability.

- Revenue Enhancement Tools
- Cost Containment Programs
- Staff Training and Team Building
- Practice Marketing
- Technology Solutions

Price Integrity

We have relationships at the local, regional, and national level—including all major GPOs—to help manage operating costs and drive value for your practice. Web-enabled reporting tools guarantee your pricing is accurate and consistent.

Guaranteed Fast Delivery

All orders received by 5 PM local time are shipped the same day. Your products arrive within 24–48 hours.

Vaccines, Pharmaceuticals, and Diagnostics

We supply every adult and pediatric vaccine recommended by the ACIP, advancing prevention and wellness in your patient community. We meet all of your flu vaccination needs with our Worry-Free Guarantee Influenza Program.

Our pharmaceutical and diagnostic solutions include brand names as well as a complete line of generics and Henry Schein Brand options for added value.





WHAT'S NEW?



HENRY SCHEIN® SURGICAL SKIN MARKERS

These surgical marking pens provide an accurate and reliable way to mark the incision site prior to surgery.

- Prep-resistant
- Nonabrasive tip
- Large ink reservoir

- Skin Markers
Non-Sterile, Fine Taper Tip, Ruler and up to 9 labels
(900-7018)100/box
- Sterile, Fine Taper Tip, Ruler and up to 6 labels
(900-7019)100/box
- Mini (WriteSite Jr.), Dark Blue
(900-7020)50/box

HENRY SCHEIN® SUTURE REMOVAL KIT



All items are sterile and latex-free.
(112-5822)ea

Contains: 1 small-loop, 4.5" straight metal iris scissor; 1 serrated, 4.75" metal forceps, 1-3" x 3" gauze pad; 1 large alcohol prep pad.

HENRY SCHEIN® SUTURE REMOVAL KIT-PLASTIC



Sterile, single-use kit. Plastic-handled Littauer scissors with plastic thumb forceps. Flatter, thinner-blade scissors provide easier suture removal by allowing blade to slip under suture.
(900-4687)ea

Contains: 1 pair of Littauer (hook-tipped) scissors, 1 pair of PosiGrip™ plastic forceps & 1-3" x 3" gauze sponge.

HENRY SCHEIN® DISPOSABLE SKIN STAPLER -35 WIDE



Our disposable skin stapler is specially designed for surgical skin suturing. With novel structure and application security, it has the advantages of quick speed of sewing up during the operation, little histological reaction, and good wound matching and healing. Leaves a small surgical scar; no pain during removal. Enables you to not only work with higher efficiency but also assists in avoiding cross-infection effectively.

- Skin Stapler
(900-7010)ea
- Skin Staple Remover
(900-7011)ea



HENRY SCHEIN® POVIDONE IODINE 10%

- Prep Pad, Medium
(112-5637)100/box
- Swab Stick 1's
(112-5638)50/box
- Swab Stick 3's
(112-5639)25/box

HENRY SCHEIN® GRADUATED PLASTIC MEDICINE CUPS



Plastic medicine cups with easy-to-read gradations.
1 oz
(112-5555)100/bag

HENRY SCHEIN® CRITERION® N200



Nitrile Powder-Free Exam

- Eliminates type-1 allergic reactions associated with natural latex
- Innovative formulation emulates the fit and feel of latex with added strength and puncture resistance
- Eco-friendly, 200-count box increases your storage by 50%, reduces cardboard waste, and still fits your glove box holders. You save money while reducing your environmental impact
- "Real-feel" grip provides optimal tactile sensitivity. Textured fingers provide enhanced, superior grip in wet and dry applications
- Easy donning: very smooth inside, so glove slides onto your hand with ease. No need to fight with it when you are in a hurry

- Gloves200/box
- Specify:
- X-small **(900-7437)**
- Small **(900-7438)**
- Medium **(900-7439)**
- Large **(900-7440)**
- X-large (180/box) **(900-7441)**

HENRY SCHEIN® URINE CONTAINER



This 4.5-oz container has a blue, tamper-evident, leak-resistant top, and is suitable for routine urine or specimen collection and transport.

Sterile
(900-4803)100/case

Non-Sterile
(900-7432)500/Cs

HENRY SCHEIN® EMESIS BASIN



Plastic, kidney shape, 16-oz capacity.
8.5"
(112-5809)25/box

GLOVE CLUB®



NYTRILE™ X200

Nitrile Powder-Free Exam

- Excellent value when compared to similar nitrile gloves
- Soft nitrile formulation contours to hand for comfort
- Great tactile sensitivity
- Textured fingers to provide superior grip in wet and dry conditions
- Eco-friendly (200/box) packaging

Gloves200/box

Specify:

X-small **(112-5658)**

Small **(112-5659)**

Medium **(112-5660)**

Large **(112-5661)**

X-large 180/box **(112-5662)**



HENRY SCHEIN® UNDERPADS

Disposable underpads made of waterproof blue poly backing, absorbent core of tissue layers and a white, hypoallergenic spunbond facing sheet.

17" x 24"
(112-5633)300/case

23" x 24"
(112-5634)200/case



HENRY SCHEIN® TOP HAT SPECIMEN COLLECTOR

800 cc
(112-5810)25/box

MEDICAL WEB SITE



Quick & Easy Ordering!

- 1 Login
- 2 Order
- 3 Checkout

What dynamic features does the Henry Schein Medical Web site offer your practice?

- 24 hour/365 days availability to supply ordering
- Track orders and shipments and process returns
- Generate comprehensive utilization reports via our dynamic report wizard
- Set budgets and electronically approve orders for greater security and accuracy
- Access purchase history and real-time pricing
- Custom shopping list capability for easy ordering
- Paperless billing statement option
- Shop by procedure, chronic conditions, or simple everyday essentials
- Access free MSDS documents for quick OSHA compliance
- Retrieve electronic statements and invoices for paperless environment
- 6 different specialty Web sites including Family Practice, Internal Medicine, and Pediatrics
- Access to promotions created uniquely for your specialty



www.henryschein.com/hsibrand



HENRY SCHEIN® HAND SANITIZER

- Alcohol-based waterless hand gel
- Kills 99.9% of most common disease-causing germs spread by the hands
- Rinse-free; no water required
- Fresh clean scent and soothing emollients leave hands soft and clean

8-oz Pump Bottle
(900-4443)ea

16-oz Pump Bottle
(900-4442)ea

HENRY SCHEIN® FOAMING HAND SANITIZER

Alcohol-free, waterless foaming hand sanitizer. Special formulation contains emollients and a fresh pleasant Mountain Meadow fragrance, leaving hands feeling soft and clean after use, without any stickiness or drying out of the skin.

Rinse-free: no water required. Kills 99.9% of the most common disease-causing germs spread by hands. Contains benzalkonium chloride (active ingredient) to kill germs.



- Alcohol-free formulation that is nondrying to the skin
 - Kills 99.9% of most common disease-causing germs spread by the hands
 - Rinse-free, no water required
 - Fresh clean scent and soothing emollients leave hands soft and clean
- 2-oz Pump Bottle
(900-4445)ea
- 18.6-oz Pump Bottle
(900-4444)ea

HENRY SCHEIN® HAND LOTION

- Luxurious hand cream
- Soothing to skin
- Contains emollients to leave hands feeling soft
- Pleasant Cucumber and Aloe fragrance

8-oz Pump Bottle
(900-4446)ea



HENRY SCHEIN® FOAMING ANTIBACTERIAL HAND SOAP

- Luxurious foaming soap
- Provides a rich, dense lather
- Rinses easily using minimal water
- Leaves hands feeling soft and smooth
- Pleasant Spring Water Clean fragrance
- Active ingredient is 0.2% triclosan

18.6-oz Pump Bottle
(900-4441)ea



HENRY SCHEIN® ANTIBACTERIAL LIQUID HAND SOAP

Luxurious Liquid Hand Soap for Sinks and Dispensers

Provides a luxurious feel, leaving hands smooth and soft. Mild, gentle, and moisturizing. Pleasant Cucumber and Aloe fragrance. Active ingredient is 0.2% triclosan.

8-oz Pump Bottle
(900-4439)ea

16-oz Pump Bottle
(900-4438)ea

1-gal Refill Bottle
(900-4440)ea



HENRY SCHEIN® MAXICLENS® ANTIMICROBIAL/ANTISEPTIC SKIN CLEANSER

Formulated with detergents, moisturizers, skin protectors, and humectants that cleanse without drying. Fast-acting formula that is proven effective at 15 seconds against a broad spectrum of bacteria. Effective for up to 6 hours after use.

- Clear, viscous liquid, 4% chlorhexidine gluconate
- Formulated with detergents, moisturizers, skin protectors, and humectants that cleanse without drying
- Fast acting, proven effective at 15 seconds against a broad spectrum of bacteria
- Persistent; effective for up to 6 hours after use; high level residual action

1 Quart
(431-0032)ea

1 Gallon
(431-0099)ea

Hand Pump for Quart Bottle
(431-1478)ea



HENRY SCHEIN® QUICK WIPES

Nonallergenic, premoistened, ready-to-use hand towelettes. Pleasant fragrance with lanolin, alcohol, and aloe. Refreshing aid for the office, home, or car. Use for easy, gentle cleansing of hands and face.

- Soft, strong towelettes: allow easy and gentle cleansing of hands and face
- Effective and gentle formulation that contains aloe: formulated for fast drying and gentleness
- Convenient packaging: towelettes easily dispensed through orifice; each towelette is saturated with cleansing liquid

Pop-Up Container, 9" x 5½"
(104-4176) 150/cont



HENRY SCHEIN® MAXISPRAY® PLUS

Hospital-level, general-purpose disinfectant/cleaner proven to kill TB in 5 minutes; and HBV, HIV, and HCV in 3 minutes. It is also EPA-registered as a broad-spectrum disinfectant for both surface and immersion use. Bactericidal, virucidal, tuberculocidal, and fungicidal. Ready to use in cleaning and disinfecting equipment surfaces and noncritical areas.

- Easy, ready-to-use formula: no mixing, measuring, or rinsing required
- Compatible with most medical/dental device material
- Safe to use in neonatal and infant-care areas
- Contains no phenols
- Nontoxic and nonirritating
- Satisfies OSHA's Bloodborne Pathogens Standard

24-oz Trigger-Spray Bottle
(101-4114)ea

1-gal Refill
(101-0285)ea



HENRY SCHEIN® MAXIWIPE™ GERMICIDAL CLOTHS

Cleans and disinfects common hard, nonporous surfaces. Bactericidal, tuberculocidal, and virucidal. Effective against 13 microorganisms including TB, HBV, MRSA, VRE, and RSV. Contains 55% alcohol. Meets CDC and OSHA regulations.

Large
(106-6794) 160/cont

X-Large
(112-5305) 65/cont





HENRY SCHEIN® MAXIZYME®

Dual Enzymatic Detergent
Nonfoaming, fast-working, dual-enzymatic detergent for use as a presoak, manual detergent, ultrasonic detergent, and evacuator cleaner. Dissolves blood, tissue, mucus, and other protein-rich body fluids from instruments. Noncorrosive, nonflammable, and rust inhibiting. Concentrated and economical.

- Contains 2 gold-standard proteolytic enzymes for more effective cleaning
- Safe for instruments
- Less corrosive than distilled water
- Dissolves and releases blood and tissue: minimizes handling and scrubbing of contaminated items
- Nonfoaming: compatible with evacuators, ultrasonic units, and presoaks



(101-9031) 1-Gallonea

HENRY SCHEIN® MAXICIDE®

Can be used as a sterilant and high-level disinfectant for up to 28 days. 90 minutes at 77°F (25°C) for high-level disinfection; 10 hours at 77°F (25°C) for sterilization. Contains a surfactant to enhance the disinfection process and a rust inhibitor to reduce corrosion. No dilution required.

- Sterilization in 10 hours
- 28-day reuse life after activation
- Clear color before activation; turns blue after activation
- Activated solution contains anticorrosion and antioxidant ingredients
- Rapid and effective broad-spectrum kill
- Pleasant pine scent



1-qt Bottle
(102-6003)ea
1-gal Bottle
(102-7394)ea

HENRY SCHEIN® SELF-SEAL STERILIZATION POUCH WITH DUAL INDICATOR

Superior quality pouch with dual indicators that comply with the recommendations of the CDC for effective sterilization. Preprinted internal and external indicators change color when processed with steam and EtO. Made with medical-grade paper and an easy-to-open thumb notch. Lead-free and latex-free.

- (112-4856) 2 3/4" x 9"200 per Box
- (112-4853) 3 1/2" x 9"200 per Box
- (112-4854) 5 1/4" x 10"200 per Box
- (112-4858) 7 1/2" x 13"200 per Box
- (112-4862) 10 1/2" x 17"100 per Box
- (112-4865) 13" x 20"100 per Box



HENRY SCHEIN® MAXICIDE® ACTIVATED

Can be used as a sterilant and high-level disinfectant for up to 30 days. Tuberculocidal in 45 minutes. Ready to use; has a fresh, pleasant scent.

1-gal Bottle
(104-4681)ea



HENRY SCHEIN® INDICATOR STRIPS

For steam and chemical sterilizers. Indicator ink changes from white to black, indicating sterilization process has been completed.

(102-6332)250/box



HENRY SCHEIN®

MAXIMA® MULTIPARAMETER SELF-SEAL STERILIZATION POUCHES WITH SURE-CHECK® TECHNOLOGY

Confirms that all 3 parameters for sterilization have been met: time, temperature, and presence of steam. Internal and external multiparameter indicators are printed directly on the pouch. Eliminates the extra step and cost of adding separate multiparameter strips to each pouch. Medical-grade paper with triple seals; chevron safety seal at base for added protection; thumb notch for easy opening. Latex-free. Made in the USA.

- 2 3/4" x 9"200/box
(900-4602)
- 3 1/2" x 5 1/4"200/box
(900-4604)
- 3 1/2" x 9"200/box
(900-4601)
- 5 1/4" x 10"200/box
(900-4603)
- 7 1/2" x 13"200/box
(900-4782)
- 12" x 18"100/box
(900-4783)



HENRY SCHEIN® CSR WRAP

For general-purpose sterilization, providing a barrier against airborne and waterborne bacteria. Conforms easily around various sizes and shapes. Provides a barrier, yet allows excellent penetration for various methods of sterilization, including steam, EtO, radiation, and E-beam. Autoclavable.

- (100-3347) 12" x 12"100/case
- (102-5193) 15" x 15"100/case
- (100-4683) 18" x 18"100/case
- (100-6950) 24" x 24"100/case
- (100-8187) 30" x 30"100/case
- (100-2016) 36" x 36"100/case



HENRY SCHEIN® MAXITEST®

Mail-in service includes 3 bacterial spore test strips (2 test strips and 1 control strip) per test package. Postage is prepaid and culturing service is included. Monitors all sterilizers including steam, chemical vapor, dry heat, and EtO. Third-party documentation serves as tangible, legal document where proof of sterility testing must be provided. Immediate and confidential notification of test failures. Test results available electronically.

- Easy to use: place MaxiTest strips within an average load; process according to normal procedure; when cycle is complete, remove load from sterilizer and remove test strips; forward to bacteriological laboratory for culturing
- Suitable for use with all sterilizers: dual-organism spore strips contain both *Bacillus stearothermophilus* and *Bacillus subtilis* so they can be used to monitor all sterilizers, including unsaturated chemical vapor
- 2 test strips and 1 control strip: provide increased safety to patients and staff; 2 test strips placed in separate areas of the sterilizer provide greater accuracy of test results
- Biological Sterilization Monitoring Service: MaxiTest includes test sets, culturing service, laboratory reports, all postage and handling, and any device concerning sterilization procedures; the lab cultures are completed in 7 days; however, if any test strip cultures positive (spore growth indicating sterilization failure), the office will be telephoned immediately; if the test strips do not culture positive, the office will receive a laboratory report in approximately 2 weeks, indicating sterilization was achieved

12-Test Package
(101-5830)ea

Contains: 12 packets, ea containing 2 test strips & 1 control strip.

48-Test Package
(101-1253)ea

Contains: 48 packets, ea containing 2 test strips & 1 control strip.

Not for export.





HENRY SCHEIN® CRITERION® PC POWDER-FREE LATEX EXAM GLOVES



Non-sterile, Polymer Coated, White
Textured • Nonchlorinated

- Certified Protein Claim: None
- Sets the standard in quality, powder-free gloves. Nonchlorinated, synthetic polymer coating eases donning. Textured finish for enhanced grip. Excellent tactile sensitivity.
- Textured finish for enhanced grip
- Excellent tactile sensitivity
- Powder-free gloves eliminate potential aerosolization of glove powder and powder-associated transfer of chemical allergens
- Multilayered and fusion bonded for enhanced strength and durability
- Microroughened texturing provides reliable wet or dry gripping while maintaining tactile sensitivity

Gloves100/box

Specify:

| | |
|------------------------|------------|
| X-small | (102-5418) |
| Small | (102-5419) |
| Medium | (102-5421) |
| Large | (102-6730) |
| X-large (90/box) | (102-5422) |

HENRY SCHEIN® CRITERION® POWDER-FREE NITRILE EXAM GLOVES



Non-sterile, Blue
Fully textured nitrile glove that offers excellent tactile sensitivity. Soft and stretchy, this nonchlorinated glove is the ideal alternative to latex.

Gloves100/box

Specify:

| | |
|------------------------|------------|
| X-small | (111-8533) |
| Small | (111-8535) |
| Medium | (111-8536) |
| Large | (111-8537) |
| X-large (90/box) | (111-8539) |



HENRY SCHEIN® CRITERION® POWDER-FREE VINYL EXAM GLOVES

Non-sterile
Contains no natural latex, providing an excellent solution for those at risk of latex sensitivities.

- Supple touch for maximum tactile sensitivity
 - Powder-free to eliminate glove-powder aerosolization
- Smooth
Gloves100/box

Specify:

| | |
|------------------------|------------|
| X-small | (102-5721) |
| Small | (102-6721) |
| Medium | (102-5338) |
| Large | (102-7319) |
| X-large (90/box) | (102-5377) |



HENRY SCHEIN® CRITERION® CL POWDER-FREE LATEX EXAM GLOVES

Non-sterile
Textured • Chlorinated

- Certified Protein Claim: None
- An on-line chlorinated glove, textured for superior grip, comfortable fit, and low protein levels.
- Produced using state-of-the-art on-line chlorination technology for better consistency
- Textured for superior grip

Gloves100/box

Specify:

| | |
|---------------|------------|
| X-small | (102-5354) |
| Small | (102-8296) |
| Medium | (102-5355) |
| Large | (102-7763) |
| X-large | (102-5357) |



GLOVE CLUB® *Tomorrow's Gloves. Today!*

CYBER-GEL® Powder-Free Latex Exam Gloves

Non-sterile
Textured • Chlorinated • Certified Protein
Claim: None

- Gel coated with enhanced donning effect
- Textured for optimal grip
- Excellent tactile sensitivity
- Soft latex provides excellent fit to ensure ultimate comfort
- Chlorinated

Gloves100/box

Specify:

| | |
|---------------|------------|
| X-small | (103-5123) |
| Small | (103-5133) |
| Medium | (103-5165) |
| Large | (103-5196) |
| X-large | (103-5121) |



HENRY SCHEIN® CRITERION® LP LIGHTLY POWDERED LATEX EXAM GLOVES

Non-sterile, Micro-Textured
Textured
Choose Criterion Latex for superior value. These low-powder gloves offer the grip and durability you need at a low price.

- Lightly powdered to assist in donning
- Bisque finish with excellent tactile sensitivity
- Maximum tensile strength to prevent tearing

Gloves100/box

Specify:

| | |
|------------------------|------------|
| X-small | (102-5425) |
| Small | (102-7597) |
| Medium | (102-5426) |
| Large | (102-5427) |
| X-large (90/box) | (102-5428) |



HENRY SCHEIN® CRITERION® POWDER-FREE LATEX SURGEONS' GLOVES

Compare to Micro-Touch® PF
Latex Surgeons' Gloves
Sterile

Textured • Chlorinated • Certified Protein
Claim: None
Henry Schein's Criterion brand of latex surgical gloves helps eliminate powder-related complications and irritation, yet still provides outstanding quality and performance you would expect from a premium surgical glove.

Gloves50 pair/box

Specify:

| | |
|----------|------------|
| 5½ | (104-6125) |
| 6 | (104-5073) |
| 6½ | (104-9218) |
| 7 | (104-3502) |
| 7½ | (104-3503) |
| 8 | (104-3563) |
| 8½ | (104-7315) |
| 9 | (104-7227) |



HENRY SCHEIN® CRITERION® LIGHTLY POWDERED LATEX SURGEONS' GLOVES

Compare to Triflex® Sterile Powdered Surgeons'
Latex Gloves
Sterile

Henry Schein's Criterion brand of low-powder latex surgical gloves provides optimal balance between barrier protection and tactile sensitivity.

Smooth • Chlorinated • Certified Protein
Claim: None
Gloves50 pair/box

Specify:

| | |
|----------|------------|
| 5½ | (104-7721) |
| 6 | (104-3500) |
| 6½ | (104-9242) |
| 7 | (104-3501) |
| 7½ | (104-5920) |
| 8 | (104-4696) |
| 8½ | (104-5787) |
| 9 | (104-9192) |



**HENRY SCHEIN®
ULTRASOFT PROCEDURAL MASKS**
Latex-Free

Fluid-resistant spunbond outer layer/inner white layer made from extremely soft hypoallergenic cellulose. Will not lint, tear, or shred. Meets ASTM standards as a low (primary) fluid-resistant barrier @ 80mm Hg; PFE 99.75% @ 1 micron; BFE 99.3% @ 3.1 microns; Delta P-3.4. Made in the USA.
 UltraSoft.....50/box
 Specify:
 Blue(104-8072)
 Pink(104-8073)



**HENRY SCHEIN®
EAR LOOP PROCEDURAL FACE MASKS**

Made from spunbond outer/inner layers with filter media and flat ear loops for greater comfort. Meets ASTM standards as a low-barrier (fluid resistant) face mask. Made in the USA. Latex-free.
 Masks.....50/box
 Specify:
 Blue(104-7799)
 Lavender.....(104-8071)
 Pink(104-9200)
 White with White Embossed Hearts
 (112-5052).....50/box



**HENRY SCHEIN®
PLEATED MASK**

Triple-layer, fluid-resistant, and high filtration rate. Soft, lightweight, comfortable design. Manufactured to exceed all current ASTM F2100-03a standards. Cloth-like ties.
 (104-6809).....50/box



**HENRY SCHEIN®
EAR LOOP PROCEDURAL FACE MASKS**

Light, comfortable masks for extended wear when necessary. Ear loops are made from special stretch yarn: no irritating rubber or plastic touches the skin. Three-ply, glass-free filter media is surrounded by soft, nonirritating inner and outer facing. Full-width flexible nosepiece provides fit and security. Manufactured to exceed all current ASTM F2100-07 standards. Latex-free.
 • Light, comfortable masks for extended wear when necessary
 • Ear loops are made from special stretch yarn: no irritating rubber or plastic touches the skin
 • Three-ply, glass-free filter media
 • Soft, nonirritating inner and outer facing
 • Full-width flexible nosepiece provides fit and security
 • Manufactured to exceed all current ASTM F2100-04 standards
 • Latex-free
 Masks.....50/box
 Specify:
 Blue(104-3809)
 White(104-6611)
 Lavender.....(104-8600)
 Pink(104-3730)
 Yellow(104-2849)



**HENRY SCHEIN®
ECONOMY MOLDED FACE MASK**

Blue, preformed with aluminum nose guard.
 #MK-6829, Blue
 (101-6829).....50/box



**HENRY SCHEIN®
DISPOSABLE FACE SHIELDS**

Optically clear, see-through visor provides eye protection from splash and splatter. "Soft-strap" has an impenetrable foam barrier providing comfortable, wide coverage without compromising safety. Each shield is packaged individually in a plastic bag intended for single-patient use to prevent cross contamination.
 Standard, 13"W x 7½"L
 (101-2254).....24/box
 Bulk Packaging
 (112-5656).....72/box



**HENRY SCHEIN®
SAFETY GLASSES**

Clear plastic lenses with side shields for eye protection in the operating room or lab. Unisex.
 (100-7382).....ea



**HENRY SCHEIN®
EAR LOOP FACE MASK WITH EYE SHIELD**

Soft, comfortable loops and absorbable inner facing make this mask very comfortable for extended wear. Manufactured to exceed all current ASTM F2100-03a standards. Superior fluid protection; passed the Standard Fluid Penetration Resistance Test. The optically clear shield is treated on both sides with antifog material to ensure that vision is not compromised. Latex-free.
 Teal
 (104-7321).....25/box



HENRY SCHEIN® PREMIUM TABLE PAPER ROLLS

Designed and manufactured to provide strength, protection, and value.



- Crepe, 18" x 125' (100-9865)12/case
- Crepe, 21" x 125' (100-1657)12/case
- Smooth, 18" x 225' (100-9491)12/case
- Smooth, 21" x 225' (100-0177)12/case
- Smooth, Heavy-duty, 14" x 225' (104-4501)12/case
- Smooth, Heavy-duty, 18" x 225' (104-6092)12/case
- Smooth, Heavy-duty, 21" x 225' (104-5273)12/case

HENRY SCHEIN® ESSENTIALS TABLE PAPER

Quiet, absorbent barrier. Designed and manufactured to provide strength, protection, and value.

- White, Smooth, 225' 18" (900-4978)12 rolls/case
- 21" (900-4979)12 rolls/case
- White Crepe, 125' 18" (979-9865)12 rolls/case
- 21" (979-1657)12 rolls/case



HENRY SCHEIN® EXAM GOWNS

3-Ply Tissue, 30" x 42"

Soft tissue provides patient comfort, while the embossed pattern provides extra strength. Plastic ties attached.

- White (100-4471)50/case
- Blue (100-4975)50/case
- Mauve (100-7115)50/case



HENRY SCHEIN® EXAM CAPES

3-Ply Tissue, 30" x 21"

- White (100-8821)100/case
- Blue (100-2393)100/case
- Mauve (100-2010)100/case



HENRY SCHEIN® COTTON BALLS

- Uniform size and shape: hold together even when saturated
- Highest quality, 100% pure-white cotton: soft and absorbent
- Tamperproof package: safe
- Machine manufactured for exact control of weight and density: consistent shape

- Sterile, Medium (100-4752)500/box
- Non-sterile, Medium (100-0468)4000/case
- Non-sterile, Large (100-1491)2000/case

HENRY SCHEIN® FACIAL TISSUES

Compare to Kleenex®
2 Ply, White, 8" x 8.33"

Soft, absorbent, elegantly embossed white facial tissue. Your patients and employees will notice the softness and quality.
(900-4214)100/box



HENRY SCHEIN® PREMIUM EXAM DRAPES

2-Ply, All-Virgin Tissue

Made from 2-ply tissue to provide soft and strong coverage for most exams. They are designed to blanket the patient in light-weight comfort and are available in an array of colors to enhance office decor.



- White, 36" x 48" (101-3658)100/case
- White, 40" x 48" (100-1797)100/case
- Blue, 40" x 48" (100-3360)100/case
- Mauve, 40" x 48" (100-6428)100/case
- Peach, 40" x 48" (100-3845)100/case
- White, 40" x 60" (101-8521)100/case



HENRY SCHEIN® PILLOWCASE

Tissue/Poly/Tissue

Provides linen-like softness for patient comfort and an inner poly barrier for extra protection. Strong and durable.

- White, 21" x 30" (104-9206)100/case

HENRY SCHEIN® PREMIUM GOWNS

Tissue/Poly/Tissue, 30" x 42"

Layer of poly provides strength; layers of tissue absorb moisture. Opens front or back with attached elastic ties.

- Gowns50/case
- Specify:
- White (100-2392)
- Blue (100-5728)
- Mauve (100-9652)
- Peach (100-6999)



HENRY SCHEIN® DELUXE PREMIUM EXAM DRAPE

3-Ply, All-Virgin Tissue

Soft, opaque, and heavyweight, this embossed tissue drape with dependable wet strength is ideal for short to moderate exams. It conforms to the body and offers the ultimate in absorbency.

- White, 40" x 48" (101-7084)100/case





**HENRY SCHEIN®
ONESTEP+ PREGNANCY
STRIP TEST**



- Urine test
 - CLIA waived
 - 25 mIU/mL detection level
 - Room-temperature storage
 - Read time 3 minutes or less
- (900-4072)**25/box

Contains: 25 pouched tests per kit & 1 package insert.

CPT® Code: Urine 81025QW†

**HENRY SCHEIN®
ONESTEP+ PREGNANCY
CASSETTE TESTS**



- Urine test
 - CLIA waived
 - 25 mIU/mL detection level
 - Room-temperature storage
- 25-Patient Kit
(900-4073)ea

Contains: 25 pouched tests per kit & 1 package insert.

- 100-Patient Kit
(900-4074)ea

Contains: 100 pouched tests per kit & 1 package insert.

CPT® Code: Urine 81025QW†



**HENRY SCHEIN®
ALCOHOL PREP PADS**

Medium, 1 1/2" x 2 3/4"
Pad is saturated with 70% isopropyl alcohol and used for topical cleansing prior to injection or venipuncture.

- Non-sterile
(104-6352)200/box
- Sterile
(104-8298)200/box

**HENRY SCHEIN®
DISPOSABLE NEEDLES**

- 30 ga x 1" **(900-4474)**100/box
- 30 ga x 3/4" **(900-4467)**100/box
- 25 ga x 3/8" **(900-4468)**100/box
- 25 ga x 1" **(900-4470)**100/box
- 25 ga x 1 1/2" **(900-4472)**100/box
- 23 ga x 1" **(900-4473)**100/box
- 18 ga x 1" **(900-4469)**100/box
- 18 ga x 1 1/2" **(900-4471)**100/box

**HENRY SCHEIN®
1-CC TB SYRINGE WITH NEEDLE**

- 25 ga x 3/8" **(900-4459)**100/box

**HENRY SCHEIN®
LUER SLIP 1-CC DISPOSABLE SYRINGE
WITH NEEDLE**

- 27 ga x 1/2" **(900-4463)**100/box

**HENRY SCHEIN®
LUER LOCK 3-CC DISPOSABLE
SYRINGES WITH NEEDLES**

- Syringes100/box
- Specify:
- 25 ga x 3/8"**(900-4458)**
- 25 ga x 1"**(900-4461)**
- 23 ga x 1"**(900-4460)**
- 21 ga x 1"**(900-4462)**



**HENRY SCHEIN®
SHARPS COLLECTORS**

Antidefluxion wing sizes prevent outflow of waste. Maximum filling level indicator. No emission of any harmful gases during the incineration process. Simple and practical snap-type closure. Easy maneuvering for temporary opening with just one hand. Wide opening for the disposal of large amounts of sharps waste. Mechanism for the removal of needles on some units.

- Redea

- Specify:
- 1 qt**(900-4980)**
- 1.5 qt**(900-4981)**
- 4 qt**(900-4982)**
- 5 qt**(900-4983)**
- 2 gal**(900-4984)**
- 3 gal**(900-4985)**
- 8 gal**(900-4986)**

- Red, 18 gal
(900-7032)7/box
- Chemotherapy–Yellowea

- Specify:
- 3 gal**(900-4988)**
- 8 gal**(900-4989)**

- Chemotherapy–Yellow, 18 gal
(900-7033)7/box

**HENRY SCHEIN®
SYRINGES WITHOUT NEEDLES**

- 1 cc **(900-4480)**100/box
- 3 cc **(900-4475)**100/box
- 5 cc **(900-4477)**100/box
- 10 cc **(900-4476)**100/box
- 20 cc **(900-4478)**50/box
- 30 cc **(900-4481)**50/box
- 60 cc **(900-4479)**50/box



**HENRY SCHEIN®
INFECTIOUS WASTE COLLECTION BAGS**

Constructed of high-grade polyethylene for safe collection, handling, and disposal of infectious wastes destined for incineration or landfill. Red and labeled with a biohazard symbol. Twist ties provided.

- Small, 10 gal, 15" x 9" x 24"
(100-9305)100/box
- Medium, 18 gal, 15" x 9" x 31"
(100-8878)100/box
- Large, 33 gal, 23" x 8" x 41"
(100-3416)100/box



Superior Point: Needles are honed with a unique multiedged point to provide greater ease of penetration. Superior attachment: Needles are drilled for a seamless needle-suture juncture, providing a smoother transition with less tissue disruption. Superior packaging: Unique double-sterile barrier packaging ends costly waste.

- Precision-honed needle: offers surgeon easier penetration and precise tracking
- Finely crafted swaged needles: tissue trauma is reduced because swaged needles are used only once and then discarded; dull needles are eliminated
- Sterile packaging: ensures sterility
- Unique, see-through packaging: permits positive identification without the necessity to touch or otherwise contaminate suture
- Wide assortment of suture material (Silk–Black Braided, Gut–Plain, Gut–Chromic, Polyester–Green Braided): practitioner can select from several suture materials



HENRY SCHEIN® **CHROMIC GUT ABSORBABLE SUTURES**

- Reverse Cutting, $\frac{1}{2}$ Circle
Needle Code: C-6 (FS-2*), Length: 27",
Size: 4-0, Mfr. #635-CG
(101-2787) 12/box
- Needle Code: C-6 (FS-2*), Length: 27",
Size: 3-0, Mfr. #636-CG
(101-1582) 12/box
- Needle Code: C-3 (P-3*), Length: 18",
Size: 5-0, Mfr. #687-CG
(101-8824) 12/box
- * Ethicon needle equivalent

HENRY SCHEIN® **PLAIN GUT ABSORBABLE SUTURES**

- Reverse Cutting, $\frac{1}{2}$ Circle; Needle Code:
C-6 (FS-2*), Length: 18", Size: 4-0, Mfr.
#7771-PG
(100-7771) 12/box
- Reverse Cutting, $\frac{1}{2}$ Circle; Needle Code:
C-6 (FS-2*), Length: 27", Size: 4-0, Mfr.
#H821-PG
(101-0668) 12/box
- Reverse Cutting, $\frac{1}{2}$ Circle; Needle Code:
C-6 (FS-2*), Length: 27", Size: 3-0, Mfr.
#H822-PG
(101-2976) 12/box
- Reverse Cutting, $\frac{1}{2}$ Circle; Needle Code:
C9, (X-1*), Length: 18", Size: 3-0, Mfr.
#612-PG
(100-6120) 12/box
- Precision Reverse Cutting, $\frac{1}{2}$ Circle;
Needle Code: C-3 (P-3*), Length: 18",
Size: 5-0, Mfr. #686-PG
(101-3192) 12/box
- *Ethicon needle equivalent

HENRY SCHEIN® **NYLON BLACK MONOFILAMENT NONABSORBABLE SURGICAL SUTURES**

- 18" Length
Reverse Cutting,
 $\frac{1}{2}$ Circle 12/box
- Specify:
Needle Code: C22 (FS-3*),
Size: 6-0, Mfr. #660-NY **(101-6053)**
Needle Code: C6 (FS-2*),
Size: 5-0, Mfr. #661-NY **(101-7137)**
Needle Code: C6 (FS-2*),
Size: 4-0, Mfr. #662-NY **(101-1209)**
Needle Code: C7 (FS-1*),
Size: 3-0, Mfr. #663-NY **(101-3728)**
- *Ethicon needle equivalent

HENRY SCHEIN® **POLYPROPYLENE BLUE MONOFILAMENT NONABSORBABLE SURGICAL SUTURES**

- 18" Length
Reverse Cutting,
 $\frac{1}{2}$ Circle 12/box
- Specify:
Needle Code: C22 (FS-3*),
Size: 6-0, Mfr. #8660-PP **(101-4378)**
Needle Code: C6 (FS-2*),
Size: 5-0, Mfr. #8661-PP **(101-6409)**
Needle Code: C6 (FS-2*),
Size: 4-0, Mfr. #8683-PP **(101-7477)**
Needle Code: C7 (FS-1*),
Size: 3-0, Mfr. #8684-PP **(101-7700)**
- *Ethicon needle equivalent

HENRY SCHEIN® **SILK BLACK BRAIDED NONABSORBABLE SUTURES**

- 18" Length
Reverse Cutting, $\frac{1}{2}$ Circle; Needle Code:
C-22 (FS-3*), Size: 6-0, Mfr. #681-SK
(101-2636) 12/box
- Reverse Cutting, $\frac{1}{2}$ Circle; Needle Code:
C-6 (FS-2*), Size: 5-0, Mfr. #682-SK
(101-3207) 12/box
- Reverse Cutting, $\frac{1}{2}$ Circle; Needle Code:
C-6 (FS-2*), Size: 4-0, Mfr. #683-SK
(100-6830) 12/box
- Reverse Cutting, $\frac{1}{2}$ Circle; Needle Code:
C-7 (FS-1*), Size: 3-0, Mfr. #684-SK
(101-2049) 12/box

HENRY SCHEIN® **PGA (POLYGLYCOLIC ACID) ABSORBABLE SUTURES**

Compare to Vicryl
3/8 Circle, 18"
Synthetic absorbable braided sutures have a smooth, low-friction coating that reduces tissue drag. They retain approximately 50% of their initial strength at 2 weeks and 20% at 3 weeks. PGA is ideal for skin and mucosa closings and can be used in most areas of the body.

- Reverse Cutting, $\frac{1}{2}$ Circle; Needle Code:
C-6 (FS-2*), Size: 3-0, Mfr. #7772-SK
(100-7772) 12/box
- Reverse Cutting, $\frac{1}{2}$ Circle; Needle Code:
C-9 (X-1*), Size: 3-0, Mfr. #632-SK
(100-6320) 12/box
- Reverse Cutting, $\frac{1}{2}$ Circle; Needle Code:
C-43, Size: 3-0, Mfr. #7762-SK
(100-7762) 12/box
- Precision Reverse Cutting, $\frac{1}{2}$ Circle; Needle
Code: C-3 (P-3*), Size: 5-0, Mfr. #640-SK
(101-4148) 12/box
- Precision Reverse Cutting, $\frac{1}{2}$ Circle; Needle
Code: C-3 (P-3*), Size: 4-0, Mfr. #641-SK
(101-2590) 12/box
- *Ethicon needle equivalent

- Precision Reverse Cutting, Undyed Braided
C3, 5/0
(900-4417) 12/box
- C3, 4/0
(900-4418) 12/box
- PC31, 4/0
(900-4419) 12/box
- PC31, 3/0
(900-4421) 12/box
- Reverse Cutting, Violet Braided
C6, 4/0
(900-4420) 12/box



HENRY SCHEIN®
BACKHAUS TOWEL CLAMP 3½"
(101-4934).....ea



HENRY SCHEIN®
LORNA NONPERFORATING TOWEL CLAMP, 3½"
(104-8244).....ea



HENRY SCHEIN®
SPENCER STITCH SCISSORS 3½"
Economy (Value)
(100-9827).....ea



HENRY SCHEIN®
LISTER BANDAGE SCISSORS 5½"
Extra Fine
(100-1087).....ea



HENRY SCHEIN®
IRIS SCISSORS 4½"
German Stainless Steel, Straight
(100-8118).....ea



HENRY SCHEIN®
DISPOSABLE BIOPSY PUNCHES
Our disposable skin biopsy punches have beneficial properties not found with generic products. Features include:
• Stainless steel ground blade
• Seamless punch
• Ergonomic blue handle
Available in 9 different sizes from 1.5 mm to 8.0 mm. Sterile, individually pouched with a tamper-evident seal.
Biopsy Punchesea
Specify:
1.5 mm(900-4863)
2 mm(900-4864)
2.5 mm(900-4865)
3 mm(900-4866)
3.5 mm(900-4867)
4 mm(900-4868)
5 mm(900-4869)
6 mm(900-4870)
8 mm(900-4871)



HENRY SCHEIN®
ADSON DRESSING FORCEPS 4¾"
German, Serrated, Delicate Pattern
(101-6279).....ea



HENRY SCHEIN®
SPENCER STITCH SCISSORS 3½"
German Stainless Steel
(100-6024).....ea



HENRY SCHEIN®
IRIS SCISSORS 4½"
Economy (Value)
Straight
(100-2767).....ea
Curved
(100-8820).....ea




HENRY SCHEIN®
IRIS SCISSORS 4"
German, Curved
(101-8111).....ea



HENRY SCHEIN®
GRADLE EYE SUTURE SCISSORS 3¾"
German, Slightly Curved, Sharp Points
(101-7676).....ea



HENRY SCHEIN®
IRIS SCISSORS 4½"
German Stainless Steel, Curved
(100-5880).....ea



HENRY SCHEIN®
SURGICAL HANDLES
Economy, #3
(100-7520).....ea
#3 is for blades #10, #11, #12B, #15 & #15C.



HENRY SCHEIN®
SURGICAL BLADES
Individual foil packs.
Computer-controlled grinding ensures a consistently sharp edge.
Stainless Steel100/box
Specify:
#10(100-8976)
#11(100-5794)
#15(100-0249)
Carbon Steel #15
(100-1259).....10/box
#15C
(100-4013)10/Pkg



HENRY SCHEIN®
DISPOSABLE SCALPELS
Individually peel-backed for convenience. Plastic handle with millimeter scale and durable stainless steel blade.
• European raw material: high-quality steel
• Sharp: clean, true, consistent cut
• Strong: reduced breakage
• Packaging: superior protection and handling
• Testing: consistent high quality; assured quality to international standards
Scalpels.....10/box
Specify:
#10(100-8037)
#11(100-5622)
#15(100-5255)



HENRY SCHEIN® ALL-GAUZE SPONGES

100% Woven Cotton, Non-sterile
Soft, pliable, and absorbent. Available in 3 sizes in both 8 ply and 12 ply.

- Precision sizes and folding: consistency within functions
- Soft and pliable: nonabrasive touch
- Highly absorbent: soaks up unwanted fluid
- Can be used for general cleaning, debriding, and prepping

| | |
|----------------------------|---------|
| 8 Ply, 2" x 2" (100-2524) | 200/pkg |
| 8 Ply, 3" x 3" (100-6458) | 200/pkg |
| 8 Ply, 4" x 4" (100-5952) | 200/pkg |
| 12 Ply, 2" x 2" (100-1570) | 200/pkg |
| 12 Ply, 3" x 3" (100-5004) | 200/pkg |
| 12 Ply, 4" x 4" (100-1444) | 200/pkg |



HENRY SCHEIN® PREMIUM ALL-GAUZE SPONGES

Non-sterile

| | |
|---|---------|
| 8 Ply, 2" x 2" (104-8197) | 200/pkg |
| 8 Ply, 4" x 4" (104-8194) | 200/pkg |
| 12 Ply, 2" x 2" (104-8588) | 200/pkg |
| 12 Ply, 3" x 3" (104-8195) | 200/pkg |
| 12 Ply, 4" x 4" (104-8198) | 200/pkg |
| 8 Ply, 2" x 2" (101-4947) | 100/pkg |
| Contains: 2 sponges per envelope, 50 envelopes per package. | |
| 8 Ply, 4" x 4" (101-5519) | 50/pkg |
| Contains: 2 sponges per envelope, 25 envelopes per package. | |
| 12 Ply, 3" x 3" (101-1768) | 80/pkg |
| Contains: 2 sponges per envelope, 40 envelopes per package. | |
| 12 Ply, 4" x 4" (101-0005) | 50/pkg |
| Contains: 2 sponges per envelope, 25 envelopes per package. | |



HENRY SCHEIN® CRITERION™ NON-WOVEN SPONGES

Compare to NuGauze®
Non-sterile, 4 Ply
• High quality
• Excellent construction
• Rayon/poly-blend is soft and gentle to the touch
• Thicker and more absorbent than standard sponges
• Multi-use sponge can be used for all applications of wound-care treatments
• Easy-to-open, peel-back packaging

| | |
|--------------------|---------|
| 2" x 2" (104-2930) | 200/pkg |
| 3" x 3" (104-2410) | 200/pkg |
| 4" x 4" (104-2411) | 200/pkg |



HENRY SCHEIN® GAUZE PADS

100% Woven Cotton, Sterile 1's in Peel-Back Packages, 12 Ply

- Highly absorbent
- 100% woven cotton
- Can be used for general cleaning, debriding, and prepping
- Box is perforated at bottom for easier gauze removal

| | |
|--------------------|---------|
| 2" x 2" (100-4380) | 100/pkg |
| 3" x 3" (100-3963) | 100/pkg |
| 4" x 4" (100-2634) | 100/pkg |



HENRY SCHEIN® CRITERION® NON-WOVEN SPONGES

Sterile 2's, 4 Ply
High-quality, excellently constructed rayon/poly-blend sponge that is soft and gentle to the touch and very comfortable against skin. Thicker than the standard sponge and, with its open-weave construction, offers more absorbency than other sponges available. Can be used for all applications in wound-care treatment, making it an excellent multiuse sponge. Easy-to-open, peel-back packaging.

- High quality
 - Excellently constructed
 - Rayon/poly-blend that is soft and gentle to the touch
 - Thicker and more absorbent than standard sponges
 - Multiuse sponge that can be used for all applications of wound-care treatments
 - Easy-to-open, peel-back packaging
- | | |
|--------------------|--------|
| 2" x 2" (104-7897) | 50/pkg |
| 3" x 3" (104-4132) | 50/pkg |
| 4" x 4" (104-8665) | 50/pkg |



HENRY SCHEIN® CONFORMING STRETCH GAUZE ROLLS

- 4.1-1-1/2 Rolls
- Economical knitted bandage in various widths
 - High-quality materials consist of nylon, polyester, and cotton; provides light compression as it cushions and protects
 - Offers softness, breathability, and alleviation of irritation
 - Low lint content reduces contamination to wound site
 - Improved packaging is resealable, durable, easy to store, and remains clean



HENRY SCHEIN® NON-WOVEN COTTON SPONGES

Non-sterile, 4 Ply
Highly absorbent 100% cotton. Packaged in a see-through bag with a convenient sterilization indicator.

- 4-ply gauze sponges: thick, soft gauze; highly absorbent
 - Water bleach: no discoloration after sterilization
 - See-through hygienic package imprinted with a sterilization indicator: easy identification
 - Tear notch: convenient bag opening
- | | |
|--------------------|---------|
| 2" x 2" (100-1090) | 100/pkg |
| 3" x 3" (100-3835) | 100/pkg |
| 4" x 4" (100-4420) | 100/pkg |



HENRY SCHEIN® BULK ABSORBENT COTTON ROLLS, USP

1-lb Bulk Rolls
Sterile (100-0145)ea
Non-Sterile (100-0147)ea



| | |
|----------------------------|--------|
| Non-Sterile, 1" (104-2704) | 12/bag |
| Non-Sterile, 2" (104-8825) | 12/bag |
| Non-Sterile, 3" (104-8997) | 12/bag |
| Non-Sterile, 4" (104-2728) | 12/bag |
| Non-Sterile, 6" (104-9424) | 6/bag |
| Sterile, 2" (104-2739) | 12/bag |
| Sterile, 3" (104-9724) | 12/bag |
| Sterile, 4" (104-6708) | 12/bag |
| Sterile, 6" (104-4293) | 6/bag |



HENRY SCHEIN®
WOUND CLOSURE STRIPS

Latex-free
Adhesive wound closures with an exceptionally porous, spun-bond fiber backing include bilateral nylon filaments for high tensile strength. Hypoallergenic, medical-grade adhesive affords excellent adhesion, yet remains gentle to the skin. To be used in addition to or in place of sutures. Each size is in a color-coded package for quick and easy size identification. Latex-free.

- **Highly aggressive adhesive properties:** Resistant to moisture, with proven superior longevity, on the patient, to popular brands
- **Bilateral nylon filaments:** Provide resistance to stretching while maintaining closure and high tensile strength
- **Hypoallergenic:** Designed for patient comfort and longer-term exposure
- **Porous:** Allows breathing for advanced healing and moisture control

5/8" x 3"
(900-4317)250/box

Contains: 5 strips per envelope, 50 envelopes per box.

1/4" x 3"
(900-4315)150/box

Contains: 3 strips per envelope, 50 envelopes per box.

1/4" x 4"
(900-4318)500/box

Contains: 10 strips per envelope, 50 envelopes per box.

1/4" x 1 1/2"
(900-4316)300/box

Contains: 6 strips per envelope.

1/2" x 4"
(900-4314)300/box

Contains: 6 strips per envelope, 50 envelopes per box.

1" x 5"
(900-4319)100/box

Contains: 4 strips per envelope, 25 envelopes per box.



HENRY SCHEIN®
CLOTH ADHESIVE TAPE

10-yd Rolls
Strong, high-quality cloth tape that is ideal for support strapping and patient positioning. Available in white or tan in a variety of sizes. Manufactured with zinc-oxide adhesive to help prevent irritation.

Tanea

Specify:
1", 12 rolls/box**(100-7220)**
1 1/2", 8 rolls/box**(100-2028)**
2", 6 rolls/box**(100-7004)**
3", 4 rolls/box**(100-7874)**

White, Regularea

Specify:
1/2", 24 rolls/box**(100-5334)**
1", 12 rolls/box**(100-5387)**
1 1/2", 8 rolls/box**(100-4058)**
2", 6 rolls/box**(100-5929)**
3", 4 rolls/box**(100-6296)**



HENRY SCHEIN®
PREMIUM SURGICAL HYPALLERGENIC PAPER TAPE

Latex-Free, 10-yd Rolls
• Excellent adhesion
• Pinked, nonfraying edges
• Strong woven backing
• Gentle removal

Tapebox

Specify:
1/2", 24/box**(101-4305)**
1", 12/box**(101-9673)**
2", 6/box**(101-1332)**



HENRY SCHEIN®
PREMIUM HYPALLERGENIC TRANSPARENT SURGICAL TAPE

Latex-Free, 10-yd Rolls
Excellent adhesion. Porous, with flexible backing. Tears easily in both directions.

Tapebox

Specify:
1/2", 24/box**(101-6655)**
1", 12/box**(101-4645)**
2", 6/box**(101-9226)**



HENRY SCHEIN®
PACKING STRIPS

5-yd Roll
• 100% cotton
• Nonfraying
• Nonlinting
• Selvage edge
• Absorbs fluid; less chance of leaving foreign matter in the wound

1/4"
(104-7976)ea

1/2"
(104-6304)ea

1"
(104-8748)ea

2"
(104-6302)ea



HENRY SCHEIN®
PREMIUM HYPALLERGENIC CLOTH SURGICAL TAPE

Latex-Free, 10-yd Rolls
• Reliable adhesion
• Highly porous
• Thin, flexible backing
• Gentle removal

Tapebox

Specify:
1/2", 24/box**(101-3025)**
1", 12/box**(101-2323)**
2", 6/box**(101-5979)**

HENRY SCHEIN®
FLEXIBLE FABRIC ELASTIC STRIPS

Stretch and conform to permit free movement. Will not stick to wound.

3/4" x 3"
(900-4499)100/box

1" x 3"
(900-4500)100/box

1 1/2" x 3" Knuckle
(900-4496)100/box



HENRY SCHEIN®
SHEER STRIPS

Sterile
Individually wrapped adhesive bandages for the protection of minor cuts, punctures, and abrasions. Highly absorbent pads won't stick to wound. Air vents let air in and keep moisture out.

3/4" x 3"
(900-4503)100/box


Single Box
1" x 3"
(900-4504)100/box

Sterile, 7/8"
(900-4505)100/box

Professional Tray Pack, 3/4" x 3"
(900-4501)1500/case

Professional Tray Pack, 1" x 3"
(900-4502)1500/case






Call: **1-800-P-SCHEIN (1-800-772-4346)**
8am-9pm (et)

Fax: **1-800-329-9109** 24 Hours

www.henryschein.com/hsibrand

Due to manufacturer price increases, prices are subject to change without notification. Offer valid on our full array of products, excluding flu vaccine and Synagis®. Upon any request by Medicare, Medicaid, Tricare or other payer, you must report the value of any benefit received from a discount program (e.g., points redemptions, gifts or other special awards). ©2012 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors. No adjustments from prior sales. Offer not to be combined with other promotions and/or special contract pricing. Void where prohibited. Local regulatory requirements may apply to use or installation of certain products. It is your responsibility to understand and comply with any such requirements prior to purchase, use or installation of products. †For informational purposes only. Source: ReimbursementCodes.com and/or www.cms.gov. Customer is responsible for verification of billing/coding in accordance with applicable specific circumstances.



12MS6111B



Henry Schein supplies standard, as well as custom surgical procedure packs, for health care practices across the United States. Henry Schein works directly with each of our customers to design a comprehensive surgical supply solution for their individual practices.

SurgiTeam takes pride in being an integral part of solutions that deliver increased cost savings, greater office efficiencies, higher profits and optimal patient results.

- ❖ SurgiTeam™ Packs are fully CUSTOMIZED to meet your exact surgical requirements.
- ❖ SurgiTeam™ Packs SAVE your facility money by reducing direct supply costs and surgical set-up time.
- ❖ Design a FREE custom surgical pack for your facility today.

For more information contact your Sales Consultant.