



PRODUCT CATALOG

A full range of PPE options just
right for health care professionals.

FROM HERE



TO HERE

WEAR IT WITH
CONFIDENCE.



CRITERION[®] by  HENRY SCHEIN[®]

CRITERION[®] provides a smart selection of PPE products that let you outfit your practice from head-to-toe with quality you can rely on.

Contents

GLOVES

8	Nitrile Exam N100 Gloves
8	Nitrile Exam N200 Gloves
9	Nitrile Exam N250 Gloves
9	Nitrile Exam N300 Gloves
9	Nitrile Exam Low Derma Gloves
9	Nitrile Exam Aloe + Vitamin E Gloves
10	Nitrile Exam Colloidal Oatmeal Gloves
10	Chloroprene Exam Gloves
10	Chloroprene Exam Ultrathin Gloves
10	Latex Exam Polymer Coated Gloves
11	Latex Exam Aloe + Vitamin E Gloves
11	Latex Exam Chlorinated Gloves
11	Vinyl Exam Gloves
12	Latex Surgical Gloves
12	Polychloroprene Surgical Gloves
12	Polyisoprene Surgical Gloves

MASKS

13	Earloop Masks Level 1
14	Sensitive Skin Earloop Masks Level 1
14	Fitted Earloop Masks Level 1
14	Extra Breathing Room Earloop Masks Level 1
14	Pink Ribbon Earloop Masks Level 1
14	Surgical Tie-On Masks Level 1
15	Earloop Masks Level 2
15	Earloop Masks Level 3
15	Earloop Masks Level 3 with Eye Shield
15	Pediatric Earloop Masks

DISPOSABLE GARMENTS

16	SMS Jackets
17	SMS Lab Coats
18	SMS Isolation Gowns
18	Polypropylene Chemo Gowns
19	SMS Cover Gowns
19	Polypropylene Cover Gowns
19	Polycoated Polypropylene Aprons
19	Bouffant Caps
19	Shoe Covers

PROTECTIVE EYEWEAR

20	Protective Eyewear
20	Safety Glasses
20	Temple Grip Protective Eyewear
21	Foam Padded Protective Eyewear
21	Disposable Face Shields

GLOVE ACCESSORIES

22	Acrylic Glove Box Holders
22	Wire Glove Box Holders

Product Branding Updates

WHAT TO EXPECT FROM HEAD TO TOE

All your favorite CRITERION® gloves will start showing up in new packaging over the next several months. In addition, Maxima® Masks and Maxi-Gard™ Masks, Jackets, Lab Coats, and Protective Eyewear will transition to the CRITERION® brand. While the name and packaging will look different, the product inside the box will remain the same.

CURRENT NAME	NEW NAME	PAGE NUMBER
CRITERION® GLOVES		
Criterion® Nitrile Exam Glove	CRITERION® Nitrile Exam N100 Gloves	8
Criterion® N200 Nitrile Exam Glove	CRITERION® Nitrile Exam N200 Gloves	8
Criterion® N200 Sky Blue Nitrile Exam Glove	CRITERION® Nitrile Exam N200 Gloves	8
Criterion® N250 Nitrile Exam Glove	CRITERION® Nitrile Exam N250 Gloves	9
Criterion® N300 Ultra Nitrile Exam Glove	CRITERION® Nitrile Exam N300 Gloves	9
Criterion® Pure Freedom Nitrile Exam Glove	CRITERION® Nitrile Exam Low Derma Gloves	9
Criterion® Radiant Razz Nitrile Exam Glove	CRITERION® Nitrile Exam Aloe + Vitamin E Gloves	9
Criterion® Coats Nitrile Exam Glove	CRITERION® Nitrile Exam Colloidal Oatmeal Gloves	10
Criterion® CR Chloroprene Exam Glove	CRITERION® Chloroprene Exam Gloves	10
Criterion® CR4G Chloroprene Exam Glove	CRITERION® Chloroprene Exam Ultrathin Gloves	10
Criterion® PC Latex Exam Glove	CRITERION® Latex Exam Polymer Coated Gloves	10
Criterion® Aloe Advanced Latex Exam Glove	CRITERION® Latex Exam Aloe + Vitamin E Gloves	11
Criterion® CL Latex Exam Glove	CRITERION® Latex Exam Chlorinated Gloves	11
Criterion® PF Vinyl Glove	CRITERION® Vinyl Exam Gloves	11
Criterion® Latex Surgeons Glove	CRITERION® Latex Surgical Gloves	12
Criterion® CR Surgeons Glove	CRITERION® Polychloroprene Surgical Gloves	12
Criterion® PI Surgical Glove	CRITERION® Polyisoprene Surgical Gloves	12
CRITERION® MASKS		
Procedure Masks Level 1	CRITERION® Earloop Masks Level 1	13
Ultrasoft Procedure Masks Level 1	CRITERION® Sensitive Skin Earloop Masks Level 1	14

Product Branding Updates (CONT.)

CURRENT NAME	NEW NAME	PAGE NUMBER
MaxiGard™ SFT Masks Level 1	CRITERION® Fitted Earloop Masks Level 1	14
Maxi-Gard™ Air Procedure Masks With Air Space™ Technology Level 1	CRITERION® Extra Breathing Room Earloop Masks Level 1	14
Procedure Pink Ribbon Mask Level 1	CRITERION® Pink Ribbon Earloop Masks Level 1	14
Pleated Surgical Masks Level 1	CRITERION® Surgical Tie-On Masks Level 1	14
Procedure Masks Level 2	CRITERION® Earloop Masks Level 2	15
Maxima® Procedure Masks Level 3	CRITERION® Earloop Masks Level 3	15
Earloop/Eye Shield Procedure Masks Level 3	CRITERION® Earloop Masks Level 3 with Eye Shield	15
Pediatric Face Mask	CRITERION® Pediatric Earloop Masks	15

CRITERION® DISPOSABLE GARMENTS

Maxi-Gard® Jackets	CRITERION® SMS Jackets	16
Maxi-Gard® Lab Coats	CRITERION® SMS Lab Coats	17
Isolation Gowns	CRITERION® SMS Isolation Gowns	18
Chemo Gown	CRITERION® Polypropylene Chemo Gowns	18
Maxi-Gard® Cover Gowns	CRITERION® SMS Cover Gowns	19
Cover Gowns	CRITERION® Polypropylene Cover Gowns	19
Polycoated Aprons	CRITERION® Polycoated Polypropylene Aprons	19
Bouffant Caps	CRITERION® Bouffant Caps	19
Shoe Covers	CRITERION® Shoe Covers	19

CRITERION® PROTECTIVE EYEWEAR

Protective Eyewear	CRITERION® Protective Eyewear	20
Safety Glasses	CRITERION® Safety Glasses	20
Maxi-Gard® Protective Eyewear	CRITERION® Temple Grip Protective Eyewear	20
Maxi-Gard® Protective Eyewear	CRITERION® Foam Padded Protective Eyewear	21
Disposable Face Shields	CRITERION® Disposable Face Shields	21

CRITERION® GLOVE ACCESSORIES

Maxi-Dispenz® Acrylic Glove Box Dispensers	CRITERION® Acrylic Glove Box Holders	22
Maxi-Dispenz® Wire Glove Box Dispensers	CRITERION® Wire Glove Box Holders	22

The New CRITERION®

Quality head-to-toe PPE, now in one reliable product line



1. UPDATED PACKAGING

makes it easy to know what's inside the box with attractive imagery and bold design elements.



2. INFO ZONE

appears in the bottom right of all packaging. The icons indicate key product information: PPE standards, size, color, and quantity.



3. REORDER PANEL

helps to simplify the reordering process by clearly stating the full product name and reorder item number.

CRITERION® NEW LOOK



You may receive either package during this transition



original packaging



original packaging



original packaging



original packaging

Gloves

Comparison Chart

GLOVE NAME	SIZES	COLOR	MIL*	KEY FEATURES
NITRILE				
CRITERION® Nitrile Exam N100 Gloves	XS–XL	Blue	3.15	Textured fingers. Chemotherapy tested to ASTM D6978 standard.
CRITERION® Nitrile Exam N200 Gloves	XS–XL	Periwinkle Blue	2.36	Textured fingers. Chemotherapy tested to ASTM D6978 standard.
CRITERION® Nitrile Exam N200 Gloves	S–XL	Sky Blue	2.36	Textured fingers. Chemotherapy tested to ASTM D6978 standard.
CRITERION® Nitrile Exam N250 Gloves	S–XL	Lilac	2.00	Textured fingers. Chemotherapy tested to ASTM D6978 standard.
CRITERION® Nitrile Exam N300 Gloves	XS–XL	Ice Blue	2.00	Eco-Friendly Packaging. Textured fingers. Chemotherapy tested to ASTM D6978 standard.
CRITERION® Nitrile Exam Low Derma Gloves	XS–XL	Slate Blue	2.00	Low Dermatitis Potential. Textured fingers. Chemotherapy tested to ASTM D6978 standard.
CRITERION® Nitrile Exam Aloe + Vitamin E Gloves	XS–L	Pink	2.00	Aloe + Vitamin E Coated. Non-Chlorinated. Textured fingers.
CRITERION® Nitrile Exam Colloidal Oatmeal Gloves	XS–XL	White	2.00	Colloidal Oatmeal Coated. Gluten Free. Textured fingers.
CHLOROPRENE				
CRITERION® Chloroprene Exam Gloves	XS–XL	Green	3.15	Textured fingers.
CRITERION® Chloroprene Exam Ultrathin Gloves	XS–XL	Lime	2.00	Textured fingers.
LATEX				
CRITERION® Latex Exam Polymer Coated Gloves	XS–XL	White	3.15	Polymer Coated. Non-Chlorinated. Textured fingers and palm.
CRITERION® Latex Exam Aloe + Vitamin E Gloves	XS–XL	White	3.15	Aloe + Vitamin E Coated. Non-Chlorinated. Textured fingers and palm.
CRITERION® Latex Exam Chlorinated Gloves	XS–XL	White	4.33	Chlorinated. Textured fingers and palm.
VINYL				
CRITERION® Vinyl Exam Gloves	XS–XL	Clear	3.15	Smooth. Accelerator Free.
SURGICAL				
CRITERION® Latex Surgical Gloves	5.5–9.0	White	6.30	Polymer Coated. Microtextured. Extended Cuff. Sterile.
CRITERION® Polychloroprene Surgical Gloves	5.5–9.0	White	5.91	Polymer Coated. Palm and Finger Textured. Extended Cuff. Sterile.
CRITERION® Polyisoprene Surgical Gloves	5.5–9.0	White		Skin-friendly accelerators. Microtextured. Extended Cuff. Sterile.

*Thickness at Palm

Gloves



CRITERION® NITRILE EXAM N100 GLOVES

NITRILE

BLUE COLOR

- Durable nitrile construction
- Chlorinated for reduced friction and easy donning
- Textured fingers provide enhanced grip
- Chemotherapy tested to ASTM D6978 standard
- Beaded cuff. 3.15-mil thickness at palm

X-Small	100/box	111-8533
Small	100/box	111-8535
Medium	100/box	111-8536
Large	100/box	111-8537
X-Large	90/box	111-8539

CRITERION® NITRILE EXAM N200 GLOVES

NITRILE

PERIWINKLE BLUE AND SKY BLUE COLORS

- Durable nitrile construction, slightly thinner than the N100
- Chlorinated for reduced friction and easy donning
- Textured fingers provide enhanced grip
- Chemotherapy tested to ASTM D6978 standard
- Beaded cuff. 2.36-mil thickness at palm

PERIWINKLE BLUE

X-Small	200/box	900-7437
Small	200/box	900-7438
Medium	200/box	900-7439
Large	200/box	900-7440
X-Large	180/box	900-7441

SKY BLUE

Small	200/box	570-2098
Medium	200/box	570-2099
Large	200/box	570-2100
X-Large	180/box	570-2101

CRITERION® Gloves (CONT.)

CRITERION® NITRILE EXAM N250 GLOVES

NITRILE

LILAC COLOR

- Ultrathin formulation for tactile sensitivity and dexterity
- Chlorinated for reduced friction and easy donning
- Textured fingers provide enhanced grip
- Chemotherapy tested to ASTM D6978 standard
- Beaded cuff. 2.0-mil thickness at palm

Small	250/box	570-1947
Medium	250/box	570-1948
Large	250/box	570-1949
X-Large	230/box	570-1950

CRITERION® NITRILE EXAM N300 GLOVES

NITRILE

ICE BLUE COLOR

- Ultrathin formulation for tactile sensitivity and dexterity
- Chlorinated for reduced friction and easy donning
- Textured fingers provide enhanced grip
- Chemotherapy tested to ASTM D6978 standard
- Beaded cuff. 2.0-mil thickness at palm
- Eco-friendly packaging, with 300 gloves in a box, reduces waste and saves space

X-Small	300/box	112-6829
Small	300/box	112-6830
Medium	300/box	112-6831
Large	300/box	112-6832
X-Large	250/box	112-6833

CRITERION® NITRILE EXAM LOW DERMA GLOVES

NITRILE

SLATE BLUE COLOR

- Made from a soft and durable nitrile formulation without the use of chemical accelerators or sulfur, preventing type-IV allergic reactions and reduces the potential for dermatitis
- Chlorinated for reduced friction and easy donning
- Textured fingers provide enhanced grip
- Chemotherapy tested to ASTM D6978 standard
- Beaded cuff. 2.0-mil thickness at palm

X-Small	200/box	570-0630
Small	200/box	570-0631
Medium	200/box	570-0632
Large	200/box	570-0633
X-Large	180/box	570-0634

CRITERION® NITRILE EXAM ALOE + VITAMIN E GLOVES

NITRILE

PINK COLOR

- Enhanced inner coating of Aloe + Vitamin E to soothe, moisturize, and rejuvenate hands
- Non-chlorinated, fast and easy donning
- Textured fingers provide enhanced grip
- Beaded cuff. 2.0-mil thickness at palm

X-Small	200/box	112-7143
Small	200/box	112-7144
Medium	200/box	112-7146
Large	200/box	112-7174

CRITERION® Gloves (CONT.)

CRITERION® NITRILE EXAM COLLOIDAL OATMEAL GLOVES

NITRILE

WHITE COLOR

- Made from 100% all natural colloidal oatmeal
- Provides continuous protection from prolonged glove use and hand washing
- Gluten-free
- Chlorinated for reduced friction and easy donning
- Textured fingers provide enhanced grip
- Beaded cuff. 2.0-mil thickness at palm

X-Small	200/box	112-6220
Small	200/box	112-6221
Medium	200/box	112-6222
Large	200/box	112-6223
X-Large	200/box	112-6224

CRITERION® CHLOROPRENE EXAM GLOVES

CHLOROPRENE

GREEN COLOR

- Provides the comfort and fit of latex, but with the synthetic strength of nitrile
- Puncture and tear resistant for high-risk procedures
- Textured fingertips for enhanced grip
- Proprietary polymer coating for "glide-on" donning
- Beaded cuff. 3.15-mil thickness at palm

X-Small	100/box	112-6857
Small	100/box	112-6858
Medium	100/box	112-6859
Large	100/box	112-6860
X-Large	90/box	112-6861

CRITERION® CHLOROPRENE EXAM ULTRATHIN GLOVES

CHLOROPRENE

LIME COLOR

- Ultrathin polychloroprene glove provides greater dexterity for precision procedures and extended-wear comfort
- Puncture and tear resistant for high-risk procedures
- Textured fingertips for enhanced grip
- Proprietary polymer coating for "glide-on" donning
- Beaded cuff. 2.0-mil thickness at palm

X-Small	100/box	112-7306
Small	100/box	112-7307
Medium	100/box	112-7308
Large	100/box	112-7309
X-Large	90/box	112-7310

CRITERION® LATEX EXAM POLYMER COATED GLOVES

LATEX

WHITE COLOR

- Polymer coating provides a smooth interior finish for fast and easy donning
- Comfortable fit and feel provide superior range of motion
- Textured finish for enhanced grip
- Beaded cuff. 3.15-mil thickness at palm

X-Small	100/box	102-5418
Small	100/box	102-5419
Medium	100/box	102-5421
Large	100/box	102-6730
X-Large	90/box	102-5422

CRITERION® Gloves (CONT.)

CRITERION® LATEX EXAM ALOE + VITAMIN E GLOVES

LATEX

WHITE COLOR

- Enhanced inner coating of Aloe + Vitamin E to soothe and moisturize
- Fast and easy donning
- Textured finish for enhanced grip
- Beaded cuff. 3.15-mil thickness at palm

X-Small	100/box	112-5840
Small	100/box	112-5841
Medium	100/box	112-5842
Large	100/box	112-5843
X-Large	90/box	112-5844

CRITERION® LATEX EXAM CHLORINATED GLOVES

LATEX

WHITE COLOR

- Chlorinated for reduced friction and lower protein levels
- Textured finish for enhanced grip
- Beaded cuff. 4.33-mil thickness at palm

X-Small	100/box	102-5354
Small	100/box	102-8296
Medium	100/box	102-5355
Large	100/box	102-7763
X-Large	90/box	102-5357

CRITERION® VINYL EXAM GLOVES

VINYL

CLEAR COLOR

- Non-latex formulation is an ideal solution for those at risk of latex sensitivity
- Accelerator-free and non-chlorinated
- Smooth finish for maximum tactile sensitivity
- Beaded cuff. 3.15-mil thickness at palm

X-Small	100/box	570-0087
Small	100/box	570-0088
Medium	100/box	570-0089
Large	100/box	570-0090
X-Large	90/box	570-0091



CRITERION® Gloves (CONT.)

CRITERION® LATEX SURGICAL GLOVES

SURGICAL

WHITE COLOR

- Advanced polymer coating on the inner surface allows for easy donning when hands are wet, damp, or dry
- Microtextured surface provides optimal grip in all conditions
- Supersoft latex film and anatomical design reduces hand fatigue during long procedures
- Minimum cuff length of 275mm
- Beaded cuff. 6.3-mil thickness at palm
- Sterile

Size 5.5	50 pairs/box	570-2102
Size 6.0	50 pairs/box	570-2103
Size 6.5	50 pairs/box	570-2104
Size 7.0	50 pairs/box	570-2105
Size 7.5	50 pairs/box	570-2106
Size 8.0	50 pairs/box	570-2107
Size 8.5	50 pairs/box	570-2108
Size 9.0	50 pairs/box	570-2109

CRITERION® POLYCHLOROPRENE SURGICAL GLOVES

SURGICAL

WHITE COLOR

- Proprietary anatomical former design for a more natural fit
- Polymer-coating finish that allows for damp or dry-hand donning
- Textured fingers and palm for enhanced grip and control
- Ultra-low modulus for enhanced comfort and reduced hand fatigue
- Minimum cuff length of 300mm for extended coverage and protection
- Beaded cuff. 5.91-mil thickness at palm
- Sterile

Size 5.5	50 pairs/box	112-7079
Size 6.0	50 pairs/box	112-7080
Size 6.5	50 pairs/box	112-7081
Size 7.0	50 pairs/box	112-7082
Size 7.5	50 pairs/box	112-7083
Size 8.0	50 pairs/box	112-7085
Size 8.5	50 pairs/box	112-7086
Size 9.0	50 pairs/box	112-7088

CRITERION® POLYISOPRENE SURGICAL GLOVES

SURGICAL

WHITE COLOR

- Similar structure to latex
- Fully anatomical shape reduces hand fatigue
- Skin-friendly accelerators reduce the potential for allergic reaction
- Micro-rough texture for enhanced grip
- Easy donning
- Extended and beaded cuff
- Sterile

Size 5.5	50 pairs/box	570-1119
Size 6.0	50 pairs/box	570-1120
Size 6.5	50 pairs/box	570-1121
Size 7.0	50 pairs/box	570-1122
Size 7.5	50 pairs/box	570-1123
Size 8.0	50 pairs/box	570-1124
Size 8.5	50 pairs/box	570-1125
Size 9.0	50 pairs/box	570-1126

Masks



ASTM Chart

Not all masks are created equal.

Know the key facts about ASTM-rated mask levels so you can ensure you have the right protection.

LEVEL 1

- Low barrier protection
- For general use
- Not used for aerosols, spray, or fluids

LEVEL 2

- Moderate barrier protection
- Use for low to moderate levels of aerosols, spray, and/or fluids

LEVEL 3

- Maximum/high barrier protection
- Use for high risk of fluid, spray, and/or fluids

ASTM REQUIREMENTS	LEVEL 1	LEVEL 2	LEVEL 3
Particle Filtration Efficiency (PFE)	≥ 95% @ 0.1 micron	≥ 98% @ 0.1 micron	≥ 98% @ 0.1 micron
Bacteria Filtration Efficiency (BFE)	≥ 95% @ 3.0 microns	≥ 98% @ 3.0 microns	≥ 98% @ 3.0 microns
Breathability (Delta P)	<5 mm H ₂ O/cm ²	<6 mm H ₂ O/cm ²	<6 mm H ₂ O/cm ²
Fluid Resistance	≥ 80 mm Hg	≥ 120 mm Hg	≥ 160 mm Hg
Flame Spread	Class 1	Class 1	Class 1

CRITERION® EARLOOP MASKS LEVEL 1

LEVEL 1

- Reliable, fluid-resistant mask meets ASTM level 1 Performance Standards
- Lightweight and comfortable for extended wear
- Three-ply design with soft inner and outer facing
- Extra-long adjustable nosepiece
- Not made with natural rubber latex

Blue	50/box	104-3809
Blue (U.S. Made)	50/box	570-2443
Lavender	50/box	104-8600
Pink	50/box	104-3730
White	50/box	104-6611
Yellow	50/box	104-2849

CRITERION® Masks (CONT.)

CRITERION® SENSITIVE SKIN EARLOOP MASKS LEVEL 1

LEVEL 1

- Meets ASTM Level 1 Performance Standards
- Hypoallergenic, free of chemicals, inks, and dyes
- Fluid-resistant mask will not tear, shred, or lint
- Spunbound outer layer with white inner facing
- Not made with natural rubber latex

Blue	50/box	104-8072
Pink	50/box	104-8073

CRITERION® FITTED EARLOOP MASKS LEVEL 1

LEVEL 1

- Meets ASTM Level 1 Performance Standards
- Extra-long adjustable nose and chin piece secures a tight custom fit
- Fluid resistant, nonirritating soft inner and outer facing
- Not made with natural rubber latex

Blue	50/box	900-7451
Pink	50/box	900-7452

CRITERION® EXTRA BREATHING ROOM EARLOOP MASKS LEVEL 1

LEVEL 1

- Meets ASTM Level 1 Performance Standards
- Contoured design enhances airflow space and provides complete coverage and fit
- Extra breathing room allows for comfort and extended wear
- Extra-long adjustable nosepiece

Blue	50/box	570-0722
Lavender	50/box	570-0724
White	50/box	570-0723

CRITERION® PINK RIBBON EARLOOP MASKS LEVEL 1

LEVEL 1

- White Mask with pink-ribbon design supporting breast cancer awareness
- Meets ASTM level 1 Performance Standards
- Lightweight and comfortable for extended wear with fluid resistance
- Three-ply, soft inner and outer facing with extra long adjustable nosepiece
- Not made with natural rubber latex

Pink Ribbon	50/box	570-2829
-------------	--------	----------

CRITERION® SURGICAL TIE-ON MASKS LEVEL 1

LEVEL 1

- Meets ASTM Level 1 Performance Standards
- Pleated with cloth-like ties behind the head for a secure fit
- Fluid-resistant mask
- Not made with natural rubber latex

Blue	50/box	104-6809
------	--------	----------

CRITERION® Masks (CONT.)

CRITERION® EARLOOP MASKS LEVEL 2

LEVEL 2

- Meets ASTM Level 2 Performance Standards
- Fluid resistant, light, and comfortable fit for extended use
- Three-ply, breathable spunbound outer layer with white inner facing
- Extra long adjustable nosepiece
- Not made with natural rubber latex

Blue	50/box	104-7799
Lavender	50/box	104-8071
Pink	50/box	104-9200

CRITERION® EARLOOP MASKS LEVEL 3

LEVEL 3

- Meets ASTM Level 3 Performance Standards
- Comfortable and lightweight for extended wear
- Three-ply soft inner and outer facing with extra long adjustable nosepiece
- Not made with natural rubber latex

Blue	50/box	570-2044
Blue (U.S. Made)	50/box	570-2448
Lavender	50/box	570-2047
Pink	50/box	570-2045
Teal	50/box	570-2042
White	50/box	570-2043
Yellow	50/box	570-2046

CRITERION® EARLOOP MASKS LEVEL 3 WITH EYE SHIELD

LEVEL 3

- Meets ASTM Level 3 Performance Standards
- Highest filtration rate and fluid-penetration resistance
- Stretchable, soft, earloops provide a comfortable fit
- Antifog material, optically clear shield is treated on both sides
- Not made with natural rubber latex

Teal	25/box	104-7321
------	--------	----------

CRITERION® PEDIATRIC EARLOOP MASKS

- Designed to fit small faces
- Soft, non-woven material for enhanced breathability and comfort
- Optimal respiratory protection for children ages 4–12
- Adjustable nosepiece provides custom fit
- Smooth, stretchy, ear loops prevent irritation and allow easy removal
- Not made with natural rubber latex

Blue	50/box	570-2268
------	--------	----------

Disposable Garments



CRITERION® SMS JACKETS

Stylish hip-length jacket made from high-quality SMS fabric meets OSHA fluid-resistance requirements for protection. Lightweight, breathable material. 3-pocket design and snap front closure. Not made with natural rubber latex.

SMALL 10/pack

● Black	570-2288
● Blueberry	570-2289
● Ceil Blue	570-2290
● Light Pink	570-2291
● Purple	570-2292
● Raspberry	570-2293
● Sky Blue	570-2294
● Teal	570-2295
● White	570-2296

LARGE 10/pack

● Black	570-2327
● Blueberry	570-2328
● Ceil Blue	570-2329
● Light Pink	570-2297
● Purple	570-2278
● Raspberry	570-2330
● Sky Blue	570-2331
● Teal	570-2332
● White	570-2333

2XL 10/pack

● Ceil Blue	570-2307
-------------	----------

MEDIUM 10/pack

● Black	570-2298
● Blueberry	570-2299
● Ceil Blue	570-2300
● Light Pink	570-2301
● Purple	570-2302
● Raspberry	570-2303
● Sky Blue	570-2304
● Teal	570-2305
● White	570-2306

X-LARGE 10/pack

● Black	570-2318
● Blueberry	570-2319
● Ceil Blue	570-2320
● Light Pink	570-2321
● Purple	570-2322
● Raspberry	570-2323
● Sky Blue	570-2324
● Teal	570-2325
● White	570-2326

3XL 10/pack

● Ceil Blue	570-2308
-------------	----------

CRITERION® Garments (CONT.)

CRITERION® SMS LAB COATS

- Stylish knee-length jacket made from high-quality SMS fabric meets OSHA fluid-resistance requirements for protection
- Lightweight, breathable material
- 3-pocket design and snap front closure
- Not made with natural rubber latex

SMALL 10/pack

Black	570-2334
Blueberry	570-2335
Ceil Blue	570-2336
Light Pink	570-2337
Purple	570-2338
Raspberry	570-2339
Sky Blue	570-2340
Teal	570-2341
White	570-2342

LARGE 10/pack

Black	570-2279
Blueberry	570-2280
Ceil Blue	570-2281
Light Pink	570-2282
Purple	570-2283
Raspberry	570-2284
Sky Blue	570-2285
Teal	570-2286
White	570-2287

Ceil Blue	570-2343
-----------	----------

MEDIUM 10/pack

Black	570-2346
Blueberry	570-2347
Ceil Blue	570-2348
Light Pink	570-2349
Purple	570-2350
Raspberry	570-2351
Sky Blue	570-2352
Teal	570-2353
White	570-2354

X-LARGE 10/pack

Black	570-2309
Blueberry	570-2310
Ceil Blue	570-2311
Light Pink	570-2312
Purple	570-2313
Raspberry	570-2314
Sky Blue	570-2315
Teal	570-2316
White	570-2317

3XL 10/pack

Ceil Blue	570-2344
-----------	----------



Black



Blueberry



Ceil Blue



Light Pink



Purple



Raspberry



Sky Blue



Teal



White

CRITERION® Garments (CONT.)

AAMI Chart

CHOOSE THE RIGHT ISOLATION GOWN FOR THE RIGHT PROCEDURE

The FDA recognizes the Association of the Advancement of Medical Instrumentation (AAMI) for protective apparel. The standard assists end-users in determining the types of protective products most appropriate for a particular task or situation.

<https://www.fda.gov/medical-devices/personal-protective-equipment-infection-control/medical-gowns#g1>

LEVEL 1	Minimal Risk Procedures
LEVEL 2	Low Risk Procedures
LEVEL 3	Moderate Risk Procedures

CRITERION® SMS ISOLATION GOWNS

Single-use, trilayered isolation gowns made from SMS fabric. Antistatic properties, elastic cuffs, and waist ties. AAMI Level 1–3.

AAMI LEVEL 1

Medium	Blue	10/bag	570-2357	Medium	Yellow	10/bag	570-2358
Large	Blue	10/bag	570-2359	Large	Yellow	10/bag	570-2360
X-Large	Blue	10/bag	570-2361	X-Large	Yellow	10/bag	570-2362

AAMI LEVEL 2

Medium	Blue	10/bag	570-2375	Medium	Yellow	10/bag	570-2376
Large	Blue	10/bag	570-2377	Large	Yellow	10/bag	570-2378
X-Large	Blue	10/bag	570-2363	X-Large	Yellow	10/bag	570-2364

AAMI LEVEL 3

Medium	Blue	10/bag	570-2371	Medium	Yellow	10/bag	570-2372
Large	Blue	10/bag	570-2373	Large	Yellow	10/bag	570-2374
X-Large	Blue	10/bag	570-2355	X-Large	Yellow	10/bag	570-2356

CRITERION® POLYPROPYLENE CHEMO GOWNS

- Protects against cross contamination, dirt, and grime
- Disposable non-sterile, durable gown neck ties for an adjustable fit

Large	Blue	10/pack	570-2211
X-Large	Blue	10/pack	570-2212

CRITERION® Garments (CONT.)

CRITERION® SMS COVER GOWNS

- Knee-length cover gown with complete front coverage
- Dual-fabric design
- Sleeves and front of gown are made from SMS fabric that meets OSHA fluid resistant requirements
- Back of gown is made of natural lightweight, breathable material
- Adjustable knitted cuffs and open-back design with neck and waist ties
- Not made with natural rubber latex

Regular	Blue	10/pack	112-6245
---------	------	---------	----------

X-Large	Blue	10/pack	112-6246
---------	------	---------	----------

CRITERION® POLYPROPYLENE COVER GOWNS

- Full-length gown made of fluid-resistant polypropylene material

Small	Pink	10/pack	100-6784
-------	------	---------	----------

Medium	Blue	10/pack	100-5870
--------	------	---------	----------

Large	Pink	10/pack	100-9392
-------	------	---------	----------

X-Large	Blue	10/pack	100-5806
---------	------	---------	----------

CRITERION® POLYCOATED POLYPROPYLENE APRONS

- Full-length aprons are excellent for hygiene and protection
- Made of polycoated polypropylene material

White	25/pack	100-7483
-------	---------	----------

CRITERION® BOUFFANT CAPS

- Comfortable, lightweight and breathable bouffant cap provides total hair coverage
- Secure fit prevents hair loss from contaminating sensitive area or site
- Elastic not made with natural rubber latex

21"	Blue	100/box	104-6032
-----	------	---------	----------

21"	White	100/box	102-4151
-----	-------	---------	----------

24"	Blue	100/box	104-6058
-----	------	---------	----------

24"	White	100/box	102-6614
-----	-------	---------	----------

CRITERION® SHOE COVERS

- Non-sterile, nonconductive, with a nonskid design
- Easy-to-use, one size fits all
- Not made with natural rubber latex

One Size	Blue	50 pairs/box	104-9452
----------	------	--------------	----------

One Size	Blue	150 pairs/box	104-6350
----------	------	---------------	----------

Protective Eyewear



CRITERION® PROTECTIVE EYEWEAR

- Lightweight protective eyewear
- Wraparound 10-base curve for complete coverage
- Gray lens and temples
- Ideal for patient and staff use

Gray Lens

570-0435

CRITERION® SAFETY GLASSES

- Clear plastic lens with side shields for eye protection

Clear Lens

6/box

100-7382

CRITERION® TEMPLE GRIP PROTECTIVE EYEWEAR

- Rubber temples provide extra grip and comfort
- Clear, shatter-resistant polycarbonate lens with scratch-resistant coating
- Wraparound 10-base curve allows for complete coverage
- Meets ANSI Z87+ standards*

Black

6/box

900-4436

Blue

6/box

900-4488

Pink

6/box

900-4779

Red

6/box

900-4781

CRITERION® Protective Eyewear (CONT.)

CRITERION® FOAM PADDED PROTECTIVE EYEWEAR

- Foam inserts offer comfort and protection
- Clear, shatter-resistant polycarbonate lens with scratch-resistant and antifog coating
- Wraparound 10-base curve allows for complete coverage
- Meets ANSI Z87+ standards*

Black	570-0434
Purple	570-0432
Pink	570-0431
Blue	570-0433

CRITERION® DISPOSABLE FACE SHIELDS

- Optically clear, full-length visor provides increased protection against splash, splatter, and aerosols
- Features a 1.5" thick, 9" long foam headband that keeps the face shield off the face
- Ultrasonically welded, 13" elastic headband for added strength and increased wearer comfort
- Special fog-free coating on inside and outside visor
- Not made with natural rubber latex

24/box	101-2254
--------	----------

CRITERION®
COMFORT IN PROTECTION™

A Smart Selection
of PPE Products for
Head-to-Toe Protection

by **HENRY SCHEIN®**



*OSHA RECOGNIZES THE AMERICAN NATIONAL STANDARDS INSTITUTE (ANSI) Z87.1 AND Z87+ STANDARDS FOR EYE PROTECTION.
<https://www.osha.gov/eye-face-protection/standards>

Glove Accessories



CRITERION® ACRYLIC GLOVE BOX HOLDERS

- Clear acrylic glove box holders
- Available in 4 sizes to meet any office or facility's needs
- Accommodates most brands and sizes
- Option to wall mount vertically or horizontally with predrilled holes (screws not included)

1-Box Holder	102-6467
2-Box Holder	102-6468
3-Box Holder	102-8344
4-Box Holder	102-7201

CRITERION® WIRE GLOVE BOX HOLDERS

- Durable steel wire frame glove box holders with white-powdered finish
- Selection of 5 different choices meets any office or facility's needs
- Hardware is included

Dual Box Holder	112-5471
2-Box Horizontal Holder	112-5473
2-Box Vertical Holder	112-5472
3-Box Horizontal Holder	112-5475
3-Box Vertical Holder	112-5474

Terms of Sale

THE HENRY SCHEIN PRICE POLICY:

We endeavor to maintain prices for the duration of a catalog, but we reserve the right to make price adjustments in response to manufacturers' price increases or extraordinary circumstances. Prices are subject to change without notice.

Henry Schein, Inc. ("Henry Schein") and customer agree that the terms and conditions hereinafter set forth shall govern the relationship between Henry Schein and the customer to the extent that the parties do not have a written agreement in effect that conflicts with such terms and conditions. Customer acknowledges and accepts all such terms and conditions by placing an order for goods with Henry Schein, and upon Henry Schein's delivery of the order to the customer.

Choose Your Payment Method

Reduce the cost and administration of paying Henry Schein—

Pay electronically (ACH Debit) or set up AutoPay. Please call Customer Service for details.

For your convenience, we provide several payment alternatives. Orders billed to your account may be paid by ACH Debit, Check by Phone, or Check. If you prefer, you may use your Henry Schein Credit Card, American Express, Visa, MasterCard or Discover Card when placing your order. All sales are subject to our normal terms and conditions.

Unless otherwise instructed, check payments must be mailed to: Henry Schein, Inc.

• Dept. CH 10241 • Palatine, IL 60055-0241

All sales are subject to credit approval. Invoices are payable within agreed terms of sale. Henry Schein and its affiliates reserve the right to offer a dual pricing system for payment by check, Henry Schein-approved ACH and/or Henry Schein-approved EFT vs. any other payment method and, to the extent permitted by law, reserves the right to surcharge for payment by credit card and/or other methods incurring a fee to Henry Schein.

Open Accounts Receivable:

All unpaid accounts receivable past due are subject to a 1.5% finance charge.

DELIVERY TERMS:

Delivery dates are not guaranteed. Purchase orders are subject to product availability. Customers' sole and exclusive remedy for failure to deliver shall be refund of monies paid for the products at issue. Added to Delivery Terms section in first column. Unless otherwise agreed, freight terms are FOB Shipper's Dock ("Ex Works" outside North America). Except as noted below, title passes at the time the shipment is loaded at the shipper's dock.

California: For all shipments of goods to customers located within California, title will pass upon receipt of goods by California customers.

Continental U.S.: All orders will be subject to a handling charge. This charge includes freight, except for additional carrier charges related to special delivery services and hazardous material shipments. Special orders are subject to additional freight charges.

Alaska, Hawaii & Pacific Protectorates: Standard shipping methods provide direct, reduced cost, expedited air delivery service to all accounts in Alaska and Hawaii. Customers in the Pacific Protectorates are offered direct surface transport, or postal services for reliable delivery. No additional surcharges apply, except when special services are requested. Low-level hazardous items (dangerous goods in accepted quantities and Consumer Commodity ID 8000) are now available via UPS 2nd-day air.

Guam, Puerto Rico, U.S. Trust Territories & Virgin Islands: • All orders will be subject to a handling charge. This charge includes freight through the United States Postal Service (USPS). • Special delivery orders and hazardous material shipments can be shipped via United Parcel Service (UPS) for an additional charge. No minimum order amount or weight applies. Speak to your International Representative for details.

Outside U.S. (50 States): If your order is being shipped outside the U.S. (50 states), please refer to the International Terms & Conditions at henryschein.com. Unless otherwise agreed, freight terms are FOB Shipper's Dock ("Ex Works" outside North America). Title passes at the time the shipment is loaded at the shipper's dock. Customer is responsible for compliance with any applicable import requirements.

RX PRODUCTS & CONTROLLED SUBSTANCES:

• Regulations require us to limit the sale of Rx and controlled substances only to registered, licensed healthcare professionals. If you are a new customer or have recently moved, please furnish us with a copy of your updated state and federal registrations verifying your shipping address. Please note that all orders for controlled substances are subject to a due diligence review process. Schedule II controlled substances can be ordered electronically or by mail. For more information on our Controlled Substance Ordering System please visit www.henryschein.com/e222; if you prefer to continue using Federal 222 Forms to order Schedule II controlled substances, mail the form to: Henry Schein, Inc. 5315 West 74th Street • Indianapolis, IN 46268

• Henry Schein restricts the sale or other transfer of medications to prisons/correctional facilities for use in lethal injections, based on our manufacturer agreements. The goods Henry Schein sells are intended to be used for their label-approved purposes or applicable standards of care, which do not include human lethal injection.

REGULATORY REQUIREMENTS:

Local regulatory requirements may apply to use or installation of certain products. Be sure to understand and comply with any such requirements prior to purchase, use, or installation of products.

THE DRUG SUPPLY CHAIN SECURITY ACT (DSCSA):

(MN, DM, WH, M2) The Drug Supply Chain Security Act (DSCSA) information related to prescription drug products is available on our website www.HenrySchein.com/pedigree. If you have any problems accessing our website or would like to receive a copy of DSCSA documentation via fax, mail or email, please contact our customer service department.

RETURNS:

WE CANNOT ACCEPT ANY RETURNS WITHOUT PRIOR AUTHORIZATION.

To arrange for a return, simply call our Customer Service department or contact your Sales Consultant.

The following conditions must be complied with:

• All returns must be accompanied by a copy of your invoice and a reason for the return • Merchandise must be returned in its original container, unmarked, and properly packaged • Returned products must have been purchased within the previous thirty (30) days. Any returns past thirty (30) days are subject to a restocking

DISCOUNTS, REBATES AND DISCLOSURES: Invoice or statement prices may reflect or be subjected to a bundled discount or rebate pursuant to purchase offer, promotion or discount program.

You must fully and accurately report to Medicare, Medicaid, Tricare and/or any other federal or State program, upon request by such program, the discounted price(s) or net price(s) for each invoiced item, after giving effect to any applicable discounts or rebates, which price(s) may differ from the extended prices set forth on your invoice. Accordingly, you should retain your invoice and all relevant information for your records. It is your responsibility to review any agreements or other documents, including offers or promotions, applicable to the invoiced products/prices to determine if your purchase(s) are subject to a bundled discount or rebate. Any such discounts must be calculated pursuant to the terms of the applicable purchase offer, promotion or discount program. Participation in a promotional discount program is only permissible in accordance with discount program rules. By participation in such program, you agree that, to your knowledge, your practice complies with the discount program requirements.

fee • Shortages or errors in shipments must be reported within seven (7) days of invoice date to issue credit (if applicable) • Shipping charges will apply on all returns

Exceptions:

The following special, customized, or government-regulated items are not returnable:

• Opened handpieces, small equipment, and custom-ordered equipment • Special order items (products that we do not ordinarily stock) • Personalized and imprinted items • Opened computer hardware and software • Hazardous/flammable materials • Expired products • Items that cannot be returned to the manufacturer • Any item marked non-returnable • Items required to be shipped and stored frozen • Any drop-shipped products

Prescription Drug Returns:

Please note that, in order to comply with Federal and State traceability requirements, prescription drugs may be returned providing that the following key elements are met:

- 1) Returns of prescription drugs will only be accepted if Henry Schein is notified within 30 calendar days of shipment date and valid return authorization is issued by Henry Schein.
- 2) The Prescription Drug Marketing Act requires any customer returning prescription drugs to complete and return a Prescription Drug Return Authorization form. Federal law requires that the healthcare entity returning prescription drugs document that the product was kept under proper storage and handling conditions while in their possession and during the return of the product. To get a copy of the form and proper return authorization, please contact Customer Service.
- 3) In addition, traceability regulations require that the healthcare entity returning prescription drugs certifies that the product being returned is the same exact product purchased from HSI.
- 4) Henry Schein will not issue credit for any returned prescription drugs which return was not authorized as provided herein, have been tampered with or where the labeling has been altered in any way.

INSTITUTIONAL, GOVERNMENT & CORPORATE ACCOUNTS:

Terms of Sale follow the same guidelines unless denoted differently in a contract. When applicable, freight is calculated by the weight. There is no minimum order; however, there is a nominal handling charge on orders under \$125. Some offers and promotions outlined in this catalog may not apply. Requests for bids and proposals may be sent to:

Henry Schein, Inc. • Bid Department—E255 • 135 Duryea Rd. • Melville, NY 11747

INTERNATIONAL ORDERS:

We proudly serve healthcare professionals, governments, and dealers throughout the world. To place orders or for inquiries on export terms and conditions please contact the International Department (USA) at: export@henryschein.com.

ARBITRATION: All Claims related to or arising under or relating to this Agreement (except for Claims for non-payment of amounts due to Henry Schein for the sale of products) are to be exclusively and finally determined by binding arbitration in the state of New York, or another location mutually agreeable to the parties. Any and all Claims must be arbitrated on an individual basis, and there shall be no right or authority for any Claims or disputes to be arbitrated on a class action or collective basis. For avoidance of doubt, each party irrevocably waives any right to: (i) have any Claim resolved in connection with any class action or collective action, or (ii) recover any damages or relief directly or indirectly as part of any class action or collective action. The arbitration shall be conducted on a confidential basis pursuant to the Commercial Arbitration Rules of the American Arbitration Association, or if applicable, under its Procedures for Large, Complex Commercial Disputes. Any decision or award as a result of any such arbitration proceeding shall be in writing and shall provide an explanation for all conclusions of law and fact and shall include the assessment of costs, expenses and reasonable attorneys' fees. Any such arbitration shall be conducted by an arbitrator experienced in the disputed subject matter and shall include a written record of the arbitration hearing. The parties reserve the right to object to any individual who shall be or has been at any time employed by or affiliated with a competing organization or entity. An award of arbitration may be confirmed in a court of competent jurisdiction. To the extent that any Claim or dispute is determined to not be subject to arbitration, all other Claims or disputes that would otherwise be subject to arbitration must be arbitrated. As used in this Agreement, "Claims" shall mean any and all liabilities, disputes and expenses whatsoever including, without limitation, claims, adversary proceedings (whether before a court, administrative agency or any other tribunal), damages (whether compensatory, multiple, exemplary or punitive), judgments, awards, penalties, settlements, investigations, costs, responses to subpoenas or other governmental directives and reasonable attorneys' fees and disbursements with respect to any claims that may be sustained, suffered or incurred by a Party hereto.

WARRANTIES: Henry Schein will pass through to the customer, at the time of sale, any transferable product warranties, indemnities and remedies provided to Henry Schein by the applicable manufacturer.

EXCEPT AS OTHERWISE PROVIDED HEREIN, TO THE EXTENT PERMITTED BY LAW, HENRY SCHEIN PROVIDES NO WARRANTIES OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION ANY WARRANTY OF MERCHANTABILITY, FITNESS FOR PARTICULAR PURPOSE OR NON-INFRINGEMENT, AND THE CUSTOMER SHALL LOOK TO THE MANUFACTURER OF THE PRODUCT FOR ANY WARRANTY THEREON.

LIMITATION OF LIABILITY: THE CUSTOMER AGREES TO LOOK SOLELY TO THE MANUFACTURER OF THE PRODUCT FOR ANY CLAIM ARISING DUE TO LOSS, INJURY, DAMAGE OR DEATH RELATED TO THE USE OR SALE OF PRODUCTS. HENRY SCHEIN SHALL NOT BE LIABLE FOR INDIRECT, INCIDENTAL, PUNITIVE, SPECIAL OR CONSEQUENTIAL DAMAGES INCLUDING, BUT NOT LIMITED TO, LOST PROFITS AND LOSS OF GOODWILL, ARISING FROM OR RELATING TO ANY BREACH TO THIS AGREEMENT (OR ANY DUTY OF COMMON LAW, AND WHETHER OR NOT OCCASIONED BY THE NEGLIGENCE OF HENRY SCHEIN OR ITS AFFILIATES), REGARDLESS OF ANY NOTICE OF THE POSSIBILITY OF SUCH DAMAGES.

Terms of Sale may change without notice. For current information, go to www.henryschein.com/legalterms.

Den-8/16/21-7

CRITERION[®]

COMFORT IN PROTECTION[™]

23PL1013

Visit us on the Web at:
HenrySchein.com/Criterion

© 2023 Henry Schein, Inc.
No copying without permission.
Not responsible for typographical errors.

by  HENRY SCHEIN[®]

 HENRY SCHEIN[®]
MEDICAL

