

CATEGORY SPOTLIGHT

HENRY SCHEIN DENTAL

ACCLEAN Ultrasonic Scaler: Power Up & Go

Space-saving ultrasonic scaler features dual-frequency technology, an easy-to-use design, and a cruise-control option for continuous power



eBOOK

Scan the QR code to learn how the ACCLEAN 3-stage system can organize workflows for optimal patient care.



Seen one, seen them all? Not necessarily. In the case of ultrasonic scalers, sometimes subtle differences can make a world of difference.

Take, for instance, unit size and ease of use. The ACCLEAN Ultrasonic Scaler from Henry Schein Dental has both features covered, and they're 2 of the top reasons why Tara Cabral, RDH, BS, likes using the unit in her Cape Coral, FL, practice.

Smaller Footprint

Weighing in at just 1.3 pounds, the ACCLEAN Ultrasonic Scaler features a conveniently small footprint of just under 6 inches long and about 3 inches high, which means the scaler is a perfect fit for space-challenged countertops. "The nice size means that less counterspace is taken up and the unit is easier to move," Cabral shared.

The scaler features an easily detachable, autoclavable handpiece sheath and user-friendly controls, features that are helpful to hygienists who find themselves using it for longer periods of time due to increased stain and calculus accumulation on patients who've missed appointments during the pandemic.

"I find the scaler to be more comfortable and easier to use than other ultrasonic scalers," she explained.



Tara Cabral,
RDH

"I really like the size of the ACCLEAN ultrasonic scaler, and the adjustability of the unit is very convenient."



A Host of Features

The ACCLEAN Ultrasonic Scaler features dual-frequency technology and accepts both 25K and 30K inserts, automatically detecting the inserts and then quickly integrating them into existing equipment. The unit's Auto Gain Control technology maintains power throughout procedures. Additionally, a hands-free, on-demand boost mode can be activated by a foot pedal, a design feature that can provide additional force if it's needed when removing challenging stains and calculus.

Cruise Control Your Cleanings

In addition to user-friendly controls, ACCLEAN features a cruise-control option, another ergonomic benefit that allows hygienists to "select power and go" while not needing to apply continuous pressure on the foot pedal. The unit also features a 1-button purge mode that provides automatic line cleaning.

"I really like the size of the ACCLEAN Ultrasonic Scaler, and the adjustability of the unit is very convenient," Cabral concluded.

FOR MORE INFORMATION:

800.372.4346

henryschein.com/
accleanpreventives