

VOCO

Meron Plus QM

PRODUCT REVIEW SNAPSHOT

CRITERIA BASED ON AVERAGE SCORE (OUT OF 5)



Delivery system **4.9**

Consistency of the cement **4.7**

Time savings/number of steps needed **4.6**

Radiopacity **4.7**

Lack of postop sensitivity **4.7**

SECTION A
AVERAGE **4.7**

SECTION B
OVERALL SATISFACTION **4.8**

takeaways

- Paste/paste formulation is easy to mix and dispense with QuickMix syringe
- Self-adhesive formula allows for quick, secure bonding without primer or conditioner
- 2-minute working time is sufficient for single-tooth or multi-unit restorations
- Excess material is easy to remove during gel phase
- Shear bond strength to dentin is 11.3 MPa and bond strength to zirconium dioxide is 11.8 MPa

[CLICK HERE TO READ THE FULL REVIEW](#)

How was this evaluated

FINAL SCORE
AVERAGE OF SECTIONS A AND B

4.8

