

DIRECTA DENTAL

Hygovac Bio by Orsing



PRODUCT REVIEW SNAPSHOT CRITERIA BASED ON AVERAGE SCORE (OUT OF 5)

Environmental friendliness **4.8**

Provides continuous suction **4.4**

Prevents drawing in/
suctioning of tissue **3.1**

Variety of sizes—short and long **3.7**

Ergonomics—design allows for
comfortable grip and posture **3.7**

Patient comfort **3.9**

Ability to retract tissue **3.7**

SECTION A
AVERAGE **3.9**

SECTION B
OVERALL SATISFACTION **4.0**

takeaways

- Made of 100% biobased polyethylene from sugarcane
- Comes in 2 lengths, with the shorter version measuring 95 mm
- Shorter length enhances ergonomics for the clinician
- Same diameter and suction capacity as other Hygovac products
- Reduces the use of fossil resources and carbon footprint

FINAL SCORE
AVERAGE OF SECTIONS A AND B

4.0