

GC AMERICA

EXAFAST and EXAMIX



PRODUCT REVIEW SNAPSHOT CRITERIA BASED ON AVERAGE SCORE (OUT OF 5)

Working time **4.3**

Setting time **4.6**

Viscosities **4.1**

Tear strength **4.3**

Accuracy at the margin **4.6**

Smell and taste **4.3**

Fit of the final restoration **4.7**

SECTION A
AVERAGE **4.4**

SECTION B
OVERALL SATISFACTION **4.5**

takeaways

- EXA family includes 5 VPS impression materials to meet any clinical need
- EXAFAST has a 60-second working time and sets in 90 seconds
- EXAMIX presents a working time of 3 minutes and a setting time of 4 minutes
- Predictable, precise, and durable impressions
- Excellent elasticity and tear strength

FINAL SCORE
AVERAGE OF SECTIONS A AND B

4.5