

Rapid and Effective Sterilization

Up to 28,000 Cycles

Maxima™ Advanced Chamber Sterilizer

Boost efficiency with advanced sterilization technology. Built for 28,000 cycles, this large, 11" chamber sterilizer offers rapid drying, dynamic pulsed air removal, and customizable programs.



**CLOSED-DOOR
VACUUM-ASSISTED
DRYING**



**SINGLE USE
WATER MANAGEMENT**



**AUTOMATED
DIGITAL CYCLE
DATA RECORDING**



**2-YEAR
WARRANTY**



maxima®

by  HENRY SCHEIN®

Performance Handpieces & Small Equipment

Explore the New Maxima: HenryScheinDental.com/Maxima

Maxima™ Advanced Chamber Sterilizer

Highlights

- **Designed to last up to 28,000 cycles*** for long-term durability and performance
- **Closed-door vacuum-assisted drying** ensures sterilization meets stringent hygiene standards
- **Dynamic pulsed air removal** effectively removes air from the chamber, providing thorough sterilization even with porous materials and instruments
- **11" large capacity chamber** with four trays to accommodate up to 18lbs
- **Flexible cycle times for high efficiency:** 18 lbs load (unwrapped) sterilized and dried within 30 minutes! Customizable programming options allows for reduced overall cycle time
- **Single-use distilled water management** supports stringent hygiene standards using fresh water every cycle
- **No periodic cleaning cycles required**, reducing costs and downtime
- **Automated digital cycle** records data via USB stick for compliance and quality assurance
- **Automatic water filling and draining** options for convenience
- **User-friendly, color touchscreen**
- **2-Year Warranty**

*Following the instructions for use and general maintenance checks outlined in the IFU



Maxima™ Advanced Chamber Sterilizer

1 Sterilizer Unit, 1 Universal Chamber Rack, 4 Large Aluminum Trays, 1 Tray Holder, 1 Power Cord, 1 Drain Tube, 1 USB Drive, 1 Quick Start Guide

External Dimensions: 19.3" x 17.8" x 24.3" | Tray Dimensions: 0.8" x 8.5" x 14.9"
Chamber Dimensions: 11" x 18" | Reservoir Capacity: 1.7 Gal. (6.4 L)

Maxima™ Advanced Chamber Sterilizer 110V
Maxima™ Advanced Chamber Sterilizer 220V

(570-3795)
(570-3796)

ACCESSORIES

Maxima Advanced Chamber Sterilizer WiFi Dongle Key	(570-4333)
Maxima Advanced Chamber Sterilizer Autofill kit w/ Valve	(570-4340)
Maxima Advanced Chamber Sterilizer Autodrain Kit	(570-4337)
Maxima Advanced Chamber Sterilizer Autofill kit w/ Pump (110V)	(570-4341)
Maxima Advanced Chamber Sterilizer Autofill kit w/ Pump (220V)	(570-4338)

PREVENTATIVE MAINTENANCE (Buy every 1200 cycles)

Maxima Advanced Chamber Sterilizer Door Gasket	(570-4335)
Maxima Advanced HEPA Filter	(570-3820)

REPLACEMENT PARTS

Maxima Advanced Chamber Sterilizer Standard Chamber Rack	(570-4334)
Maxima Advanced Chamber Sterilizer Large Aluminum Tray	(570-4328)
Maxima Advanced Chamber Sterilizer Tray Holder	(570-4329)
Maxima Advanced Chamber Sterilizer Emergency Door Opening Tool	(570-4330)
Maxima Advanced Chamber Sterilizer Drain Tube	(570-4331)
Maxima Advanced Chamber Sterilizer Mains Cable 125V	(570-4339)
Maxima Advanced Chamber Sterilizer Mains Cable 250V	(570-4336)
Maxima Advanced Chamber Sterilizer USB Pen Drive	(570-4332)

