

MICROCOPY

NeoShine Zirconia Polishing System



PRODUCT REVIEW SNAPSHOT CRITERIA BASED ON AVERAGE SCORE (OUT OF 5)

Smoothness of zirconia surface after using step 1

4.2

Glossiness of zirconia after using step 2

4.6

Ease of access for all areas of the tooth with the mini cup and point

4.6

Durability of the polishers when polishing zirconia crown

3.5

Product color ID system—matching polisher to substrate and polishing step

4.2

Clarity of the packet label information

4.6

SECTION A
AVERAGE

4.3

SECTION B
OVERALL SATISFACTION

4.3

takeaways

- Mini cups and points used for occlusal adjustments of an already sintered or lab-fabricated zirconia restoration
- 2-step zirconia polishing system removes striations and produces a high-gloss, shiny surface
- Quickly polishes crowns to a high shine
- Single-use zirconia polishers prevent cross contamination
- Each cup or point is sterile, individually packaged, and clearly labeled

FINAL SCORE
AVERAGE OF SECTIONS A AND B

4.3