

PREMIER DENTAL
BioCoat



PRODUCT REVIEW SNAPSHOT CRITERIA BASED ON AVERAGE SCORE (OUT OF 5)

Flowability/viscosity **4.4**

Bubble-free consistency **4.4**

Ability to achieve smooth margins **4.6**

Convenience of disposable applicator tip **4.4**

Working time/curing time **4.6**

SECTION A
AVERAGE **4.5**

SECTION B
OVERALL SATISFACTION **4.7**

takeaways

- Dispensed in a 1.2 mL syringe and easily applied with disposable applicator tips
- Excellent flowability and bubble-free consistency for smooth and consistent application
- Thixotropic properties allow it to penetrate and adapt into small pits and fissures
- SmartCap Technology allows for long-term release of bioactive ions of fluoride, calcium, and phosphate
- Strength comparable to flowable composites for excellent retention

FINAL SCORE
AVERAGE OF SECTIONS A AND B

4.6