

SINGLE-USE AEROSOL EVACUATOR

A simple solution for aerosol control to protect your patients and staff.



✓ Disposable, cost-effective device for aerosol reduction in the operator.

✓ Use during hygiene, restorative, oral surgery, and other procedures.

✓ Hands-free, ease-of-use for both the hygienist and dentist.

✓ Use in conjunction with saliva ejector.



SINGLE-USE AEROSOL EVACUATOR

NEW


The Single-Use Aerosol Evacuator is a disposable device that attaches to any standard HVE connection to provide an effective level of aerosol reduction. The extraoral product placement allows for hands-free suction of aerosols at the exit of the oral cavity, delivering comfort for the patient and ease-of-use for the hygienist or dentist. This low-profile, single-piece device provides minimal interference at the procedure location, especially when compared with intraoral options.

Advantages

- Proprietary patented design.
- Universal size with malleable “wings” that comfortably adjust to patient’s mouth.
- Affordable aerosol control.
- Low-profile provides unobstructed view of the procedure area.

Specs

- Material: High-density polyethylene (HDPE)
- Color: White
- Quantity: 50/pack
- Dimensions: 4.25” x 3” (L x W)
- Use with saliva ejector



1. Internal Wings
2. External Wings
3. Lower Teeth Bumper
4. Suction Inlet
5. Valve Stem

Ordering Information

Single-Use Aerosol Evacuator
Item No.: (570-2441-WSH)
Catalog Price: \$59.99/50pk

Maxima PRO 2AR™ High-Speed
Air-Driven Handpieces
Specifically engineered for aerosol
reduction and anti-retraction.



Maxima PRO 2LAR™ Optic (570-2094)
Maxima PRO 2SAR™ Non-Optic (570-2095)

Save 10% with promo code WSH.*

*Offer valid through May 31, 2021.

For more information, please see product Instructions for Use.

Available through

HENRY SCHEIN®
DENTAL

www.henryscheindental.com

© 2021 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors. To guarantee special pricing/offers, where applicable, please include the three-digit promo code when ordering. Henry Schein Dental reserves the right to discontinue this promotion at any time [including, without limitation, if supply or demand for the promotional products is affected due to extraordinary events]. Note: Participating in a promotional discount program (e.g., points, discount redemptions or other special awards) is only permissible in accordance with discount program rules. By participation in such program, you agree that, to your knowledge, your practice complies with the program requirements. For HS promotions: This promotion offers a bundled discount. The cost of any “buy/get” (or “no charge”) items is a discount applied to your other purchases on the same invoice. The reportable value (unit price) of the “no charge” item will be identified on the invoice. To calculate the applicable discount on the other items on the same invoice, divide the total value (unit price(s)) of any no-charge item(s) by the total value (extension price(s)) of all other items purchased on the same invoice, and apply the resulting percentage equally to each item to determine net prices. For manufacturer promotions: The prices stated in the invoice or statement may reflect a discount or be subject to a rebate. You must fully and accurately report this stated discount price, or if applicable, any net pricing, after giving effect to any rebates, to Medicare, Medicaid, Tricare and any other federal or state program upon request by any such program. Accordingly, you should retain these records. It is your responsibility to review any agreements or other documents applicable to these prices to determine if they are subject to a rebate.