

VOCO

Meron Plus QM



PRODUCT REVIEW SNAPSHOT CRITERIA BASED ON AVERAGE SCORE (OUT OF 5)

Delivery system

4.9

Consistency of the cement

4.7

Time savings/number
of steps needed

4.6

Radiopacity

4.7

Lack of postop sensitivity

4.7

SECTION A
AVERAGE

4.7

SECTION B
OVERALL SATISFACTION

4.8

takeaways

- Paste/paste formulation is easy to mix and dispense with QuickMix syringe
- Self-adhesive formula allows for quick, secure bonding without primer or conditioner
- 2-minute working time is sufficient for single-tooth or multi-unit restorations
- Excess material is easy to remove during gel phase
- Shear bond strength to dentin is 11.3 MPa and bond strength to zirconium dioxide is 11.8 MPa

FINAL SCORE
AVERAGE OF SECTIONS A AND B

4.8