



LEARN MORE

To view the Henry Schein Brand Products catalog and discover additional product information, visit henryschein.com/dental.com/hsbrand.

HENRY SCHEIN BRAND

Natural Elegance Restorative Line

Easy-to-use composites, bulk fills, and bonding agents designed for optimal restoration results at a great value for your practice

From office managers to hygienists to lab techs, having a solid dental team to rely on is essential to building a successful practice. So is having a go-to line of products designed to work in tandem during the most common procedures. For that, you can rely on Henry Schein's Natural Elegance restorative products—an artful blend of composites, bulk fills, and bonding agents designed to achieve optimal results in a variety of restorative scenarios.

Meet the Natural Elegance Line Up

The Natural Elegance Premium Flow Nano Composite is a reliable solution when restoring in smaller, hard-to reach areas. Its optimized flowability makes it ideal as a liner/base and it is formulated with nanoparticles to provide strength, polishability, and wear resistance. Additionally, it is fluoride releasing to help prevent caries and offers good radiopacity.

With better mechanical properties due to improved polymerization, the Natural Elegance Premium Composite is a universal material suited for both posterior and anterior fillings. It features nonsticky handling and contains a unique



combination of photo initiators that provide total shade predictability. Its spheroidal zircon silicate filler offers enhanced brightness and wear resistance.

In the bulk fill department, Natural Elegance Flowable Bulk Fill is self-leveling and ideal for small, non-occlusal stress-bearing Class I restorations or as a base/liner in Class I and IIs. Maxima Natural Elegance Packable Bulk Fill offers a new multi-hybrid filler and innovative photo initiator technology. Available in syringes and capsules, it is best suited for core buildups and Class I and II posterior fillings. Both materials are highly radiopaque and offer a 4-mm depth of cure with minimal shrinkage.

Last, but certainly not least, Henry Schein Natural Elegance total-etch, self-etch, and universal bonding agents feature drop control dispensing technology for less waste. Natural Elegance Universal Bond has a very low film thickness and offers excellent marginal sealing to minimize postoperative sensitivity. It can be used to prime zirconia, alumina, metals, and porcelain. Try the full line of Natural Elegance restorative products for optimal results at a value your practice will love.

3 WAYS TO BOND

Henry Schein Natural Elegance bonding agents offer a solution for every technique:

Total Etch Bond: Good shear bond strength—Dentin: 22.7 MPA; Enamel: 20.9 MPA

Self Etch Bond: Excellent shear bond strength—Dentin: 25.9 MPA; Enamel: 19.5 MPA

Universal Bond: Ultimate shear bond strength—Dentin: 27.5 MPA; Enamel: 22.2 MPA

FOR MORE INFORMATION:

800.372.4346

henryschein.com/dental.com/NErestoratives