



Flexible



3M™ Medipore™ H Soft Cloth Surgical Tape

Made to move with  
your patients.



## The choice for stretchable, conformable securement.

Medipore H Soft Cloth Surgical Tape features a soft, conformable backing with bi-directional stretch to accommodate swelling and movement, while minimizing the risk for skin injury. It's now available in individually packaged single-use rolls to help reduce the risk of cross-contamination.



### Conformability

Bi-directional stretch accommodates swelling, distention and movement.



### Breathability

Highly breathable, soft cloth tape.



### Easy to use

Easy tear at perforations with no scissors necessary.



### Helps reduce risk

Individually packaged single-use length rolls can help reduce cross-contamination risk by 99% compared to unpackaged tape.\*

Helping you focus on  
patient care and safety



**74%**  
of tape rolls

In one study, 74% of partially used tape rolls from various sites within the hospital had some bacterial growth.<sup>1</sup>

\*3M's packaged single-use length rolls can help reduce exposure to bacteria and spores. In vitro testing on contaminated hard dry surfaces showed as much as 99% protection from pathogens such as MRSA, CRE, VRE and C. difficile when compared to unpackaged rolls of tape.<sup>2</sup> Individual roll packaging helps minimize exposure of the tape to environmental contaminants from surfaces and hands. Practice good infection prevention technique per facility protocols and published guidelines, including environmental cleaning and proper hand hygiene before and after opening package.

# 3M™ Medipore™ H Soft Cloth Surgical Tape

## Clinical Application: Flexible Securement

- Dressings and added pressure
- Secondary securement for:
  - ▶ IV lines and tubing
  - ▶ Central venous catheters
  - ▶ Chest tubes
  - ▶ Surgical drain tubes
- When swelling or movement is anticipated
- When repeated taping is needed

**Ideal for patients who need to maintain movement or are experiencing swelling or distension.**

**Tape material:** Perforated nonwoven soft cloth



**Stretches to accommodate swelling and movement**



**Conforms to challenging contours**

## Ordering information

Catalog Number	Size	Individually packaged single-use length rolls	Rolls	Case	HCPCS Code
2860S-1	1 in. x 2 yd. (2,5 cm x 1,8 m)	✓	72 per bag	1 bag	A4452
2860S-2	2 in. x 2 yd. (5 cm x 1,8 m)	✓	48 per bag	1 bag	A4452
2860S-4	4 in. x 2 yd. (10,1 cm x 1,8 m)	✓	24 per bag	1 bag	A4452
2860S-6	6 in. x 2 yd. (15,2 cm x 1,8 m)	✓	16 per bag	1 bag	A4452
2861	1 in. x 10 yd. (2,5 cm x 9,1 m)	—	24 per box	1 box	A4452
2862	2 in. x 10 yd. (5 cm x 9,1 m)	—	12 per box	1 box	A4452
2863	3 in. x 10 yd. (7,6 cm x 9,1 m)	—	12 per box	1 box	A4452
2864	4 in. x 10 yd. (10,1 cm x 9,1 m)	—	12 per box	1 box	A4452
2866	6 in. x 10 yd. (15,2 cm x 9,1 m)	—	12 per box	1 box	A4452
2868	8 in. x 10 yd. (20,3 cm x 9,1 m)	—	6 per box	1 box	A4452



For more information about 3M™ Medipore™ H Soft Cloth Surgical Tape and other 3M Medical Securement Solutions, visit [3M.com/MedicalSecurement](https://www.3m.com/MedicalSecurement).

**Medical Solutions Division**  
3M Health Care  
2510 Conway Ave.  
St. Paul, MN 55144 USA

Phone 1-800-228-3957  
Web [3M.com/medical](https://www.3m.com/medical)

<sup>1</sup>Lavelle BE. Reducing the risk of skin trauma related to medical adhesives, Managing Infection Control, June 2004.

<sup>2</sup>3M data on file. TECH REPORT-05-742819.

3M and Medipore are marks and/or registered marks of 3M. Unauthorized use prohibited.  
Please recycle. Printed in USA. © 3M 2021. All rights reserved. 70-2011-7899-6