

Rely on Us™



To Order: 1-800-521-1472 8am-6pm (et)
To Fax: 1-800-372-5702 24 Hours

www.henryschein.com/dialysis

SUCTION UNITS



Madavac 172BS Aspirator

- Bacteria filter prevents contaminant damage
- Vibration-free, easy-to-read gauge
- Quick-connect tubing for all connections
- Draws up to 22" Hg (560mm Hg) with a flow rate of 28 lpm
- Large 800cc collection bottle with overflow cutoff
- Carrying handle for easy portability
- Quiet yet powerful 115-V, 60-Hz, 2.9-A unit

(923-3316)ea



Schuco-Vac® 310 Aspirator

Regular assembly has been replaced with a bacteria filter, eliminating the cap-and-float assembly. Filter prevents fluids and other contaminants from damaging pump. This durable unit continues to meet the need for controlled vacuum regulation from 0" up to 22" Hg with a flow rate of 28 lpm.

Standard features include:

- Bacteria filters (3 included) prevent contaminant damage
- Vibration-free, easy-to-read regulator gauge
- Wire basket accommodates a variety of canisters
- Quick-connect tubing for all connections
- Hydrophobic filter eliminates head flooding, the most common cause of pump damage

Aspirator, with 1200-cc Canister (115 V/60 Hz)ea
(268-7332)ea

Replacement Canister, Plastic, 1200 cc (141-2411)4/case



Laerdal Suction Units

Features industry-first device diagnostics, a no-tools replaceable battery, built-in charger, and an integrated vacuum regulator and power-on knob. All this in a rugged bump-proof and splash-proof design that is very quiet.

- Oral, nasal, and tracheal suction applications for infants through adults
- Multisource power from internal battery and AC or DC sources
- Powerful—meets ISO 10079-1, the international standard for medical suction equipment
- Adjusts from 80-mm to 550-mm Hg
- Spacious and durable soft-pack carrying case holds extra supplies
- Operate unit right from the carrying case: convenient for storage

With Disposable Canister (827-2329)ea

Replacement Canister #484410 (116-9264)ea



7305 Series Suction Unit

Easy-to-use and compact design with vacuum adjustments from 80mm to 550mm Hg. The 7305 series meets the RTCA/DO160D Aircraft Standard as well as the ISO-10079 Suction Standard. Backed by a 2-year warranty.

#7305P-D, Unit with Battery 6.3 lb, 8"H x 9"W x 7"L (401-6284)ea



Hydrophobic Canister

#424410, 800cc, Red Top (645-5651)ea

OXYGEN CONCENTRATORS & THERMOMETRY



5-Liter Compact Concentrator

34% more compact than traditional concentrators. Exclusive turn-down technology reduces power consumption by 15% or more at or below 2.5 lpm. Includes oxygen sensing device; no maintenance required for 2 years. Dimensions: 24½"H x 13½"W x 12"L Weight: 36 lb
#525DSea
(110-8052)ea



NewLife Elite Oxygen Concentrator

The world's most in-demand unit for patients who require supplemental oxygen. Its ease of use, reliability, and classic styling make it the popular choice for oxygen therapy in more than 100 countries.

- 5-LPM performance with easy-to-access patient controls
- Units available with optional oxygen monitor or optional air outlet for medicated nebulizer treatments
- Low service requirements make the unit virtually maintenance-free

(919-3556)ea

THERMOMETRY



Braun Thermoscan® Pro 4000 Ear Thermometer

- Improved accuracy and repeatability through a preheated probe tip and sensor
- Automatic probe cover eject button for quick and easy disposal of used probe covers
- Exclusive ExacTemp technology provides an active user feedback light, guiding proper probe positioning and improved accuracy
- Preheated probe tip minimizes cool-down effect experienced when taking repeat temperatures, improving both accuracy and repeatability
- Probe cover detection system improves accuracy by ensuring precise infrared heat transfer
- Innovative security features reduce the occurrence of theft in the hospital environment (*only available with Base Station*)



Thermometer
(566-3430)ea
Probe Cover
#05075-005
(678-2442)200/box



Genius™ 2 Tympanic Thermometer

- "Peak Select" technology records and analyzes up to 100 temperature measurements and selects the most accurate reading in 1 to 2 seconds
 - Redesigned tapered probe tip facilitates proper ear-canal alignment, minimizing technique-related errors and increasing patient comfort
 - Rugged tip-based sensor construction improves durability and accuracy, minimizing the need for recalibration
 - No-touch probe cover application and ejection reduce potential cross contamination
 - Programmable display shows oral, rectal, axillary, and core or an absolute temperature reading in ear mode
 - Temperature range: 91.4°F to 107.6°F
- #303000, Thermometer with Base
(106-6433)ea
#303030, Probe Covers
(194-3810)96/box



Non-Contact Infrared Thermometer

Specially designed to take a human body temperature while never touching the patient. Uses an infrared detection system and displays temperature on a backlit screen.

- For use on adult and pediatric patients
- Measures body temp, ambient temp and surface temp
- Audible alarm if temp is too high
- 32-reading memory
- Celsius or Fahrenheit reading
- Instantaneous reading
- 2 AA batteries (included)

(700-0646)ea



RediScan™ Infrared Thermometer

Professional accuracy and wide-range readout: 50°F–122°F (person) and 32°F–221°F (object/liquid). Features the Nite-Glo display and automatic shut-off, 1-second readout, °F/°C, peak-temperature tone, 100.4°F fever alarm, 10-memory recall, and 5-year limited warranty.

(114-2518)ea

Contains: 1 thermometer, 1 drawstring storage pouch & 2 AAA batteries.



Adtemp™ 427

Noninvasive forehead temperature measurement using patented conductive technology delivers clinically accurate audible readings in about 6 seconds. Reduces the risk of cross contamination and avoids the need for costly probe covers. Dual scale (°F/°C), temperature range 95°F to 107.6°F ± .2°F (35°C to 42°C ± .1°C).

(114-7121)ea



Filac™ 3000 EZ Electronic Thermometer

Fastest, Accurate, Predictive-type

- Accurate reading in under 4 seconds or advanced accuracy reading in under 6 seconds
- Compact, robust, user-friendly design will stand upright on a flat surface
- Intuitive icon-based interface system with backlit LCD readout
- Unique progress bar indicates status of temperature measurement
- Last-temperature-recall function
- Push-button probe cover ejection
- Patented isolation chamber design and durable probe covers minimize cross contamination
- Operates on 4 AA batteries for 6000 measurements
- Equipped with self-calibrating electronics and auto on/off function
- 1-year full manufacturer warranty

Oral Thermometer

(110-3415)ea

Probe Covers

(110-3414)2000/box



SureTemp® Plus 690 Electronic Thermometer

SureTemp® Plus Thermometer Probe Covers

Combines the speed of ear thermometry with the accuracy and reliability of traditional oral thermometry. Takes fast, 4–6-second oral temperatures, 10-second pediatric axillary and rectal temperatures, and 15-second adult axillary temperatures.

- Accurate to within ± 0.2°F (0.1°C)
- Interchangeable, removable probe wells to reduce the risk of cross contamination
- Large LCD that can display temperatures in Fahrenheit or Celsius
- Last-temperature recall
- Mode selection and recall buttons
- Four different temperature modes
- Convenient storage for 25 probe covers
- Takes about 6000 readings on 3 AA batteries (included)

Thermometer, with Oral Probe, Well, and 4' Cord (203-3856).....ea

Thermometer, with Wall Mount, Oral Probe, Well, and 4' Cord (178-8374).....ea

Thermometer with Wall Mount, Oral Probe and 9' Cord (454-2689).....ea

Wall Holder Mount for the Thermometer (206-8985)ea

Probe Cover (566-0460)ea





HENRY SCHEIN®

Pro Plus Series Palm Aneroid

Latex-Free

This revolutionary palm-style manometer has an integral latex-free bulb and trigger-style air release valve for one-hand convenience.

Features:

- Precision-crafted, 300mm-Hg, CE-marked, no-pin stop movement tested to 50,000 cycles
- Movement housed in a zinc safety cage for unsurpassed durability
- Ambidextrous infinite control trigger-release air valve
- Nylon cuff with proprietary size guide marking system, to prevent miscuffing
- Dip-molded latex inflation bladder withstands heavy use
- Quick-release Luer-slip connector
- Complete with leatherette carrying case

Pro Plus Palm Aneroid

Adult, Black
(112-6067)ea

Pro Plus Cuff & Bladder, 1 Tube
Child
(112-6090)ea

Small Adult
(112-6087)ea

Adult
(112-6088)ea

Large Adult
(112-6092)ea

Thigh
(112-6091)ea



HENRY SCHEIN®

Pro Pocket Aneroid Sphygmomanometer

- Value-priced aneroid that meets the needs of the most demanding professionals
- Precision-crafted movement is rated to 30,000 cycles and backed by a 20-year calibration warranty
- Extended graduation marks eliminate parallax error
- White numbers on contrasting black scale with red indicator needle for easier viewing
- Proprietary Size Guide Marking system enhances measurement ease and prevent miscuffing
- Chrome-plated brass deflation valve
- Inspected, assembled, and packaged in the USA
- Leatherette zippered carrying case
- Latex-free

Child
(112-6073)ea

Small Adult
(112-6071)ea

Adult
(112-6072)ea

Large Adult
(112-6075)ea

Thigh
(112-6074)ea



Blood Pressure Unit

This high-quality BP unit features an accurate, no-needle stop gauge, a durable, easy-to-adjust nylon cuff with hook-and-loop fasteners and a standard pressure-release valve. It was created for the cost-conscious care provider who is not willing to compromise on accuracy and quality. Each unit is sold complete with matching carrying case.

Blood Pressure Unit
Infant
(499-2556)ea

Pediatric
(499-2555)ea



Adult
(499-2488)ea

Large Adult
(499-2489)ea

Thigh
(499-5520)ea

BLOOD PRESSURE MONITORING



HENRY SCHEIN®

Pro Plus Series Clockface Aneroids

- Integral ABS cuff-storage compartment
- Manometer is angled back 15° for easy viewing
- 18" base diameter

Adult, Black
Wall Aneroid, Latex-Free
(112-6069)ea

Adult, Black, Latex-Free, Mobile
(112-6070)ea



W.A. Baum Co. Inc.

Wall Aneroid

- 300mm Hg swivel-type (110°-arc) sphygmomanometer
- Meets ANSI/AAMI SP-10 for accuracy and performance
- Neutral white ABS plastic manometer with chrome-plated steel bracket
- Unique, high-contrast scale and red pointer for ease of viewing
- Calibrated V-Lok® inflation system

#0950NL, With Adult Nonlatex V-Lok® Inflation System and Wall Bracket
(603-7721)ea



Advanced Wrist Blood Pressure Monitor

- 4-user memory with averages
- WHO risk indicator with 6-color LED
- Irregular heartbeat detector
- Accurate vital health information in seconds

#1147
(635-0059)ea

WelchAllyn®

Tycos® 767-Series Mobile Aneroid

Mounted on a sturdy, quiet mobile stand for portable accuracy from room to room. Features all of the benefits of the 767-series Wall Aneroid, plus it has an adjustable-height stand and 360° swivel for clear viewing, large wire storage basket for cuffs, and 4' coiled tubing to provide great range.



(294-5567)ea



OMRON.

10 Series™ BP785N Upper Arm Blood Pressure Monitor

- Memory storage
- Irregular heartbeat detector
- Two-user mode
- ComFit™ cuff
- BP level indicator
- AC adapter (can also run on 4 AA batteries, not included)
- Tru Read™ with Advanced Averaging
- AM/PM averaging
- Calibration check system
- Cuff wrap guide
- Carrying case include

BP Monitor
(635-0278)ea



OMRON.

Automatic Blood Pressure Monitor

- Memory storage: 60
- Irregular heartbeat detector
- 2-user mode
- Advanced averaging
- Hypertension indicator
- IntelliSense® inflation technology
- Cuff Size D, Ring cuff arm 9" to 13" (optional Large cuff, 13" to 17", H-003D)
- 4 AA batteries included; AC adapter not included
- 2-year warranty

#HEM-712-C BP Cuff, Large Adult
(427-1666)ea

BLOOD PRESSURE MONITORING/ DIAGNOSTIC SYSTEMS



Elemano™ Blood Pressure Monitor

A hybrid design blood pressure monitor

- Gold standard accuracy: +/- 3mm Hg
- Minimal user errors: Ease of use and measurement standardization
- Total flexibility: Normal mode—routine blood pressure measurements; Slow mode—for patients with hypotension and weak pulse; Auscultation mode—conventional manual reading for blood pressure; Improved technology—double-cuff technology to minimize variety of blood pressure readings
- Inflation bulb and air bag contain natural rubber latex

X-small Cuff Cover (315-7179)	ea
Small Cuff Cover (315-9163)	ea
Large Cuff Cover (315-9330)	ea
X-small Cuff Cover with Air Bag and Tubing (315-5187)	ea
Small Cuff Cover with Air Bag and Tubing (315-5582)	ea
Medium Cuff Cover with Air Bag and Tubing (315-1986)	ea
Large Cuff Cover with Air Bag and Tubing (315-7422)	ea



HEM-907XL Intellisense® Professional Digital Blood Pressure Monitor and Accessories

The most recent advancement in digital, noninvasive blood-pressure monitors. Engineered for accuracy, precision, and clinical effectiveness, it is mobile, easy to use, time-saving, and versatile. Developed specially for professional use in a clinical setting, including physician's office, hospitals, or other point-of-care environments.

- Automatic inflation level setting with IntelliSense® allows for automatic cuff inflation, eliminating the need to preset inflation levels
- Average mode measures up to 3 readings and averages the total, according to current blood pressure, for a more accurate measurement
- Large, easy-to-read display
- Auto cuff pressure control for manual (*auscultatory*) measurement
- Quiet operation eliminates unnecessary patient anxiety and clinical noise
- Quick-and-easy one-button operation

Monitor (524-3175)	ea
Contains: large, medium & small cuff/bladder set, 1.0-m air tube & AC adapter.	
Stand (635-3856)	ea
Wall Mount Kit (256-6232)	ea
Pole Mounted Kit (524-0630)	ea



DIAGNOSTIC SYSTEMS

WelchAllyn® 97100-M Diagnostic Set

3.5 V (566-1543)	ea
Contains: 1—#23820 MacroView™ otoscope with throat illuminator, 1—#11710 3.5-V halogen ophthalmoscope, 1—#71000-A rechargeable 3.5-V handle, 4 reusable specula, & 1—#05259-M hard case.	



3.5-V Otoscope Head

#25020 (566-5884)	ea
Contains: illuminator complete with 4 polypropylene specula (<i>sizes 2.5 mm, 3 mm, 4 mm & 5 mm.</i>)	





Pneumatic/Operating Otoscope Heads

Ultrasealed for easy pneumatic otoscopy.

- #20201, without Tips
(729-3654)ea
- #20200, with Tips
(566-2532)ea




HENRY SCHEIN®

Otoscope Speculum Tips

These disposable tips are carefully manufactured to be free of rough edges or burrs. They are designed for a snug fit on Welch Allyn®, Propper™, and Heine™ otoscopes. They will not fit the Welch Allyn® Macroview™.

- 2.5 mm
(112-5527)1000/bag
- 4.0 mm
(112-5528)1000/bag



WelchAllyn®

Green Series™ 777 Integrated Wall System

Providing you with accurate, reliable, energy-efficient products you trust in a workflow that will help improve your practice's throughput.

- Standardize your facility with a sturdy, predrilled mounting board that mounts easily to any wall
- Choose either PanOptic™ Ophthalmoscope (11810 or 11820) or the Coaxial Ophthalmoscope (11720)
- Choose either MacroView™ Otoscope (23810 or 23820) or Pneumatic Otoscope (20200)
- Help improve workflow with contoured handles that automatically turn on/off
- Configure to meet your needs with the choice of the Wall Aneroid, ProBP™ Digital Blood Pressure Device, Spot Vital Signs® Device or Spot Vital Signs® LXi Device
- Optional KleenSpec® Disposable Diagnostic or Pneumatic Specula Dispenser
- Optional SureTemp® Plus 690 Electronic Thermometer
- Durable, easy-to-clean, high-impact ABS housing

#77510, With Specula Dispenser
(672-0015)ea

GS 777 Wall Transformer
#77791-1MPX, With PanOptic Ophthalmoscope, MacroView Otoscope, Wall Aneroid, KleenSpec Specula Dispenser, SureTemp Plus 690 Thermometer & Wall Board
(672-0016)ea

#77791-2MPX, With PanOptic, MacroView Otoscope, Wall Aneroid, KleenSpec Specula Dispenser, SureTemp Plus 690 Thermometer & Wall Board
(672-0018)ea

#77791-MX, With Coaxial Ophthalmoscope, MacroView Otoscope, Wall Aneroid, KleenSpec Specula Dispenser, SureTemp Plus 690 Thermometer & Wall Board
(118-3650)ea

#77791-MXNOBP, With Coaxial Ophthalmoscope, MacroView Otoscope, KleenSpec Specula Dispenser, SureTemp Plus 690 Thermometer & Wall Board
(672-0020)ea



ZOLL

ZOLL AED Plus®

Bringing Lifesaving Defibrillation to Everyone
The only AED with CPR feedback and an ECG waveform. By supporting CPR as well as defibrillation, this device can help every unresponsive victim. Features Real CPR Help®, a tool that is able to actually see what user is doing and provides audio and visual prompts to help them do it with confidence and clarity.

- Provides a one-piece electrode for accurate placement
- Uses consumer lithium batteries that are easy to replace, long-lasting, and inexpensive

Specify:

Fully Automatic Unit	(119-7910)
Semi-Automatic Unit.....	(866-7559)
Contains: AED, electrodes, batteries, ECG, wall-mount bracket & operator's guide.	
CPR-D padz®—Adult Defibrillator Pad (253-2611)	ea
Pedi padz® II—Pediatric Electrode Pad (203-7054)	ea
Soft Carrying Case (106-1937)	ea
3D Wall Sign (110-5384)	ea
Wall Cabinet (108-7520)	ea
AED Plus Simulator (499-7385)	ea



defibtech
Lifeline AED

Designed for speed and simplicity. Uses biphasic defibrillation—the most effective at resuscitating in ventricular fibrillation, the most common form of sudden cardiac arrest. Text, lights, and voice instructions guide the responder step by step through the lifesaving process. Electrode pads are stored in the back pocket, ready to apply. “Power On” and “Shock” buttons are oversized and brightly lit. Lightweight and durable. 4.2 lb. 5-year limited AED warranty, 3-year warranty on the battery pack, 1-year warranty on the 9-V lithium battery, and 2-year shelf life on the defibrillation pads.

Complete Package

(499-1362)ea

Contains: 1 of each: unit, battery pack, set of defibrillation pads, 9-volt status indicator battery, quick-start card & user manual.

Defibrillation and Monitoring Pads

Adult Pad Set

Preconnected, single-use, nonpolarized, and disposable. Self-adhesive electrodes with cable and connector.

(499-1367)ea

Pediatric Pad Set

(499-1368)ea

Wall Mount

(499-3395)ea



LifePak® CR® Plus

Designed for the first person to respond. This AED features an easy two-step operation, turn it on and then apply the electrode pads. The unit then takes over, giving prompts and delivering the shock. 8-year warranty.

Semi-Automatic Unit
(365-0088)ea

Fully Automatic Unit
(106-8280).....ea

Contains: 2 pair of Quik-Pak® packing/defibrillation/ECG electrodes with 1 Redi-Pak™ preconnect system, 1 Charge-Pak™ charging unit, 1 quick-reference instruction card, 1 carrying case, 1 Ambu® Res-Cue mask, 1 AED program implementation starter kit, 1 reference manual CD & 1 orientation DVD.



Infant/Child Reduced Energy Defibrillation Electrode

(985-6621)ea

Adult Defibrillation Electrodes

(920-4738)pair

Electrodes and Charge-Pak™ Replacement Kit

(601-0223)ea

Contains: 2 sets electrodes, 1 battery charger, instructions & discharger for safe disposal of used charger.

Replacement Battery and Adult Electrodes

(373-3406)ea

Contains: 1 of each: pair of adult electrodes & Charge-Pak™ charging unit.

Soft-Shell Carrying Case

(486-1462)ea

Hard-Shell, Watertight Carrying Case

(111-6660)ea

Wall Cabinet with Alarm

(499-7215)ea

SummitDoppler

LifeDop™ Dopplers and Probes

Offers great sensitivity, a recording option, and superior sound. Long-lasting battery saves money. Comfortable, easy-to-hold probes. Very durable.

- Large, easy-to-read LCD and technically advanced software display extended 50–220 bpm heart rate as well as remaining battery life, signal quality, and probe identification; available in both display and nondisplay units
- Built-in Audio Recorder LifeDop™ is the first Doppler to feature a built-in solid-state recorder; records up to 32 seconds of Doppler sounds (*nonrecording Dopplers also available*)
- Quality speaker and audio processing for superior sound
- Improved battery management design to reduce frequency of battery charges

#L150R, Nondisplay Doppler, with Recharger (118-8427)	ea
#L150A, Nondisplay Doppler, with Audio Recorder (118-8694)	ea
#L250, Display Doppler, Basic Model (118-3501)	ea
#L250AR, Display Doppler, with Recharger and Audio Recorder (118-9924)	ea
#8P800, Vascular, 8MHz Probe (118-9695)	ea



Newman Medical

DigiDop™ 300 and DigiDop™ 301 Non-Display Digital Dopplers

Feature Digitally Optimized Processing, DOP™, for high performance and durability for obstetrical or vascular applications. Tough and durable polycarbonate and stainless steel case with unique magnetic probe holder and full 2-year product warranty. Made in the USA.

With 8-MHz Vascular Probe (118-4769)	ea
DigiDop™ 301 Rechargeable Doppler with Recharger & 8-MHz Vascular Probe (119-5396)	ea



SCHILLER
A M E R I C A

Cardiovit AT-2

Low-cost ECG system with performance features previously only available in high-priced systems. Provides A4 reports within seconds at the touch of a button. Features a powerful integrated, rechargeable battery that generates up to 400 reports.

- 6-channel and 12-channel representation of all 12 simultaneously acquired leads
- Automatic and manual modes
- Full-size (8½" x 11") printouts (A4)
- Computerized ECG interpretation program
- Accurate measurements
- Single-button operation
- Portable: only 4.8 kg (11 lb)
- Battery operation

(600-4272)	ea
------------------	----

AT-2 Light ECG

Offers the economics of a low-cost unit with many of the high-end features needed to care for patients. Ideal for a small office that wants to offer more. 12-channel system offers full-size reports to fit patient records, alphanumeric keyboard with large LCD screen for easy entry of patient demographics, and the option of spirometry and interpretation.

- Battery operation of up to 6 hours
- 3-year warranty
- 8½" x 11" printouts
- Automatic printout function for 6 and 12 channels
- Adult and pediatric ECG interpretation and measurement
- 6-channel printout in real time at 5 mm/S, 25 mm/S, and 50 mm/S

(104-2998)	ea
------------------	----

SCHILLER
A M E R I C A



CP 50™ Resting Electrocardiograph



12 Lead

Get the power of a full-featured ECG in a compact, affordable, easy-to-use package. With simple, one-button operation, this device makes performing ECG tests easy for trained caregivers of all skill levels. Other advanced features: a full-color touch-screen display, high-resolution thermal printer, and a full range of EMR/EHR connectivity options.

- Reliable, full-color touch-screen display
- On-screen instructions that guide you through a test
- Programmable report formats
- High-resolution thermal printer
- External printer integration
- Programmable patient entry fields
- Optional MEANS interpretation
- Optional mobile stand and carrying case
- Flexible connectivity options:
EHR, EMR, HIS, e-mail, Web, and telemedicine



#CP50-1EN1, CP 50 Non-interpretive (566-0338)	ea
#CP50A-1EN1, CP 50 Interpretive (566-0339)	ea
#CP50P-1EN1, CP 50 Plus, Non-interpretive (566-0340)	ea
#CP50AP-1EN1, CP 50 Plus, Interpretive (566-0341)	ea

Accessories

#45008-0000, Disposable Tab Electrodes (701-6335)	1000/box
#406815, CP 50 Mobile Stand (566-0345)	ea
#715107, Carrying Case (566-0344)	ea
#406021, CP 50 Printer Paper, Z-fold (566-0342)	4 packs/case
#406020, CP 50 Printer Paper, Roll (566-0343)	4 rolls/case

**Burdick ELI™ 250c
12-Lead Resting Electrocardiograph**



- **Portable solution.** Compact and lightweight, the ELI 250c electrocardiograph provides comprehensive functionality, in an easy-to-use, portable device.
- **Best 10 technologies.** Capture the best 10 seconds of ECG data from system memory, reducing the need for repeat ECGs.
- **Wireless data acquisition.** Increase workspace flexibility by gathering 12-lead data without the need for a tethered cable.
- **Connectivity your way.** Multiple options ensure patient data transfers easily from the ELI 250c to the location of our choice.



BUR250C-WIX with WAM & 802.11 Wireless data acquisition (907-3637)	ea
---	----



HIGH-PERFORMANCE BLOOD GLUCOSE MONITORING SYSTEM

TRUEresult™

TRUERESULT™ BLOOD GLUCOSE MONITORING SYSTEM

- A no-coding system with high performance features
- Clinically proven accurate
- Lifetime warranty
- Only requires a tiny, 0.5-mL blood sample
- Provides accurate results in as fast as 4 seconds



Features audible fill detection, easy-to-read display, 500-test memory, a weekday display, 7-, 14-, and 30-day averaging, glucose control detection, 4 testing reminders, ketone testing reminder, a strip release button, auto on/off, and data management capabilities.

(112-6826)	ea
------------------	----

Contains: TRUEresult® meter, installed 3-volt battery, lancing device, 10 lancets, instructions for use, self-test log book & compact carrying case.



TRUETEST™ TEST STRIPS

For use with the TRUEresult® Blood Glucose Monitoring Systems. TRUEtest™ Test Strips use Quad-Electrode Laser Accuracy for enhanced precision and accuracy. They also use the patented TRUEfill® beveled tip, which allows for precise entry of blood sample.

(112-6827)	50/box
------------------	--------



Rely On Us for Quality, Selection and Performance.

GLUCOSE MONITORING



Freestyle Precision H



Blood Glucose Monitoring System

- Used for whole blood: arterial, venous, fingertip capillary, and neonates
- Simple 2-step testing
- 0.6- μ L sample size, 20 second test time for glucose
- Top- or end-fill provides easy sample application with visual confirmation
- Fill trigger ensures reaction starts only when enough blood applied
- Individual foil-wrapped package
- Proven accuracy with GDH-NAD chemistry
- Minimizes potential interfering effects of common medications and substances

Freestyle Precision H Meter (657-0005)	ea
Freestyle Precision H Strips (657-0004)	100/box
High/Low Control Solution (657-6290)	2 vials
TriLevel Glucose Control (657-3970)3/box



EmbracePRO™ Glucometer

Professional Blood Glucose Meter



- Easy to use
- Test strip release button
- Auto coding
- Fast—results in 5 seconds
- Approved for multipatient use
- Test with less pain
- Small sample size: 0.5 UL
- Alternate site testing: fingers, palm or forearm
- Manage your results
- 500-test memory with time/date stamp
- Event Recording Mode for more effective data management

Glucometer (240-0042)	ea
Test Strips (240-0043)	50/box
Control Solution Hi/Lo (240-0044)	ea



Hb 201+ Handheld Hemoglobin Analyzer

Offers the same important benefits as its predecessor, the B-Hemoglobin Photometer. Smaller size gives more flexibility and mobility to use the system anytime and anywhere. Performs an internal self-test. No control cuvette required.

- Produces results in less than 1 minute with lab-quality accuracy
- 3 easy steps
- Precalibrated, portable analyzer that requires minimal maintenance
- No calibration or instrument manipulation (between cuvette/reagent batches) is needed
- Safe and easy to use
- Reimbursable
- CLIA '88 waived
- 2-year warranty



CONTOUR® Blood Glucose Meter



- No coding required
- Automatic underfill detection
- Automatic control marking
- Automatic temperature control
- Small sample size (0.6 μ L)
- Quick testing time—results in 15 seconds
- Saves information for up to 240 tests with time and date

#7151, Meter (147-1405)	ea
#7099, Strips (147-9303)	50/vial
#7109, Control—Normal, 2.5 mL (147-0874)	ea
#7111, Control—High, 2.5 mL (147-8283)	ea
#7110, Control—Low, 2.5 mL (147-5579)	ea

CPT Code; 82962QW

A1CNow+® Tests

This handheld monitor provides accurate in-office A1c results, allowing physicians to make on-the-spot treatment decisions for diabetes patients.



- CLIA waived
 - In-office testing; no waiting for lab results
 - Enables A1c testing in every exam room
 - Results in just 5 minutes
 - No refrigeration necessary if used within 4 months
 - Made in the USA
- | | |
|---------------------------|--------|
| #3024
(261-0324) | 10/box |
| #3021
(261-0323) | 20/box |
- CPT Code: 83037QW†



HENRY SCHEIN®

Fingertip Pulse Oximeter

- Adult and pediatric use: ideal for all specialties
 - Range for SpO₂: 70%–100%
 - Range for pulse: 30–235 bpm
 - Dual-color LED with display modes and brightness control
 - Includes pulse bar and plethysmogram
 - Visual and audible low-battery alarm
 - 2 AAA alkaline batteries with 30 hours of continuous use
 - Electromagnetic compatibility
- (900-6368)ea



smths medical
bringing technology to life

Digit® Finger Oximeter

Combines monitor and sensor into one unit that provides fast, reliable SpO₂, pulse rate, and pulse strength measurements on patients from pediatric to adult.

- Portable, lightweight, pocket-sized monitor
- One-button keypad for ease of operation

#3420Y, Yellow
(488-1236)ea



MedSource
INTERNATIONAL

Fingertip Pulse Oximeter MS-74001

Affordability in Pulse Oximetry
This unit features a single directional LED display that includes a pulse-rate monitor, an SpO₂ monitor, and a pleth bar. The automatic switch-off feature preserves battery life. Suitable for adult, pediatric, and neonatal use.

(499-4035)ea



smths medical
bringing technology to life

BCI® 3301 Hand-held Pulse Oximeter

- Bright, easy-to-read LED numerals
- Quantitative pulse strength bar depicts actual signal strength
- Provides fast, reliable SpO₂ and pulse rate measurements on any patient from neonate to adult
- HPO software corrects background noise for enhanced artifact recognition
- Unique memory feature stores data for up to 99 different patients
- 3 C batteries (included) provide 24½ hours of continuous run time (approximately 1500 one-minute spot checks)
- Auto shutoff feature turns unit off when probe is removed from patient
- Reusable finger probe #3044 and operator's manual are included with the unit
- 2-year unlimited warranty

Low-cost, handheld pulse oximeter with an integral printer. A lightweight, highly portable device providing SpO₂ and pulse rate measurements, it has easy-to-read red LEDs, runs on 4 AA alkaline batteries, and has a rubber boot for added protection.

#3301A1, Pulse Oximeter
(117-1807)ea



MedSource
INTERNATIONAL

Fingertip Pulse Oximeter

Offers 6 display modes, including SpO₂, Pulse Rate, and Pleth Waveform. This easy-to-read, 2-color OLED display makes monitoring clear and simple. The automatic switch-off feature preserves battery life. Suitable for adults and pediatric patients.

(499-4744)ea

PULSE OXIMETERS/ VITAL MONITORS



smiths medical
bringing technology to life

Spectro2™ 10 Pulse Oximeter

- Large, bright, easy-to-read LED for at-a-glance assessment
- Unique onboard cradle to store and protect finger sensor
- Multiple power options
- Compatible with docking station, attachable printer, and BCI® spot-check sensor
- Pulse amplitude index (PI)
- Low-battery indicator

#WW1000A1EN
(499-8711)ea

VITAL MONITORS

WelchAllyn

Spot Vital Signs™ Devices

- Fast and accurate—completes BP, pulse rate, temperature, and SpO2 in approximately 30 seconds
- Multiparameter and modular so you can select precisely the vital signs you need
- Optional Masimo® or Nellcor® pulse oximetry
- Optional SureTemp® thermometry
- IR connectivity to interface with computerized patient records
- Portable—carry it, wheel it, or mount it
- Big, bright, easy-to-read LCD screen
- Offers you one of the lowest priced multiparameter instruments on the market
- Fully upgradeable to meet your changing needs
- Integrated long-life battery enhances portability and continuous operation



#4200B-E1, Spot Vital Signs with NIBP only (pulse rate and MAP)
(151-1071)ea

#420TB-E1, Spot Vital Signs with NIBP and SureTemp Thermometer
(464-7803)ea

#42N0B-E1 Spot Vital Signs with NIBP and Nellcor Pulse Oximetry
(611-4160)ea

#42NTB-E1, Spot Vital Signs with NIBP, Nellcor Pulse Oximetry and SureTemp Thermometer
(522-6109)ea

WelchAllyn

Spot Vital Signs® LXI

- SureBP™ measures blood pressure in about 15 seconds as the cuff is inflating, resulting in time savings for you and improved comfort for your patients
- Choice of SureTemp® Plus or Braun ThermoScan® PRO 4000 ear thermometer
- Choice of Masimo or Nellcor pulse oximetry
- Ability to manually enter weight, height, respiration rate, and pain level
- Calculates body mass index (BMI) or device connects to Health-o-meter or Tanita digital scales
- 50-reading memory allows retrieval of previous data for review or to download later
- Comes connectivity ready, allowing you to connect to the electronic medical record system you now have or one you will get in the future
- Bigger, brighter LCD readout provides greater visibility from farther away



#450E0-E1, with SureBP™ and Braun ThermoScan PRO 4000 Thermometer
(566-0606)ea

#45ME0-E1, with SureBP™, Masimo SpO2, and Braun ThermoScan PRO 4000 Thermometer
(566-3405)ea

#45NT0-E1, with SureBP™, Nellcor SpO2, and SureTemp® Plus Thermometer
(566-0052)ea

#45MT0-E1, with SureBP™, Masimo SpO2, and SureTemp Plus Thermometer
(566-3759)ea

#45NE0-E1, with SureBP™, Nellcor SpO2, and Braun ThermoScan PRO 4000 Thermometer
(566-7992)ea



Four-Drawer Crash Cart with Specialty Accessory Package

Features:

- Breakaway lock bar
- 100 breakaway seals included
- 24" of vertical drawer space
- Stabilizing bumper frame with extended front casters and 4 corner bumpers
- 5" ball-bearing casters, all full swivel, 2 with brake
- Predrilled threaded holes for V-Series side-mounted accessories
- Ball-bearing, full-extension drawer slides
- 12-year warranty
- Ships with wheels installed and ready to use
- Minimal assembly required at delivery
- Latex-free
- Durable, hard-baked powder-coat finish, Red
- Steel construction

Emergency package accessories include:

- Cardiac board and brackets
- Side-mounted drop shelf
- IV pole
- Oxygen tank holder (*D/E size cylinders*)
- Powder-coated steel top rail/push handle



#V24-4B-HSACC, Crash Cart (738-0008)ea



Emergency Crash Carts

- Stainless steel breakaway locking system
- Full wraparound bumper comes standard on the 6400 and 6300 Series
- Stainless steel top rail/push handle comes standard on the 6400, 6401, 6300, 6301, E24-4EMS and E30-6EMS carts
- Drop shelf comes standard on the 6401, 6301, E24-4EMS and E30-6EMS carts
- All come with a removable plastic top
- All steel construction
- Ball-bearing, full-extension drawer slides
- Powder coat finish in 12 colors
- 5-year warranty

- 44½"H x 22"D
- #6401, 6 Drawers with IV Pole and Oxygen Tank Holder, 38"W (600-4365)ea
- #6400, 6 Drawers, 38"W (280-5287)ea
- #6201, 2 Drawers with Lower Cabinet, 32"W (189-3578)ea



Treatment Carts

- Key lock with 2 keys
- Removeable plastic top
- Aluminum push handle
- Ball-bearing full extension drawer slides
- Unibody steel construction
- 40"H x 18"W x 18"D
- Powder coat finish available in 12 colors

- #3156B-EMG, 6 Drawers with IV Pole and Oxygen Tank Holder (921-2755)ea
- #3156B, 6 Drawers with Breakaway Locking System (701-6338)ea
- #3156K, 6 Drawers (306-1676)ea
- #3154K, 4 Drawers (104-8029)ea

CARTS



2-Shelf Utility Carts

Perfect for transporting equipment and heavy loads in any environment. Constructed of sturdy, structural foam, so it will not rust, dent, chip or peel. Large, nonmarking casters reduce touch-up cleaning. Rounded corners protect walls and furniture.

- Lightweight and maneuverable
- Foam construction
- Casters
- Mass capacity: 500 lb

33¼"H x 25⅞"W x 45¼"L
Black
(606-9546)ea

33¼"H x 25⅞"W x 45¼"L
Beige
(606-9547)ea



Standard-Duty Utility Carts—Stainless Steel

Made from stainless steel, they feature 3 stain-resistant shelves, 3½" casters, and a load capacity of 300 lb.

- Durable, unitized electronically welded stainless steel angle leg design
- Reinforced, hemmed front on 22-ga stainless steel shelves
- Legs are 20-ga stainless steel
- Sound-deadening panels under shelves restrict vibrations and noise
- Bumpers on legs and handles
- Cleans and sanitizes easily
- 300-lb capacity

33"H x 18⅝"W x 30¾"L
(522-7623)ea

32⅞"H x 16¼"W x 27½"L
(522-0083)ea



Brewer

Plastic All-Purpose Cart

- 200-lb capacity
 - 3" swivel casters
 - Easy to assemble
 - 3 ABS plastic shelves (27¼"L x 17"W)
- (301-7983)ea



Brewer

Multi-Purpose Cart

- Strong, chrome-plated steel frame
- 3" rubber wheels; ball-bearing casters
- Convenient height of 35"
- Coated steel basket supports; no assembly required

#63400, Cart Frame Only
(923-8191)ea

#63021, 12" Wire Basket
(452-6707)ea



Oxygen Crash Cart

- Dimensions: 36" x 13" x 10"
- Weight: 19½ lb complete

#1629
(923-3034)ea

Contains: #1603 aluminum cylinder (680 liters filled), #1442 fixed-flow regulator (6 LPM avg.), #370C adjustable chrome cart, #4000 manual resuscitator, mask, & tube.

♦ Hazardous material



Chrome Hydraulic Lift

Features:

- 450 lb weight capacity
- Lightweight construction and easy disassembly allow quick transport and setup
- 20"-64" range allows for patient to be picked up from lying position on floor
- Elongated handle design serves a broad range of operator heights and offers improved handling ability
- Offset mast and boom style provides better lift path, maximizing full travel range
- Pump handle can rotate from side to side for convenience of attendant
- Six-point swivel bar attachment easily adapts to all sling styles and positions



- #9805, Hydraulic Lift (120-9593)ea
- #9402, Standard Sling (120-8369)ea
- #9071, Sling Chains for Standard Series Slings (120-7757)ea



Invacare® Reliant Series 450 Full Body Lift

- Total pinch-point protection covers all moving parts for added safety
- Height range allows lifting from floor and high surfaces
- 450-lb capacity



The RPL450-1 Battery-powered Lift's electronic system includes 2-24-volt batteries: 1 for the lift, 1 for the battery-charging station; a manual emergency lowering device; and anti-entrapment if boom meets any resistance while lowering. Covered base provides protection against wear, dirt, and moisture. The RPL450-2 Power Lift's electronic system includes 2-24-volt batteries, a manual emergency lowering device, and anti-entrapment if boom meets any resistance while lowering. A padded swivel bar with 360° rotation and 6-point hookup allows use of Invacare chain-free slings as well as older-style slings with chains or straps. Power-opening base legs adjust easily using the hand pendant.

- #RPL450-1, Battery-Powered Lift with Low Base (246-0993)ea
- #RPL450-2, Power Lift with Power Opening Low Base (834-5641)ea

*Slings sold separately



Reliant 600 Heavy Duty Patient Lift

Designed to meet the needs of the growing bariatric population up to 600 lb, the Reliant 600 lifts provide a large lifting range that allows for patient transfers to high surfaces. Low-friction casters significantly improve maneuverability and contribute to caregiver safety. The boom is controlled through the hand pendant. The manual low base provides access to beds and furnishings that are close to the floor with a 4.5" base height. An essential device for safe patient handling in hospitals and LTC (long-term care) facilities. Total pinch-point protection covers all moving parts for added safety.

Features:

- Manual lowering for instant response to power loss
- Safety tested with CE mark for highest quality standards
- Base legs adjust easily and lock securely into open position with the padded shifter handle
- Extra-wide coated swivel bar with 360° rotation and 6-point hookup
- Unique sling fabric designed to provide the utmost in comfort and strength
- Warranty includes 3 years on lift, 1 year on electric components

- #RPL600-1, Power Lift with Low Base (108-7889)ea
- #RPL600-2, Power Lift w/Power Opening Low Base (118-4333)ea

*Slings sold separately



Reliant Digital Scale



- Weight capacity: 600 lb
 - Works with Reliant 450/600; Model 9805 may require the swivel bar kit
 - Easy-to-use, 4-button display
 - Automatic shutoff to prolong battery life
- #RLS6 (120-1238)ea

B|BRAUN

Vista® Basic Infusion Pump

Large-volume infusion pump, with Continuous, Piggyback, Ramp/Taper, and Program modes. Operates with DEHP-free and straight-line IV sets. Using prompts, it walks user through programming process. Easy-to-load pump is less expensive than cassettes and fosters quick learning. Easy to carry, lightweight (6.8 lb), and doesn't require an IV pole.

(187-1304)ea



IV Stand

Custom

(106-5072)ea



Heavy Base Four-Hook IV Pole

- Chrome-plated steel pole
- Knob height adjustment
- Detachable dual ram horn hooks with plastic end caps
- Epoxy-coated, 22"-diameter steel base
- 4 rubber-wheeled, ball-bearing casters


IV Pole
(635-4559)ea



Value IV Pole with 2 Hooks

- Chrome-plated steel pole
- Twist lock adjusts height 51½"-93"
- Detachable dual ram's horn hooks with plastic tips
- Dual wheel casters
- Epoxy-coated, 4-leg base with bumper guard (22" diameter)

(100-5982)ea



Five-Leg IV Pole with 2 Hooks

- Chrome-plated steel pole
- Knob height adjustment
- Detachable dual ram horn hooks with plastic tips
- Dual wheel casters
- Five-leg, cast-aluminum base

(635-5666)ea



3-Panel Mobile Screen

Mobile design. Each 27"W panel contains a flame-retardant, bacteriostatic, white vinyl curtain. Unit contains 4-2" twin-wheel casters. Overall height is 69".

(120-0545)ea



3-Panel Privacy Screen

- Sturdy 1" diameter anodized aluminum tubing
- Screens are 6-mil flame-retardant bacteriostatic white vinyl
- Four 3" hooded casters
- Easy to maintain and clean
- Overall height: 70½"
- Product weight: 15 lb
- Shipping weight: 20 lb
- Each hinged panel measures 56"H x 27½"W
- Limited lifetime warranty

#13508
(109-9968)ea

HENRY SCHEIN®

PowerSpin™ LX Centrifuge

Spin blood, urine, fecal materials, or semen at the best speed for each type of specimen with this powerful, variable-speed, 6-place or 8-place centrifuge.



Maintenance-free, brushless motor runs quietly at a fixed speed up to 4000 rpm, emitting only 58 dB, all with low-temperature rise. Suction-cup rubber feet hold unit firmly in place and virtually eliminate vibration. Lid features a gasket to reduced aerosol emission, a window for strobe speed check, and a safety interlock switch that cuts power to the motor if lid is raised. Angled rotor accepts 2–10-mL tubes (maximum tube size 16 mm x 110 mm).

- Speed range: 300–4000 rpm
- 30-minute timer
- Stainless steel rotor chamber
- Lid latch safety switch
- Dimensions: 10.2" x 11.8" x 7.9"
- Weight: 15 lb
- 115 V/60 Hz

#C-856, Variable Speed, 6 Place
(112-6907)ea

#C-858, Variable Speed, 8 Place
(112-6909)ea

Tuttnauer
Your Satisfaction Is Our Greatest Priority



Valueklave 1730

The speed and quality you demand with the value you deserve.

- Provides all basic sterilization needs at an economical price
- Compact size (7" chamber) makes it a perfect unit for small offices or a second autoclave for larger practices
- Fast cycle time: 9 minutes from a hot start (*standard unwrapped cycle time*)
- Overall dimensions: 12"H x 17.4"W x 17.9"D
- Chamber size: 7" x 12"

#1730
(668-2244)ea

HENRY SCHEIN®

Urispec® Plus Urinalysis Reader

Provides dependable urine status results to detect early stages of many diseases, such as diabetes, kidney disease, and urinary tract infections. The comprehensive interface options and the optimized printouts minimize risk associated with manual transcriptions. Reliable results can be obtained virtually immediately at the point of care.

- CLIA-waived
- Readings eliminate subjectivity of visual color interpretation, which makes urinalysis simpler and more reliable
- Features no-touch operation: After strip placement, measurement starts automatically, which makes urinalysis faster and more hygienic
- Easy switch between different test strips: Urispec Plus is for use with Henry Schein Brand Urispec 11-Way, 5-Way Plus, 4-Way, and GP+A urine test strips



Reads 11-Way, 5-Way Plus, 4-Way & GP+A Test Strips

(900-7507)ea

Contains: 1 Urispec Plus Urinalysis Reader & 12 bottles of Urispec 11-Way Urine Strips.

HENRY SCHEIN®

G380 Microscope

Designed for viewing comfort, this microscope provides crisp, clear imaging on a versatile, ergonomic platform. Standard 10x high-eyepoint eyepieces are engineered to reduce eyestrain for users who wear glasses. Each eye tube features independent focusing diopter adjustment to compensate for the natural vision differential between each eye. Base design adds stability to the body and allows extra room for hand access to the oversized coarse-focus and fine-focus controls and easy access to the coaxial mechanical stage controls.



- Siedentopf binocular head
- Dual diopter adjustments
- 10x high eyepoint widefield eyepieces
- 4x, 10x, 40x & 100x (oil) objectives
- Variable LED illumination
- NA 1.25 Abbe condenser

(112-6914)ea

HENRY SCHEIN®

Dry Heat Sterilizers

- All stainless steel
- 60 minutes sterilization set at 320°F or 160°C
- Low operating cost: economical
- Maintenance-free chamber sealing door: no gasket to be replaced
- Burnout-proof heating elements: maintenance-free
- No pressure and no moisture: safe; prevents rusting of instruments
- Trays can be sterilized with lids on: convenient
- Heat is even; heats up rapidly for fast operation: prevents scorching or incomplete sterilization in spots
- Precision thermostat: ensures a controlled temperature
- Will not dull sharp instruments: promotes instrument longevity



#200 Unit
Dimensions: outside: 15½"W x 10½"H x 9½"D;
inside (useful space): 12½"W x 5"H x 7½"D;
no. of trays: 2; weight: 17 lb
(100-5400)ea

#300 Unit
Dimensions: outside: 19"W x 14½"H x 11"D;
inside (useful space): 15"W x 7½" H x 9"D;
no. of trays: 3; weight: 27 lb
(100-7250)ea

SCALES

HENRY SCHEIN®

Digital Waist-High Scale

Compact waist-high digital scale has integrated mechanical height rod. Large, backlit LCD display shows lb or kg, and has capacity up to 500 lb (225 kg). Low profile for easy access—platform is only 2½" high. Wheels for easy movement. Contains BMI calculator and a clock. Includes 2-year warranty. 22" x 15" x 36"



(112-7155)ea

HENRY SCHEIN®

Digital Electronic Eye-Height Scale

Compact eye-height digital scale features integrated mechanical height rod, large backlit LCD display, and BMI calculator. Weight capacity up to 500 lb (225 kg). Low profile for easy access—platform is only 2½" high. Sta-Bilt stabilizing software ensures an accurate measure, even if the patient is moving. Autohold feature allows the weight to continue to be displayed after the patient has stepped off the scale. An integrated USB interface allows data to be transferred to a computer. Includes 2-year warranty. 22" x 15" x 53"



(112-7154)ea

seca

Digital Column Scale

High-Capacity Scale

With very high capacity and BMI function. Extremely robust, precise, and stable: with a capacity of 550 lb, a large but low platform and the integrated BMI function, this column scale has the prerequisites for the diagnosis and therapy of heavyweight patients. The energy-saving battery operation—up to 16,000 weighings with just one set of batteries—means that this scale can be used anywhere, as it is not dependent on AC power.



- Includes seca #220 measuring rod
- Capacity: 550 lb/250 kg
- Graduation: 0.2 lb/100 g
- Functions: TARE, HOLD, BMI, kg/lb/sts switch over
- Transport casters
- Optional adapter set: RS 232 seca #460

#703
(600-5578)ea

HENRY SCHEIN®

Mechanical Beam Scale With Height Rod

Standing mechanical beam scale with built-in height rod. Weight capacity of 450 lb (200 kg). 1-year warranty. 22¼" x 18¼" x 58¼"



(112-7157)ea

DETECTO

Eye-Level Beam Scale

- Capacity: 400 lb x 4 oz
- Platform size: 10½" x 14½" (267 mm x 368 mm)
- Height rod measures: 30"–78" (76 cm–198 cm)
- Scale height: 59" (1.5 m)



#438, with Wheels
Shipping weight: 39 lb (17.6 kg)
(545-8586)ea

#439, without Wheels
Shipping weight: 38 lb (17.2 kg)
(545-3857)ea

#338, with Wheels and Pound/Kilo Increments
Shipping weight: 39 lb (17.6 kg)
(545-0533)ea

#CW-30X, Counterweight
Provides an additional 100 lb of capacity.
(545-6822)ea

seca

#840 Digital Floor Scale

A timeless beauty, with perfectly executed curves, clear LC display, and AUTO-HOLD function. It embodies grace, convenience, and precision.



- Capacity: 308 lb, (140 kg)
- Graduations: 0.2 lb (100 g)
- Weight: 4½ lb
- Power supply: 2–3-V batteries (included)
- 1-year warranty

#803
(109-2860)ea

HENRY SCHEIN®
Digital Floor Scale



Compact digital floor scale has a large backlit LCD display. Weighs accurately in both lb and kg, up to 400 lb (200 kg). Features antislip feet and a low profile for added safety. Automatically turns on. Includes 1-year warranty. 14" x 11 $\frac{5}{8}$ " x 1" (112-7153)ea

Health o meter®
Digital Scale



Non-slip mat with a 1 $\frac{1}{2}$ " thin platform. 1 $\frac{1}{2}$ " easy-to-read, red LED. Instant on/off feature for immediate weight readout and battery conversion.
 • Capacity: 390 lb x $\frac{1}{2}$ lb (136 kg x 0.2 kg)
 • 5-year limited warranty
 #800KL
 (110-7427)ea

TANITA®
Digital Baby And Infant Scale With LED—Model 1584



Affordable digital infant scale that offers 1-oz (20-g) sensitivity up to 20 lb.
 • Capacity: 44 lb.
 • Graduations: 0–20 lb x 1 oz (20 g) and 22–44 lb x 2 oz (50 g)
 • $\frac{3}{4}$ " LED display
 • Functions: auto off and weight lock
 • 14" x 14" x 23" tray
 • 4 AAA batteries (included)
 #1584
 (611-0994)ea

seca
Mechanical Telescopic Stadiometer



This full-length model is wall mounted, extends to the floor, and comes with a foot plate for increased stability. Additionally, this stadiometer—made from sturdy anodized aluminum—has an additional 4" head plate to ensure further accuracy. Because it telescopes up to eye level, readings will become easier.
 • Measuring range: 2 $\frac{1}{2}$ "–90 $\frac{1}{2}$ " (6–230 cm)
 • Graduations: $\frac{1}{8}$ " x 1 mm
 #222
 (266-8199)ea

TANITA®
Certified Digital Scales



Perfect for those instances where absolute measurement of weight gain or loss is a necessity. Approved as a Class III NTEP Certifiable Legal for Trade Scale, the WB-100 series is available in weight capacities of 600 and 440 lb with automatic BMI calculations standard. Both models can be easily wall-mounted or placed on a countertop with the built-in desktop stand.
 • Platform size: 11" x 13" x 3"
 • 1" LCD display
 • Weight: 11 lb
 • Runs on AC adapter (included) and 6 AA batteries (not included).
 #WB-100A, 440 lb (200 kg) Capacity (600-3250)ea
 #WB-110A, 600 lb (270 kg) Capacity (605-6856)ea

HENRY SCHEIN®
Digital Pediatric Table Scale



Small countertop pediatric scale has a curved design to support and protect the infant. Sta-Bilt stabilizing software ensures an accurate measure even if the infant is moving. Large backlit LCD display measures in lb and oz. Weight capacity of 45 lb. Integrated tape measurer to measure size of the baby. The autohold feature allows the weight to be displayed even after the baby has been taken off the scale. An integrated USB interface allows data to be transferred to a computer. Includes a 2-year warranty. 16 $\frac{3}{4}$ " x 27 $\frac{3}{4}$ " x 4" (112-7156)ea

seca
Accu-Hite Stadiometer



From its built-in leveling bubble to its handsome design, this stadiometer has long been the stadiometer of choice. It features a headpiece that comes with a lock knob for measurement accuracy and safety. Plus, readings can be taken in feet, inches, and centimeters.
 • Measuring range: 24"–83" (61–212 cm)
 • Graduations: $\frac{1}{16}$ " (1 mm)
 #216
 (266-5011)ea

CHAMPION

Standard Medical Recliner-Manual Recline

Patient comfort, ease of use, dependability, and long life make our recliner one of the most popular for oncology, surgery, infusion therapy, dialysis, and other treatment needs.

- **Comfort** : Pressure-reducing foam conforms to patients' bodies, and patients can easily adjust the recliner to the most comfortable position, from full upright to full recline
- **Wall-hugger mechanism** : Move to a full recline position without making contact with the wall
- **Durability** : All-steel frame, premium components and quality construction give years of trouble-free service with proper use
- **Swing-away arms** : Secured with heavy-duty latches, our swing-away arms make patient lateral transfers easy and secure, while improving cleanability for improved infection control (*Series 85 and 89*)
- **Advanced technology options** : State-of-the-art components add



Call today for various chair options and configurations.



Lumex Deluxe Clinical Care Recliner

Sets a new standard of comfort for the patient and ease of use for the caregiver. Its unique, patent-pending recline mechanism enables the patient to position the back from TV position to full recline. The wide, contoured seat and back, along with an integrated footrest, provide an enhanced level of comfort. Built-in arm pads in the back give improved support for the arms while in Trendelenburg. Wall-saver armrest caps increase durability of the armrests and protect walls from marring. Higher floor-to-back height enables the caregiver to raise the patient up with less effort, reducing the chance for a back injury. Four Tente® swivel casters facilitate transport and maneuverability in confined areas. Left- and right-rear top-locking casters provide positive lock. Directional wheel lock on right-rear caster reduces "crabbing" of recliner and enables the chair to be moved down hallways with less effort.



Call today for various chair options and configurations.



Winco CareCliner

The standard in the clinical recliner line. Nonporous, blow-molded, fold-down side trays are mounted on both sides. Each chair in the series features a standard Trendelenburg release with handle accessible from either side. It also has a built-in ergonomic headrest, adjustable elastic bands on the back frame that provide exceptional comfort, an elastic, open-seat base, and nylon casters. All of these features combine to create the ultimate in comfort, durability, and value. Also available with a Dual Drop Arm to allow for easy cleaning and safe transferring of patients to and from the recliner.

- #6530, Standard (*Weight Capacity: 275 lb*)
- #6540, XL (*Weight Capacity: 450 lb*)
- #6550, Drop Arm (*Weight Capacity: 275 lb*)
- #6570, XL Drop Arm (*Weight Capacity: 450 lb*)



Call today for various chair options and configurations.



Winco Elite CareCliner

Design, Comfort, Durability, and Value

Simple concepts that other manufacturers have yet to grasp. These concepts are what sets the Winco Elite CareCliner clinical chair series apart from the rest, and make ours among the best-selling clinical chairs in the market today. Wide elastic band on back frame provides exceptional comfort. Standard Trendelenburg release on back with handle accessible from either side of the chair. Nonporous, blow-molded, fold-down side tray comes standard on both the right and left sides. Elastic open-seat base provides the ultimate in comfort and support.

Call today for various chair options and configurations.





**#204
Table Base**

Option 1

- Pass-through drawers
- Electrical receptacle (*base only*)

Option 2

- Pass-through drawers
- Drawer heater
- Pelvic tilt
- Electrical receptacle (*base only*)

(702-5876)ea



**Treatment Tables with
Adjustable Backrest**

31"H x 72"L

- Straight line table
- Solid wood legs with H-brace

#1010-24, 24"W
(173-7409)ea

#1010-27, 27"W
(122-3886)ea

#1010-30, 30"W
(701-3142)ea



**Alpha Series Treatment
Table with Drawers**

- Features all-steel frame, laminate shelf, 2 drawers with steel sides and adjustable backrest
- Paper dispenser included
- 31"H x 27"W/30"W x 72"L

#3013-27
(122-0010)ea

#3013-30
(122-0011)ea



**Walnut Overbed
Table H-Base**

- Easy spring load height adjustment
- Height 26.75" – 43.25"
- Weight capacity 50 lbs

#11610, Overbed Table
(605-6802)ea



**LUMEX®
Everyday
Overbed
Table,
Non-Tilt**

- Easy to assemble
- Height adjusts infinitely between 28" and 44"
- Twin-wheel swivel casters provide smooth-rolling action
- 2 locking and 2 nonlocking casters
- Walnut, wood-grain laminated ¾" particleboard core with plastic-edge molding
- Chrome plated "H" base and column for aesthetics and durability
- Designed to meet the requirements of HCPCS code E0274
- Supports up to 30 lb of dispersed weight
- 1-year warranty

Overbed Table
(635-0158)ea

STOOLS/ WASTE RECEPTICLES



#2102 Revolving Stools without Backs

- 18½"–25½" height range
- 14" diameter
- Chrome base
- Thread adjustment
- 4"-thick polyfoam cushioning



Stools.....ea

#2102, Chrome/Rubber Casters (122-3678)

Please specify color when ordering: Emerald, Black, Desert Tan, Dove Gray, Mulberry, Burgundy, Royal Blue, or Slate Blue.



#2153 Economy Air-Lift Stool

- Hand-operated air-spring pneumatic height adjustment
- 16"-diameter seat with 3½" thick foam
- Durable 5-leg, black composite base with 2" twin-wheel casters
- 250-lb weight capacity
- 24"D x 20–27"H
- Available in 9 Pro-Form vinyl colors



Black (143-5033).....ea

Specify Color When Ordering (143-0094).....ea



Baked Enamel Steel Step-On Cans

Each can, with epoxy polyester finish, features unitized construction, full-width hinges for long life, heavy-duty foot pedal and linkage system, self-closing lid, and bag-securing mechanism. Plastic protective molding hems bottom. "Biohazard" and "Infectious Waste" adhesive-backed labels included.

#P-24, 24 qt.....ea

Red.....(545-2145) White.....(545-3315)

#P-32, 32 qt.....ea

Red.....(545-3427) White.....(545-5228)

#P-48, 48 qt.....ea

Red.....(545-5926) White.....(545-6276)



Step-On Waste Containers

- Deodorant block holder and tight-fitting, overlapping lid help contain odors
- Smooth contours allow easy cleaning and disinfecting
- Quiet, heavy-duty pedal allows hands-free sanitary use
- 18 gal
- 19¾" W x 26½" H x 16½" D

Red (120-5445).....ea

Beige (120-5926).....ea

Shipped directly from manufacturer. Please allow extra time for delivery.



Rubbermaid® Step-On Containers

For hands-free, sanitary refuse disposal. Approved for fire safety. Tight-fitting, overlapping lid helps contain odor. All cans have rounded corners and smooth contours for easy cleaning. Tough, rugged, all-plastic construction provides years of use.

- Assists in complying with OSHA regulation 29 CFR 1910.1030, "Occupational Exposure to Bloodborne Pathogens"
- Factory Mutual (FM) and California State Fire Marshal (CSFM) approved for fire safety
- Tough, puncture-resistant, all-plastic construction
- Tight-fitting, overlapping lid to help contain odors
- Quiet, heavy-duty pedal for sanitary, hands-free use
- Optional rigid liners for containers #6143, #6144, #6145, and #6146 provide leak-proof bases to help contain liquids

#6143, 8 Gallon, 13¾" D x 15¾" H.....ea
Light Beige.....(707-0542)

Red.....(307-0320)

#6146, 23 Gallon, 19.8" x 16.1" x 32½"ea

Light Beige.....(613-0183)

Red.....(324-5203)





22" Removable-Arm Wheelchair

- 22" Sentra EC with detachable full arms
- Dual cross-brace
- Swing-away footrest (replaces STD22ECDFA)
- Carbon steel frame/black upholstery
- Weight capacity: 450 lb

#STD22ECDFA-SFea
(115-5651)



18" Fixed-Arm Sport Wheelchair

- 18" Silver Sport 2
- Silver vein fixed arm
- Swing-away footrests (replaces 905N)
- Silver vein/black upholstery
- Weight capacity: 300 lb

#SSP218FA-SFea
(116-8651)



Wheelchairs

- Highly polished chrome-plated steel
- Padded armrests
- Embossed upholstery
- Dual axle
- Die-case aluminum footrests
- Adjustable front fork on 16" and 18" only
- Weight capacity: 250 lb for 16" and 18"; 350 lb for 20", 22", and 24"
- 24" x 1" mag solid composite rear wheels with sealed bearing
- 8" x 1" solid front wheels
- 1-year warranty on frame

18" Wide, with Permanent Full Arm and Fixed Footrestsea
(635-4562)

18" Wide, with Permanent Full Arm and Elevated, Detachable Footrestsea
(103-7172)

18" Wide, with Permanent Full Arm and Swingaway Footrestsea
(103-6178)

18" Wide, with Removable Desk-Length Armrests and Swingaway, Detachable Footrestsea
(635-4311)

20" Wide, with Removable Desk-Length Armrests and Swingaway, Detachable Footrestsea
(635-9603)

22" Wide, with Removable Desk-Length Armrests and Swingaway, Detachable Footrestsea
(407-0549)



24" Wide, with Removable Desk-Length Armrests and Swingaway, Detachable Footrestsea
(635-6711)

18" Wide, with Removable Desk-Length Armrests and Elevated Footrestsea
(635-3974)

20" Wide, with Removable Desk-Length Armrests and Elevated Footrestsea
(635-4561)

24" Wide, with Removable Desk-Length Armrests and Elevated Footrestsea
(635-6736)

BLOOD DRAWING CHAIRS



Blood Drawing Chair

- One-piece, easy-to-clean, plastic chair
- Welded, tubular steel frame
- Adjustable-height arms, 27"-33"
- Laminate pivot armrest
- Black

(122-4867)ea

Item shipped directly from manufacturer.
Some assembly required.



Padded Blood Drawing Chairs

- Padded and upholstered seat, back, and armrests
- Adjustable-height arms, 28"-33"
- Double-pivot armrest
- One-piece drawer
- White welded tubular steel chair frame

Desert Tan
(122-8216)ea

Mulberry
(127-8370)ea

Slate Blue
(102-5765)ea

Specify upholstery color.
Item shipped directly from manufacturer.
Some assembly required.



DIALYSIS WEBSITE

What dynamic features does the Henry Schein Dialysis Website offer your facility?

Quick & Easy Ordering!



Login Order Checkout



- 24 hour/365 days availability to supply ordering
- Track orders and shipments and process returns
- Generate comprehensive utilization reports via our dynamic report wizard
- Set budgets and electronically approve orders for greater security and accuracy
- Access purchase history and real-time pricing
- Custom shopping list capability for easy ordering
- Paperless billing statement option
- Shop by procedure, chronic conditions, or simple everyday essentials
- Access free MSDS documents for quick OSHA compliance
- Retrieve electronic statements and invoices for paperless environment
- Access to promotions created uniquely for dialysis facilities

www.henryschein.com/dialysis

accucold
by SUMMIT APPLIANCE

Compact 1.8 Cubic-Foot All-Refrigerator

- Perfectly sized compact to fit tight spaces
- Front lock
- Automatic defrost saves you time and maintenance
- Factory-installed lock provides security you can depend on
- Cold-wall interior for easier cleanup
- 19"H x 18 $\frac{5}{8}$ "W x 17 $\frac{1}{4}$ "D

(117-7318)ea



accucold
by SUMMIT APPLIANCE

Slim 2.4 Cubic-Foot All-Refrigerator

AccuCold by Summit Appliance offers the MEDDT series all-refrigerator for use in medical, scientific, pharmaceutical, and laboratory institutions.

- Digital thermostat allows more precision in temperature control
- Factory-installed lock provides security you can count on
- Simple maintenance with auto-defrost interior
- Internal fan complete with gel packs for stabilized interior
- High temperature alarm includes external temperature readout
- Includes a hospital-grade cord with a green-dot plug
- 26"H x 18 $\frac{1}{2}$ "W x 17 $\frac{3}{4}$ "D

#FF28LWHMEDDT
(116-6631)ea



ABS American BioTech Supply

**Undercounter Refrigerator
4.5 Cubic Foot**

- Standard undercounter refrigerator (*freestanding*)
- 4.5 cubic-foot capacity
- 1 solid swing door, right hinged, not reversible
- Mechanical thermostat
- Temperature range: 0°C to 10°C
- 3 shelves (*1 fixed, 2 adjustable*)
- Keyed door lock
- White cabinet and door
- High density urethane foam cabinet and door insulation
- 2-year parts and labor warranty, plus an additional 3-year compressor parts warranty
- Exterior dimensions: 33 $\frac{1}{2}$ "H x 21 $\frac{1}{2}$ "W x 22 $\frac{1}{2}$ "D

#ABT-UCFS-0504W
(121-9890)ea



accucold
by SUMMIT APPLIANCE

**Counter-Height FF6L Refrigerator
5.5 Cubic Foot**

- Full 5.5 cubic-foot capacity inside a convenient 24" footprint
- Factory-installed lock provides security you can count on
- Automatic defrost ensures minimal user maintenance
- Includes door storage and adjustable shelves
- Keyed lock offers added security
- Fully finished white cabinet
- Reversible door
- 33.5"H x 23.63"W x 23.5"D

#FF6L
(113-1696)ea



Enhance patient safety with products
that help you meet the challenge of controlling cross-contamination

NEW



Welch Allyn CareTemp™ Touch Free Thermometer

- Infrared sensor reads temperature 4-6 cm from center of forehead—No patient skin contact is required supporting a safe and comfortable patient experience
- Uses touch free sensor head to help minimize the risk of cross-contamination
- Measurement feedback system detects proper operating distance to provide consistent accuracy

Welch Allyn FlexiPort® Disposable Blood Pressure Cuffs

- Assign a single FlexiPort Disposable Cuff to each patient to help reduce cross-contamination risks
- Utilize our complete range of color-coded sizes for fast, accurate selection of the correct cuff size to help enhance patient safety

For more information call your local Welch Allyn representative or customer service at 1.800.535.6663 opt. 3

WelchAllyn®

©2015 Welch Allyn MCI2432

HENRY SCHEIN®
DIALYSIS



To Order: **1-800-521-1472** 8am-6pm (et)
To Fax: **1-800-372-5702** 24 Hours
www.henryschein.com/dialysis

Due to manufacturer price increases, prices are subject to change without notification. Offer valid on our full array of products, excluding flu vaccine and Synagis®. Upon any request by Medicare, Medicaid, Tricare or other payer, you must report the value of any benefit received from a discount program (e.g., points redemptions, gifts or other special awards) ©2015 Henry Schein, Inc.
No copying without permission. Not responsible for typographical errors. No adjustments from prior sales.
Offer not to be combined with other promotions and/or special contract pricing.

Void where prohibited. Local regulatory requirements may apply to use or installation of certain products.
It is your responsibility to understand and comply with any such requirements prior to purchase, use or installation of products. †For informational purposes only. Source: ReimbursementCodes.com and/or www.cms.gov.
Customer is responsible for verification of billing/coding in accordance with applicable specific circumstances.

15LS3849

HENRY SCHEIN®

Henry Schein, Inc. • 135 Duryea Road • Melville, NY 11747

PRRST STD
U.S. POSTAGE

PAID

HENRY SCHEIN, INC.