



Renfert USA Fall Promo

CONTACT ZAHN DENTAL AT (800) 496-9500



USE PROMO CODE: KPE

FALL into SAVINGS!

GET CASH BACK ON ALL PURCHASES!

Receive \$100 back for every \$1,000 spent*



NEW!
SILENT powerCAM EC
Powerful extraction unit
for CAM systems

Item Code: 228-5803

~~\$2,429.00~~
SALE PRICE

\$1,998.99



*See last page for details and Rules and Regulations

Valid September 2nd through November 30th, 2019

www.renfert.com | Se habla español

making work easy

CONTACT ZAHN DENTAL AT (800) 496-9500

Renfert

For digital dentistry solutions

SILENT powerCAM EC

Powerful extraction unit for CAM systems

ADVANTAGES

- Integrates seamlessly with a wide variety of heavily used CAM units. Low maintenance EC Motor (> 5,000 operating hours).
- Powerful, bag-free unit with automatic filter cleaning.
- High suction power of 38 m/s and low operating noise of max. 54.3 dB (A).



Item Code: 228-5803

~~\$2,429.00~~

SALE PRICE

\$1998.99

NEW

Item Code: 228-4785

~~\$1,424.00~~

SALE PRICE

\$1,199.99



SILENT compactCAM

Extraction unit for compact CAM systems

ADVANTAGES

- Extraction unit designed for optimized use with small to medium size dental CAM units.
- Active health protection due to mechanical filter cleaning technology.
- Powerful, bag-free unit with automatic filter cleaning.

Renfert-Scanspray

Surface marker

ADVANTAGES

- Ideal for all optical scanner systems and also for scanning impressions.
- Achieves maximum precision thanks to a very thin layer application providing accurate edge representation in 3D models.



Item Code: 228-2324

~~\$100.00~~

SALE PRICE

\$84.99

200 ml (6.80 fl oz.)

Item Code: 228-6944

~~\$22.00~~

SALE PRICE

\$19.49



75 g

GEO Classic

Milling wax

ADVANTAGES

- Scannable wax for digital duplication of wax up.
- No chip formation due to ideal wax hardness.

making work easy

CONTACT ZAHN DENTAL AT (800) 496-9500

Renfert

For finishing of CAD/CAM restorations

Universal polishing paste beige

High-luster polishing paste

- Hard polishing paste for mirror-like surface on conventional acrylics or 3D printed resin appliances.
- Fine grained consistency for high-quality polishing results with brushes or buffs.



Item Code: 228-6325

~~\$72.00~~

SALE PRICE

\$64.99

(box of 6 bars 200 g ea.)

Item Code: 228-0517

~~\$17.00~~

SALE PRICE

\$14.79



35 g

Opal L

High-luster polishing paste

- Achieve high shine on acrylic, PMMA, and other milled polymers. Particularly suitable for light-curing acrylics.
- Optimal results with use of Bison or Goat hair brushes.

lay:art style set

Natural bristle brushes

- Top-quality, natural hair brush. Extremely durable and dimensionally stable.
- Perfect for stain & glaze, micro layering.



Item Code: 228-2434

~~\$63.00~~

SALE PRICE

\$55.79

(set of 3 brushes)

Item Codes: all-in-one: 228-5529
ZrO₂: 228-5530
LiSi: 228-5531

~~\$74.00~~ (each)

SALE PRICE

\$64.99 (each)



10 g

Renfert Polish

Diamond polishing pastes

- Diamond polishing pastes for universal application and specific use for high-gloss polishing.
- Excellent for use with lithium disilicate, zirconium oxide, ceramics and high-performance polymers.

Bison

Polishing brush

- Provides a powerful and quick high-lustre polish on alloys, porcelain and acrylic facings.
- Long service life and excellent take-up of polishing pastes.



Item Code: 228-2694

~~\$126.00~~

SALE PRICE

\$111.99

(18mm - box of 100)



~~\$21.91~~ (each)

SALE PRICE

\$19.99 (each)

Dynex Brilliant Discs

Separating discs

- Fiberglass reinforced, diamond discs created for separating, cutting and trimming veneering ceramic, zirconium oxide and lithium disilicate.
- NEW Brilliant Disc (228-5900) is a specialized grinding disc for all ceramic restorations. Low heat, extremely stable and high removal capacity. .80x20mm, 1 pc.



Item Codes: Dynex Brilliant .20x20mm, 10pcs.: 228-0217
Dynex Brilliant .25x20mm, 10pcs.: 228-2229
Dynex Brilliant .30x20mm, 10pcs.: 228-0808

~~\$54.00~~ (each)

SALE PRICE

\$44.99 (each)

making work easy

CONTACT ZAHN DENTAL AT (800) 496-9500

Renfert

For model production workflows

Twister

Vacuum mixer

- Easy-to-operate vacuum mixer that produces homogeneous mixing results even with the thickest consistency and largest volume of materials.
- Automatic bowl coupling enables convenient, single-handed operation.



Item Code: 228-5266

~~\$1,725.00~~

SALE PRICE

\$1,424.99

(stand sold separately)

Item Code: 228-0451

~~\$641.00~~

SALE PRICE

\$549.99



Vibrax

Vibrator

- Produces a wide range of vibrations to optimally process materials of all types of viscosity.
- No transfer of vibrations to the workbench and extremely quiet.

MT3 Klettfix

Wet trimmer

- High-performance, wet trimmer with a powerful 1-3/4 HP motor & compact design.
- Klettfix version includes quick change trimmer disc.
- Designed for easy use and fast cleaning.



Item Code: 228-2064

~~\$1,014.00~~

SALE PRICE

\$833.99

Item Code: 228-2627

~~\$1,213.00~~

SALE PRICE

\$998.99



Top spin

Pin drilling unit

- Unique design ensures a <0.01 mm drilling accuracy without resistance for all degrees of plaster hardness.
- Automatic start/stop function with easily removable, clear dust drawer.
- Easy, tool-free drill change.

Vortex compact 3L

Compact dust extractor

- Ideal for both wet and dry extraction at a low cost.
- Powerful, lightweight and portable for use in any size lab.



Item Code: 228-1556

~~\$968.00~~

SALE PRICE

\$819.99

making work easy

CONTACT ZAHN DENTAL AT (800) 496-9500

Renfert

For model production workflows

AUTO spin

Automatic pin-drill system

ADVANTAGES

- Automatic pin-drilling system for precise fabrication of models with plastic base plates.
- The Renfert Level Controller design helps achieve safe positioning and alignment in challenging impressions.



Item Code: 228-5804

~~\$1,672.00~~

SALE PRICE

\$1,373.99

NEW

Item Code: 228-3214

~~\$1,068.00~~

SALE PRICE

\$877.99



Millo pro

Dental arch trimmer

ADVANTAGES

- Provides the basis for quick and effortless trimming of dental arches.
- Adjustable cutter height for easy trimming of difficult tooth positions.

Power pillo

Deflasking chisel

ADVANTAGES

- Deflask restorations quickly and easily without the risk of creating stresses, cracks or splinters.
- Quick-release coupling connection.



Item Code: 228-1477

~~\$649.00~~

SALE PRICE

\$534.99

making work easy

CONTACT ZAHN DENTAL AT (800) 496-9500

Renfert

For sandblasting using precision nozzles, optimum viewing and unique mixing technology

Basic eco

Compact sandblaster

ADVANTAGES

- Compact, yet provides maximum field of vision thanks to a large viewing panel with innovative LED lighting.
- Uniform flow of abrasive guarantees consistent surface conditioning with very low consumption.
- Allows for quick switching between abrasives without interrupting work.



Item Code: 228-2682

~~\$796.00~~

SALE PRICE

\$579.99

Item Code: 228-2066

~~\$1,295.00~~

SALE PRICE

\$1,069.99



Basic master

Fine spacious sandblaster

ADVANTAGES

- Precise blasting with the special mixing chamber technology of Renfert's Venturi principle.
- Large volume blasting chamber provides maximum freedom of movement.
- With many extras including ergonomics, ventilation and durability, this unit is perfect for labs with high workloads.

Basic quattro IS

Fine sandblaster with Immediate Stop technology

ADVANTAGES

- Suitable for extremely delicate fissure blasting as a result of the IS – Immediate Stop technology.
- IS technology reduces blasting costs by up to 80%.
- Can be quickly retrofitted for up to 4 tanks, enabling several applications simultaneously.
- Quick and safe operation due to exact pressure adjustment and pressure control for every tank.



Item Code: 228-2068

~~\$2,442.00~~

SALE PRICE

\$2,012.99

(Optional 3rd or 4th tank sold separately)

Item Code: 228-4782

~~\$1,158.00~~

SALE PRICE

\$981.99



SILENT compact

Single station extraction unit

- Features intelligent switch-on function controlled by the integrated appliance.
- Can be easily built-in to work benches due to compact construction.
- Bagless dust collection with mechanical filter cleaning function.

making work easy

CONTACT ZAHN DENTAL AT (800) 496-9500

Renfert

For cleaning dentures and restorations

SYMPRO

Denture cleaning unit

ADVANTAGES

- Automated mechanical removal of surface accretion using a rotary magnetic field and specialized needles - for optimal cleaning of dentures, orthodontic appliances and other dental restorations.
- Up to 80% time savings compared to conventional methods, with virtually no labor costs.
- Easily service denture wearers and in turn, increase customer retention and loyalty.



Item Code: 228-1622

~~\$978.00~~

SALE PRICE

\$803.99

Item Code: 228-1619

~~\$68.00~~

SALE PRICE

\$59.99



2 x 2 L (2 x 0.53 gal)

SYMPROfluid Universal

Cleaning agent

ADVANTAGES

- Effectively aids in removal of deposits in conjunction with SYMPRO device.
- Removes standard accretion such as plaque, calculus, discoloration and adhesive denture creams.

EASYclean

Ultrasonic cleaning unit

ADVANTAGES

- Provides state-of-the-art ultrasonic cleaning technology for all types of ultrasonic cleaning. Includes degas function to optimize cleaning solutions.
- Features a temperature regulated ultrasonic function and automatic cleaning start at a preselected temperature.



Item Code: 228-0405

~~\$952.00~~

SALE PRICE

\$783.99

Item Code: 228-5559

~~\$52.00~~

SALE PRICE

\$45.99



2 L (0.53 gal)

GO 2011 speed

Plaster & alginate solvent

ADVANTAGES

- Efficient removal of stubborn residues such as dental stones and alginates in a cleaning bath and even faster in an ultrasonic unit.
- 50% faster cleaning of different instruments compared to conventional solvents.

Valid September 2nd through November 30th, 2019

making work easy

CONTACT ZAHN DENTAL AT (800) 496-9500

Renfert

For any step in the porcelain modeling & waxing-up process

Brushes

Synthetic & natural

SYNTHETIC BRISTLE BRUSHES

Takanishi set

Item Code: 228-6020 (set of 6 brushes)

- Quality synthetic brush at an affordable price.



~~\$68.00~~
SALE PRICE
\$58.99

Ceramicus set

Item Code: 228-0434 (set of 6 brushes)

- Durable, special microhair structure with elegant, wood-like handle.



~~\$105.00~~
SALE PRICE
\$93.99

NATURAL BRISTLE BRUSHES

Genius set

Item Code: 228-8090 (set of 4 brushes)

- High quality natural bristle with protective handle and replaceable brush inserts.



~~\$157.00~~
SALE PRICE
\$140.99

Kolinsky set

Item Code: 228-3116 (set of 6 brushes)

- Ideal combination of stability and elasticity for porcelain processing with multi-functional handle.



~~\$80.00~~
SALE PRICE
\$70.99

Trays

Mixing trays

lay:art oasis

Item Code: 228-5720



~~\$86.00~~
SALE PRICE
\$75.99

- Continually moistens ceramic materials, keeping them processable for long periods of time.
- Gentle on the brush thanks to a completely glazed surface.

lay:art crystal aqua xl

Item Code: 228-2373



~~\$216.00~~
SALE PRICE
\$193.99

*small version also available

- Gentle on the brush and completely abrasion free using float glass.
- Homogeneous moistening and consistency of the porcelain thanks to dimensionally stable and durable moistening strips.

Waxes

Modeling waxes

GEO Expert Functional

Item Code: 228-5907



~~\$99.00~~
SALE PRICE
\$87.99

- Wax set with step-by-step guide for the production of functional models and restorations based on occlusal compass & NAT technique.
- Easy implementation of functional tooth structures by following wax-up instructions.

GEO Expert A. Bruguera

Item Code: 124-9246



~~\$118.00~~
SALE PRICE
\$104.99

- Esthetic wax-up system w/ step-by-step manual helps patient and dentist envision and evaluate the final restoration.
- Esthetic, true to nature, hyper realistic waxups – Gingival Kit for tissue also available in natural red and pink.

Rules and Regulations:

Renfert equipment, instrument and material purchases must be made between September 2nd through November 30th, 2019. Does not include taxes and shipping costs. Purchases may be made through any participating Renfert USA dealer in the US ONLY.

To redeem your gift card, please download the redemption form at www.renfert.com/usafallpromo2019. You are required to completely fill out and submit ONE redemption form between December 1st-16th, 2019. Invoices must be dated 9/2/19 - 11/30/19. NO EXCEPTIONS. CANNOT BE COMBINED WITH ANY OTHER OFFERS. Please allow 4-6 weeks from the date we confirm receipt of your submitted invoices to receive your gift card. For more info or any questions please call 1-800-336-7422 or email: promos@renfertusa.com.