

LABORATORY GUIDE

2013-2014 EDITION

Henry Schein introduces the 2013/2014 **Physician Office Lab Supply Guide**. This guide is intended to serve as an overview of our most popular diagnostic testing options and supplies needed at the point of care. Don't see what you are looking for, ask your Sales Consultant. From rapid test kits to large analyzers, from CLIA waived to highly complex, we have POL solutions for all your needs.

SIEMENS

TERUMO®

Alere™

BD

UNICO®

OraSure Technologies, Inc.

RPS Right Diagnosis
Right Treatment
Right Now

globe
SCIENTIFIC INC.

ABAXIS

**BECKMAN
COULTER**
We're better together

ALLERAVISION

Alere™ Toxicology

Bayer HealthCare
Diabetes Care

ENTERIX®
A Quest Diagnostics Company

HORIBA
Medical

HEMOSURE® iFOB TEST
-Immunochemical Fecal Occult Blood Test (FIT)

Roche

HealthLink®
A Family of Medical Products
Joining Forces to Serve the Healthcare Community

**CLOROX
HEALTHCARE®**

esa

MAGELLAN BIOSCIENCES

osom

ALFA WASSERMANN
Diagnostic Technologies, LLC

If you have questions regarding diagnostic, pharmaceutical, or vaccine items used in your practice, the DxRx Solutions Hotline is here to help! The DxRx Solutions Hotline can provide CPT Code, package inserts, CLIA product status, product storage requirements, and brand-to-generic comparisons. Let us help you today!

DxRx Solutions Hotline
1-877-523-SHOT 8:30am – 5:00pm, ET
Dxrxsolutions@henryschein.com

TABLE OF CONTENTS

Respiratory	5
Strep A	9-11
Diabetes Management.....	12-18
Kidney Function	19
Urinalysis	20-25
Fecal Occult	26-28
Fecal Occult Immunochemical	29-30
Women's Health	31-34
Infectious Disease Testing.....	35-39
Drugs of Abuse Testing	40-41
Specialty Testing	42-46
Cardiac.....	47

General Chemistry and Lipids	48
Coagulation.....	49
Blood Collection	50-53
Lab Equipment	54
Lab Supplies	55-57
MOHs Lab Solutions.....	58
Cytology/Histology Supplies	60-61
Stains & Reagents	62-63
Plated Media	65
Clinical Chemistry	66-81
Hematology	82-85
Immunoassay	86-111



small package

BIG REVOLUTION

*It will change
the way you view
rapid testing...*



New BD Veritor™ System

CLIA-waived for Rapid Detection of Flu A+B



Helping all people
live healthy lives

- **Proven Performance vs PCR** — PCR is currently the highest sensitivity standard available
- **Reproducibility You Can Count On** — Unique Adaptive Read Technology helps reduce false-positive results
- **Results You Can Trust** — Easy-to-read digital display eliminates subjective interpretation
- **Convenient Workflow** — Testing requires minimal hands-on time

Medical technology that captures the trusted answer

www.bd.com/ds/veritorsystem

BD, BD Logo and Veritor are property of Becton, Dickinson and Company ©2011 BD

BD Diagnostics
7 Loveton Circle
Sparks, MD 21152-0999
www.bd.com/ds

HENRY SCHEIN®
BRAND PRODUCTS

Rely on Us

BRAND PROMISE

QUALITY SELECTION PERFORMANCE



THE ONESTEP LINE OF DIAGNOSTIC TESTS FEATURES SUPERIOR CLINICAL PERFORMANCE, ACCURACY, AND QUICK RESULTS

LOOK WHAT'S BEEN ADDED TO THE HENRY SCHEIN ONESTEP+ PRODUCT LINE
HENRY SCHEIN ONESTEP+ INFLUENZA A & B TEST



OneStep+ Influenza A & B Test


This test provides a simple, rapid method to aid in the differential diagnosis of influenza A & B using nasal swab specimens.

- CLIA Waived
- 10-minute results
- Sensitivity (A/B): 93.8%/77.4%
- Specificity (A/B): 95.8%/98%

(112-6219).....25/box

CPT® Codes: Influenza A: 87804QW†; Influenza B: 87804-59†

Serving you Quality, Selection and Performance.
It's our promise to you.



BD Veritor™ System

A rapid chromatographic immunoassay for use with symptomatic patients. This CLIA-waived kit is to be used with the BD Veritor System Reader to test nasal and nasopharyngeal swab samples that are freshly collected, processed and dispensed directly onto an assay test device.

- Proven performance vs. PCR
- Unique Adaptive Read Technology provides reproducibility you can count on
- Easy-to-read digital display eliminates subjective interpretation and delivers the "Trusted Answer"
- Convenient workflow with minimal hands-on time

#256055, BD Veritor™ System Reader

(992-0003)ea
#256045, CLIA-waived BD Veritor™ System Flu A+B Kit

(992-0002)30 tests/pkg
CPT® Codes: Influenza A: 87804QW† / Influenza B: 87804-59†; Using the Objective and Differentiated Reader System: 87999†



BinaxNOW® Influenza A & B Tests

Rapid test for identification and differentiation of influenza A and B, using specimen collected by nasal washes/aspirates, nasal swabs and nasopharyngeal swabs.

- CLIA waived
- Differentiates types A and B on the same card
- Sensitivity (A/B): 70–89%/50–69%
- Specificity (A/B): 90–99%/94–100%
- Easy to use: no additional reagents needed
- On-board extraction
- 15-minute results
- Same specimen can be used with BinaxNOW® RSV test
- Made in the USA

#416-022
(153-9781)22/box

#416-110
(348-7480)10/box

CPT® Codes: Influenza A: 87804QW†/Influenza B: 87804-59†

HENRY SCHEIN®

ONESTEP+ INFLUENZA A & B TEST

This test provides a simple, rapid method to aid in the differential diagnosis of influenza A & B using nasal swab specimens.

- CLIA Waived
- 10-minute results
- Sensitivity (A/B): 93.8%/77.4%
- Specificity (A/B): 95.8%/98%

#1126219, Influenza A & B Test (112-6219)25/box

CPT® Codes: 87804QW† (Influenza A); 87804-59† (Influenza B)



NEW!



Quidel®

QuickVue® Influenza Tests

An immunoassay for detecting influenza types A and B. Three simple steps with results in 10 minutes or less. Color-coded procedure card helps demonstrate each step of test. Kit contains positive and negative controls, and stores at room temperature.

- Detects both A and B types of influenza
- CLIA waived
- Sensitivity (A/B): 94%/70%
- Specificity (A/B): 90%/97%
- Color-coded procedure card helps demonstrate each step of the test
- No instrumentation needed
- Results in 10 minutes or less, right in the exam room
- Room-temperature storage with up to a 24-month shelf life from date of manufacture
- Made in the USA

#20183, A+B Test, Differentiates Type A or Type B (227-5611)25/box



#00317, Non-differential (227-9012)25/box

#20103, Swab Pack* with 25 Individually Pouched Sterile Foam Swabs (986-3022)25/pkg

#20171, Transport Tubes* with 25 Sterile Foam Swabs Inside Plastic Transport Tubes (227-3550)25/box

CPT® Codes: Influenza A, 87804QW†/Influenza B, 87804-59†



OSOM® Influenza A & B Tests

In vitro diagnostic test for the qualitative detection of influenza A and B antigens.

- Differentiates between influenza A and B
- Sensitivity (A/B): 73.8%/60%
- Specificity (A/B): 96.4%/96.4%
- Easy to read: 3 distinct lines
- Nasal swab specimens
- Less than 1 minute hands-on time; results in 10 minutes
- Room-temperature storage
- Moderately complex
- 2 extra test sticks included for external QC testing
- Made in the USA



#190 (106-3126)27/box

#191, Influenza A & B Control Kit (106-3132)ea


CPT® Codes: Influenza A and B, 87804-59†



The #1 Rapid Tests Used in Hospitals and Physician Offices Nationwide¹

BinaxNOW[®] Influenza A & B

- Packaged with NP Swab, a CDC preferred specimen type²
- Validated for NP Swab, Nasal Swab and Nasal Wash/Aspirate
- Compatible with a dozen common viral transport medias³
- Easy Procedure with On Board extraction
- Enclosed card format reduces pathogen exposure
- 24 month shelf life⁴



CLIA
Waived

BinaxNOW[®] RSV

- Results in 15 minutes with no hands-on time
- Onboard extraction: no extra steps
- Nasopharyngeal swab and nasal wash
- 24-month shelf life



Stock up today!
Contact your sales representative!

Rapid Results, Smarter Decisions

alere.com

1. Data on file 2. CDC website. 3. BinaxNOW Influenza A&B Package Insert 4. From date of manufacture. BinaxNOW and the Alere Logo are trademarks of the Alere group of companies. ©2010 Alere. All rights reserved. 10-LF-0291vB

**BinaxNOW® RSV Tests**

- CLIA waived
- Easy to use: no additional reagents needed
- On-board extraction
- 15-minute results
- Built-in control to ensure accuracy
- Same specimen can be used with BinaxNOW Influenza A & B Test
- **Made in the USA**

#430-100, Test
(301-8899)10/box

#430-122, Test
(512-5059)22/box
CPT® Code: 87807QW†

Clearview® RSV Test

A 15-minute, four-step test for the qualitative detection of RSV. For use with nasopharyngeal swabs and nasal aspirate specimens (approved claim: in children less than 6 years and adults more than 60).

- CLIA waived
- Sensitivity: 93.9%
- Specificity: 97.7%
- Room-temperature storage
- Built-in control to ensure accuracy
- Shelf life up to 12 months from date of manufacture
- **Made in the USA**



#135060
(512-8933)20/box


Contains: 20 test devices, 2 mL extraction reagent solution, 1 mL each positive & negative controls & 40 disposable pipettes.
CPT Code: 87807QW

**BinaxNOW® S. pneumoniae Urinary Antigen Test**

An *in vitro* rapid immunochromatographic assay for the qualitative detection of *S. pneumoniae* antigen in the urine of patients with pneumonia and in the cerebral spinal fluid (CSF) of patients with meningitis. In conjunction with culture and other methods, it is intended to aid in the diagnosis of both pneumococcal pneumonia and pneumococcal meningitis.

- Dependable accuracy
- Accuracy/sensitivity: urine: 86% sensitivity/94% specificity; CSF: 97% sensitivity/99% specificity
- No results negatively affected by previous antibiotic use
- Self-contained, swab-based device—no tubes or pipettes
- Easy to use: insert urine- or CSF-saturated swab into test card, add reagent, and read results in 15 minutes
- Both built-in and external positive and negative controls ensure test viability
- **Made in the USA**

#710-000, 22-Test Kit
(512-0000)ea
CPT® Code: 87449†



QUIDEL®
CORPORATION

QuickVue® RSV Tests

- CLIA waived
- Sensitivity: 92%
- Specificity: 92%
- Easy-to-use dipstick format
- 3-step procedure
- 2-color results in 15 minutes
- Room-temperature storage
- Built-in internal controls
- Kit includes external controls
- **Made in the USA**

#20193
(227-0581)20/box
CPT Code: 87807QW†

**BD Veritor™ System**

A rapid chromatographic immunoassay for use with symptomatic patients. This CLIA-waived kit is to be used with the BD Veritor System Reader to test nasal and nasopharyngeal swab samples that are freshly collected, processed and dispensed directly onto an assay test device.


- Proven performance vs. PCR
- Unique Adaptive Read Technology provides reproducibility you can count on
- Easy-to-read digital display eliminates subjective interpretation and delivers the "Trusted Answer"
- Convenient workflow with minimal hands-on time

#256055, BD Veritor™ System Reader
(992-0003)ea



The BD Veritor™ System RSV Kit is a chromatographic immunoassay with an instrumented read for the direct and qualitative detection of RSV fusion protein from nasopharyngeal wash, aspirate and swab in transport media samples from patients suspected of having a viral respiratory infection. This test is intended for *in vitro* diagnostic use to aid in the diagnosis of RSV infections in infants and pediatric patients under the age of 20 years. Moderately complex.

#256042, BD Veritor™ System RSV Kit
(992-0005)30/box
CPT® Code: 87807†



Start here.

Questions? We have answers. Today's healthcare environment is evolving and often confusing. We understand that you are constantly challenged by this rapidly changing system, which is why we developed **Answers for Healthcare**.

The **OSOM**® rapid diagnostics team will strive to help you meet these goals:

- Improve outcomes
- Create better patient experiences
- Lower costs
- Coordinate care
- Prevention

Answers for Healthcare. Awesome...yes, OSOM®.

osom®

To learn more visit www.osomtests.com or contact your Territory Manager at 800-332-1042.





osom



OSOM® Ultra Strep A Tests

Ultrasensitive strep A test with the ease of use and rapid results your busy practice needs. Studies show the sensitivity of OSOM Ultra Strep A test and the standard single-swab culture method are not statistically different.

- CLIA waived
- Easy to use
- Easy-to-read, 2-color results in minutes
- Sensitivity: 92.6%
- Specificity: 92.8%
- Room-temperature storage
- Includes positive and negative controls
- **Made in the USA**

#147
(392-7389)25/box
#149 (Includes 2 additional tests for external QC testing.)
(102-7617)52/box
CPT® Code: 87880QW†



osom



OSOM® Strep A Tests

- CLIA waived
- Easy-to-read 2-color results
- Results in 5 minutes
- Sensitivity: 96%
- Specificity: 98%
- Room-temperature storage
- Includes positive and negative controls
- Includes 2 additional tests for external QC testing
- **Made in the USA**

#141
(224-7320)52/box
CPT® Code: 87880QW†

HENRY SCHEIN®
BRAND PRODUCTS

Rely on Us

BRAND PROMISE
QUALITY SELECTION PERFORMANCE

THE ONESTEP LINE OF DIAGNOSTIC TESTS FEATURES SUPERIOR CLINICAL PERFORMANCE, ACCURACY, AND QUICK RESULTS

ONESTEP + STREP A CASSETTE KIT

- CLIA waived
- Easy-to-use flip design; all-in-one extraction, incubation, and test
- Results in 5–10 minutes
- Sensitivity of 96.2% and specificity of 98.7%
- External negative and positive controls included

(900-7570)25 tests/kit
CPT® Code: 87880QW†



NEW!

ONESTEP+ STREP A TEST STRIPS

- CLIA waived
- Room temperature storage
- Canister packaging for ease of multiple testing
- Walk-away procedure

(107-4050)25/box
CPT® Code: 87880QW†



Serving you Quality, Selection and Performance.
It's our promise to you.

BRAND PROMISE

QUALITY SELECTION PERFORMANCE



NOW, MORE THAN EVER, RELY ON HENRY SCHEIN BRAND PRODUCTS

Save without compromising on excellence with our own high-quality Henry Schein Brand.

The **One Step** line

of diagnostic tests

features superior

clinical performance,

accuracy, and

quick results

Item#	Description	Size
900-4071	One Step+ H Pylori Test	20/box
900-4070	One Step+ Mono Test	20/box
900-4074	One Step+ Pregnancy hCG Urine Cassette Test	100/box
900-4073	One Step+ Pregnancy hCG Urine Cassette Test	25/box
900-4075	One Step+ Pregnancy hCG Combo Test	25/box
900-4072	One Step+ Pregnancy hCG Urine Strip Dip	25/box
107-4050	One Step+ Strep-A Test Dip Stick	25/box
112-6219	One Step+ Influenza A & B Test	25/box

Serving you Quality, Selection and Performance.

It's our promise to you.



ICON® SC Strep A Kits

- Easy-to-use cassette with on-board extraction.
- CLIA waived
- Easy-to-read, 5-minute results
- Accuracy: 98%
- Sensitivity: 96.2%
- Specificity: 98.7%
- Made in the USA**

#395161A
(331-6584)25/box

Contains: 25 test cassettes, sterile rayon swabs, vial of reagent A, vial of reagent B, positive control & negative control.
CPT® Code: 87880QW†



QuickVue In-Line® One-Step Strep A Tests

True one-step test with a unique in-line extraction and rapid procedure that allows testing to occur right at the point of care. Two-color results give you clear readability, as well as proven accuracy and dependable results. Kit stores at room temperature and includes positive and negative external controls.

- Built-in controls for CLIA compliance
- Easy to use and easy to read
- 2-color, rapid results in 5 minutes
- Sensitivity: 92%
- Specificity: 99%
- Room-temperature storage
- Made in the USA**



#00343
(227-4207)25/box
CPT® Code: 87880QW†



ICON® DS Strep A Kits

- Dipstick-format test.
- CLIA waived
- Easy-to-read, 5-minute results
- Accuracy: 98%
- Sensitivity: 96.2%
- Specificity: 98.7%
- Made in the USA**

#395098A
(331-1419)25/kit

Contains: 25 test strips, extraction reagent A, extraction reagent B, positive control, negative control, extraction tubes & throat swabs.
CPT® Code: 87880QW†



QuickVue® Dipstick Strep A Tests

- CLIA waived
- Results in 5 minutes or less
- 2-color endpoint: results are easy to read and interpret
- Sensitivity: 92%
- Specificity: 98%
- Overall agreement: 95%
- Room-temperature storage
- Built-in positive and negative controls
- Made in the USA**

#20108
(227-0611)50/box
CPT® Code: 87880QW†



QuickVue+® Strep A Tests

- Built-in controls for CLIA compliance
- Easy test procedure
- ± endpoint
- Built-in workstation
- Sensitivity: 95%
- Specificity: 98%
- Overall agreement: 98%
- Room-temperature storage
- Moderately complex
- Made in the USA**

#20122
(489-0075)25/box
CPT® Code: 87880†



culture plates & media



Taxo™ Differentiation Discs

- #231552, A Discs
For differentiating group A streptococci.
(224-1056)50/box
- #231042, A Discs
For differentiating group A streptococci.
(693-1886)50/box
- #231562, 1-mg Kanamycin Discs
For differentiating gram-negative anaerobic bacilli.
(987-1965)50/box
- #231041, A Disc Cartridges
For differentiating group A streptococci.
(693-8964)6/box



Selective Strep Culture Plates

#221935
(120-4787)10/pkg



SSA Selective Strep AGAR Plate

#221779
(603-4472)20/pkg



Plated Media

- #1146, Strep Select Agar SSA
Plates, 100 mm
(437-3477)10/pkg
- #1160, TSA 5% Sheep Blood, 100 mm
(437-5777)10/pkg



Bacitracin Differentiation Disc Strep A

#2851HL, Disc
(437-0018)ea



Taxo A Differentiation Disc

#2852, Disc
(146-6565)ea



HENRY SCHEIN®

BRAND PRODUCTS



Rely on Us

BRAND PROMISE

QUALITY SELECTION PERFORMANCE



TRUEtest™

TRUETEST™ TEST STRIPS

For use with the TRUEresult® Blood Glucose Monitoring Systems. TRUEtest™ Test Strips use Quad-Electrode Laser Accuracy for enhanced precision and accuracy. They also use the patented TRUEfill® beveled tip, which allows for precise entry of blood sample.

(112-6827)50/box



Coming Soon!



TRUEresult®

TRUERESULT® BLOOD GLUCOSE MONITORING SYSTEM

- A no-coding system with high performance features
- Clinically proven accurate
- Lifetime warranty
- Only requires a tiny, 0.5-mL blood sample
- Provides accurate results in as fast as 4 seconds

Features audible fill detection, easy-to-read display, 500-test memory, a weekday display, 7-, 14-, and 30-day averaging, glucose control detection, 4 testing reminders, ketone testing reminder, a strip release button, auto on/off, and data management capabilities.

(112-6826)ea

Contains: TRUEresult® meter, installed 3-volt battery, lancing device, 10 lancets, 10 TRUEtest™ strips, instructions for use, self-test log book & compact carrying case.

NOW, MORE THAN EVER, RELY ON HENRY SCHEIN BRAND PRODUCTS

Serving you Quality, Selection and Performance.

It's our promise to you.



CONTOUR® Blood Glucose Meter

Cleared for multiple-site testing, which means more choice and flexibility about where you test: finger, forearm, palm, abdomen, or thigh—it's up to you. However, be sure to talk with your health care professional first to ensure that alternate site testing is right for you.

- No coding required
- Automatic underfill detection
- Automatic control marking
- Automatic temperature control
- Small sample size (0.6 µL)
- Quick testing time—results in 15 seconds
- Saves information for up to 240 tests with time and date



#7151, Meter
(147-1405)ea

#7099, Strips
(147-9303)50/vial

#7109, Control-Normal, 2.5 mL
(147-0874)ea

#7111, Control-High, 2.5 mL
(147-8283)ea

#7110, Control-Low, 2.5 mL
(147-5579)ea

CPT® Code: 82962QW



BREEZE®2 Meter

This meter makes testing simple. Its 10-test disc eliminates the struggle with small, individual test strips, and the use of No Coding™ technology ensures easy accuracy every time you test.

- Results in just 5 seconds
- Small sample size
- Large, easy-to-read display
- Ergonomic design is easy to handle
- Single-function buttons for easy handling
- Improved 420-test memory and 1-, 7-, 14-, and 30-day averages

#1440, Meter Kit

(795-0077)ea

Contains: 1 meter, 1–10-strip cartridge, 1 lancing device, 10 lancets & 1 carrying case.

BREEZE® was the first meter to receive the "Ease of Use" commendation from the Arthritis Foundation.



BREEZE®2 Test Strip Disc Cartridges

- Less chance of dropping or losing test strips
- Simply load and go

#1467

(795-7966)50/pkg

Contains: 10 strips per cartridge, 5 cartridges per box.



Henry Schein DxRx Hotline

The Henry Schein DxRx Hotline is a resource to assist medical professionals with various types of questions on diagnostics, pharmaceuticals, and vaccines that are purchased from Henry Schein. The hotline provides information on the following topics for diagnostics, pharmaceuticals, and vaccines:

- Indications and usage as provided on a package insert
- Package inserts
- Brand to generic comparisons
- CPT Codes and J Codes*
- ACIP recommendations for vaccinations
- State Law Requirements for School and Day Care Entry**
- Latest updates on supply issues and much more

The Henry Schein DxRx Hotline can be reached Monday – Friday, 8:30AM to 5:00PM EST at **1-877-523-SHOT** (7468), via E-mail at DxRxSolutions@henryschein.com, or via live chat with a representative on our Web site, www.henryschein.com/dxrxsolutions.



*For informational purposes only.

Source: www.ReimbursementCodes.com and/or www.cms.gov

Customer is responsible for verification of billing/coding in accordance with applicable specific circumstances.

**For informational purposes only.

Does not constitute legal advice.

Source: www.immunize.org



OneTouch® Ultra® 2 Blood Glucose Monitor

- Very small blood sample
- More choices for testing: 3 sites, including finger and arm
- FastDraw™ design automatically draws blood into test strip
- Fastest test time available: just 5 seconds

#02109800

(567-0390).....ea



OneTouch® Ultra® Blue Strips

Calibrated for use with a normal control solution.

#020994

(567-6191).....25/box

02289601

(567-0091).....50/box

#02289501

(567-0092).....100/box

CPT® Code: 82962QW†

OneTouch® SureSoft™ Lancing Devices

Single-use lancing devices are electropolished, which results in more blood than an acid-washed needle. This means you can use a smaller gauge needle while still ensuring that a sufficient sample will be drawn for point-of-care testing.

- Electropolished needle means a smaller gauge and more comfort for your patients
- Meets OSHA guidelines: nonreusable, retractable needle
- Distinctive wing design prevents reloading and reusing
- Preset gauge for no guessing, no maneuvering and easy training
- On-demand training through Web and CD-ROM programs

Gentle, 28-ga Needle

(567-0000).....200/box

Regular, 21-ga x 2-mm Needle

(115-5267).....200/box



AgaMatrix PRESTO™ PRO Blood Glucose Monitoring System

Key features:

- No code
- 0.5 µL sample size
- Test results in As Fast As 1-2-3®
- Backlight
- Alternate site testing
- 300-test memory
- 14-, 30-, and 90-day averages
- On-screen graphing
- 6 user-settable alarms
- Hypo/hyperglycemic warning alarms
- **Made in the USA**

PRESTO™ PRO System Kit (Includes built-in test strip ejector for safely removing used test strips.)

#8000-04088 Meter

(476-0003).....ea

Contains: 1 of each: meter, carrying case & owner's guide.

PRESTO™ System Kit

#8000-02649

(476-6255).....ea

Contains: 1 meter, 1 carrying case, 1 owner's guide, 1 lancing device, 10 lancets, 10 test strips & 1 bottle of Control solution.

PRESTO™ Test Strips

#8000-03329

(476-7212).....50/btl

#8000-03337

(476-2624).....100/btl



Call for more options!

TruTol® Glucose Tolerance

Screen for hyperglycemia (diabetes mellitus and gestational diabetes) and hypoglycemia (low blood sugar). Volume graduations are imprinted on the label to aid in dispensing. Pediatric dose chart lists precalculated dosage based on body weight.

Orange, 50 g

(122-0183).....ea

Fruit Punch, 50 g

(917-1244).....24/case

Orange, 100 g

(122-0184).....ea

Lemon-Lime, 100 g

(275-4299).....24/case

Lemon-Lime, 50 g

(486-4583).....24/case

TRUEbalance Test Kit

- Patented TRUEfill design that promotes first-test success and accurate results
- High performance, low acquisition cost and no-coding technology with a small blood sample size and fast test results

H4H0181, Starter Kit

(890-0056).....ea

H3H01-81, Test Strips

(890-0055).....50/box



NIPRO DIAGNOSTICS™



TRUE2go™ Blood Glucose Meter Kit

- Easy to use: simply twist the TRUE2go on the top of a new vial of strips
- Test results within 4 seconds
- No coding: eliminates a step that can lead to inaccurate results
- Convenient: the smallest meter in the world

The world's smallest meter system makes blood glucose testing easy and affordable. Just twist TRUE2go™ on to the top of a new vial of TRUEtest strips and you're ready to go.

Glucose Meter Kit

(890-0020).....ea

TRUEtest™ Test Strips

(890-0034).....50/box



NIPRO DIAGNOSTICS™



TRUEresult™ Blood Glucose Meter Kit

Proven performance and accuracy at a great value.

- No coding: eliminates a step that can lead to inaccurate results
- Results within 4 seconds: means virtually no wait time
- 0.5-mL sample size: less blood means less painful testing
- 500-test memory: includes time/date and 7-day, 14-day, and 30-day averaging
- Affordable: no-coding convenience at an affordable price

Glucose Meter Kit

(890-0019).....ea

TRUEtest™ Test Strips

(890-0034).....50/box



Glucose 201 Analyzer

The HemoCue Glucose Systems are based on a glucose-dehydrogenase method and consist of a small, dedicated analyzer and a unique, disposable microcuvette. They will produce lab-accurate results for diagnosing or monitoring diabetes, as well as monitoring neonatal blood glucose levels.

- More accurate and 68% higher reimbursement than glucose meters
- FDA-compliant unit when screening and diagnosing
- Analyzer: 0–444 mg/dL (0–24.6 mmol/L)
- Only 5 µL of capillary, venous, or arterial blood is needed
- No inherent limitations related to hematocrit or known drug interferences
- No recalibration needed between cuvette lots
- 2-year warranty

#120706, Analyzerea
(104-2290)

#110706, Microcuvette, Glucose (439-9158)100/box

#110402, Microcuvette, Glucose (609-2005)100/case



Hb 201+ Handheld Hemoglobin Analyzer

Offers the same important benefits as its predecessor, the B-Hemoglobin Photometer. Smaller size gives more flexibility and mobility to use the system anytime and anywhere. Performs an internal self-test. No control cuvette required.

- Produces results in less than 1 minute with lab-quality accuracy
- 3 easy steps
- Precalibrated, portable analyzer that requires minimal maintenance
- No calibration or instrument manipulation (between cuvette/reagent batches) is needed
- Safe and easy to use
- Reimbursable
- CLIA '88 waived
- 2-year warranty

Analyzer (345-2344)ea
Contains: analyzer, AC adapter, manual & training CD-ROM.

#3201 Microcuvettes (110-9167)100/box
Contains: 4 packs of 25.

(345-8733)200/box
Contains: 4 packs of 50.



Accutrend® Plus Meter Kit

- CLIA waived
- Stores 100 results for each parameter (glucose and cholesterol)
- Out-of-meter dosing
- Onboard QC
- Positive strip lot ID
- Fast, reliable results
- Portable
- Room-temperature storage
- Made in the USA

Complete Kit #53467544160 (891-1038)ea

Contains: 1 meter, 4–1.5-V AAA alkaline batteries, 1 user's manual, 1 quick-reference guide, 1 carrying case & 1 warranty card.



Cholesterol Test Strips #5213312160
2.9% low sensitivity.
173–224 mg/dL range.
(891-1036)25/pkg
Contains: 25 strips per vial & 1 code strip.



Glucose Test Strips #11447475160 (891-0008)25/pkg
Contains: 25 strips per vial & 1 code strip.

Cholesterol Control 2 Level #5219957001 (891-1037)ea
Contains: 1 of ea 4-mL cholesterol control level: low & high.



Glucose Control 2 Level #5213231160 (891-1039)ea
Contains: 1 of ea 4-mL glucose control level: 2.5% low & 2.4% high.
CPT® Code: 82465QW†



CardioChek™ P+ A

Perform on-the-spot testing for cholesterol and lipid levels with 1 convenient, portable, and accurate whole-blood test system.

- CLIA waived
- Easy to use
- Single-testing and panel-testing capability
- Lipid profile results from a single finger-stick whole blood sample

CardioChek P+ A (133-2551)ea

Glucose Test Strips (261-5308)25/box

Total Cholesterol Test Strips (412-5021)25 tests

HDL Test Strips (107-4803)25 tests

HDL Cholesterol Controls, High & Low (261-5338)ea

Triglyceride Test Strips (872-5107)25 tests

Lipid Panel Test Strips (261-7537)15/box

Multi Chemistry Control (261-8361)ea



Alere™ HemoPoint® H2 System

Using a single drop of blood, the HemoPoint® H2 System delivers accurate hemoglobin and hematocrit results in less than a minute. The touch screen is easy to learn and operate, minimizing training time. Rechargeable battery can last 100 hours between charges. Recalls the last 4000 test results with the touch of a button.

#55118, HemoPoint® H2 Meter (117-9599)ea

Contains: HemoPoint® H2 Photometer, user guide, quick reference guide, optics cleaner, control cuvette & power adapter.

#80549, HemoPoint® H2 Microcuvettes (117-9605)100/box

Contains: 2 containers of 50 cuvettes each.

#88772, Hemoglobin Controls (117-9656)ea

Contains: 1–1.5 mL Low & 1–1.5 mL High.

#52219, HemoPoint™ H2 Optics Cleaner (118-5057)ea

Contains: 5 cleaning units.



Chemcard™ Glucose Test

A simple, effective screening test that is used to identify abnormally high fasting plasma glucose (blood sugar) levels which are often a warning sign of diabetes.

- 3-minute whole-blood screening test
- No instrumentation required
- Room temperature storage, not to exceed 86°F (30°C)

#47024, Glucose Test (844-1247)24/box

Chemcard™ Cholesterol Test

- 3-minute whole-blood screening test
- Color interpretation from 150 mg/dL to 300 mg/dL in increments of 25
- No instrumentation required
- Room-temperature storage, not to exceed 86°F (30°C).

(100-7110)24/box





Afinion™

Simple ■ Fast ■ Reliable

Afinion™ HbA1c Test

- CLIA-waived - no special training required
- Easy sample collection and operation
 - Sample size: 1.5 µL of blood
 - Sample types: capillary and venous whole blood
- 3 minute test time
 - Increased efficiency
- Reliable, lab quality results
 - Precision: <3% CV
 - NGSP certified
 - No interferences from hemoglobin variants

Now available through Henry Schein!

For more information, call your Henry Schein Sales Consultant at **800.772.4346**



1.877.441.7440 | alere.com

Connectivity

The Afinion™ Data Connectivity Converter (ADCC) allows for simple transfer of patient and control results from the Afinion™ AS100 Analyzer to LIS/HIS system of choice.

- Results are transferred via Ethernet connection
- Available standard formats of HL7 or ASTM

Other Available Tests

- Afinion™ ACR





Alere

Afinion™ Test System

The Afinion™ AS100 Analyzer provides valuable near-patient testing for today's health-care professionals. It utilizes the latest technology and provides a simple, fast, and reliable method to monitor patients' disease progression.

- Multi-assay analyzer system
- 3-minute HbA1c, CLIA-waived
- 5-minute albumin/creatinine ratio (ACR)
- Short test time provides the opportunity for more timely treatment changes
- Easy sample collection and operation, minimal training required
- Reliable, lab-quality results
- No manual calibration required
- Unique integrated fail-safe system
- Allows for input of patient and operator IDs

AS100 Analyzer
(119-8787)ea

Cleaning Kit
(119-8810)50/box

Universal Printer
(119-8807)ea

Thermal Printer Labels
(119-8813)ea

Barcode Scanner
(119-8889)ea



Alere

Afinion™ HbA1c Test

Test for quantitative determination of glycated hemoglobin (HbA1c) in human whole blood, used to monitor metabolic control in patients with diabetes.

- 3-minute test time
- Excellent precision (<3% CV)
- No interference from Hb variants
- NGSP and IFCC certified
- CLIA-waived

HbA1c
(119-8788)15/box

HbA1c Controls
(119-8791)ea



Alere

Afinion™ ACR Test

Test for quantitative determination of albumin, creatinine, and albumin/creatinine ratio (ACR) in human urine, used for early identification of renal disease in patients with diabetes and/or hypertension.

- Lab-quality results in 5 minutes
- 3.5 µL of random/spot urine sample
- Quantitative results
- Increased accuracy by measuring creatinine

ACR
(119-8789)15/box

ACR Controls
(119-8793)ea



SIEMENS

DCA™ Vantage Analyzer

Point-of-care immunoassay analyzer for diabetes management. Provides precise and accurate quantitative test results for HbA1c in whole blood during a patient's visit to help improve decision making. With a single urine sample, the analyzer also provides quantitative test results for low levels of albumin, creatinine, and A:C ratio in the physician's office or clinic.

- #5075US, Analyzerea
(110-0743)
 - #5035C, HbA1c Reagent Kit
(147-8045)10/pkg
 - #5068A, Hemoglobin Like Control Kit
(147-6067)ea
- Contains: 2 vials of each control: high & low.

Microalbumin/Creatinine Reagent Kits
Detects and monitors microalbuminuria, an early indication of kidney disease. Achieve test results in just 7 minutes from a urine specimen, using the DCA™ Vantage Analyzer. Existing DCA analyzers may be updated with a new program card.

- **Made in the USA.**
 - #6011A
(147-8917)10/pkg
- CPT® Codes: 83036QW†/82043†/82570†

DCA Microalbumin/Creatinine Control Kit

- #6012A
(147-8477)ea
- Contains: 2 vials of each control (high & low) & reconstitution fluid.

Can be used with Bayer A1C or Siemens DCA

Liquid Whole Blood hbA1C Controls

- #NOD13122-100
4-Vial Kit, 1 mL
2 Level 1 and 2 Level 2
(111-9856)ea

Liquid Whole Blood hbA1C Controls

- #NOD13111-100
2-Vial Kit, 1 mL
1 Level 1 and 1 Level 2
(111-9857)ea



Monitor diabetes treatment with unparalleled precision.

The DCA Vantage™ analyzer provides HBA1c results in minutes.

Answers for life.






SIEMENS

Multistix® 10LS Pro Reagent Strips

More reliable tool for detecting kidney disease. Test for glucose, ketone, high protein, low protein, specific gravity, creatinine, pH, blood, nitrite, and leukocytes.

- Creatinine ratio to correct for varying urine concentrations
- Decreases false positive protein and false negative protein rates
- For over-65 patients and those with diabetes, CVD, or UTI
- Indication of a dilute requires patient to submit specimen, triggering another urine sample
- Automatic identification of strip configuration for Clinitek® 50 users
- **Made in the USA**



#1566
(147-8028)25/btl

#1554
(147-2360)100/btl

CPT® Codes:
81000†/81001†/81002†/81003†/82570QW†



Chemstrip® Micral Test Strips

Test for microalbumin. Fast, easy test to detect subclinical neuropathy.

- 1-minute test
- In office
- Accurate compared to radioimmunoassay

#417146
(891-6088)30/btl
CPT® Code: #82044QW†




SIEMENS

Clinitek® Microalbumin Test Strips

Test for semiquantitative determination of microalbumin/creatinine in urine. Semiquantitative tests for protein (2–4 mg/dl) and creatinine. Includes an A:C ratio result.

- CLIA waived on CT status analyzer
- Dip-and-read test for albumin and creatinine in urine
- Determines an A:C ratio for greatest accuracy
- Use a random urine sample—reduces need for a 24-hour collection
- Timeframe: 60 seconds
- **Made in the USA**

#2083
(147-4712)25/btl
CPT® Codes: 82044†/82570†



Microalbumin/Creatinine Reagent Kits

Detects and monitors microalbuminuria, an early indication of kidney disease. Achieve test results in just 7 minutes from a urine specimen, using the DCA™ Vantage Analyzer. Existing DCA analyzers may be updated with a new program card.

• **Made in the USA.**
#6011A
(147-8917)10/pkg
CPT® Codes: 83036QW†/82043†/82570†



HemoCue® Albumin 201 System

Can be used for quantitative determination of low levels of albumin in urine for the purpose of screening, diagnosis, monitoring, and supporting clinical evidence in the treatment of microalbuminuria. The system consists of a small, specially designed analyzer and individually packaged microcuvettes containing freeze-dried reagent. The test is performed at the point-of-care, using a microcuvette filled by capillary action. To obtain a quantitative result, fill the microcuvette, place it in the analyzer, and read the result. The system is based on an immunoturbidimetric reaction using antihuman albumin antibodies. The antigen-antibody complex creates turbidity, which is measured photometrically at 610 nm. The albumin concentration is proportional to the turbidity, and the result is shown on the display in mg/L.

Analyzer Kit with 1 Box of Microcuvettes
(111-5593)ea
Albumin Microcuvettes
(111-5591)50/box

Roche

Urisys 1100™ Urine Analyzer

This CLIA-waived urine analyzer makes urine testing easier. With the simple press of a button, analyzer automatically evaluates urine test strips and prints a patient report. This easy-to-use system is an ideal choice for physician offices looking to standardize test results and improve office efficiency through automated urinalysis testing. Uses Chemstrip® 7, 10 MD, and 5 OB strips.

#3617556, Analyzer
(891-7039)ea



Roche

Chemstrip® 10 MD Test Strips

Intended to be used for visual testing or for automated testing on the Chemstrip® 101 urine analyzer and Criterion II urine analyzer. Measures blood, protein, leukocytes, nitrate, pH, glucose, ketones, bilirubin, urobilinogen, and specific gravity.

#3260763
(891-4205)100/btl

Roche

Chemstrip® 10 SG Test Strips

For visual or automated testing on the Chemstrip 101 urine analyzer and Criterion II urine analyzer. Measures leukocytes, nitrite, pH, protein, glucose, ketone, urobilinogen, bilirubin, blood, and specific gravity.

#417145
(891-8002)100/btl



Roche

Chemstrip® 9 Test Strips

Test for leukocytes, nitrite, pH, protein, glucose, ketone, urobilinogen, bilirubin, and blood.

(891-4564)100/btl



Roche

Chemstrip® 5 OB Urine Test Strips

#1893467
(891-4294)100/btl



Roche

Chemstrip® GP Test Strips

Test for glucose and protein.

#200743
(891-5552)100/btl

Roche

Chemstrip® 7 Test Strips

Test for leukocytes, pH, protein, glucose, ketone, nitrite, and blood.

#1008552
(891-1914)100/btl





URISCAN

Urine Test Strips

- DOCTORS' CHOICE SINCE 1995
- CLINICALLY PROVEN QUALITY
- TEST FOR LESS



HS ITEM NO.

URISCAN URINALYSIS TEST STRIPS

106-9405	1 KETONE STRIP (Ketones)
106-9418	1 GLUCOSE STRIP (Glucose)
257-5782	2 GP STRIP (Glucose, Protein)
106-9421	3 GPH STRIP (Glucose, Protein, pH)
106-9425	4 HemoGPH STRIP (Blood, Glucose, Protein, pH)
106-9428	5 HemoKetoGPH STRIP (Blood, Ketones, Glucose, Protein, pH)
106-9432	6L STRIP (Blood, Protein, Nitrite, Glucose, Leucocytes, pH)
257-0142	10 SGL STRIP (Blood, Bilirubin, Urobilinogen, Ketones, Protein, Nitrite, Glucose, pH, Specific Gravity, Leucocytes)

Manufactured For:



Los Angeles, CA
1.800.365.6146
www.biosyslabs.com

Uri-Check® Plus Test Kits

A semiquantitative dip-slide culture system for the screening and early detection of urinary tract infections.

More than 100,000 bacteria per mL indicates infection.

1½" diameter with patient label affixed to vial. Refer to Uri-Check reference chart to help identify bacterial growth present. 3–5-month shelf life.

Nonreturnable. Requires Boekel's Micro Incubator #260700 (100-6303).

#3000UC, CLED Agar/EMB Agar

(224-9604)10/box



Uricult Test Kits

A simple, self-contained urine test method for the detection of bacteria that commonly cause urinary tract infections. Uricult urine-culture paddles are capable of detecting bacteriuria concentrations as low as 1000 CFU/mL of urine.

#1000, CLED/EMB

(120-9110)10/box

Contains: paddles, patient ID stickers & package insert.



TROY
Biologicals

Uri-Check® Test Kits

A semiquantitative dip-slide culture test for the detection of urinary tract infections. Kit includes 1" vials and caps, paddles slightly smaller than 1", and patient labels affixed to the outside of the vial.

#1000UC, CLED Agar/EMB Agar

(224-7373)10/box



Accutest® Uriscree™ Tests

- Rapid urinary-tract infection screen test for bacteriuria and presence of somatic cells in urine
- Detects true infections with high sensitivity and high NPV
- Positive screening result indicates presence of bacteria and/or somatic cells

(106-9335)20/box

CPT® Code: 87899QW†





Precision Uroflowmetry System by

HENRY SCHEIN®
MEDICAL

What is the Flo-Med 1100 by Henry Schein?

It is a simple, non-invasive diagnostic screening device used to calculate the flow rate of urine over time. This aids in diagnosis of bladder and sphincter function as well as prostate conditions.

- Easy to use (hygienic & convenient)
- Height adjustable
- Battery operated
- Wireless communication to laptop iPad, or other PDA (hardware capable)
- CPT 51741†
- Compact/portable
- One-touch technology
- Easy calibration
- PDF can easily be stored into an EHR

System Requirements:

Windows 7, XP, Vista
Hardwire compatible
with PCL5 Printer

(112-5811)

Flo-Med Precision Uroflowmetry System,
1 Unit

(112-5812)

Flo-Med Cup & Funnel Disposables,
24/Sets



Other accessories such as Laptops
and printers are available separately.



†For informational purposes only. Source: ReimbursementCodes.com and/or www.cms.gov. Customer is responsible for verification of billing/coding in accordance with applicable specific circumstances.

View Demo Video at:

<http://vimeopro.com/user8610595/flo-med-1100-demo>



Call Your Sales Consultant Today For More Information!



**Reads
11-Way,
5-Way Plus,
4-Way & GP+A
Test Strips**

URISPEC® PLUS URINALYSIS READER & URISPEC® 11-WAY URINE STRIP PROMOTION

Provides dependable urine status results to detect early stages of many diseases, such as diabetes, kidney disease, and urinary tract infections. The comprehensive interface options and the optimized printouts minimize risk associated with manual transcriptions. Reliable results can be obtained virtually immediately at the point of care.

- CLIA-waived
- Readings eliminate subjectivity of visual color interpretation, which makes urinalysis simpler and more reliable
- Features no-touch operation: After strip placement, measurement starts automatically, which makes urinalysis faster and more hygienic
- Easy switch between different test strips: Urispec Plus is for use with Henry Schein Brand Urispec 11-Way, 5-Way Plus, 4-Way, and GP+A urine test strips

(900-7507)ea

Contains: 1 Urispec Plus Urinalysis Reader & 12 bottles of Urispec 11-Way Urine Strips.



URISPEC® 5-WAY PLUS URINE REAGENT TEST STRIPS

- Tests for glucose, blood, protein, nitrite, leukocytes, and urobilinogen
- Strips for visual read or for use with the Henry Schein Urispec Plus Urinalysis Reader only

(900-7508)100/btl

HENRY SCHEIN®

BRAND PRODUCTS

URISPEC™ 11-WAY URINE REAGENT TEST STRIPS

Urine Reagent Test Strips

- Unique "read all" feature at 60 seconds
- Similar to Multistix® 10 SG Plus Ascorbic Acid Test and all brands of 10 SG products
- Glucose, bilirubin, ketone, specific gravity, blood, pH, protein, urobilinogen, nitrite, and leukocytes
- Reader available
- Ascorbic acid pad to ensure uncompromised results
- Room-temperature storage

(102-2285)100/btl

CPT® Codes: 81000†/81001†/81002†/81003†



URISPEC™ GP+A URINE REAGENT TEST STRIPS

Test for glucose and protein, with an ascorbic acid test for accuracy.

- Protein and glucose
- Ideal for metabolic, kidney, and diabetic monitoring
- Dip and read in 60 seconds

(100-6567)100/btl



URISPEC™ 4-WAY URINE REAGENT STRIPS

- Protein, nitrite, glucose, and leukocytes
- Ideal for OB/GYN
- Infection, diabetes, and kidney function
- Dip and read in 60 seconds

(107-7157)100/btl

**Urispec Dipper (112-6841)
& Dropper (112-6842)
Controls now available!**



URISPEC™ PRO 2 GP URINE STRIPS

Strips for Visual Read

Tests for glucose and protein.

(900-4996)100/btl



URISPEC™ 5 OB URINE STRIPS

Strips for Visual Read

Tests for glucose, protein, blood, nitrites, and leukocytes.

(900-4995)100/btl



URISPEC™ 10SG URINE STRIPS

Strips for Visual Read

Tests for glucose, bilirubin, ketone, urobilinogen, specific gravity, blood, pH, protein, nitrites, and leukocytes.

(900-4994)100/btl



SIEMENS

Siemens Adds Value— Every Step of the Way.



Innovation

Bundled connectivity featuring RAPIDComm® Data Management and Connectivity Solution, DCA Vantage® and CLINITEK Status® Connect Analyzers, minimal preparation for tests, self-contained reagents, walk-away convenience, market leading position, and first to market innovation with connectivity

Quality

Product exchange warranty, 24/7 call center, technical support, extended service agreements, internal electronic checks and built in error messages in analyzers, QC management and operator management for security, audit trail data



Training

Proficiency testing, Personalized Education Plan (PEP) web site, telephone installation and training, webinars, patient education material, key opinion leaders, field support, video training, and procedure manuals

Service

Technical support, reimbursement information, 24/7 call center, customer service specialists, government affairs team, and territory sales representatives in market

Contact your local sales representative or visit us at:
www.usa.siemens.com/diagnostics

Order No. A91DX-POC-110347-UC1-4A00 | 03-2012 | All rights reserved
© 2012 Siemens Healthcare Diagnostics Inc.

Siemens Adds Value

www.usa.siemens.com/diagnostics

Answers for life.



SIEMENS

Clinitek Status+ Urine Analyzer

Provides convenient point-of-care testing with improved data integrity for more confidence in the results. Simple, fast, and accurate testing produces a result in about 60 seconds. Status analyzer can process basic urine strips, ratio tests, and hCG cassette tests. Full alphanumeric keyboard, easy-to-use touch-screen display, and automatic strip identification make this unit convenient to use. Eliminates need for visual reading, timing errors, and transcription errors. Memory to store 950 patient records. Future-ready for EMR when you are.

#1780
(795-0055)ea



SIEMENS

Clinitek® Advantus™ Analyzer

With Bar-Code Reader

For laboratories or practices requiring medium to high throughput (up to 500 tests per hour). Respond to demands for higher productivity and high quality. Streamline your workflow with flexible operation. Immediate start-up. Automatic calibration, network ready, and offers the widest range of test options. Added QC features.

- Easy-to-use color touch-screen
- Ability to consolidate microscopy results onto one patient record
- QC lockout available or automatic prompting if regimen is not followed
- Testing flexibility with simple, one-touch switching between different strip configurations
- Three sieve functions—confirmatory and microscopic—to flag abnormal results for quick review
- Convenient data retrieval for up to 500 patients by patient ID
- 200 control results
- Streamlined operation with immediate start-up
- Improve efficiency with nonpacing processing
- Improve accuracy with automatic calibration for every strip
- Continued productivity boost with throughput of 7 seconds per sample
- Minimize errors with simple, bar-code identification of sample, color, and clarity values
- Flexibility in reporting options, on-screen display, printout and/or uploading to laboratory information system

(108-6007)ea

it's
back!

SIEMENS

Clinitek® hCG Pregnancy Tests

hCG test cassette for use with the CLINITEK Status® Analyzer. Quick, reliable, instrument-read results eliminate interpretation and provide accurate results for your patient.

- CLIA waived
- 99% accuracy
- Sensitivity equal to or greater than 25 mIU/mL
- Time frame: 5 minutes
- Room-temperature storage

(186-4896)25/box

CPT® Code: 84703QW†



SIEMENS

Multistix® 10LS Pro Reagent Strips

More reliable tool for detecting kidney disease. Test for glucose, ketone, high protein, low protein, specific gravity, creatinine, pH, blood, nitrite, and leukocytes.

- Creatinine ratio to correct for varying urine concentrations
- Decreases false positive protein and false negative protein rates
- For over-65 patients and those with diabetes, CVD, or UTI
- Indication of a dilute requires patient to submit specimen, triggering another urine sample
- Automatic identification of strip configuration for Clinitek® 50 users
- **Made in the USA**

#1554
(147-2360)100/btl

CPT® Codes:
81000†/81001†/81002†/81003†/82570QW†

Multistix® Reagent Strips

Test for pH, protein, glucose, ketone, bilirubin, blood, and urobilinogen in urine.

#2179
(147-5635)100/btl

SIEMENS

Multistix® 5 Reagent Strips

- Includes glucose, protein, blood nitrite, and leukocyte
- Easy, accurate, and reliable results in 60 seconds

#2309
(147-0768)100/btl

Multistix® 9 SG Reagent Strips

Test for glucose, bilirubin, ketone, specific gravity, blood, pH, protein, nitrite, and leukocytes in urine.

#2163
(147-5776)100/btl



SIEMENS

Clinitek® Microalbumin Test Strips

Test for semiquantitative determination of microalbumin/creatinine in urine. Semiquantitative tests for protein (2–4 mg/dl) and creatinine. Includes an A:C ratio result.

- CLIA waived on CT status analyzer
- Dip-and-read test for albumin and creatinine in urine
- Determines an A:C ratio for greatest accuracy
- Use a random urine sample—reduces need for a 24-hour collection
- Timeframe: 60 seconds
- **Made in the USA**

#2083
(147-4712)25/btl
CPT® Codes: 82044†/82570†



SIEMENS

Multistix® 10 SG Reagent Strips

Test for glucose, bilirubin, ketone, specific gravity, blood, pH, protein, urobilinogen, nitrite, and leukocytes in urine.

#2161
(147-5913)100/box
CPT® Codes: 81000†/81001†/81002†/81003†



Albustix® Strips

Test for proteinuria.

(147-1625)100/btl

Combistix® Reagent Strips

Test for pH, glucose, and protein in urine.

(147-9259)100/btl

Hemastix®

Test for blood in urine.

(147-1005)50/btl

Labstix®

Test for pH, protein, glucose, ketone, and occult blood.

#2181
(147-6828)100/btl
CPT® Code: Urinalysis #81000†/81002†



SIEMENS

Uristix® 4 Reagent Strips

Test for glucose, protein, nitrite, and leukocytes in urine.

#2166
(147-4912)100/btl
CPT® Codes: Urinalysis #81000†/81002†

Uristix®

Test for glucose and protein in urine.

#2184
(147-9199)100/btl




When **quality** is at the **core** of every decision you make...

Trust Hemoccult ICT—the only Immunochemical fecal occult blood test clinically tested over an extended duration and proven to work. Our dry-slide process has provided excellent proven results for over 40 years.

Excellent 14-day stability—developed with dry-slide technology, which increases stability of the sample—allowing for multiple-day collection and necessary transit/development time.

Higher sensitivity and specificity—clinically proven to have a relative sensitivity of 100% and apparent specificity of 99.1%.



No drug or dietary restrictions—Hemoccult ICT is specific to intact human hemoglobin—eliminating interactions from red meat, Vitamin C and other factors that lead to false positives.

The right choice—meeting the American Cancer Society and the U.S. Task Force for Colorectal Cancer (CRC) Foundation recommendation of multiple-day testing for CRC.

For information on developing a colorectal cancer screening program in your office, contact your Henry Schein representative today.

Blood Banking
Immunoassay

Centrifugation
Information Systems

Chemistry
Lab Automation

Flow Cytometry
Molecular Diagnostics

Hematology
Diagnostics

Hemostasis
Rapid Diagnostics



Beckman Coulter, the stylized logo and Hemoccult are trademarks of Beckman Coulter, Inc. and are registered in the USPTO.
© Copyright 2011 Beckman Coulter, Inc.



HENRY SCHEIN®

BRAND PRODUCTS



ONESTEP™ OCCULT BLOOD SINGLE SLIDES

- Detects hemoglobin through chemical reaction
 - Compares to high-sensitivity products
 - Room-temperature storage
 - 100–2-hole slides and 100 sticks
 - **Made in the USA**
- (147-1246)100/box



ONESTEP™ TRIPLE SLIDE PATIENT MAILERS

Effective and rapid method for detecting occult blood in stool. Aids in diagnosis of asymptomatic gastrointestinal problems such as colorectal cancer, ulcers, polyps, anemia, and diverticulosis.

- Detects hemoglobin through chemical reaction
 - Compares to high-sensitivity products
 - Room-temperature storage
 - **Made in the USA**
- (147-7182)50 kits/box

Contains: 50 triple 2-hole slides, 150 sticks, developer & 50 US postal mailers.

CPT® Code: #82270



Hemocult II® Dispensapak™

#61100, 100-Patient Kit
(254-5061)ea

Contains: 100 patient kits, 6–15 mL bottles developer & 1 product instruction.
Kit contains: 1 Hemocult II triple-test card, 3 applicators & 1 illustrated diet/test instructions.
Barrier mailing pouch may be purchased separately.
CPT® Co2de: 82272QW



Hemocult II® SENA® Dispensapak™ Plus

The fecal occult blood test that provides improved readability versus Hemocult and other standard FOBT products.

- Easier readability
- Higher clinical sensitivity
- Detects more true positives
- **Made in the USA**

#64130A, 40-Patient Kit
(254-4454)ea

Contains: 40 patient kits, 2–15-mL bottles developer & 1 product instruction manual.
Patient kit: 1 Hemocult II® SENA triple-test card, 3 flushable sample collection tissues, 3 applicators, 1 illustrated diet & test instructions, 1 patient pamphlet & 1 Barrier® mailing pouch.

CPT® Code: 82272QW

#61200A, Hemocult II Slides
(224-9819)102/box



Hemocult® Single-Test Card

#60151A
(254-9052)ea

Contains: 100 Hemocult single-test cards, 2–15-mL bottles of developer, 100 applicators & 1 product instruction manual.
CPT® Code: 82270QW



Hemocult II® Dispensapak™ Plus

Most clinically proven and widely used FOBT. Trusted for over 30 years.

#61130A, 40-Patient Kit
(254-5748)ea

Contains: 40 patient kits, 2–15-mL bottles developer & 1 product instruction manual.
Kit contains: 1 Hemocult II triple-test card, 3 flushable sample collection tissues, 3 applicators, 1 illustrated diet & test instructions, 1 patient pamphlet & 1 Barrier mailing pouch.
CPT® Code: 82270QW



Hemocult II® SENA® Elite Dispensapak™

- Easier diet and patient instructions
- Easier sample collection
- Easier processing
- **Made in the USA**

#395035A
(331-2723)40/box
CPT® Code: 82270QW



Hemocult® Tape Dispenser Kit

#63202
(254-8235)2/box

Contains: Hemocult tape dispenser (100 tests/dispenser), 1–15-mL bottle developer & product instructions.

Hemocult® Barrier Mailing Pouch

Designed to meet US Postal Service regulations.

#62200
(254-6204)100/box

Hemocult® Developer

Most clinically proven and widely used FOBT. Trusted for over 30 years.

#62115, 15-mL Bottle
(254-2360)ea



Hemocult® SENA® Single-Test Card

#64151A
(254-3896)ea

Contains: 100 Hemocult SENA single-test cards, 2–15-mL bottles developer, 100 applicators & 1 product instruction manual.
CPT® Code: 82270QW

Hemocult® SENA® Developer

#64115, 15 mL
(254-8308)ea



Colon Cancer Kills.

Are You Screening Your Patients?

Choose the Leader in Colorectal Cancer Screening.



HEMOSURE[®], INC.

Immunochemical FOB test (FIT)

- Clinically Proven for Hospitals and Physicians
- High Patient Compliance
- Cost Effective
- Higher Sensitivity & Specificity
- No Drug or Diet Restrictions
- Greatly Reduces False Positives & False Negatives.



HEMOSURE[®] iFOB Test (T1-CK30)*



Hemosure[®] is a Worldwide Recognized Brand known for Quality, Accuracy and Reliability. We are dedicated to saving lives and educating people on the importance of getting screened for Colorectal Cancer, one of America's top cancer killers.

Hemosure[®] GPO Contracts:

- Amerinet: Contract #VL10050
- Yankee Alliance: Contract #YA-LA-012
- Novation: Contract #LB01228
- HPG: Contract #HPG-5829

(patient pamphlet)



(screening poster)



FREE CRC Awareness Materials

CRC Awareness Materials includes:
2 CRC Screening Posters
20 Patient Pamphlets

Request your CRC Materials today!
Visit www.hemosure.com
or e-mail us at
sales@hemosure.com
and use promo code:

HSPOL2013

Place your order today by contacting your Henry Schein representative or call 888-HEMOSURE (436-6787) or visit our website at www.hemosure.com

Henry Schein Product Number	Hemosure [®] Product Number	Description
108-5284	T1-CK30T	IFOB 30 TESTS AND 30 SPECIMEN COLLECTION TUBES
792-8076	T1-CK30*	IFOB 30 TESTS AND 50 COLLECTION TUBES/HOME KITS
109-3145	T1-CK50	IFOB 50 TESTS AND 50 SPECIMEN COLLECTION TUBES



QuickVue® iFOB Test

- Results in 5 to 10 minutes
- Analytical sensitivity:
50 µg hHb/mL buffer or
50 µg hHb/g feces
- Specificity: specific to human hemoglobin
- Physician office lab (POL) analytical accuracy: 98%
- Specimen collection tube doubles as convenient, sanitary dispenser
- Self-contained test cassette
- Bold burgundy end points for ease in reading and interpreting results
- Room-temperature storage with 24-month shelf life from date of manufacture
- External controls available separately
- **Made in the USA**



#20194, Complete Test Kit
(950-2866)20/box

#20201, Tray Packs
(227-2458)50/box

#20197, Control
(223-4668)ea

CPT® Codes: 82274QW†/G0328QW†



Hemosure® iFOB Tests

- Best-value package (792-8076): 40% more collection tubes and home kits
- 5X higher reimbursement (ROI)
- No dietary or drug restrictions
- Specific to human hemoglobin
- Higher patient compliance, due to ease of use
- Easy-to-read results, similar to pregnancy test
- User-friendly and more sanitary
- Complete with home kit
- CLIA waived

#T1-CK10, Complete with Home Kit
(792-5630)10/box

#T1-CK30, Complete with Home Kit
(792-8076)30/box



#T1-CK30T, Complete Kit
(108-5284)30/box

#T1-CK50, Complete Kit
(109-3145)50/box

#T1-TC01, Control
(108-1373)ea

CPT® Codes: 82274QW†/G0328W†



Hemocult® ICT Tests

Next-generation fecal occult blood test (FOBT) for colorectal cancer screening. Specific for human hemoglobin, this immunochemical test eliminates all diet and drug restrictions. Specific for lower GI bleeding, resulting in fewer false positives and a reduction in unnecessary invasive procedures.

- CLIA waived; easy test procedure
- **Made in the USA**
- Screening reimbursement code: G0328QW

#395065A, Collection Card Set
(834-7610)ea

Contains: 100 cards, 100 applicator sticks, physician instructions & sample collection instructions.

#395066A, Screening Kit
(834-3249)ea

Contains: 1 box of 40 patient screening kits. Each patient screening kit contains: 1 dispensing envelope with patient instructions, 3 Hemocult ICT collection cards, 3 flushable collection tissues, 3 applicator sticks & 1 mailing pouch.

#395067A, Test Kit
(834-6389)ea

Contains: 20 test devices & 1-1-mL bottle Hemocult ICT buffer (to be used with Hemocult ICT collection cards & patient screening kits).

#395068A, Controls
(834-5620)4/box

#395261A, 2-Day Patient Screening Kit
(237-0085)50/box

Contains: 1 dispensing envelope, 2 Hemocult ICT collection cards, 2 flushable collection tissues, 2 applicator sticks & 1 mailing pouch.



HENRY SCHEIN® BRAND PRODUCTS

ONESTEP® iFOBT OCCULT BLOOD MONOCLONAL ANTIBODY KITS

Innovative dry-base collection method that utilizes 1 sampling card for the collection of 2 fecal specimens, and requires only 1 immunochemical specimen preparation tube/test cassette system for analysis.

- Meets the American Cancer Society and the US Task Force for Colorectal Cancer (CRC) Foundation recommendation of multiple-day testing
- The sampling tab technology replaces unpleasant liquid buffer tubes and eliminates human error due to patient oversampling
- No dietary or drug restrictions
- 30-day stability/24-month shelf life from date of manufacture
- CLIA waived
- **Made in the USA**



Cassettes
(900-7458)30/box
Contains: 30 test cassettes & 30 buffer tubes.

25-Test Complete Kit
(900-4413)25 tests/kit
Contains: 25 test cassettes, 25 buffer tubes & 25 patient return mailers.

20 Patient Mailer Kits
(900-4416)20/pkg

Contains: 20 postal-approved envelopes, each containing 1 sampling card slide, 2 applicators & patient instructions.

25-Test Complete Kit
CPT® Codes: 82274QW†/G0328QW†

Making FIT Happen

The Inconvenience and Mess of Fecal Occult Blood Tests Are Over

InSure® FIT™ Fecal Immunochemical Test for Early Detection

InSure® FIT™ is the sensible take-home multiple sample method made simple, quick, and easy for you and your patients. It features the only patented, water-based collection method with two samples on one test card meeting American Cancer Society Guidelines* for Colorectal Cancer Screening.

Unlike other FOBTs that require fecal handling and test-card smears, InSure® FIT™ is the one test with no mess. With a gentle brushing of the stool in water, InSure® FIT™ provides 87% sensitivity and 98% specificity for colorectal cancer. What's more, this unique, convenient method promotes high patient compliance.

Help Your Patients Get FIT and Stay Healthy

For a limited time, InSure® FIT™ Sample and Educational Packages are available for you to trial FREE. To learn more about how you can make FIT happen, contact your InSure® FIT™ Representative, visit www.insuretest.com, or call 1.800.531.3681.

*The American Cancer Society recommends for people over 50 a multiple sample method for take home FOBT or FIT screening test; FOBT or FIT testing done during a digital rectal exam in the doctor's office is not recommended.

Manufactured by Enterix Inc., Edison, NJ 08837 USA, a Quest Diagnostics Company. Quest, Quest Diagnostics, the associated logo and all associated Quest Diagnostics marks are the trademarks of Quest Diagnostics. Copyright © 2013, Enterix Inc. All rights reserved. InSure® FIT™ is a trademark of Enterix Inc. For in vitro diagnostic use. ENTAd01/2013

InSure® FIT™ Fecal Immunochemical Tests (FIT)

Test provides high sensitivity and ease of use with a simple, water-based collection method designed to be more patient-friendly. The patented brush sampling replaces unpleasant fecal collection or smearing. No dietary or medicinal restrictions, enhanced specificity and sensitivity for more productive colonoscopies.

- 87% sensitivity for colorectal cancer*
- 98% specificity*
- 66% increase in patient compliance vs. a guaiac-based fecal occult blood test (gFOBT) was found in one study*
- No handling/smearing of dry, solid fecal materials*

• Made in the USA

#90030, Combo Kit 25/25
(450-0009)ea

Contains: 25 test strips, 1–9-mL bottle of run buffer & 25 patient collection kits.

#90035, Combo Kit 25/35
(450-0008)ea

Contains: 25 test strips, 1–9-mL bottle of run buffer & 35 patient collection kits.

#90010, Patient Collection, 10-Pack
(450-0003)ea

Contains: Each envelope contains: 1–2-day sample card, 2 brushes & specimen collection instructions.



#90025, Patient Collection, 25-Pack
(450-0002)ea

Contains: Each envelope contains: 1–2-day sample card, 2 brushes & specimen collection instructions.

#80025, Developer Kit
(450-0011)ea

Contains: 25 test strips & 1–9-mL bottle of run buffer.

#11018, QC Pack, 50 Test Cards
(450-0012)ea

Contains: 50–2-day sample test cards.

#16800, FOBT Controls
(450-0010)ea

Contains: 2–1.5-mL bottles, 1 of each: positive & negative.

*Data on file.
CPT® Code: 82274QW†;
Screening: G0328QW†



HENRY SCHEIN®
ONESTEP+ PREGNANCY STRIP TEST

- Urine test
 - CLIA waived
 - 25 mIU/mL detection level
 - 99% sensitivity
 - Room-temperature storage
 - Read time 3 minutes or less
- (900-4072)25/box
CPT® Code: 81025QW† (Urine)

HENRY SCHEIN®
BRAND PRODUCTS



HENRY SCHEIN®
ONESTEP+ PREGNANCY COMBO CASSETTE TEST

- Urine or serum test
 - CLIA waived for urine only
 - 25 mIU/mL detection level
 - 99% sensitivity
 - Room-temperature storage
- (900-4075)25/box
CPT® Codes: 81025QW† (Urine); 84703† (Serum)



HENRY SCHEIN®
ONESTEP+ PREGNANCY CASSETTE TESTS

- Urine test
 - CLIA waived
 - 25 mIU/mL detection level
 - 99% sensitivity
 - Room-temperature storage
- 25-Patient Kit
(900-4073)ea
100-Patient Kit
(900-4074)ea
CPT® Code: 81025QW† (Urine)



osom

OSOM® hCG Urine Test

- CLIA waived
- Dipstick: simple, 1-step procedure
- Fast results: as soon as 1 minute
- Sensitivity: 25 mIU/mL
- 2-color results
- Room-temperature storage
- More than 99% accurate in clinical studies
- **Made in the USA**

#101, Test
(236-2470)50/box
CPT® Code: 81025QW†



osom

OSOM® hCG Combo Test

- Urine/Serum Pregnancy Test
 - CLIA waived for urine only
 - Easy to read: black-and-white results
 - Fast: one step; no reagents
 - Sensitivity: 10 mIU/mL (serum), 20 mIU/mL (urine)
 - Room-temperature storage; 24-month shelf life
 - Easy to run
 - **Made in the USA**
- #124, Combo Test
(327-9522)25/box
CPT® Codes: 81025QW† (Urine); 84703† (Serum)

osom

OSOM® Card hCG Pregnancy Test

- Traditional physician-office hCG test with the added benefit of easy-to-read, 2-color results.
 - Add 3 drops of urine and read result.
 - CLIA waived
 - Easy to run: no-drop counting pipette
 - Results in as soon as 1 minute
 - Sensitivity: 25 mIU/mL (urine)
 - More than 99% accurate in clinical studies
 - **Made in the USA**
- #102, Test
(327-0743)25/box
CPT® Code: 81025QW†
#134, Urine Control
(327-6476)ea
CPT® Code: 81025†



Alere™ hCG Combo Cassette Test (20/10 mIU/mL)

#92215, hCG Combo Cassette Test
(512-0014)40/box
CPT® Code: hCG 81025QW† (Urine); 84703† (Serum hCG)
Combo tests are classified as moderate when testing serum samples.



- CLIA-waived (for urine only), easy-to-use test
- Cut-off level: 20 mIU/mL in urine; 10 mIU/mL in serum
- Single-step test provides positive results in as little as 3 minutes for urine and 5 minutes for serum
- Highly sensitive assay with built-in control ensures accuracy
- Test can be used with either urine or serum samples
- Sensitivity: >99%
- Specificity: >99%
- Shelf life: 24 months from date of manufacture



Alere™ hCG Dipstick Test



- CLIA-waived, easy-to-use lateral flow test
 - Cut-off level: 25 mIU/mL
 - Single-step test provides positive results in as little as 3 minutes
 - Sensitivity: >99%
 - Specificity: >99%
 - Shelf life: 24 months from date of manufacture
- #92211, hCG Dipstick Test
(512-0012)50/box
CPT® Code: 81025QW†



Alere™ hCG Cassette Test (25mIU/mL)

- CLIA-waived, easy-to-use lateral flow urine test
 - Cut-off level: 25 mIU/mL
 - Single-step test provides positive results in as little as 3 minutes
 - Sensitivity: >99%
 - Specificity: >99%
 - Shelf life: 24 months from date of manufacture
- #92217, hCG Cassette Test
(512-0013)40/box
CPT® Code: 80125QW†



RAPID DIAGNOSTICS

www.sekisuidiagnostics.com



*Which one
is easier to
read?*

(600X MAGNIFICATION)



OSOM® Trichomonas Rapid Dipstick
(ACTUAL SIZE)

No More Wet Mounts...

OSOM® Trichomonas Rapid Test

Are you still relying on wet mount to
diagnose Trichomoniasis?

The OSOM® Trichomonas Rapid Test is a rapid antigen
test with 83% sensitivity versus culture and wet mount
combined. The easy-to-read, two-color results appear
within 10 minutes, and it's the only CLIA-waived rapid
test for trichomonas available.

For more information, call us at 800 332 1042 or visit us
at www.sekisuidiagnostics.com. No more wet mounts,
awesome, yes, OSOM.



SEKISUI
DIAGNOSTICS

THE AMERICAS
Sekisui Diagnostics, LLC
Phone: 800 332 1042
Fax: 800 762 6311
Email: questions@sekisuidiagnostics.com
ISO 13485:2003 Certified

INTERNATIONAL
Sekisui Diagnostics (UK) Limited
Phone: +44 (0) 1732 220022
Fax: +44 (0) 1732 220024
Email: info@sekisuidiagnostics.com
ISO 9001:2000 Certified, ISO 13485:2003 Certified

OSOM® is a registered trademark of Sekisui Diagnostics, LLC

1265-2 06/11



ICON® 25 hCG Pregnancy Test

A one-step, rapid, combination serum/urine test designed for the early detection of pregnancy.

- CLIA waived for urine
- 3-minute results for urine, 5-minute results for serum
- Sensitivity: 25 mIU/mL
- Specificity: >99%
- Requires just 3 drops of specimen

#395097A
(254-2726)25/box



ICON® 20 hCG

A rapid test for detecting hCG qualitatively.

#395097A
(331-0431)25/box



QuickVue® One-Step hCG Urine Test

- CLIA waived for urine only
- Results in 3 minutes or less
- 25 mIU/mL sensitivity
- Accuracy: 99%
- Sensitivity: 99%
- Specificity: 99%
- Easy to use
- Built-in controls for CLIA-waived compliance
- Made in the USA

#20109
(227-6481)25/box
CPT® Code: 81025QW†



QuickVue+® One-Step hCG Combo Tests

- 20 mIU sensitivity in urine, 10 mIU sensitivity in serum
- Results in 5 minutes for serum, 3 minutes for urine (CLIA waived urine only)
- 99% sensitivity, 99% specificity
- Room-temperature storage
- Built-in controls for CLIA compliance
- One-step ± endpoint

#0178
(371-6739)30/box

#0179
(371-1811)90/box

CPT Codes: 81025† (Urine); 84703† (Serum)
(Combo tests are classified as moderate when testing serum samples.)



QuickVue® One-Step hCG Combo Tests and Control Sets

- CLIA waived for urine only
- Results in 5 minutes (serum), 3 minutes (urine)
- 25 mIU/mL sensitivity (urine), 10 mIU sensitivity (serum)
- Sensitivity: 99%
- Specificity: 99%
- Built-in controls for CLIA compliance
- Made in the USA

#20110 Combo Test
(227-7244)50/box

#00272, Control Set (Urine)
(227-4849)ea

#00281, Control Set (Serum)
(227-5659)ea

CPT® Codes: 81025† (Urine hCG); 84703† (Serum hCG)
(Combo tests are classified as moderate when testing serum samples.)



OSOM® BVBLUE® Test

For *Gardnerella vaginalis*

The only CLIA-waived rapid test for detecting bacterial vaginosis, including gardnerella, bacteroides, prevotella, and mobilincus.

- CLIA waived
- Instant color change provides clear, easy-to-read results
- Sensitivity: 92.8% vs. Gram stain
- More sensitive and specific than wet mount
- One-minute hands-on time
- Room-temperature storage
- Made in the USA

#183, Test
(414-4170)25/box

#184, Control
(327-4117)ea

CPT® Code: 87905QW†



OSOM® Trichomonas Rapid Test

The only CLIA-waived rapid test for trichomonas infection. Detects antigen directly from vaginal swab.

- More sensitive than wet mount
- Easy-to-run dipstick
- Easy-to-read, 2-color results
- Results in 10 minutes or less
- Sensitivity: 83% vs. culture
- Agreement: 95% with culture and wet mount combined
- Included in CDC and ACOG guidelines
- Includes QC controls
- Made in the USA

#181, Test
(406-7631)25/box

#182, Control
(611-1157)ea

CPT® Code: 87808QW†



VS-SENSE™ Test

A rapid acidity test for the evaluation of vaginitis (bacterial vaginosis and trichomonas). This 3-second, user-friendly test accurately detects mixed vaginal infections, which allows for immediate treatment, resulting in health care provider and patient satisfaction.

- Easy to use and operate
- Immediate results
- No need for a color scale
- No need to send samples to a lab and wait for results
- Sensitivity: 86.32%
- Specificity: 93.43%
- 3–5 second timeframe
- Room-temperature storage

#FPHSCSVS25
(411-0001)25 swabs/box
CPT Code: 83986QW†



BD Affirm™ VPIII Microbial Identification System



The Problem:

Mixed vaginal infections.

One of every four vaginal infections are caused by two or more organisms: *Candida* species, *Gardnerella vaginalis* or *Trichomonas vaginalis*.

Failure to diagnose these common causes of vaginitis can carry significant risk for your patients. Trichomoniasis and bacterial vaginosis (BV) are associated with increased risk of HIV infection, delivering preterm, low birth weight babies and, if untreated, PID.^{1,2}

¹ Trichomoniasis, CDC, Div. STD. Fact Sheets Sept 2000.

² Bacterial Vaginosis, CDC, Div. STD. Fact Sheets Sept 2000.

³ Data from package insert.

BD, BD Logo and Affirm are trademarks of Becton, Dickinson and Company. ©2009 BD

The Solution:

1 molecular test that detects & differentiates 3 organisms.

The BD Affirm™ VPIII system quickly identifies and differentiates three pathogens from one single test. The Affirm VPIII technology provides RNA accuracy to enable targeted therapies.

- Results in < 45 minutes
- Enables proper therapy, the first visit
- Unaffected by difficult specimens: menses, sperm or creams from self medication.³

Learn more about Affirm VPIII testing in your laboratory, call 800.553.3470.



Helping all people
live healthy lives

BD Diagnostics
800.638.8663
www.bd.com/ds

CYTOLOGY



Affirm™ VPIII Kit

Easy-to-perform, moderately complex DNA probe technology. Conclusive identification of *Candida* species, *Gardnerella vaginalis*, and *Trichomonas vaginalis* from a single vaginal sample.

- Accuracy up to 30% higher than slide microscopy
 - Less than 2 minutes hands-on time required
 - Easy-to-read results
 - Built-in procedural controls for each sample
- (987-7212)24/box



pHizatest™ Nitrazine Indicator Paper

Compare
to NitraTest™

- Economical
- Instant, easy-to-read results
- Reliable: manufacturer has made nitrazine-indicator paper for more than 50 years
- Made in the USA

(122-0181)ea
CPT® Code 83986†



Cytology Brush

#60317, 2.3 x 240 x 2.0
(889-0701)10/case

Cytology Wood Scraper

Non-sterile
#2173
(113-5208)5000/case

Spray-Cyte® Cytological Fixative

1½ oz
(168-9697)ea

Aerocell™

3½-oz Cytology Fixative Spray
(100-0370)ea



Histobrush®

Cytology Brush
(100-8560)100/bag



Reagents Solutions

#400420, Acetic Acid 3%, 16 oz
(740-8432)ea
#400450, Acetic Acid 5%, 16 oz
(666-5674)ea
#400567, Trichloroacetic Acid 80%, 15 mL
(877-9174)ea



Histology/Cytology Solutions and Stains

#400631, Cytology Fixative, 4-oz Pump Bottle
(437-7476)ea
#400600, Eosin Y Stain 1% (Alcoholic), 16-oz Bottle
(113-4902)ea
#400604, Formalin 10% Buffered, 1-gal Bottle
(108-2287)ea
#400605, Formalin 10% Buffered, 5-gal Bottle
(108-2289)ea





herpes

biokitHSV-2 Rapid Test

For Herpes Simplex Virus 2

A type-specific glycoprotein G (gG-2) immunoassay kit—a leading aid in the diagnosis of genital herpes. Two-reagent system.

- Results in less than 10 minutes
- Finger-stick blood or serum specimens
- Built-in procedural codes
- **Made in the USA**



#3000-28001
(630-3080)20/box
CPT Code: 86696†

syphilis



RPR Card Tests for Syphilis

RPR (rapid plasma reagin) macroscopic nontreponemal flocculation test kits for the qualitative and semiquantitative detection of reagin antibodies in human serum. The macroscopic RPR carbon antigen enhances the visual discrimination between reactive and nonreactive results. Kits have a shelf life of up to 24 months from date of manufacture.

- Easy-to-interpret results available in 8 minutes
- Includes controls
- **Made in the USA**

Syphilis Test Cards

(653-4966)100/box
(117-8951)500/box
(115-5717)5000/box

Contains: Each kit contains RPR carbon antigen, reactive control, weak reactive control, nonreactive control, a dropping bottle, a 20-g dispensing needle, 10- or 30-well test cards & disposable stirrer pipettes.

c. difficile



OSOM® C. Difficile Toxin A/B Test

- Detects both toxins: A & B
- Sensitivity: 88.2%
- Specificity: 96.8%
- Just 4 simple steps
- Includes 2 additional test sticks for external QC testing
- Less than 1 minute hands-on time; results in 20 minutes or less
- **Made in the USA**




#173
(309-0085)27/box
CPT® Code: 87324†

chlamydia



QuickVue® Chlamydia Test

- Results in 12 minutes or less
- Simple extraction procedure requires only 2 reagents
- Everything needed is in the kit
- Sensitivity: 92%
- Specificity: 99%
- Room-temperature storage
- Built-in controls, external controls are included
- **Made in the USA**



#B006
(227-0527)25/box
CPT® Code: 87810†



Impact RPR® Tests

For Syphilis

Rapid detection of syphilis in serum or plasma.

#60C1, Tests
(512-1386)50/box

#61C2, Dry Control Cards
(512-4000)10/box

#61C3, Liquid Controls
(512-3595)3/box
CPT® Code: 86592†

UPDATED: USPSTF "B" RATED GUIDELINE –

recommends screening for hepatitis C virus in persons at high risk for infection and a one-time screening for ALL adults born between 1945 and 1965 (baby boomers)*.

About **4.1 million** people have **hepatitis C** and **3 in 4** don't know it

and About **1.1 million** people have **HIV** and **1 in 5** don't know it

Change a Life in 20 Minutes...

Only early diagnosis and linkage to care can improve your patient's health outcomes.

The CDC and the new USPSTF "A" rated guidelines recommend HIV screening of ALL adolescents and adults age 15 to 65. Additionally, both groups urge testing for ALL persons at risk for hepatitis C.

Deliver results in 20 minutes at point of care

Learn more at: OraSure.com/RapidTesting

OraQuick[®] Rapid Antibody Test
HCV

*The **first** and **only** FDA-approved, CLIA-waived test for hepatitis C (HCV) antibodies*

- Rapid Results, Lab Accurate, Easy to Use



OraQuick[®] Rapid Antibody Test
ADVANCE[®] HIV-1/2

*The **only** CLIA-waived oral fluid rapid HIV-1/2 antibody test*

- Simple, Accurate, Versatile
- Proven Patient Acceptance and Increased Test Result Delivery



OraSure Technologies

220 East First Street
Bethlehem, PA 18015
1-800-672-7873
www.OraSure.com

*Moyer, Screening for Hepatitis C Virus Infection in Adults: U.S. Preventive Services Task Force Recommendation Statement, Annals of Internal Medicine, www.annals.org 07/09/2013.

©2013 OraSure Technologies, Inc. OraQuick[®] and OraQuick ADVANCE[®] are registered trademarks of OraSure Technologies, Inc. IFD0022 (rev. 07/13)

To order, contact your Henry Schein sales consultant at 1-800-772-4346



hiv


Uni-Gold™ Recombigen® HIV-1 Test

Single-use rapid immunoassay for qualitative detection of antibodies to HIV-1 in serum, plasma, and whole blood (venipuncture and finger stick). Intended for use in point-of-care settings as an aid in diagnosis of infection with HIV-1.

- CLIA Waived
- Results in 10 minutes or less
- 100% sensitivity

#1206506, Test
(104-3399)20/box
#1206530, Test Control
(174-3798)ea
CPT® Code: 86703QW†

Alere™
Clearview® HIV 1/2 Stat-Pak®

A simple, cassette-based format for the detection of antibodies to HIV-1 and HIV-2 in human whole blood, serum, and plasma.

- CLIA waived for fingerstick and venipuncture whole blood
- Results in 15 minutes
- Sensitivity: 99.7%
- Specificity: 99.9%
- 24-month shelf life from date of manufacture
- **Made in the USA**



#92110
(512-9386)20/box
CPT Code: 86703QW†

Alere™
Clearview® Complete HIV 1/2 Test

A unique 3-in-1, self-contained barrel format that handles the collection, processing, and analysis of whole blood, serum, or plasma for detecting HIV-1 and HIV-2 antibodies.

- CLIA waived for fingerstick and venipuncture whole blood
- Includes a bandage and lancet with each device
- Results in 15 minutes
- Sensitivity: 99.7%
- Specificity: 99.9%
- 24-month shelf life from date of manufacture
- **Made in the USA**



#92111
(512-1454)25/box
CPT® Code: 86703QW†

hepatitis c


OraQuick™ HCV Rapid Antibody Test Kit

- The first and only FDA-approved CLIA-waived HCV rapid antibody test
- Lab accuracy (greater than 98% accurate)
- Tests for multiple HCV genotypes
- Easy to use, simple platform
- Rapid results in 20 minutes
- Fingerstick and venipuncture whole blood collection
- **Made in the USA**

#1001-0181
(810-0012)25/kit
#1001-0180
(810-0013)100/kit
#1001-0182, Oraquick HCV Rapid Controls
(810-0011)ea
Contains: 1 bottle of each: positive & negative.
CPT® code: 86803QW†


OraQuick® HIV Rapid Tests

- Rapid human immunodeficiency virus (HIV) types 1 and 2 (antibody test is a single-use, qualitative immunoassay)
- Detects antibodies to HIV
- Only rapid HIV test for use with specimens for oral fluids, venipuncture whole blood, fingerstick whole blood, and plasma
- CLIA waived for specimens for oral fluids, venipuncture whole blood, and fingerstick; whole blood only
- Simple, 2-step, oral-fluid procedure
- Reliable results in just 20 minutes
- Greater than 99% accuracy
- **Made in the USA**



#1001-0079, Test
(810-3221)25/box
#1001-0078, Test
(810-6864)100/box
#1001-0077, Control
(810-3390)ea
CPT® Code: 86703QW†



Make one call or e-mail to our diagnostics, pharmaceutical, and vaccine experts and let us do the research.

Henry Schein's DxRx Hotline is a single point of contact for your organization for answers to various types of questions related to diagnostics, pharmaceuticals, and vaccines.

FEATURES:

- Latest updates on supply issues
- Indications and usage as provided on package inserts
- Brand-to-generic comparisons
- CPT® and J codes*
- CLIA product status
- Product storage requirements
- ACIP recommendations for vaccinations
- Flu vaccine and flu-related topics

**For informational purposes only.*

*Source: www.reimbursementcodes.com or www.cms.gov.
Customer is responsible for verification of billing/coding
in accordance with applicable specific circumstance.*



HENRY SCHEIN®
DxRxSolutions™
DIAGNOSE • TREAT • PREVENT

BENEFITS:

- Expedite diagnosis and treatment at the point of care
- Free, on-demand resource for immediate needs
- Easy-to-access and use resource
- One phone call to answer pharmaceutical, diagnostic, and vaccination questions
- Live support with no automation

1-877-523-SHOT (7468)

E-mail: DxRxSolutions@henryschein.com

Visit: www.henryschein.com/drxrsolutions

 HENRY SCHEIN®
MEDICAL

 **Rely on Us**



mononucleosis



HENRY SCHEIN® ONESTEP+ MONO CASSETTE TEST

- CLIA waived
- Finger stick or collection tube
- Room-temperature storage
- Includes controls

(900-4070)20/box

Contains: 20 pouched test devices with disposable sample droppers; disposable, heparinized capillary tubes with dispensing bulb; 1 of each control (positive & negative); 1 sample buffer; 1 procedure card & 1 package insert.
CPT® Code: 86308QW†

HENRY SCHEIN® BRAND PRODUCTS

HENRY SCHEIN® ONE STEP PRO+ MONONUCLEOSIS TEST KIT

- CLIA waived for whole blood
- No age restrictions
- Measures heterophile antibodies in whole blood, serum, and plasma
- Uses best technology of bovine derivatives, which eliminate cross reactivity
- Room temperature or refrigerated storage
- Two internal controls exist to avoid additional cost and ensure good laboratory practice
- Only 25 µL blood needed
- **Made in the USA**

**No Age
Indication!**



Mono Test Kit (112-5547)20/box

Contains: 20 test devices, 1 developer solution, 1 negative control, 1 positive control, 40 sample transfer pipettes (20 for serum/plasma & 20 for whole blood), 1 procedure card & 1 package insert.

CPT® Codes: 86308QW† (Whole Blood) / 86308† (Serum/Plasma)



ICON® Mono Test Kits

Approved for testing on all ages. Combines fingerstick whole-blood convenience with EIA membrane technology for fast, accurate results. Uses direct solid-phase immunoassay technology for early qualitative detection of IM heterophile antibodies in whole blood, serum, or plasma. Serum/plasma is moderately complex.

- Whole-blood-test waived
- Relative sensitivity of serum/plasma: 99.9%
- Relative specificity of serum/plasma: 98.8%
- Relative sensitivity of whole-blood assay: 98.3%
- **Made in the USA**



#395100A (111-7123)30/box

Contains: 30 test cassettes, 1 developer solution, 1 of each control (negative & positive), 1 bulb & capillary tubes.

CPT® Code: 86308†



QuickVue® Infectious Mononucleosis Test Kits

- IgM-specific test
- CLIA waived for whole blood
- Easy-to-read ± test result
- Results in 5 minutes
- Sensitivity: 99%
- Specificity: 99%
- Flexible specimen-collection options
- **Made in the USA**

#20121 (227-0024)20/box

Contains: 20 test cassettes, 1 developer, 1 of each control (negative & positive), 20 sample pipettes & 20 capillary tubes.

CPT® Codes: 86308QW† (Whole blood); 86308† (Serum)



Acceava® Mono Cassette

The high-performance mononucleosis diagnostic test for whole-blood, serum, or plasma samples.

- CLIA waived for whole-blood procedure only
- FDA cleared with no age limits
- Positive results in 8 minutes
- Superior sensitivity for earlier IM detection: 99.9%
- Specificity: 98.8%
- External positive and negative controls provided with your kit—no additional purchase is necessary
- **Made in the USA**

#92407
(512-0004)25/pkg



OSOM® Mono Tests

Simple, in-office test for the detection of mononucleosis.

- CLIA waived for whole blood
- Easy to use
- Clinical studies included children 4 and up
- Sensitivity: 100%
- Specificity: 96%
- Two-color results in 5 minutes
- Includes 2 additional tests for external QC testing
- **Made in the USA**



#145
(327-3297)27/box
CPT® Code: 86308QW†

Alere Toxicology

iScreen™ DX Drug Tests

- CLIA waived
- Simple dip-and-read procedure
- Results in minutes
- Laboratory accurate
- **Made in the USA**

Single Dip Devices25/box

Specify:

#IS1 AMP-DX,
Amphetamine (AMP).....(395-0031)

#IS1 BAR-DX,
Barbiturate (BAR).....(395-0036)

#IS1 BZO-DX,
Benzodiazepine (BZO).....(395-0033)

#IS1 BUP-DX,
Buprenorphine (BUP).....(395-0166)

#IS1 COC-DX,
Cocaine (COC)(395-0035)

#IS1 MET-DX,
Methamphetamine
(mAMP).....(395-0037)

#IS1 MTD-DX,
Methadone (MTD).....(395-0038)

#IS1 OPI-DX,
Opiate (OPI).....(395-0039)

#IS1 OXY-DX,
Oxycodone (OXY).....(459-0011)

#IS1 MDM-DX,
Ecstasy (MDMA)(395-0043)

#IS1-PCP-DX,
Phencyclidine (PCP).....(395-0040)

#IS1 TCA-DX, Tricyclic
Antidepressants (TCA).....(395-0041)

#IS1 THC-DX,
Marijuana (THC)(395-0042)

Dip Cards

#IS5M-DX: COC, THC, OPI, AMP, mAMP
(395-0044)25/box



#IS10-22015-DX, COC, THC, OPI, AMP,
mAMP, PCP, BZO, BAR, MTD, MDMA
(395-0045)25/box

#IS10-DX: AMP, BAR, BZO, COC, mAMP,
MOP, MTD, OXY, TCA, THC
(395-0136)25/box

Medicare/Medicaid CPT® Code: G0434QW†
Private Pay CPT® Code: 80104†

ID CODE TRANSLATIONS:

AMP: Amphetamine
BAR: Barbiturates
BUP: Buprenorphine
BZO: Benzodiazepines
COC: Cocaine
THC: Marijuana
MTD: Methadone
mAMP: Methamphetamine
MDMA (Ecstasy): Methylenedioxymethamphetamine
MOP: Opiate 300
OPI: Opiate 2000
OXY: Oxycodone
PCP: Phencyclidine
TCA: Tricyclic Antidepressants

*ID Code Translations
also valid for items
on page 41*



Alere Toxicology

iCup™ DX Drug Screen Cups

- CLIA waived
- Optional sample validity tests
- Simple process
- All inclusive
- No urine exposure
- Tester activates test when ready
- 5-minute results
- **Made in the USA**

I-DXA-1127-023, 12 Panels

Tests for up to 12 drugs: COC, THC, OPI,
AMP, mAMP, PCP, BZO, BAR, OXY, MTD,
TCA, and MDMA.

(113-2722)25/box

I-DXA-1107-1-142, 10 Panels

Tests for up to 10 drugs: COC, THC, OPI,
AMP, mAMP, BZO, BAR, OXY, MTD,
and MDMA.

(395-0024)25/box

I-DXA-187-016, 8 Panels

Tests up to 8 drugs: COC, THC, OPI, AMP,
mAMP, BZO, OXY, and MTD.

(852-0015)25/box

I-DXA-167-013, 6 Panels

Tests for up to 6 drugs: COC, THC, OPI,
mAMP, BZO, OXY; plus SVT
adulteration strips.

(459-0014)25/box

Point of Care Testing + Laboratory Confirmation

Substance abuse testing with more substance.



From rapid screening devices to full lab confirmation testing, we offer a complete solution for all of your drug detection or monitoring needs.

Introducing our new and improved CLIA-waived devices.

iCup® Dx 14 Drug Test
Cup with Buprenorphine

iScreen® Dx Dip
Card 10 Panel

BUP Single
Dip Device



Alere Toxicology

aleretoxicology.com

© 2013 Alere. All rights reserved. Alere is a trademark of the Alere group of companies. 810118 REV1 7/13

Alere Toxicology

iCassette® Dx Test

CLIA waived, easy to use, and screens for 11 or 12 drugs on one device.

#I-DCB-1115-011, 11-Test Panel: COC, THC, OPI, AMP, PCP, BZO, BAR, MDMA, PPX, OXY, and TCA
(459-0007).....25/box



Alere Toxicology

iCup® AD

All-Inclusive Drug-Screening Test

- Zero exposure
- Built-in validity test
- 3–13 drug configurations*
- Results can be photocopied
- Collect and read
- 24-month shelf life



#I-DUA-157-013, 5 Panel (CLIA Waived): COC, THC, OPI, AMP, and mAMP
(195-4184).....25/box

#I-DUA-157-023, 5 Panel (CLIA Waived): COC, THC, OPI, AMP, and PCP
(673-2588).....25/box

#I-DUA-157-034, 5 Panel (CLIA Waived): COC, THC, OPI, mAMP, and PCP
(204-9814).....25/box

#I-DUA-167-022, 6 Panel (CLIA Waived): COC, THC, OPI, AMP, mAMP, PCP (OX, SG, PH)
(604-7353).....25/box

Alere Toxicology

iScreen® OFD–Cotinine

The iScreen® OFD–Cotinine device is a rapid, oral fluid screening test that can be performed any time without the use of an instrument. The test detects elevated levels of cotinine (nicotine metabolite) in oral fluid for accurate determination of smoking status.

- Accurately screen for cotinine (primary metabolite of nicotine)
- Detects up to 1-2 days after nicotine use
- Qualitative detection of cotinine in oral fluids at a cutoff of 30 ng/mL
- Simple collection and testing procedure
- Results in as little as 10 minutes
- Integrated internal procedural control

#I-DCT-B702, Nicotine Test
(395-0149).....25/box

#I-DUA-167-022, 6 Panel (Moderate): COC, THC, OPI, AMP, mAMP, and BZO
(104-2333).....25/box

#I-DUA-167-291, 6 Panel (Moderate): COC, THC, OPI, AMP300, OX4, MDMA (OX, SG, and PH)
(874-7298).....25/box

#I-DUD-187-013, 8 Panel (Moderate): COC, THC, OPI, AMP, mAMP, PCP, BZO, and BAR (OX, SG, PH, NI, GL, and CR)
(303-6739).....25/box

#I-DUE-187-071, 8 Panel (Moderate): COC, THC, MOP, AMP, mAMP, PCP, BZO, OXY (OX, CR, and PH)
(113-5160).....25/box

#I-DUE-1107-141, 10 Panel (Moderate): COC, THC, OPI, AMP, mAMP, BZO, BAR, OXY, MTD, PPX (OX, CR, and PH)
(113-2737).....25/box

#I-DUE-1127-022, 12 Panel (Moderate): COC, THC, OPI, AMP, mAMP, PCP, BZO, BAR, MTD, TCA, OXY, and PPX
(110-5436).....25/box

*Not all configurations are CLIA waived.

Medicare/Medicaid CPT® Code: G0434QW†
Private Pay CPT® Code: 80104†

Alere Toxicology

iCassette®

One-Step Drug-Screening Test

- Minimal sample volume required
- 1–11 drug configurations
- Results in about 5 minutes
- 24-month shelf life
- Simple procedure



#I-DOA-2155, 5 Panel, CLIA Waived
(447-2471).....25/box

#I-DOA-3155, 5 Panel, CLIA Waived
(724-9019).....25/box

#I-DOA-1165, 6 Panel, CLIA Waived
(922-8619).....25/box

Medicare/Medicaid CPT® Code: G0434QW†
Private Pay CPT® Code: 80104†

Q.E.D.® A150 Alcohol Saliva Screen

- DOT approved
- CLIA waived
- Quantitative results
- Saliva test
- Fast, reliable results
- Requires no calibration or controls
- Made in the USA

(102-5216).....10/box
CPT Code: 82055†



Buprenorphine (BUP) Drug-Screen Dip

Single-Dip Device (Moderate) Buprenorphine is used in the treatment of opioid addiction. Therapeutically, buprenorphine is used as a substitution treatment for opioid addicts. It is as effective as methadone but demonstrates a lower level of physical dependence.

#DBU-114
(106-8526).....25/box
CPT® Code: G0434†
Private Pay CPT® Code: 80104†



Alere Toxicology

AlcoMate Premium

- CLIA status: Moderate
- 4-digit display (clear LCD)
- DOT approved
- Low-battery indicator/auto power-off
- Internal mouthpiece storage slot (x2)
- Cumulative test counter
- No calibration service required

#ALF7000, Kit
(114-2179).....ea

Contains: 1 AlcoMate, 50 mouthpieces, 1 extra sensor, 1 hard carry case & 1 car charger.

#ATM-Mouthpiece
(114-3049).....50/box
CPT® Code: 82075†



Alere Toxicology

iScreen™ One-Step Drug Card

Efficiently detects the presence of drugs in minutes.

- On-site testing
- Saves you time and money
- Simple dip-and-read procedure
- Room-temperature storage
- Laboratory accurate
- Results in minutes

#IS4A, 4-Test Panel: COC, THC, OPI, and AMP (Moderate)
(108-1504).....25/box

#IS5A, 5-Test Panel: COC, THC, OPI, AMP, and PCP (Moderate)
(827-2314).....25/box

#IS5MP, 5-Test Panel: COC, THC, OPI, mAMP, and PCP
(395-8161).....25/box

#DOA-164-551, 6-Test Panel: COC, THC, OPI, AMP, mAMP, and OXY (Moderate)
(110-8078).....25/box

#I-DOA-364-0BC, 6-Test Panel: COC, THC, OPI, AMP, BZO, and MTD (Moderate)
(680-1961).....25/box

#IS9, 9-Test Panel: COC, THC, OPI, AMP, mAMP, PCP, BZO, BAR, and MTD (Moderate)
(172-1463).....25/box

#IS-10, 10-Test Panel: COC, THC, OPI, AMP, mAMP, PCP, BZO, BAR, MTD and TCA (Moderate)
(392-1003).....25/box

#DOA-1124-011T, 12-Test Panel (Moderate): COC, THC, OPI, AMP, mAMP, PCP, BZO, BAR, MTD, MDMA, OXY, and PPX
(106-5613).....25/box





HENRY SCHEIN® ONESTEP+ H. PYLORI CASSETTE TEST

- Whole blood, serum, or plasma
- CLIA waived for whole blood only
- Simplified
- Results in 10 minutes or less
- Sensitivity: 88%
- Specificity: 88%
- Allows for rapid treatment decisions
- Controls are included

Complete Kit
(900-4071)20/box

Contains: 20 pouched test devices with disposable sample droppers; disposable, heparinized capillary tubes & dispensing bulb; 1 of each control: (positive & negative); 1 sample buffer; 1 procedure card & 1 package insert.
CPT® Codes: 86677 (Serum/Plasma); 86318QW† (Whole Blood)

conjunctivitis

RPS AdenoPlus™ for Pink Eye



10-minute, in-office test used for conjunctivitis (pink eye). Working from an easy-to-collect eye-fluid sample, AdenoPlus™ detects all known serotypes of adenovirus. Knowing the underlying infection source helps clinicians determine the most appropriate treatment.

- CLIA waived
- Accuracy: 90% sensitivity, 96% specificity
- Cost-effective: no additional equipment needed
- Easy to use: 4 simple steps

#RPS-AD, Test
(770-0028)10/box

Contains: 2–2-mL control vials (1 positive & 1 negative.)

#RPS-ADSTD, Control
(108-7896)ea
CPT® Code: 87809QW†

h. pylori

Kimberly-Clark H. pylori Rapid Test MicroCOUNT Lite Meter and Accessories



- Up to 96% sensitivity and 96% specificity
- Active infection determination versus the presence of antibodies, as measured by serology
- Noninvasive, easy-to-swallow gelatin capsules
- Onsite test administration and analysis that can be completed in less than 20 minutes
- Testing not regulated by CLIA
- Reimbursement by Medicare and many other carriers
- **Made in the USA**

#60495, Lite Liquid,
Single-Sample Liquid Scintillation Counter
(108-1365)ea

#60496, Lite Printer
(108-1366)ea

#60442, PYtest Capsules
(370-9222)10/pkg

#60443, PYtest Capsules
(106-0064)100/pkg

#60439, PYtest Startup Package
(106-6419)ea



OSOM® H. pylori Test

For use with whole blood, serum, or plasma.

- Easy to use: single-reagent procedure fits easily into any practice or lab
- Fast: results in 10 minutes or less
- Easy to read: 1 pink line indicates a negative result; 2 pink lines indicates a positive
- Small sample size: test only requires 10 µL of serum or plasma; 25 µL of whole blood
- Accurate sensitivity: 95.9%
- Room-temperature storage
- **Made in the USA**

#175, Test
(309-0040)25/box

Control Kit
(309-0074)ea

Contains: 2–2-mL control vials (1 of each: positive & negative).

CPT® Codes: 86677QW† (Serum/Plasma); 86318QW† (Whole Blood)



Clearview® H. pylori Rapid Test

For the visual qualitative detection of antibodies specific to *Helicobacter pylori* in human serum, plasma, or whole blood. Simple, two-step procedure.



- CLIA waived for whole blood only
- Sensitivity: 88%, Specificity: 88%
- Results in 10 minutes or less
- Room-temperature storage for up to 18 months
- Built-in procedure controls to ensure accuracy
- **Made in the USA**

#4581505020
(512-4160)30/box

CPT® Codes: 86318QW†/86318†

COUNTLESS EYES AT RISK



AdenoPlus® is the only CLIA-waived, point-of-care test to help identify the viral form of conjunctivitis, commonly known as **pink eye**. Commit to making the right diagnosis and the right treatment during the office visit.

- CLIA-Waived
- Easy to use
- 90% sensitivity¹
- 96% specificity¹
- Reimbursement 87809QW
- Results in 10-minutes
- No additional equipment needed



Reference [1] FDA Section 510k number (K110722)
for RPS Adeno Detector Plus™; March 15, 2011.

info@RPSdetectors.com
RPSdetectors.com
1.941.556.1850





thyroid

ThyroChek® Thyroid Rapid-Test Kit


Elevated TSH is the first clinical indicator of hypothyroidism. This test allows you to know in minutes if your patient's TSH is greater than 5 mIU/L. Requires 1 drop of blood from a fingerstick or venous draw.

#TSH-120-A, ThyroChek® Test Kit
(818-0440)ea

Contains: 20 tests, pipettes, buffer, & instructions.

#TSH-120-N, Control Set
(818-0001)2 vials

CPT® Code: 84443QW†



dtm

TROY
Biologicals

DTM Test Agar Dermatophyte Test

Medium Vial

Selective isolation of fungi from cutaneous (skin) samples.

#7051

(224-9011)20/box

Item stored under refrigeration. May be shipped separately.

dvt



Clearview

Clearview® Simplify D-dimer Test

10-Minute Rapid Test

Rapid lateral flow test for qualitative detection of D-dimer to aid in the diagnosis of deep vein thromboses (DVT), pulmonary embolism (PE), and disseminated intravascular coagulation (DIC).

- Fast: 10-minute results
- Easy: 2 easy steps
- Flexible: uses whole blood or plasma
- Reliable: built-in control
- Made in the USA

(891-9725)10/box

CPT® Code: 85380†

Microhematocrit Capillary Tubes

Red Tip, Heparinized

(130-3223)200/vial

Blue Tip, Plain

(227-5359)200/vial

sedi-rate



globe
SCIENTIFIC INC.

Sedi-Rate™ Westergren ESR System and Accessories

The erythrocyte sedimentation rate (ESR) system is a nonspecific test that alerts physicians to the presence of disease or abnormality in the body. Sedi-Rate's pipette features a unique self-sealing aerosol barrier plug designed to seal upon contact with blood, eliminating exposure to blood aerosols. Its polypropylene vial is prefilled with 0.25 mL of 3.8% trisodium-citrate diluent and features a self-sealing stopper. Test results are ready in 60 minutes and can be read directly from the imprinted graduation marks on the pipette.

#3467, Starter Kit
(217-0930)ea

Contains: 100 Sedi-Rate pipettes & vials, 1-10-place rack with pipette support clamps, 1 leveling platform, 100 piercing funnels & 100 ESR transfer pipettes.

#3469, System
(368-5132)100/box

Contains: 100 tests (100 pipettes & 100 vials).

Individual Kit Components

#3459, Acrylic Leveling Platform

(691-4303)ea

#3452, Rack with 10 Pipette

Support Clamps, Holds 10 Tests

(496-3393)ea

#3479, Pipettes

(600-7387)100/box



Dispette 2 Sedimentation Rate Kit

With improved safety features plus greater efficiency. Modified Westergren sedimentation rate. Includes 100 ESR disposable plastic pipette tubes and 100 prefilled, pierceable cap reservoirs.

(668-9703)ea

bladder cancer



lifeSign

Status BTA Test

A Rapid Test for the Qualitative Detection of Bladder-Tumor-Associated Antigen in Human Urine

- CLIA-waived
- Uses monoclonal antibodies to detect the presence of bladder-tumor-associated antigen in urine
- Accurate test results with superior sensitivity and specificity
- Detects early stage recurrence cancer that cytology may miss
- Adjunct test to cystoscopy
- Simple 1-step procedure
- Results available in 5 minutes
- Made in the USA

(641-0001)25/box

AllerVision Offers an “Ahh” Inspiring Allergy Solution For Your Practice.

With over 50% of Americans sensitized to at least one allergen, the World Allergy Organization urges practices to “find the cause and prevent symptoms and disease progression, rather than just rely on medications to suppress the symptoms.” AllerVision can help your practice establish a comprehensive allergy program. We’ll train your entire staff to use the most effective diagnosis and treatment options available, so you can help your patients and your practice.

To learn more, call your Henry Schein Consultant or (800)772-4346 today.

ALLER **VISION**

©2013 AllerVision LLC, Scottsdale, AZ. All rights reserved.

**Call Your
Henry Schein
Sales Consultant
for more
information
& order details.**





HITACHI

Inspire the Next

CLA-1™ Luminometer

This compact luminometer automatically prints a complete report that lists the tested allergens and the severity of the patient's response to each. Easy to use, the system is the simplest, quickest, and most convenient allergy testing method available today. It can analyze up to 180 allergens in 8 minutes.

#95013
(680-1614)ea

HITACHI

Inspire the Next

CLA® Allergy Tests

Eastern, Combines IgE Foods and Inhalants
(406-5680)20/box
Contains: 20 CLA® Pettes & reagents.

Moderate Food Panel Test
(286-2678)20/box
Contains: 20 CLA® Pettes & reagents.

Northeastern Inhalant
(406-9608)20/box
Contains: 20 CLA® Pettes & reagents.

Northwestern Inhalant
(618-7516)20/box
Contains: 20 CLA® Pettes & reagents.

Pediatric Comprehensive, Combines IgE Foods and Inhalants
(905-0111)20/box
Contains: 20 CLA® Pettes & reagents.

Southeastern Inhalant
(406-3895)20/box
Contains: 20 CLA® Pettes & reagents.

Southern, Combines IgE Foods and Inhalants
(259-4912)20/box
Contains: 20 CLA® Pettes & reagents.

Southwestern Inhalant
(406-5585)20/box
Contains: 20 CLA® Pettes & reagents.

U.S. 15
(110-0538)20/box
Contains: 20 CLA Pettes & Reagents.

Western Inhalant
(168-6399)20/box
Contains: 20 CLA® Pettes & reagents.

Serum Vials, Blue Cap
(407-0810)100/box

Thermal Printer Paper
(613-3347)ea

#93026, Positive Control Serum
(835-1985)ea

#94173, Negative Control Serum
(610-7817)ea



Hb 201+ Handheld Hemoglobin Analyzer

Offers the same important benefits as its predecessor, the B-Hemoglobin Photometer. Smaller size gives more flexibility and mobility to use the system anytime and anywhere. Performs an internal self-test. No control cuvette required.

- Produces results in less than 1 minute with lab-quality accuracy
- 3 easy steps
- Precalibrated, portable analyzer that requires minimal maintenance
- No calibration or instrument manipulation (between cuvette/reagent batches) is needed
- Safe and easy to use
- Reimbursable
- CLIA '88 waived
- 2-year warranty

#121721, Analyzer
(345-2344)ea
Contains: analyzer, AC adapter, manual & training CD-ROM.



#111715, Microcuvettes
(110-9167)100/box
Contains: 4 packs of 25.

#111716, Microcuvettes
(345-8733)200/box
Contains: 4 packs of 50.



Alere™

Alere™ HemoPoint® H2 System

Using a single drop of blood, the HemoPoint® H2 System delivers accurate hemoglobin and hematocrit results in less than a minute. The touch screen is easy to learn and operate, minimizing training time. Rechargeable battery can last 100 hours between charges. Recalls the last 4000 test results with the touch of a button.

#55118, HemoPoint® H2 Meter
(117-9599)ea

Contains: HemoPoint® H2 Photometer, user guide, quick reference guide, optics cleaner, control cuvette & power adapter.

#80549, HemoPoint® H2 Microcuvettes
(117-9605)100/box
Contains: 2 containers of 50 cuvettes each.

#88772, Hemoglobin Controls
(117-9656)ea
Contains: 1-1.5 mL Low & 1-1.5 mL High.

#52219, HemoPoint™ H2 Optics Cleaner
(118-5057)ea
Contains: 5 cleaning units.



MAGELLAN BIOSCIENCES

LeadCare® II Analyzer

Accurate, On-the-Spot Lead Testing

- CLIA-waived, reimbursable test
- Perform blood-lead measurements with only 2 drops of blood
- Accuracy and sensitivity of a reference lab test (from capillary or venous blood sample)
- Two-drop sample size requires only a finger stick; allows combination of LeadCare II with other tests
- In-office immediate results in 3 minutes

Analyzer
(106-4441)ea

Test Kit
(106-5736)48/box

#70-3440, Lead Control
(707-3252)ea

#70-3447, Printer Lead Analyzer
(106-7328)ea
CPT® Code: 83655†



9:15 Conduct diagnostic tests
9:30 Discuss the results
9:45 Back on her feet

Accelerate Decision Making at the Point of Care

Smarter diagnostic tools with rapid results empower decisions that can get your patients back on their feet, back to their loved ones and back to life.

- Alere Cholestech LDX® System for Lipid, Glucose, ALT and AST testing
- Alere INRatio®2 PT/INR Monitoring System for Anticoagulation Monitoring
- Alere Triage® Meter with BNP for the Assessment of Heart Failure



Empower Your Practice
alere.com | 877.441.7440



© 2011 Alere. All rights reserved. The Alere Logo and Alere are trademarks of the Alere group of companies. 1000363-01 11/11



Alere™ Triage® MeterPro™ Meter

High-performance, comprehensive cardiovascular and drugs-of-abuse testing in the laboratory or at the point of care (POC). Designed to provide diagnostic results quickly and easily. Quantitative diagnostic panel results are displayed and printed in easy-to-read numbers and highlighted when the result is above the cutoff. Qualitative diagnostic panel results are displayed as positive (POS) or negative (NEG), eliminating the need for interpretation of lines that often hamper qualitative diagnostics.

- 16 immunoassays available
- Results in about 15 minutes
- Whole blood, plasma, or urine
- Customizable test panels
- Onsite training and implementation support
- Results printed, displayed, and interfaceable
- Full connectivity support

#55070MP, Triage Meter Plus
(496-3405)ea

#52031, Paper for Triage Meter and MeterPro, 1 Roll
(111-6321)ea



Alere™ Triage® Profiler S.O.B. (Shortness of Breath) Panel

Rapid fluorescence immunoassay of the quantitative determination of creatine kinase MB (CK-MB), myoglobin, troponin I, B-type natriuretic peptide (BNP) and D-dimer in whole blood or plasma. Test devices include 2 levels of built-in controls. Measurement of 5 markers:

- Troponin I
- CK-MB
- Myoglobin
- BNP
- D-dimer
- **Made in the USA**

#97000BD, Panel
(241-0108)25/box




Alere™ Triage® Cardiac Panel

Measures 3 key cardiac markers: myoglobin, CK-MB, and troponin I.

- Whole blood or plasma
- **Made in the USA**


#97000HS
(131-1964)25/box



Alere™ Triage® Cardio Profiler Panel

Measurement of 4 cardiovascular markers: troponin, CR-MB, myoglobin, BNP.

#97100CP, Panel
(706-7343)25/box



Alere™ Triage® D-dimer Test

Rapid, quantitative immunoassay to aid in the assessment and evaluation of patients suspected of having disseminated intravascular coagulation (DIC) and thromboembolic events, including pulmonary embolism (PE) and deep-vein thrombosis (DVT).

- **Made in the USA**

#98100, Test
(131-0268)25/box

#98113, Control Set, Control 1 (5 x 1 Level)
(176-9303)5/box

#98114, Control Set, Control 2 (5 x 1 Level)
(241-0075)3/box



Alere™ Triage® Control 5 Calibration Verification

2 x 5 Levels

#88755
(111-9927)5/box

Alere™ Triage® Total Control 5, Control 2

For Cardio Profiler, S.O.B. and Cardiac Panel
5 x 0.25 mL

#88754
(111-9843)5/box

Alere™ Triage® Total Control 5, Control 1

For Cardio Profiler, S.O.B. and Cardiac Panel
5 x 0.25 mL

#88753
(111-9841)5/box



Cardiac Reader Analyzer/D Dimer System

Cardiac Reader Analyzer
(106-5804)ea

Cardiac T Quantitative Assay
(106-5803)10/box

Cardiac M Assay
(106-5802)20/box



Alere™ Triage® BNP Test

The only CLIA-waived blood test measuring B-type natriuretic peptide (BNP). CLIA waived for use with whole-blood samples.

Rapid, objective point-of-care immunoassay is used:

- As an aid in the diagnosis and severity of CHF
- For the risk stratification of patients with ACS and HF
- **Made in the USA**

#98010XR, Triage BNP Test
(131-0086)15/box

#98000XR, BNP Test
(131-7507)25/box

#98013XR, BNP Control Set, Level 1 (5 Level 1)
(219-6429)ea

#98014XR, BNP Control Set, Level 2 (5 Level 2)
(219-6475)ea

#98015XR, BNP Calibration Verification Control Kit (1 per Level, 5 Levels)
(131-9669)ea

#98200, BNP Reagent for Beckman Coulter Immunoanalyzers
(104-4059)100/box

#98201, BNP QC Controls for Beckman Coulter Immunoanalyzers
(478-8086)ea

#98202, BNP Calibrator for Beckman Coulter Immunoanalyzers
(640-1104)ea



Alere™ Triage® TOX Drug-Screen

Simultaneously detects 11 drugs in one rapid, easy-to-use test that is read on the Triage® meter.

- Acetaminophen
- Amphetamines
- Methamphetamines
- Barbiturates • Benzodiazepines
- Cocaine • Methadone
- Opiates • PCP • THC
- Tricyclic antidepressants
- **Made in the USA**

#94400, Test
(131-0824)25/pkg

#94413, 0.25-mL Control Level 1
(131-7103)5/pkg

#94414, 0.25-mL Control Level 2
(131-0969)5/pkg





ABAXIS

Piccolo® Point-of-Care Chemistry Analyzer Accessories

The first broad-menu clinical chemistry analyzer designed for near-patient testing. *Basic Panel* contains BUN, Ca, Cl-, CRE, GLU, K+, Na+, and tCO₂.

Comprehensive Panel contains ALB, ALP, ALT, AST, BUN, Ca, Cl-, CRE, GLU, K+, Na+, TBIL, tCO₂, and TP. *Lipid Panel* contains CHOL, CHOL/HDL (calculated), HDL, LDL (calculated), TRIG, and VLDL (calculated).

- General chemistry, liver, electrolyte, BMP, CMP, lipid, renal, and hepatic panels available
- No daily controls, calibrators, or maintenance
- Intelligent Quality Control (iQCT) verifies correct performance
- No special training required
- 0.1cc whole blood, serum, or plasma
- Single-use reagent discs
- Compact and portable
- Built-in QC and calibration
- Lipid panel is CLIA waived
- Whole blood, serum, or plasma
- Less than 2 minutes hands-on
- Results in under 13 minutes

#1100-0000, Accessory Kit
(106-6368)ea

Contains: 1-0.1cc Minipette, 1-96-pack of disposable Minipette tips, 1 operator's manual, training CD, 1 power cord, 1 power supply, 10 lithium-heparin tubes, 1 screwdriver & 2-12-packs of results cards.

#400-0024, Basic Metabolic Panel
(137-3987)10/box

#400-0028, Comprehensive Metabolic Panel
(174-6978)10/box

#400-0025, Lipid Panel
(146-0189)10/box



Alere™ Cholestech LDX® System

Complete Lipid Profile, Glucose

- CLIA waived**
- Procedure saves time and enhances productivity
- Small, tough, and portable
- Eliminates anxiety of venipuncture; requires only a fingerstick
- On-the-spot adjustment of therapy and improved patient compliance
- Accuracy comparable to laboratory tests
- Saves time: fewer follow-up visits required
- Enhances compliance through teachable moments between healthcare providers and patients

#10-959, Analyzer
(738-2917)ea

Contains: 1 of each: analyzer, printer, power supply, starter pack, optics check cassette, printer cable, user's manual, procedure manual & training video. Starter pack contains 200 lancets, 200 capillary tubes, 200 capillary plungers, 1 Minipette pipette, 50 pipette tips, 1 accessory tray & 1 fingerstick tip pad.

#11-785, Cholestech Thermal Printer Labels
White labels for use with Cholestech Thermal Label Printer.
(803-6378)400/roll

#10-989 Test for Lipid Profile
(116-9995)10/box

#10-986, Test for TC
(420-7892)10/box

#10-988, Test for TC and Glucose
(956-0149)10/box

#10-987, Test for TC and HDL Panel
(600-5779)10/box



#10-990, Test for TC, HDL, and Glucose Panel
(600-5781)10/box

#10-991, Test for Lipid Profile and Glucose
(694-9360)10/box

#12-712, Multi-Analyte Controls, 2 mL
(116-4232)ea
Contains: 1 vial of each level (1-2).

#11-255, Calibration Verification Material (Lipids/Glucose), 2 mL
(116-7301)ea
Contains: 1 vial of each level (1-4).

#11-010, Disposable Pipette Tips
For use with 35µL, 40µL, and 50µL MiniPet Pipettes.
(267-4596)50/box

#10-311, Capillary Plungers
For use with all types of Cholestech capillary tubes.
(960-1587)50/box

CPT® Codes for lipid profile: #80061QW†; for lipid profile/ glucose: 80061QW†/839471QW†; for ALT/AST: 84460QW†

Printer Cable

Connects the Cholestech LDX to the Thermal Label Printer.

#11-787
(803-3213)ea



i-STAT® 1 Handheld System

Includes these important features:

- Easy to use: Cartridges are self-contained and the system is simple to learn and operate
- Lab-accurate results: Advanced biosensor technology and real-time quality assurance provide extremely accurate testing
- Small sample size: Just a few drops of blood per cartridge
- Ready when you need it: Test results in as little as 2 minutes
- Made in the USA



i-STAT Analyzer
(111-3324)ea

Cartridges
Chem8+ (waived)
(111-3408)25/box

Creatine (waived)
(111-3394)25/box

CG8+ (waived)
(115-6662)25/box

PT
(111-3385)24/box

BNP
(111-3158)25/box

i-STAT Printer Paper
(111-3398)6/box



CoaguChek XS® System

PT/INR Monitoring Solution

The next generation in point-of-care prothrombotin (PT) and international normalized ratio (INR) testing, this new system uses exclusive smart technology to help ensure accurate results and offers the flexibility and control that health care professionals need to confidently manage a wide range of anticoagulated patients. It is the only system to perform onboard quality control and determine patient results in a single-test chamber.

- CLIA waived
- Neutralizes therapeutic levels of heparin and LMWH
- INR corrected for hematocrit within specified range; manage a broad array of patients with a single device
- Small sample size (10 µL) for easy dosing
- Strip allows dosing from either top or side for simple blood application
- Small, battery-powered, handheld meter for portability and efficiency
- Onboard controls; no external QC necessary
- Auto-on with strip insertion and auto-off capabilities
- 18-month strip shelf life; no refrigeration needed
- **Made in the USA**

XS® System

(891-0581).....ea

Strips For XS and XS Plus Meters

(891-7593).....48/box

Contains: 24 strips per bottle, 2 bottles per box.



Controls, Hi & Lo

(109-9989).....ea

3-Box Starter Kit

(108-1861).....ea

Contains: XS Analyzer, 3 boxes of CoaguChek XS PT strips (144 total) & CoaguChek Lancets.

5-Box Starter Kit

(108-7757).....ea

Contains: XS Analyzer, 5 boxes of CoaguChek XS PT test strips (240 total) & CoaguChek Lancets.

10-Box Starter Kit

(108-2228).....ea

Contains: XS Analyzer, 10 boxes of CoaguChek XS PT test strips (480 total) & CoaguChek Lancets.

CoaguChek XS® Plus System

Battery

(109-9303).....ea

10-Box Starter Kit

(109-3345).....ea

Contains: CoaguChek® XS Plus Analyzer, 10 boxes (480 total) of CoaguChek XS PT test strips & CoaguChek lancets.



VerifyNow™ Aspirin Assay

A qualitative assay to aid in the detection of platelet dysfunction due to aspirin ingestion incitrated whole blood for the point-of-care or laboratory setting.

- Simple: 1-step, closed-tube sampling
- Rapid: results in 5 minutes
- Accurate: >99% specificity; 91.4% sensitivity
- Cost-effective: eligible for Medicare and private-payer reimbursement with CPT® Code 85576QW†; platelet aggregometry–arachidonic acid (agonist)

(113-2708).....ea



ProTime® Monitor

- Accepts and stores patient and/or operator ID for reduced transcription errors
- Automatic: sends results quickly and easily to printer, computer, or LIS
- Only system that offers flexibility of both onboard and optional external quality controls
- Thromboplastin reagent of ISI=1.0 (AHA, CAP, and WHO recommended)
- Off-meter sampling provides safety against disease transmission
- Bar-coded cuvette: no coding necessary
- Longer open-pouch stability: up to 16 hours

Meter

(107-0307).....ea

#PRO3-25, Disposable Cuvette, 3 Channels

(487-7457).....25/box

#PRO5-25, Disposable Cuvette, 5 Channels

(487-2270).....25/box

CPT® Code: 85610QW†



Alere™ INRatio2® Monitor

Allows you to obtain a prothrombin time (PT) and international normalized ratio (INR) result using 1 drop of blood from a fingerstick.

- CLIA waived
- Easy to use
- Single-button operation
- 1-step auto turn-on with strip insertion
- Larger LCD screen and easy-to-read digits
- Small, 1-drop (15 µL), sample size
- Wider test-strip area for easier sample application
- Test results in under 1 minute
- Test strips do not require refrigeration
- 2-channel, on-board quality control runs on every test
- No need for external liquid quality control or extra strips
- Test strips are individually wrapped and stored at room temperature

INRatio2 Professional Monitor Kit

(111-3837).....ea

#INR2QS1 Quick Start #1

(111-3840).....ea

Contains: analyzer, 1 box of 48 INRatio test strips, 50 Unistik® single-use lancets, 100 alcohol pads & 250 gauze brick pads.



#INR2QS3 Quick Start #3

(111-3844).....ea

Contains: analyzer, 3 boxes of 48 INRatio test strips (144 total), 150 Unistik® single-use lancets, 200 alcohol pads & 250 gauze brick pads.

#INR2QS5 QuickStart #5

(111-6578).....ea

Contains: analyzer, 5 boxes of 48 INRatio test strips (240 total), 250 Unistik® single-use lancets, 300 alcohol pads, 5 boxes of capillary tubes & 500 gauze brick pads.

#INR2-6PK 6-Box Promo

(112-3230).....ea

Contains: analyzer, 6 boxes of 48 INRatio test strips (288 total), 10 single-use lancets & 10 capillary tubes.

CPT® Code: #85610QW†

BD Diagnostics – Preamalytical Systems Technology Built On Trust

The BD Vacutainer® Push Button Blood Collection Set

- The only blood collection system with push-button technology
- Split-second activation for improved safety
- In-vein activation – assures retraction of the needle directly from the vein
- Intuitive design – reduces training



HS Reference No.	BD Reference No.	Product Description
9879975	367326	PBBCS 12" w/o Luer 21 x.75
9872637	367324	PBBCS 12" w/o Luer 23 x.75
9873303	367323	PBBCS 12" w/o Luer 25 x.75
2170920	367338	PBBCS 7" with Luer 21 x.75
2489708	367336	PBBCS 7" with Luer 23 x.75
9878340	367335	PBBCS 7" with Luer 25 x.75
9873358	367344	PBBCS 12" with Luer 21 x.75
9194247	367342	PBBCS 12" with Luer 23 x.75
9879168	367341	PBBCS 12" with Luer 25 x.75

BD Microtainer® Contact-Activated Lancets

- Designed with a Positive Patient Experience in Mind
- Activates on Contact
- Improves Puncture Accuracy
- Simplifies the Procedure



HS Reference No.	BD Reference No.	Dimensions (Width x Depth)	Blood Volume
1065174	366592	30 G x 1.5 mm, Purple	Low Flow
1065210	366593	21 G x 1.8 mm, Pink	Med. Flow
1065209	366594	1.5 mm x 2.0 mm, Blue	High Flow

View the product
demonstration video!

Please visit:
www.bd.com/vacutainer



BD Diagnostics
Preamalytical Systems
1 Becton Drive
Franklin Lakes, NJ 07417
www.bd.com/vacutainer

The value of any rebates, discounts, incentives, or product provided to customer may constitute a "discount or other reduction in price" under Section 1128B(b)(3)(A) of the Social Security Act [42 U.S.C. Sec. 1320a-7b(b)(3)(A)]. Customer shall satisfy any and all requirements imposed on buyers relating to discounts or reductions in price, including, when required by law, to disclose all discounts or other reductions in price received from BD and to accurately report under any state or federal health care program the net cost actually paid by customer.

BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2009 BD DIST072



HENRY SCHEIN®
BLOOD LANCETS
Sterile, Disposable
(100-0178)200/box



HENRY SCHEIN®
SAFETY LANCETS

- Finger activation
- Needle retracts automatically after use
- Cannot be reloaded after use
- Safe for patient and user

28 ga, Fine
(900-4433)200/box

23 ga, Medium
(900-4434)200/box

17 ga, Coarse
(900-7569)200/box



Tenderfoot®
Lancets
Infant heel
blood sample.
#TF50
(107-9563)50/box



OneTouch® SureSoft™ Lancing Devices

Single-use lancing devices are electropolished, which results in more blood than an acid-washed needle. This means you can use a smaller gauge needle while still ensuring that a sufficient sample will be drawn for point-of-care testing.

- Electropolished needle means a smaller gauge and more comfort for your patients
- Meets OSHA guidelines: nonreusable, retractable needle
- Distinctive wing design prevents reloading and reusing
- Preset gauge for no guessing, no maneuvering and easy training
- On-demand training through Web and CD-ROM programs

Gentle, 28-ga Needle
(567-0000)200/box

Regular, 21-ga x 2-mm Needle
(115-5267)200/box



Unistik® 2
Safety
Lancets
21 ga

Normal, for Standard Tests, Yellow, 2.4 mm
(668-3836)50/box
(668-5603)200/box

Extra, for Larger Samples, Orange, 3 mm
(159-3956)100/box
(668-6699)200/box



Latex-free Tourniquet

Blue, 1" x 18"

- Used for constricting or compression
- Feels like latex
- Hypoallergenic
- Superior strength

(933-0168)100/box




**BD™ Vacutainer®
Eclipse™ Blood-
Collection Needles**

Thin Wall, 1 1/4"

Safety-engineered, multisample needle offers a simple, effective way to collect blood while it reduces the possibility of needlestick injuries. It features a patented shielding that allows for one-handed activation to cover needle immediately upon withdrawal from vein and confirm proper activation with an audible click. New needles with preattached holders are ready to use with no assembly required.

#368608, 22 ga	
(987-2645)	48/box
#368607, 21 ga	
(987-6956)	48/box
#368651, with Preattached Holder, 22 ga	
(987-6134)	100/case
#368650, with Preattached Holder, 21 ga	
(987-5021)	100/case


**BD™ Vacutainer®
Safety-Lok™ Blood
Collection Sets**

Safety-engineered winged sets for both blood collection and infusion. Each set features a translucent, integrated, protective shield that allows for one-handed activation immediately after use to help minimize the risk of needlestick injuries and to offer clear visibility of blood flashback.

#367294, 7" Tube with Luer, 25 ga x 3/4"	
(987-5250)	50/box
#367292, 7" Tube with Luer, 23 ga x 3/4"	
(987-0529)	ea
#367287, 7" Tube with Luer, 21 ga x 3/4"	
(958-5669)	200/case
#367285, 12" Tube with Luer, 25 ga x 3/4"	
(840-3447)	ea
#367283, 12" Tube with Luer, 23 ga x 3/4"	
(987-1301)	50/box
#367281, 12" Tube with Luer, 21 ga x 3/4"	
(987-8393)	50/box
#367298, 12" Tube without Luer, 25 ga x 3/4"	
(623-4670)	50/box
#367297, 12" Tube without Luer, 23 ga x 3/4"	
(987-4692)	50/box
#367296, 12" Tube without Luer, 21 ga x 3/4"	
(987-9967)	50/box
#368653, 12" Tube with Preattached Holder, 23 ga	
(987-8940)	25/box
#368652, 12" Tube with Preattached Holder, 21 ga	
(111-7212)	200/case


**BD™ Microtainer®
Contact-Activated
Lancets**

Designed with a positive patient experience in mind. The contact activation method facilitates a consistent puncture depth and minimizes the likelihood of having to repeat the puncture. It covers only a small area at the contact point, resulting in improved visibility of the puncture site for the clinician and greater accuracy of lancet positioning when performing the puncture. Its innovative ergonomic design allows for a more comfortable grip. The lancet automatically retracts into the device, which prevents the lancet from being reused. In addition, the lot number is laser-etched on each lancet for easier tracking.

#366592, Low Blood Flow (Single Drop), Purple, 30-ga x 1.5-mm-D Needle	
(106-5174)	200/box
#366593, Medium Blood Flow, Pink, 21-ga x 1.8-mm-D Needle	
(106-5210)	200/box
#366594, High Blood Flow, Blue, 1.5mm-W x 2-mm-D Blade	
(106-5209)	200/box


**BD™ Vacutainer®
Holder/Adapters**

These adapters and holder provide a quality, low-cost, single-use product that can help customers to comply with OSHA regulations. Compatible with BD Vacutainer Eclipse™ Blood Collection Needles, BD Vacutainer Safety-Lok™ Blood Collection Sets, BD Vacutainer Push-Button Blood Collection Sets, and BD Vacutainer Multi-Sample Luer Adapter.

#364815, One-Use Holder	
(987-0051)	250/bag
Complete Kit	
#364597, Ribbed Pediatric Tube Adapter	
(991-3989)	10/box


**BD™ Vacutainer® Plus Plastic Blood-
Collection Tubes**

Offer a safer method of blood collection. Reduce risk of tube breakage and specimen spillage. Safe, simple disposal in accordance with EPA guidelines.

#367981, 3.5-mL SST™ Tube with Clot Activator/Gel for Serum Separation, Red/Gray Conventional Top, Paper Label, Silicone-coated Interior, 13 mm x 75 mm	
(987-1419)	100/box
#367856, EDTA Tube/K2EDTA, Lavender Hemogard™ Top, 3 mL, Paper Label, 13 mm x 75 mm	
(127-4412)	100/box
#367861, EDTA Tube/K2EDTA, Lavender Hemogard™ Top, 4 mL, Paper Label, 13 mm x 75 mm	
(987-7504)	100/box


**BD Vacutainer®
Push-Button
Blood Collection
Sets**

Set offers split-second protection for that single moment that could potentially change your life. Push-button safety mechanism instantly helps protect you against needlestick injury. BD offers the only blood-collection system with push-button technology. Set is designed to be simple, easy to use, ergonomic and, most important, helpful in protecting you from accidental needlestick injuries.

Collection Sets.....50/box

Specify:

#367323, without Luer Adapter, 12" Tubing, 25 ga x 3/4"	
(987-3303)	
#367324, Light Blue without Luer, 12" Tubing, 23 ga x 3/4"	
(987-2637)	
#367326, Green without Luer, 12" Tubing, 21 ga x 3/4"	
(987-9975)	
#367341, Royal Blue with Luer, 12" Tubing, 25 ga x 3/4"	
(987-9168)	
#367336, Light Blue with Luer, 7" Tubing, 23 ga x 3/4"	
(987-6287)	
#367338, Green with Luer, 7" Tubing, 21 ga x 3/4"	
(987-0431)	
#367342, with Luer, 12" Tubing, 23 ga x 3/4"	
(919-4247)	50/box
#367344, with Luer Adapter, 12" Tubing, 21 ga x 3/4"	
(987-3358)	50/box
#367344, with Luer, Green, 12" Tubing, 21 ga x 3/4"	
(405-4902)	200/case
#367341, with Luer, Royal Blue, 12" Tubing, 25 ga x 3/4"	
(679-6358)	200/case



#366643, K2EDTA, 10 mL, 16 mm x 100 mm	
(286-9036)	1000/case
#367988, 8.5-mL SST™ Tube with Clot Activator/Gel for Serum Separation, Red/Gray Conventional Top, Paper Label, Silicone-coated Interior, 16 mm x 100 mm	
(987-2786)	100/box
#364300, 2-mL Glass Tube/K ₂ EDTA, Lavender Conventional Top, Paper Label, 10.25 mm x 64 mm	
(908-1970)	1000/case



Surshield™ Safety Wing Blood-Collection Sets

One-handed, easy-to-activate safety device. Requires no change in needle insertion technique. Built-in safety design is easy to use. Latex-free.

Long, 21 ga x 3/4"50/box
(315-7504)

Long, 23 ga x 3/4"50/box
(315-0039)

12" Tubing
Long, 25 ga x 3/4"50/box
(315-0434)

7" Tubing
Short, 21 ga x 3/4"500/case
(315-5499)

Short, 23 ga x 3/4"50/box
(315-6596)



Capiject®

- Extra-wide, built-in collection lip for enhanced specimen flow
- Single-piece, tapered design for easy collection and handling
- Color-coded fill lines and caps for volume accuracy and easy identification
- Lot number and expiration date printed on each tube for easy identification
- Silicone-coated anticoagulant tubes designed to reduce clotting
- Olefin-based gel for clean separations
- Tube nests in cap, allowing tube to stand alone
- Certified for lead testing
- Latex-free

#T-MG, Cherry Red, Gel Barrier with Clot Activator (Silica Particles and Inert Gel)250/box
(315-9272)

#T-MQK, Lavender, K2-EDTA (0.78 mg)250/box
(315-0280)

#T-M, Brick Red, No Additive/Clot Tube1000/case
(335-9969)



Capiject® Safety Lancets—Blue Blade

1.5 mm x 2.0 mm, Blue200/box
(110-8852)



#T-MLH, Green, Lithium Heparin (12.5 IU)250/box
(315-0873)

#T-MLHG, Mint Green, Lithium Heparin with Gel (12.5 IU)250/box
(315-1116)

#T-MPS, Gray, Sodium Fluoride (0.72 mg), Potassium Oxalate (0.58 mg)1000/case
(476-8940)

#T-MGA, Gold-Amber Tube with Gel Barrier and Clot Activator (Silica Particles and Inert Gel)1000/case
(609-2050)

#T-MLHGA, Mint Green-Amber Tube with Lithium Heparin and Gel Barrier (12.5 IU)1000/case
(315-7531)

1.5mm x 1.5mm, Green200/box
(110-8853)

1.5 mm x 1.0 mm, Pink200/box
(110-8854)

28 ga, 1.25 mm, Purple200/box
(110-8541)

23 ga, 2.25 mm, Orange200/box
(110-8539)

Making Blood Collection Safe, Precise, Simple

Capiject® Safety Lancets

High quality, single-use lancets offer the form, fit and function you've come to rely on.

- Never see or handle the sharp
- One-step activation
- No special staging required
- Color-coded sizes



Schein No.	Catalog No.
1108854	200101
1108853	200102
1108852	200103
1108539	200104
1108541	200105

Surshield® Safety Blood Collection Sets

Help shield your staff against accidental needlesticks with Surshield.

- Simple, one-handed activation
- Audible click for confirmation that needle is securely locked
- Available in 21G and 23G in 7" and 12" tubing, and 25G with 12" tubing.



Schein No.	Catalog No.
3155499	1MN*SVS21B18
3157504	1MN*SVS21B30
3156596	1MN*SVS23B18
3150039	1MN*SVS23B30
3150434	1MN*SVS25B30

CAPJECT is a trademark of Terumo Medical Corporation, Somerset, New Jersey. SURSHIELD, TERUMO and Terumo are registered trademarks of Terumo Corporation, Tokyo, Japan. © 2013 Terumo Medical Corporation 8/13. All rights reserved. 08052013-2786.



TERUMO MEDICAL PRODUCTS



UNICO® Tube-Cubes®

Available in 3 configurations to accommodate practically all of your needs. Each cube measures 2 5/8" H x 2 5/8" W x 2 5/8" D. The 9-place version holds 8-mm to 16-mm/10-mL to 15-mL tubes. The 4-place version holds 26-mm to 30-mm/50-mL tubes. The single place version holds 66-mm cups, alcohol dispensers, small flasks, and more.

- Invert and pour: tubes will stay in place
- Great for phlebotomy trays
- Organize and store samples in refrigerator or freezer

#54900, 9 Place, for 8-mm–16-mm Tubes (112-4212)4/pkg
#54400, 4 Place, for 26mm–30-mm Tubes (112-4210)4/pkg
#54100, Single Place, for Up to 66-mm Tubes (112-4214)4/pkg



UNICO® Tube Rack

Starting with a single 50-place T-Rack, you can attach as many racks as space allows. These patented racks attach together quickly and easily with their unique built-in T-connector, which allows the flexibility to attach as many racks as needed.

- Each rack holds up to 50 tubes of varying sizes in two 25-hole sections
- Top shelf has cross-referencing alphanumeric designation in order to track exact location of specific tubes in rack
- Attach additional T-Racks together; make a single large rack as big as needed
- Made of durable, autoclavable polypropylene; will not break or become brittle in low-temperature refrigerators
- Easy to clean and disinfect
- Arrives disassembled and can be snapped together in seconds
- Dimensions: 8 3/4" x 3 3/4" x 5 1/4"


#54500, 50 Place (Expandable) (111-7761)ea



UNICO® Tube-Cube® Carrier

Durable acrylic, and designed to hold any combination of Tube-Cubes. Available 4-place configurations. Provides a safe, convenient way to carry patient samples to the lab.

#54040 (112-4213)ea



SmoothRack Test Tube Racks

Features include:

- 72-tube capacity
- 9.7" x 5.0" x 2.6"
- Autoclavable polypropylene
- Accommodates test tubes 10 mm to 16 mm in diameter

Blue (720-0000)ea



BRAND PROMISE

QUALITY SELECTION PERFORMANCE

Rely on Us

**SAVE WITHOUT COMPROMISING ON EXCELLENCE,
WITH OUR OWN HIGH-QUALITY HENRY SCHEIN BRAND**



HENRY SCHEIN® LENS PAPER

For fine lenses and optical instruments.

- Special high wet strength
- Soft and dust-free
- Will not scratch glass
- Ideal if liquid cleaner is used

4" x 6" (100-8362)50/pkg

HENRY SCHEIN® SPECIMEN CONTAINERS

These 4.5-oz specimen containers have blue, tamper-evident, leak-resistant tops and are suitable for routine urine or specimen collection and transport. Container includes patient information label.

Sterile, with Tab, Unwrapped (900-7433)100/case
Sterile, Individually Wrapped, 4 oz (900-7431)100/case
Non-Sterile, Bulk (900-7432)500/case



HENRY SCHEIN® GLASS MICROSCOPE SLIDES

22 mm x 22 mm; thickness: 0.17 mm–0.25 mm

Plain Slide (102-7277)72/box
Frosted End Slide (102-5605)72/box
Cover Glass #2, 1-oz Package (102-6803)ea

**Serving you Quality, Selection and Performance.
It's our promise to you.**

UNICO® PowerSpin™ FX Centrifuge



Features a 6-place angled rotor designed for 2–10-mL tubes (maximum tube size up to 16 mm x 110 mm) and a fixed speed of 3400 rpm, generating 1350 RCF (g) for superior specimen separation. Set the 0–30-minute timer to desired run length.

- Maintenance-free, brushless motors run quietly at a fixed speed of 3400 rpm, emitting only 58dB, all with a low temperature rise
- Suction-cup rubber feet hold unit firmly in place while virtually eliminating vibration
- Lid features a gasket to reduce aerosol emissions, a window for strobe speed checks, and a safety interlock switch that cuts power to motor if lid is raised
- Accessories available on request: 8-place rotor and pediatric tube adapters
- 1-year manufacturer's warranty
- Dimensions: 10.2" x 11.8" x 7.9"
- Weight: 15 lb
- 115 V/60 Hz

#C806, Fixed Speed, 6 Place
(715-5981)ea



UNICO® PowerSpin™ LX Centrifuge

Spin blood, urine, fecal materials, or semen at the best speed for each type of specimen with this powerful, variable-speed, 6-place centrifuge. Maintenance-free, brushless motor runs quietly at a fixed speed up to 4000 rpm, emitting only 58dB, all with low-temperature rise. Suction-cup rubber feet hold unit firmly in place and virtually eliminate vibration. Lid features a gasket to reduced aerosol emission, a window for strobe speed check, and a safety interlock switch that cuts power to the motor if lid is raised. Angled rotor accepts 2–10-mL tubes (maximum tube size 16 mm x 110 mm). Accessories available upon request: 8-place rotor and pediatric tube adapters.

- Speed range: 300–4000 rpm, generating up to 1850 RCF (g)
- 30-minute timer
- 1-year warranty
- Dimensions: 10.2" x 11.8" x 7.9"
- Weight: 15 lb
- 115 V/60 Hz

#C-856, Variable Speed, 6 Place
(715-2656)ea

UNICO® Micro-Hematocrit Centrifuge



This powerful centrifuge, with all-metal housing, spins up to 24 microhematocrit tubes at 12,000 rpm. Tubes are held in place with a screw-down metal lid and replaceable rubber gasket/ring. The 15-minute mechanical timer has a bell tone to signal a completed run, while the electronic brake can stop the unit on demand. Metal lid features a safety-lid latch that cuts engine power if lid is lifted. Unit comes complete with a metal hematocrit reader tray. 1-year warranty.

#C-MH30
(715-0877)ea

UNICO® Vortex Mini Mixer



This economical vortex mixer is rugged and sturdy for the most demanding applications. Dual-mode switch lets you vortex continuously or intermittently using the touch mode: simply press your tube to the soft rubber mixing cup. Features speed control to adjust the vortexing velocity to your requirements.

- Powerful and quiet
- Variable speed up to 3000 rpm
- Switch from Auto to Touch mode
- Power On LED indicator
- 1-year limited warranty

#L-VM1000
(110-0228)ea

UNICO® Rock-IT™ Jr. Blood-Tube Rocker



A versatile rocker for hematology (blood) tubes prior to counting; also excellent for vials and gentle mixing of injectables. Nearly silent while gently completing 25 mixing cycles per minute. Platform is designed for 8 tubes and has a reversible pad that can be autoclaved. 1-year warranty.

Dimensions: 8.6" x 4.7" x 5"

#LTTR100, 8 Place
(159-0730)ea

UNICO® Rock-IT™ Blood Tube Rocker



A versatile rocker for hematology (blood) tubes prior to counting; also excellent for vials and gentle mixing of injectables. Nearly silent while gently completing 25 mixing cycles per minute. Platform is designed for 15 tubes and has a reversible pad that can be autoclaved. 1-year warranty.

Dimensions: 13.7" x 4.7" x 5"

#L-TTR200, 15 Place
(715-2650)ea

UNICO® IP730 Series Microscopes



Designed to provide crisp, clear imaging on a versatile and ergonomic platform. Available with 30° inclined binocular or trinocular heads. Features Seidentopf optical design with 55mm–75mm interpupillary distance adjustment. Both heads rotate 360° and lock in place with a thumbscrew. Eyepieces are standard with 10X 22mm extra-wide field eyepieces for viewing comfort and clarity.

#IP730, Advanced Binocular Microscope, Achromat
Objectives: 4X, 10X, 40X, 100X Achromat
(106-9122)ea

#IP730PL, Advanced Binocular Microscope, Plan
Objectives: 4X, 10X, 40X, 100X Plan
Illumination: 6 V/20 W Halogen
(109-9806)ea

UNICO® G380 Microscope



Designed for viewing comfort, this microscope provides crisp, clear imaging on a versatile, ergonomic platform. Available in a variety of optics choices and body configurations. Standard 10X High Eye-Point eyepieces are engineered to reduce eyestrain for users who wear glasses. Each eye tube features independent focusing diopter adjustment to compensate for the natural vision differential between each eye. Base design adds stability to the body and allows extra room for hand access to the oversized coarse-focus and fine-focus controls and easy access to the coaxial mechanical stage controls.

(281-4869)ea
Contains: dust cover & 2 extra fuses.

BOEKEL® Boekel Scientific #132000 Incubator



An excellent incubator for all in-office culture needs. Large capacity (0.8 cu ft) for busiest physicians' offices.

- Interior dimensions: 11½"W x 10½"H x 12"D; exterior dimensions: 13"W x 16¼"H x 14"D
 - See-through acrylic plastic door
 - Magnetic door latch
 - Pilot light and on/off switch
 - Adjustable nontip shelves
 - Maximum temperature: 60°C
- (100-8790)ea

Diamond™ White Glass Microscope Slides

Clearly Superior.

Diamond™ White Glass Microscope Slides are super high quality slides that compare to the top brands . . . but at substantial savings. Our slides are produced from water white glass sheets and are so clear . . . *you will see the sample not the slide.*

Diamond™ White Glass Microscope Slides are ideal for all routine applications, automated instruments, long-term storage, fluorescence studies and where accuracy and the highest optical quality are of paramount importance.

Available in plain, frosted, color coded, charged and instrument specific versions. All slides are double cellophane packaged to protect against moisture and foreign particles.




Producers of Exceptional Quality Laboratory Supplies

610 Winters Avenue • Paramus, NJ 07652 • USA
www.globescientific.com



For more information or to request FREE samples, contact your Henry Schein representative



Microscope Cover Glass

Vacuum Pack, 10 oz

- #1401-10, 18 mm x 18 mm, #1 Thickness (424-0016)10/box
- #1401-20, 18 mm x 18 mm, #2 Thickness (424-0017)10/box
- #1404-10, 22 mm x 22 mm, #1 Thickness (424-0020)10/box
- #1415-10, 24 mm x 50 mm, #1 Thickness (424-0034)10/box
- #1419-10, 24 mm x 60 mm, #1 Thickness (424-0039)10/box



Microscope Slides, Glass

25 mm x 75 mm, 90° Ground Edges, Safety Corners

Premium-quality microscope slides are precleaned and ready to use. They feature finely-ground edges and slightly clipped corners for comfortable handling. Each box of slides contains a moisture-absorbent sheet to prevent sticking.

Plain (Colorless)

#1321, Unfrosted72/box

#1324, Frosted

(803-6390)72/box

Frosted, Colors—Color on 1 End/1 Side

#1324W, White (868-6492)72/box



#1358W, White, Adhesive (868-0134)72/box

#1324N, Orange (424-0005)1440/case

#1324Y, Yellow (424-0007)1440/case

#1324P, Pink (424-0006)1440/case

#1324B, Blue (424-0003)1440/case

#1328, 2 Sides, Frosted (868-0123)72/box

Ground Edges, 25 mm x 75 mm x 1 mm

#1301, Plain (424-0000)1440/case

#1304, Frosted (1 End, 1 Side) (424-0002)1440/case



Microscope Slides—Glass

25 mm x 75 mm

Frosted/colored, 1 end, both sides, beveled edges, clipped corners.

#1331, Plain (424-0008)1440/case
Contains: 72/box, 20 boxes/case (10 gross).



#1334, Frost (424-0009)1440/case
Contains: 20 boxes/case (10 gross)

#1334B, Blue (424-0010)1440/case
Contains: 20 boxes/case (10 gross)

#1334G, Green (424-0011)1440/case
Contains: 20 boxes/case (10 gross)

#1334N, Orange (424-0012)1440/case
Contains: 20 boxes/case (10 gross)

#1334P, Pink (424-0013)1440/case
Contains: 20 boxes/case (10 gross)

#1334Y, Yellow (424-0014)1440/case
Contains: 20 boxes/case (10 gross)

#1334W, White (109-9262)1440/case
Contains: 20 boxes/case (10 gross)

#1338, Double Frost (424-0015)1440/case
Contains: 20 boxes/case (10 gross)



Research-Grade Centrifuge Tubes

- Acrylic spin rate is 5,000 rpm
- Polypropylene spin rate is 10,000 rpm
- All caps are made from high-density polyethylene (HDPE)

#6267, Acrylic, Non-sterile, Bulk
(102-6540)500/case

#6264, Polypropylene, Non-sterile, Bulk
(280-5289)500/case


#6270, Acrylic, Sterile
(185-8503)500/case
Contains: 25 per bag.



#6273, Acrylic, Sterile
(424-7963)500/case
Contains: 50 per rack.

#6269, Polypropylene, Sterile
(424-6254)500/case
Contains: 25 per bag.

#6271, Polypropylene, Sterile
(243-5517)500/case
Contains: 50 per rack.



Urine Centrifuge Tubes

PS, Graduated, 12 mL

These 12-mL flared-top urine centrifuge tubes are made of high-quality polystyrene (PS) and are available in Ames-type, Kova-type and Whale-type styles. Blue and yellow snap caps are also available.

- Tube dimensions: 20 mm x 115 mm
- Polystyrene spin rate: 1500 rpm

#112017, Ames-Type
(120-9178)500/pkg

#112015, Kova Type
(424-0502)2000/case

#112030, Sediment Bulb, Whale-Type
(424-0054)1500/case

#113136, Blue Snap Caps for 12-mL Flared-Top Centrifuge Tubes
(120-6724)500/pkg

#113137, Yellow Snap Caps with Sanitary Grip for 12-mL Flared-Top Centrifuge Tubes
(120-8588)500/pkg



Specimen Transport Bag

These specimen bags are available with a zipper closure or tamper-evident glue seal. The bags are printed with a biohazard symbol and feature an attached document pouch that keeps important paperwork safe.

#4918, 6" x 9" Bag with Document Pouch
(106-5227)1000/case



Research-Grade Centrifuge Tubes

50 mL, 30-mm x 116-mm Conical Centrifuge Tubes with Screw Caps

- Acrylic spin rate is 5000 rpm
- Polypropylene spin rate is 10,000 rpm
- All caps are made from high-density polyethylene (HDPE)

#6244, Polypropylene, Non-sterile, Bulk
(424-5618)500/case

#6242, Polypropylene, Non-sterile
(424-9912)500/case
Contains: 25 per rack.

#6240, Polypropylene, Non-sterile
(424-8018)500/case
Contains: 25 per bag.

#6243, Polypropylene, Sterile
(424-8999)500/case
Contains: 25 per rack.

#6241, Polypropylene, Sterile
(424-5191)500/case
Contains: 25 per bag.



Transport Tubes, Polypropylene

Safely store and transport biological materials with self-standing polypropylene tubes. The polyethylene leak-resistant caps can be removed or sealed in a single turn. Molded graduation marks measure 1.5 mL to 14 mL.

- Tubes are self-standing with conical bottoms
- Available with caps attached or packaged separately

#6101, 5 mL, Non-sterile, with Separate Screw Cap
(110-6203)1000/case

#6101C, 5 mL, Non-sterile, with Attached Screw Cap
(356-3969)1000/case

#6101S, 5 mL, Sterile, with Attached Screw Cap
(424-6233)500/case

#6102, 10 mL, Non-sterile, with Screw Cap
(120-1971)1000/case



Transfer Pipettes

These plastic transfer pipettes are made from high-quality, flexible polyethylene and are completely unbreakable. They are offered in a variety of sizes and tip styles and are perfect for all liquid transfer applications.

- Packaged in convenient dispenser boxes
- Available in sterile and non-sterile versions
- Nontoxic and inert

#137038, 5 mL, Sterile, Graduated to 1 mL
(606-9661)400/case

#135038, 7 mL, Sterile, Graduated to 3 mL
(990-1310)400/case

#137035, 3 mL, Non-sterile, Graduated to 1 mL, Dispenser Box
(185-8544)500/box

#136030, 4.5 mL, Non-sterile, Narrow Stem, Nongraduated, Dispenser Box
(904-6934)500/box

#137030, 5 mL, Non-sterile, Graduated to 1 mL, Dispenser Box
(120-9438)500/box

#135030, 7 mL, Non-sterile, Graduated to 3 mL, Dispenser Box
(120-8751)500/box



#6102C, 10 mL, Non-sterile, with Attached Screw Cap
(120-6506)1000/case

#6102S, 10 mL, Sterile, with Attached Screw Cap
(945-1766)500/case

#6255, 50 mL, Non-sterile, with Screw Cap, Bulk
(424-6071)500/case

#6256, 50 mL, Sterile, with Attached Screw Cap, Bulk
(999-6888)500/case



globe
SCIENTIFIC INC.

Borosilicate Glass Culture Tubes

Globe Scientific introduces a new line of high-quality borosilicate glass culture tubes. These tubes offer exceptional quality at great prices. They are neatly packaged in heavy-duty boxes of 250 and are shrink wrapped to keep them particle-free. The boxes are compact and designed to fit into most test-tube dispensers.

3 mL, 10 mm x 75 mm (424-0058)	250/box
5 mL, 12 mm x 75 mm (424-0059)	250/box
7 mL, 13 mm x 100 mm (116-7973)	1000/case
10 mL, 16 mm x 100 mm (424-0060)	250/box
12 mL, 16 mm x 125 mm (424-0062)	1000/case
15 mL, 16 mm x 150 mm (424-0061)	1000/case

globe
SCIENTIFIC INC.

Formalin Filled Containers

Save time and reduce formaldehyde exposure and spills with prefilled formalin containers designed for the safe collection, preservation, transportation and storage of specimens. Produced from polypropylene, choose from highly leak-resistant Tite-Rite™ containers or leak-resistant large containers, both prefilled with 10% neutral buffered formalin.

Tite-Rite™ Containers

#6518FL, 20 mL (0.66 oz) (396-2441)	100/case
#6520FL, 40 mL (1.32 oz) (120-7084)	100/case
#6522FL, 60 mL (2 oz) (120-7083)	100/case



#6525FL, 90 mL (3 oz), Wide Mouth (120-6449)	100/case
---	----------

Large Containers

#6387FL, 240 mL (8 oz) (970-0586)	90/case
#6397FL, 360 mL (12 oz) (177-5659)	48/case
#6399FL, 480 mL (16 oz) (120-1863)	48/case

globe
SCIENTIFIC INC.

Urine Collection Containers and Cups

The 4-oz (120-mL) containers are available in both sterile and non-sterile versions. Non-sterile containers have a frosted writing area for patient information. Sterile containers are individually wrapped and have a trilingual patient ID label. They also have easy-to-turn, leak-resistant screw caps and molded graduation marks. The 6.5-oz (190-mL) cups come in both polypropylene and polystyrene versions and are available with or without paper lids.



#5912, Trilingual Patient ID Label, Sterile, Individually Wrapped (606-6028)	100/case
#5913, 4-oz Container, Sterile, Individually Wrapped (868-0112)	100/case
#5915, 4-oz Container, Non-sterile, Etched Writing Space, Bulk (642-6596)	500/case
#5916, 4-oz Container, Non-sterile, with Affixed Thermometer Strip (413-4578)	500/case
#5920, 6.5-oz Specimen Collection Cup, PP (868-8877)	25/pack

*Producers of
exceptional quality
laboratory supplies.*

EXTENSIVE SELECTION - ONE TRUSTED SOURCE

Analyzer Consumables	Microhematocrit Tubes
Beakers, Bottles, Cylinders	Microscope Slides & Accessories
Biohazard Specimen Bags	Pathology Containers
Blood Collection Products	Pipettors & Pipette Tips
Centrifuge Tubes	Serological Pipettes
Collection Cups	Specimen Containers
Cytology Products	Test Tube Racks
Drug Test Containers	Test Tubes, Vials & Caps
ESR Systems	Transfer Pipets
Formalin Filled Containers	Transport Tubes
Funnels	Urinalysis Products
Histology Cassettes	Weighing Dishes

For more information or to request FREE samples,
contact your Henry Schein representative



610 Winters Avenue • Paramus, NJ 07652 • USA
www.globescientific.com

EMBEDDING MEDIA:

(102-3734)	Cryo O.C.T.	4 oz.
(116-1233)	Paraplast® Plus	8x1kg /Case
(106-4135)	Tissue-Tek® O.C.T. Compound	4 oz.

(Frozen Tissue Matrix)

TISSUE FREEZING PREPARATIONS:

(104-3248)	Cytocool® II	11 oz.
(406-1248)	Flash Freeze Spray	11 oz.
(106-3127)	Neg-50 Clear	4oz.
(758-1368)	Neg-50 Blue	4 oz. x 2/Box
(106-7137)	Tissue Freezing Media Clear	4 oz. x 12/Cs.

MOUNTING MEDIA:

(334-7213)	Cytoseal® 60 Low Viscosity	4 oz.
------------	----------------------------	-------

(Toluene-Based For Use with Cover Glass)

(122-6901)	Permount®	100 mL
------------	-----------	--------

(Toluene-Based For Hand Coverslips)

(106-7123)	SHUR/Mount™	12 x 4 oz./Case
------------	-------------	-----------------

(Xylene-Based For Automated Coverslips)

FIXATIVES:

(108-3646)	Acetic Acid 10%	16 oz.
(120-1795)	Glacial Acetic Acid	16 oz.

(Also Known As Ethanoic Acid)

(120-0443)	Acetone	16 oz.
(120-2827)	Acetone	Gallon
(120-8724)	Ethanol Alcohol 95% 200 Proof	Gallon

(Also Known As Ethanol, Denatured Alcohol, Algarin Denatured, Ethyl Hydrate Denatured, Grain Alcohol Denatured)

(364-6848)	Ethanol 100%	16 oz.
(600-4289)	Formalin Aceto-Alcohol	Gallon

(FAA Solution or Carnoy's Solution)

(120-8911)	Formalin, 10% Buffered	10 Liter
------------	------------------------	----------

(NBF or Neutral Buffered Formalin)

(110-5855)	Pen Fix	4 Gallon/Case
(106-3078)	SafeFix II	Gallon

(Formalin Substitute)

(234-6266)	Zinfix® Zinc Formalin Fixative	Gallon
------------	--------------------------------	--------

DEHYDRATING REAGENTS:

(113-7936)	Flex 100% Alcohol	4 Gallons/Case
------------	-------------------	----------------

(Premixed Blend of Methyl and Isopropyl Alcohol For H&E Staining I)

(120-9026)	Ammonium Hydroxide	500 mL
------------	--------------------	--------

(For H&E Staining)

(120-7599)	Deionized Water (Di)	20 Liters/Box
------------	----------------------	---------------

BLUING AGENTS:

(608-5631)	Substitute for Scott's Tap Water	Gallon
(106-2876)	Bluing Reagent	Gallon

(For H&E Staining)

HISTOLOGY STAINS:

(106-3048)	Clarifier 2 Reagent	Gallon
(799-1780)	Clarifier 2 Reagent	Gallon

(For Improved Contrast & Increased Cell Transparency)

(831-6922)	Diff-Quik Stain Set	3 x 500 mL
------------	---------------------	------------

(Also known as Romanowski Stain)

(113-4902)	Eosin-Y Stain Solution 1%	16 oz.
(906-6070)	Eosin-Y Working Solution	Gallon

(Also Eosin Y & Alcohol)

(652-8573)	Eosin-Y Stain Solution 1%	32 oz.
------------	---------------------------	--------

(Also Eosin Y & Alcohol)

(424-6155)	Harris Hematoxylin Stain	500 mL
------------	--------------------------	--------

(Mercury-Free)

(398-4749)	Gills I Hematoxylin	Gallon
(106-2893)	Gills II Hematoxylin2	4 Pints/Cs.
(487-7361)	Gills III Hematoxylin	33.8 oz.
(230-4808)	Gills III Hematoxylin	Gallon
(106-7191)	Tissue Marking Dye, 6 Colors	18/Kit
(925-7517)	Toluidine Blue 0.1%	8 oz.
(925-4103)	Toluidine Blue 1%	16 oz.

(Also Called Diff-Quik or Romanowski Stain)

DECALCIFYING SOLUTION:

(120-1248)	HCl-EDTA Rapid Decal	Gallon
------------	----------------------	--------

CLEARING AGENTS:

(111-8126)	Clear-Rite 3™ Solvent	4 Gallons/Case
------------	-----------------------	----------------

(Can Be Used In All Tissue Processors and Automatic Stainers)

(357-0333)	d-Limonene, Purified	4 Gallons/Case
------------	----------------------	----------------

(Citrus-Based Xylene Substitute)

(120 -3130)	Safe-Clear	Gallon
(114-7751)	Histo-Clear™	Gallon
(120-4629)	Xylene	Gallon

TRANSPORT & LAB SUPPLIES:

(369-4647)	Petri Dish 4" (100mm)	25/Pkg.
(122-7441)	Petri Dish 4" (100mm) Glass	12/Pkg.
(115-1911)	Prefilled Formalin Jar	96 x 20 mL/Case
(115-1912)	Prefilled Formalin Jar	96 x 40 mL/Case

We stock over 100,000 products, including a full line of Histology and Cytology Containers: Tissue and Biopsy Cassettes, Slide Storage/Staining/Mailing, Sample containers, and more. If you don't see what you need for your practice, call your Henry Schein Sales Consultant or visit: www.henryschein.com/mohs




Simport

Histosette® I Tissue Processing/Embedding Cassette

Disposable plastic cassettes hold tissue specimens during the embedding process, as well as in a storage file. Molded from a special high-density polymer, these cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage. The one-piece integral lid eliminates the need for separate steel lids. Just snap apart and conveniently lock the lid into the base of the cassette. Cassettes can be opened or closed as often as necessary and always relock securely without danger of specimen loss. Large labeling areas are located on three sides of the cassette. Anterior writing area at a 45° angle. The more acute angle makes cassettes more suitable to be used with some models of cassette-labeling instruments.

Cassettes 1500/case



Specify:

#M498-10, Lilac	(115-1678)
#M498-11, Orange.....	(115-1679)
#M498-12, Aqua	(115-1680)
#M498-2, White.....	(115-1681)
#M498-3, Pink.....	(115-1682)
#M498-4, Green.....	(115-1683)
#M498-5, Yellow	(115-1684)
#M498-6, Blue.....	(111-6103)
#M498-7, Peach.....	(115-1685)
#M498-8, Tan.....	(115-1686)
#M498-9, Gray.....	(115-1687)



Simport

Histosette® II Tissue Processing/Embedding Cassettes

Disposable plastic cassettes hold tissue specimens during processing and embedding, as well as in storage. Molded from a special high-density polymer, these patented cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of histological solvents. The efficient flowthrough slots maximize fluid exchange and ensure proper drainage. The one-piece integral lid eliminates the need for separate steel lids. Cassettes can be opened and closed as often as necessary and will always relock securely without danger of specimen loss. Large labeling areas are located on three sides of the cassettes for your convenience.

Cassettes 1500/case

Specify:

#M492-10, Lilac	(115-1649)
#M492-11, Orange.....	(115-1650)
#M492-12, Aqua	(115-1652)
#M492-2, White.....	(115-1653)
#M492-3, Pink.....	(115-1654)
#M492-4, Green.....	(115-1655)
#M492-5, Yellow	(115-1656)
#M492-6, Blue.....	(114-8888)
#M492-7, Peach.....	(115-1657)
#M492-8, Tan.....	(115-1659)
#M492-9, Gray.....	(115-1660)

Simport

Histosette® I Tissue Processing/Embedding Cassettes

Disposable plastic cassettes hold tissue specimens during the embedding process, as well as in a storage file. Molded from a special high density polymer, these cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage. The one-piece integral lid eliminates the need for separate steel lids. Just snap apart and conveniently lock the lid into the base of the cassette. They can be opened or closed as often as necessary and they always relock securely without danger of specimen loss. Large labeling areas are located on three sides of the cassette for your convenience. Anterior writing area at a 30° angle. Not suitable for automated printers.

Cassettes 1500/case



Specify:

#M490-10, Lilac	(641-3848)
#M490-11, Orange.....	(630-5451)
#M490-12, Aqua	(642-0225)
#M490-2, White.....	(120-5194)
#M490-3, Pink.....	(600-9085)
#M490-4, Green.....	(600-2992)
#M490-5, Yellow	(340-4586)
#M490-6, Blue.....	(104-5652)
#M490-7, Peach.....	(115-1636)
#M490-8, Tan.....	(115-1637)
#M490-9, Gray.....	(115-1639)

Simport


Histosette® II Biopsy Processing/Embedding Cassette

Similar to Histosette® II Tissue Processing Cassette but specially designed to hold biopsy specimens during the embedding process, as well as in a storage cabinet. Anterior writing area is at a 45° angle to make the cassette more suitable to be used with certain types of cassette labeling instruments.

Cassettes 1500/case

Specify:

#M493-10, Lilac	(115-1661)
#M493-11, Orange.....	(115-1662)
#M493-12, Aqua	(115-1663)
#M493-2, White.....	(115-1666)
#M493-3, Pink.....	(115-1667)
#M493-4, Green.....	(115-1670)
#M493-5, Yellow	(115-1671)
#M493-6, Blue.....	(115-1672)
#M493-7, Peach.....	(115-1673)
#M493-8, Tan.....	(115-1674)
#M493-9, Gray.....	(115-1675)



Simport

Histosette® II Tissue Processing/Embedding Cassette

Most convenient for cassette-labeling instruments since covers are already removed from cassettes and are packaged separately in the case. Disposable plastic cassettes hold tissue specimens during processing and embedding, as well as in storage. Molded from a special high-density polymer, these patented cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage.

Cassettes1000/case

Specify:

#M485-10, Lilac	(115-1613)
#M485-11, Orange.....	(115-1614)
#M485-12, Aqua	(115-1615)
#M485-2, White.....	(115-1616)
#M485-3, Pink.....	(115-1617)
#M485-4, Green.....	(115-1618)
#M485-5, Yellow	(115-1619)
#M485-6, Blue.....	(115-1620)
#M485-7, Peach.....	(115-1621)
#M485-8, Tan.....	(115-1622)
#M485-9, Gray	(115-1623)



Simport


Unisetite™ Biopsy Processing/Embedding Cassette

Molded from a special high-density polymer, these cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage. The one-piece integral lid eliminates the need for separate steel lids. The snap-latch and hinge-lock design of the Unisetite™ prevents early separation of base and lid and allows for one-hand operation. Lids can be opened and closed as often as necessary and always relock securely without danger of specimen loss. Large labeling areas are located on two sides of the cassettes. Anterior writing area is at a 35° angle. Specially designed to hold biopsy specimens during the embedding process.

Cassettes1500/case

Specify:

#M506-10, Lilac	(115-1722)
#M506-11, Orange.....	(115-1723)
#M506-12, Aqua	(115-1724)
#M506-2, White.....	(115-1725)
#M506-3, Pink.....	(115-1726)
#M506-4, Green.....	(115-1727)
#M506-5, Yellow	(115-1728)
#M506-6, Blue.....	(115-1729)
#M506-7, Peach.....	(115-1730)
#M506-8, Tan.....	(115-1731)
#M506-9, Gray	(115-1732)



Simport

Unisetite™ Tissue Processing/Embedding Cassette

Molded from a special high-density polymer, these cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage. The one-piece integral lid eliminates the need for separate steel lids. The snap-latch and hinge-lock design of the Unisetite™ prevents early separation of base and lid and allows for one-hand operation. Lids can be opened and closed as often as necessary and always relock securely without danger of specimen loss. Large labeling areas are located on two sides of the cassettes. Anterior writing area is at a 35° angle.

Cassettes1500/case

Specify:

#M505-10, Lilac	(115-1709)
#M505-11, Orange.....	(115-1710)
#M505-12, Aqua	(115-1711)
#M505-2, White.....	(115-1712)
#M505-3, Pink.....	(115-1713)
#M505-4, Green.....	(115-1714)
#M505-5, Yellow	(115-1717)
#M505-6, Blue.....	(115-1718)
#M505-7, Peach.....	(115-1719)
#M505-8, Tan.....	(115-1720)
#M505-9, Gray	(115-1721)



Simport

Histosette® II Biopsy Processing/Embedding Cassette

Most convenient for cassette-labeling instruments but specially designed to hold biopsy specimens during the embedding process, as well as in a storage file. The covers are already removed from cassettes and are packaged separately in the case.

Cassettes1000/case

Specify:

#M486-10, Lilac	(115-1624)
#M486-11, Orange.....	(115-1625)
#M486-12, Aqua	(115-1627)
#M486-2, White.....	(115-1628)
#M486-3, Pink.....	(115-1629)
#M486-4, Green.....	(115-1630)
#M486-5, Yellow	(115-1631)
#M486-6, Blue.....	(115-1632)
#M486-7, Peach.....	(115-1633)
#M486-8, Tan.....	(115-1634)
#M486-9, Gray	(115-1635)

Simport

Micromesh™ Biopsy Processing/Embedding Cassettes

This version of the Micromesh™ offers 1676 square openings (0.38mm), allowing for a greatly improved fluid exchange without having to use biopsy pads. Large anterior and posterior slots in both cassette and cover ensure that the Micromesh™ biopsy cassette will sink rapidly. A large square compartment measuring 27mm is perfect even for needle biopsies. The cover does not protrude above the cassette, a great space-saving feature allowing more cassettes to be stacked in automatic labeling machines and tissue processors.

Cassettes1000/case

Specify:

#M507-10, Lilac	(115-1733)
#M507-11, Orange.....	(115-1734)
#M507-12, Aqua	(115-1735)
#M507-2, White.....	(115-1736)
#M507-3, Pink.....	(115-1737)
#M507-4, Green.....	(115-1738)
#M507-5, Yellow	(115-1739)
#M507-6, Blue.....	(109-7544)
#M507-7, Peach.....	(115-1740)
#M507-8, Tan.....	(115-1741)
#M507-9, Gray	(115-1742)



Simport

Embedding Cassette

Disposable plastic tissue cassettes are suitable for holding and identifying tissue samples in processing, embedding, and sectioning procedures. The cassettes fit securely in microtome chuck adapters. They are molded from a high density polymer that is totally resistant to the chemical action of histological solvents. These cassettes are designed to accept standard metal lids and will keep specimens in complete safety during processing. The slanted writing surface accepts markings easily, permitting sample identification throughout all stages of embedding and long afterwards when in archives.

Cassettes1500/case

Specify:

#M480-10, Lilac	(487-7323)
#M480-11, Orange.....	(114-9995)
#M480-12, Aqua	(114-9996)
#M480-2, White.....	(281-4837)
#M480-3, Pink.....	(769-5292)
#M480-4, Green.....	(120-0469)
#M480-5, Yellow	(114-9997)
#M480-6, Blue.....	(630-2264)
#M480-7, Peach.....	(114-9998)
#M480-8, Tan.....	(114-9999)
#M480-9, Gray	(115-1611)





Simport

Biopsy Foam Pads

Used to hold biopsies in place and prevent them from being lost during processing. Specially formulated foam is always verified for consistency throughout to achieve optimum solvent flow. Biopsy samples are sandwiched between two foam pads and placed either in tissue capsules or cassettes with metal or plastic lids. Resists temperatures from -40°C to 121°C. For use with cassettes.

#M476-1
30.2 mm x 25.4 mm x 2 mm
(120-6138)1000/case




Simport

Disposable Base Molds

Disposable base molds offer ease and convenience of operation. They are inexpensive enough to be discarded after use, yet strong enough to be reused. Specially formulated plastic material provides excellent thermal exchange. Smooth interior finish and rounded corners facilitate specimen removal. Available in the same variety of sizes as metal molds and can be used with the same styles or types of cassettes and embedding rings.

Base Molds1000/case

Specify:
#M475-1
7 mm x 7 mm x 5 mm.....(114-9988)
#M475-2
15 mm x 15 mm x 5 mm.....(114-9990)
#M475-3
24 mm x 24 mm x 5 mm.....(114-9991)
#M475-4
30 mm x 24 mm x 5 mm.....(114-9992)
#M475-5
37 mm x 24 mm x 5 mm.....(114-9993)



Simport

Modular Storage Drawer

This drawer provides permanent storage and identification of up to 165 embedding rings or 250 cassettes per drawer. It is stackable to any convenient height, thanks to interlocking ridges on top and bottom. Made of high-impact resistant plastic. Identification labels included.

15 7/8"W x 9 1/8"D x 2 7/8"H
#495-6, Storage Drawer
(115-1677)6/case




Simport

SlideFolder™ Microscope Slide Holder

The SlideFolder™ will hold up to 20 standard microscope slides 75 mm x 25 mm (3" x 1") and is made of 2 parts: a base holding the slides horizontally, offering numbered spaces for easy identification, and transparent doors that can either cover the slides or be swung behind the SlideFolder™ for space-saving purposes. The 2 transparent hinged covers offer a full view of each slide without having to remove it from its position in the SlideFolder™, and allow easy reading of ID labels with or without an optical bar-code reader. All units are stackable and take minimal space on laboratory tables or shelves.

Slide Holders.....10/case

Specify:
#M750-20AS
Assorted, 2 of Each Color...(115-1875)
#M750-20B, Blue.....(115-1876)
#M750-20G, Green.....(115-1877)
#M750-20GY, Gray.....(115-1878)
#M750-20W, White.....(115-1879)
#M750-20Y, Yellow.....(115-1880)



Simport

StainTray™ Slide Staining System

The StainTray™ has a black base made of tough ABS plastic that can withstand a wide range of chemicals (avoid chlorinated hydrocarbons). It will accept up to 20 slides on four plastic rails covered with a polymer strip to hold slides perfectly even if tray is held at an angle. When humidity is needed, wells between rails will hold up to one mL of water securely without splashing. Middle wells will hold up to 2 mL each. Rails are raised not only to avoid water touching the slides but to make them more easily retrievable. The base will also hold excess stain solution dripping from the slides. Four rubber feet ensure greater base stability. Units are stackable for space-saving purposes.



Slide-Staining Systemea

Specify:
#M920-1, Base with Clear
Lid.....(115-1904)
#M920-2, Base with Black
Lid.....(115-1905)
#M921-1, Clear Lid Only.....(115-1906)
#M921-2, Black Lid Only.....(115-1907)

Simport

EasyDip™ Slide Staining System

The system has two components: a square staining jar and a 12-position vertical slide rack. Made of resistant acetal plastic, the staining jar will not break like most glass jars do. It will resist most staining agents, including alcohol and xylene (but not phenol, iodine or ferric chloride). The wide stable base offers greater stability, while the inside is recessed, allowing for a smaller reagent volume of only 80 mL. Easy to clean and no metals to corrode. Ideal for special stains, frozen sections and special processes.

Slide-Staining Systems.....6/case



Specify:
#M900-12AS
Assorted (1 Jar of Each
Color and 1 Rack).....(115-1889)
#M900-12B, Blue.....(115-1892)
#M900-12G, Green.....(115-1893)
#M900-12P, Pink.....(115-1894)
#M900-12W, White.....(115-1895)
#M900-12Y, Yellow.....(115-1896)



Gram Stain Kits and Reagents

Used to stain for microorganisms from cultures or specimens by the differential Gram method.

- #212539, Gram Stain Kit
(106-2056)4/pkg
Contains: 1-250-mL bottle ea of Gram crystal violet, Gram iodine (stabilized), Gram decolorizer, & Gram safranin.
- #212524, Gram Stain Kit
(645-2264)4/pkg
Contains: 1-250-mL bottle of ea: Gram crystal violet, Gram iodine diluent with iodine concentrate, Gram isopropyl decolorizer & Gram safranin.
- #212525, Gram Crystal Violet
(116-6788)4/pkg
- #212527, Gram Decolorizer
(224-5022)4/pkg
- #212531, Gram Safranin
(343-3274)4/pkg
- #212542, Gram Iodine (Stabilized)
(613-3242)4/pkg
- #212529, Gram Iodine (Unstabilized)
(264-9099)4/pkg
- #212536, Acridine Orange Stain
(610-4654)ea

Isopropyl Alcohol

- 99% Purified
1 gal
(120-8443)ea



Sedi-Stain™ Concentrated Reagent Stain

Stabilized modification of the Sternheimer-Malbin stain for use in the microscopic examination of urine sediment, which is generally recognized as a valuable diagnostic technique to permit rapid and accurate identification of intrinsic renal disease.

- #421570, 12.5-mL Bottle with Dropper
(168-9062)ea
Contains: Crystal violet 0.10%, safranin O 0.25%, ammonium oxalate 0.03%, ethyl alcohol (SD-3A) 10%, & water & stabilizers 89.62%.



Erado-Sol® Stain Remover

A fast and effective biological and chemical stain remover that removes most medical and histological stains from hands, plastic, glassware, lab coats, floors, and laboratory equipment and surfaces. Phosphate-free, biodegradable, and pH neutral (7-8). Easy to apply—simply rub in gently and stain is gone in a matter of seconds. Rinse with water.

- #207-32, 32 oz
(154-6143)ea



Alcohol—Anhydrous

- ACS Reagent
#1200-1GL, 1-gal Bottle
(120-8724)ea



Monse's (Rx) Solution

Formula to help control bleeding during biopsy, LLETZ, and laser procedures. Widely known as ferric-subsulfate solution, it is less irritating than ferric-sulfate solution, owing to the smaller proportion of sulfuric acid. Convenient 8-mL, light-resistant amber glass vials for single usage. Includes applicators.

- 8-mL Unidose Vial
(812-0005)12/box



Lugol's (Rx) Solution

This powerful anti-infective agent contains an aqueous solution of iodine and is used in cervical cancer screening. Convenient 8-mL, light-resistant amber glass vials for single usage. Includes applicators.

- Used for staining cervical and vaginal mucosa during colposcopic evaluations of abnormal Pap smears
- Less painful or irritating than alcoholic solutions when applied to cuts
- Potassium iodide added to increase the solubility of the iodine

- #0111-8, 8-mL Unidose Vials
(812-0004)12/box



ORDERING OPTIONS



Streamline your ordering process in a world where there is no "one size fits all".

Henry Schein offers a wide range of ordering tools and options that you can access at your convenience.



Web



Phone



Fax



Procurement
Software



Account
Representative

FEATURES:

- Quick Ordering Methods • Online Order Management • Online Order Tools • Easy Returns
- Training • Reporting • Dedicated Specialty Web Sites • Convenient Payment Options
- e-222 for Controlled Substances

BENEFITS:

- Convenient, flexible ways to order • Easy-to-use Web site tracks your most frequently purchased items • Mobile tools, bar code scanning, and iPad
- E-commerce training available by request
- Customer service support line for troubleshooting

E-mail: hsmmarketing@henryschein.com

Visit: www.henryschein.com/mdservices

HENRY SCHEIN®
MEDICAL

Rely on Us



Hematology Solutions and Stains

#400110, QuickLink I Wright's Stain, 16 oz (104-3950)	ea
#400120, QuickLink I Wright's Stain, 32 oz (104-3898)	ea
#400129, QuickLink I Wright-Giemsa Stain, 3 oz (103-0691)	ea
#400171, QuickLink III Fixative, 1 gal (996-4989)	ea
#400160, QuickLink III Set, 16 oz (3) (833-1331)	ea
#400191, Eosinophil Stain, 16 oz Compares to: Hansel's® Stain (751-9989)	ea
#400190, Eosinophil Stain, 4 oz Compares to: Hansel's® Stain (106-9583)	ea

#400220, Neat Stain® 150 Hematology (104-3900)	6/box
#400199, Phosphate Buffer, 7.2 pH, 100 mL (108-2296)	ea
#400202, Phosphate Buffer, 7.2 pH, 500 mL (108-2297)	ea
#400131, QuickLink I Wright-Giemsa Stain, 16 oz (108-8896)	ea
#400130, QuickLink I Wright-Giemsa Stain, 32 oz (253-7496)	ea
#400170, QuickLink III Fixative, 16 oz (833-3723)	ea
#400140, QuickLink III Solution I, 16 oz (833-6589)	ea

#400141, QuickLink III Solution I, 1 gal (914-2544)	ea
#400180, QuickLink III Solution II, 16 oz (833-4052)	ea
#400181, QuickLink III Solution II, 1 gal (104-9688)	ea
#400150, QuickLink III Stat-Pak, 3 oz (3) (437-8080)	ea
#400201, Sternheimer Malbin Stain, 12.5 mL Compares to Sedistain® (268-7372)	ea
#400200, Sternheimer Malbin Stain, 60 mL Compares to Sedistain® (108-2292)	ea



Microbiology Solutions and Stains

#400308, Chlorazol Black E, 2 oz (792-4876)	ea
#400327, Decolorizer Replacement, 8 oz (108-9431)	ea
#400331, Decolorizer, 16 oz (108-9429)	ea
#400328, Decolorizer Gram Stain, 1 gal (108-9430)	ea

#400310, Gram Stain Set, 2 oz (922-5498)	4/set
#400340, Gram Stain Set, Stabilized IOD, 8 oz (107-4395)	4/set
#400344, Iodine Stabilized, 8 oz (108-5887)	ea
#400351, Lugol's Solution, 2-oz Dropper (106-0788)	ea
#400350, Lugol's Iodine, 16 oz (108-2277)	ea

#400355, Lugol's Iodine Concentrate, 2-oz Dropper (108-2283)	ea
#400353, Lugol's Iodine Concentrate, 16 oz (108-2293)	ea
#400488, Methylene Blue Loeffler Formula, 4 oz (253-4317)	ea
#400334, Safranin Replacement, 8 oz (110-0522)	ea




Reagents Solutions

#400420, Acetic Acid 3%, 16 oz (740-8432)	ea
#400422, Acetic Acid 3%, 1 gal (108-0471)	ea
#400450, Acetic Acid 5%, 16 oz (666-5674)	ea
#400458, Acetone, 4 oz (108-0473)	ea
#400460, Acetone, 16 oz (309-2688)	ea
#400461, Acetone, 1 gal (108-0480)	ea
#400466, Alcohol 99% Isopropyl, 16 oz (108-0474)	ea
#400449, Alcohol 99% Isopropyl, 1 gal (437-9659)	ea
#400455, Aluminum Chloride, 6 Normal, 4 oz (437-5264)	ea
#400463, Aluminum Chloride 20%, 4 oz (300-6116)	ea
#400464, Aluminum Chloride 20%, 8 oz (199-2408)	ea
#400471, Aluminum Chloride 25%, 16 oz (355-4417)	ea
#400459, Aluminum Chloride 35%, 16 oz (914-8949)	ea
#400472, Aluminum Chloride 35%, 4 oz (811-9287)	ea
#400474, Aluminum Chloride 40%, 1 oz (990-7613)	ea

#400454, Aluminum Chloride 50%, 16 oz (163-8578)	ea
#400590, Dichloroacetic Acid, 15 mL Same as Bichloroacetic Acid® (991-3996)	ea
Ethanol Solution 95%, 1 gal (437-6392)	ea
#400453 Glycolic Acid 70%, 2 oz (484-5537)	ea
#400525, Jessner's Solution, 4 oz (647-4464)	ea
#400491, Monsels Solution Dropper, Non-Reg., 2 oz (108-2286)	ea
#400500, Monsel's Solution, 16 oz (106-5092)	ea
#400512, Phenol 6%, 4 oz (600-8610)	ea
#400507, Phenol, Std, 16 oz (108-2294)	ea
#400524, Phenol, 20%, 2 oz (108-0988)	ea
#400511, Potassium Hydroxide 10%, 2-oz Dropper (315-2772)	ea
#400530, Potassium Hydroxide 10%, 4 oz (958-8858)	ea
#400536, Potassium Hydroxide 20%, 2-oz Dropper (106-6572)	ea
#400526, Potassium Hydroxide 20% with DMSO, 30 mL (248-9676)	ea

#400553, Trichloroacetic Acid 15%, 16 oz (291-0494)	ea
#400555, Trichloroacetic Acid, 20%, 4 oz (483-9095)	ea
#400557, Trichloroacetic Acid 25%, 4 oz (747-5379)	ea
#400566, Trichloroacetic Acid 30%, 16 oz (472-4346)	ea
#400564, Trichloroacetic Acid 35%, 4 oz (938-4830)	ea
#400562, Trichloroacetic Acid 50%, 4 oz (437-8054)	ea
#400567, Trichloroacetic Acid 80%, 15 mL (877-9174)	ea
#400569, Trichloroacetic Acid 80%, 4 oz (600-4768)	ea
#400568, Trichloroacetic Acid 80%, 16 oz (240-3626)	ea
#400572, Trichloroacetic Acid 85%, 4 oz (907-5638)	ea
#400573, Trichloroacetic Acid 85%, 16 oz (437-8961)	ea
#400577, Trichloroacetic Acid 100%, 4 oz (437-7453)	ea
#400444, Monsel's High Viscosity, 1-oz Multidose Vial of Prethickened Monsel's Solution (104-9691)	ea
#1815, KOH 10% Dropper Bottle (437-6045)	ea

DON'T FORGET ABOUT HEALTHLINK TO HELP YOU WITH YOUR MICROBIOLOGY NEEDS!



HS ITEM #	PRODUCT DESCRIPTION
(437-3477)	Strep Select AGAR SSA Plates 10/Pk
(437-5777)	TSA 5% Sheep Blood 10/Pk, 10 PK/CA
(437-8270)	Bi Plate TSA/Macconkey 10/sl 10/Pk
(437-8576)	TSA 5% Sheep Blood/EMB 10/Pk
(437-3198)	Macconkey Agar 100mmx15mm 10/Pk

To order call Henry Schein: 800-772-4346

Don't forget these other great products to help fight the spread of germs.

HealthLink® offers a wide variety of diagnostics, medical devices and infection-control products to meet all your customer needs.





HealthLink®
A Family of Medical Products



CLOROX
HEALTHCARE™

Joining Forces to Serve the Healthcare Community

PH: (800) 638-2625 www.hlk.cc

2013 Clorox Professional Products Company, 1221 Broadway, Oakland, CA 94612.
Always use products as directed.



BBL™ Prepared Plated Media

#221353, Columbia CNA Agar with
5% Blood
(101-3786)100/box
(101-3782)100/box
(168-8219)100/case
(987-3908)100/case



Culture Blood Agar 5% Media

Culture Blood Agar 5% Media
#221239
(224-2775)20/box



MacConkey Agar II

MacConkey Agar Plated
#221172
(224-5653)20/box



Mueller-Hinton Agar Plates, 150 mm

#221800
(224-0871)24/box



Selective Strep Culture Plates

#221935
(120-4787)10/pkg



Sulfamethoxazole with Trimethoprim

#231536, 23.75/1.25 mcg
(100-6160)50/box



Trypticase–Soy Agar with Blood

#221291
(423-9799)100/box



Trypticase Soy Broth Medium

#221815, 2 mL
(101-3799)100/pkg



Plated Media

#1146, Strep Select Agar SSA
Plates, 100 mm
(437-3477)10/pkg

#1728, DTM, ¾ oz
(437-6907)20/box

#1103, Mueller Hinton Agar, 100 mm
(437-2229)10/pkg

#1106, Mueller Hinton Agar, 150 mm
(437-9509)8/pkg

#1087, MacConkey Agar, Monoplate,
100 mm
(437-3198)10/pkg

#1280, TSA 5%/MacConkey, Biplate
(437-8270)10/pkg

#1160, TSA 5% Sheep Blood, 100 mm
(437-5777)10/pkg

#1272, TSA 5% Sheep Blood/EMB, Biplate
(437-8576)10/pkg

#1442, TSA 5% Sheep Blood/EMB/XLD
Citrate UR Quad Plate
(437-3280)10/pkg

#1220, CNA 5%/MacConkey,
Biplate, 100 mm
(437-8450)10/pkg

#1043, EMB Agar, Monoplate
(437-8946)10/pkg

#1502, Bullseye Urine Plate
(991-7194)10/pkg

#1685, URI-ID CLED/EMB
(932-7436)10/pkg



Sensi Disc

Sensi Disc.....50/box

Specify:

Ampicillin AM-10, 10 mcg ..(615-5378)
Amoxicillin, 20/10 mcg.....(100-5150)
Bacitracin 10U.....(693-5386)
Cephalothin CF-30,
30 mcg(693-6368)
Cephazolin, CZ-30(693-7594)
Ciprofloxacin(224-4350)
Erythromycin, 15 mcg(693-2139)
Furadantin, 300 mcg.....(693-2821)
Gentamycin 120 mcg.....(104-6129)
Penicillin, 10 U(693-8689)
Sulfamethoxazole,
23.75/1.25 mcg.....(224-6637)
Tetracycline, 30 mcg.....(693-1167)



Taxo™ Differentiation Discs

#231552, A Discs
For differentiating group A streptococci.
(224-1056)50/box

#231042, A Discs
For differentiating group A streptococci.
(693-1886)50/box

#231562, 1-mg Kanamycin Discs
For differentiating gram-negative
anaerobic bacilli.
(987-1965)50/box

#231041, A Disc Cartridges
For differentiating group A streptococci.
(693-8964)6/box



ALFA WASSERMANN

Diagnostic Technologies, LLC

ACE Alera® Clinical Chemistry System

This high performance, low maintenance chemistry analyzer is capable of improving patient care and your bottom line. Designed to meet the needs of the POL, it has been cleared by the Food & Drug Administration specifically for use in Physician Office Labs.

Instrument features:

- Processes up to 165 photometric tests per hour and up to 120 potentiometric (ISE) tests per hour
- Patented, closed-tube sampling assures safety by effectively minimizing the risk of bloodborne pathogens, eliminating need to pour off samples
- Refrigerated compartment maintains specimen and reagent integrity as well as eliminating overnight reagent transfer



- Up to 40 test reagents can be stored on-board at any time
- Software identifies and tracks samples, patients, and reagent details, including expiration date and remaining reagent
- Stores up to 1,000 patient records and 5,000 requisitions
- Offers a unique stat interrupt process for single or panel testing

- Ensures fast, accurate reporting of abnormal test results by automatically diluting and repeating those tests without operator intervention while offering a unique stat test interrupt process for single tests

#402900-1
(168-9581)ea

ACE Alera and ACE Axcel Reagents/Calibrators/Controls

Σ = Not available for Analyzer

Item Code	Manuf Part No.	Item Description	Size	Ace Alera	ACE Axcel
120-6765	SA2001	ACE® Albumin Reagent (400 Tests/Kit)	Ea	•	•
120-3292	RX2002	ACE® Alkaline Phosphatase Reagent (200 Tests/Kit)	Ea	•	•
120-3571	SA2002	ACE® ALP Reagent (750 Tests/Kit)	Ea	•	•
642-0022	SA1046	ACE® ALT Reagent (450 Tests/Kit)	Ea	•	•
120-1477	RX1047	ACE® AST Reagent (200 Tests/Kit)	6/Bt	•	•
985-0010	SA1047	ACE® AST Reagent (450 Tests/Kit)	6/Bt	•	•
120-5468	SA1004	ACE® Amylase Reagent (300 Tests/Kit)	Ea	•	•
120-2959	SA2024	ACE® BUN Reagent (600 Tests/Kit)	6/Bx	•	•
120-9447	SA1009	ACE® Calcium Reagent (600 Tests/Kit)	6/BX	•	•
120-5266	RX1010	ACE® Cholesterol Reagent (267 Tests/Kit)	Ea	•	•
120-4820	SA1010	ACE® Cholesterol Reagent (800 Tests/Kit)	Ea	•	•
120-2320	SA2011	ACE® CK Reagent (750 Tests/Kit)	480/Bx	•	•
116-2176	SA2039	ACE® CO2 Reagent (500 Tests/Kit)	1/Bx	•	•
153-9698	RX2011	ACE® Creatine Kinase Reagent (300 Tests/Kit)	Ea	•	•
120-0363	RX1012	ACE® Creatinine Reagent (188 Tests/Kit)	Ea	•	•
120-3232	SA1012	ACE® Creatinine Reagent (560 Tests/Kit)	560/Test	•	•
255-3471	SA1007	ACE® Direct Bilirubin Reagent (120 Tests/Kit)	3/Bx	•	•
399-7821	RX1026	ACE® Direct LDL Reagent (30 Tests/Kit)	Ea	•	•
330-2574	SA1040	ACE® Direct LDL Reagent (100 Tests/Kit)	Ea	•	•
114-5616	ACI-32	ACE® Direct TIBC Reagent (120 Tetsts/Kit)	1/Kt	•	•
114-8840	S2-35	ACE® Ferritin Calibrator	Ea	•	•
114-8841	ACI-35	ACE® Ferritin Reagent (90 Tests/Kit)	1/Kt	•	•
120-3607	SA2013	ACE® GGT Reagent (480 Tests/Kit)	Kit	•	•
120-5177	SA1014	ACE® Glucose Reagent (600 Tests/Kit)	Ea	•	•
109-3655	RX2038	ACE® HDL Reagent (30 Tests/Kit)	1/Bx	•	•
209-3448	SA2038	ACE® HDL Reagent (100 Tests/Kit)	1/Bx	•	•
436-0000	SA1043	ACE® HDL-C Reagent (300 Tests/Kit)	Ea	•	•
271-0011	SA1044	ACE® Hemoglobin A1C Reagent (100 Tests/Kit)	1/Bx	•	•
363-0915	RX1016	ACE® Inorganic Phosphorous Reagent (300 Tests/Kit)	Ea	•	•
120-6989	SA1016	ACE® Inorganic Phosphorous Reagent (900 Tests/Kit)	6/Bx	•	•
271-1374	RX1017	ACE® Iron Reagent (80 Tests/Kit)	Ea	•	•
120-2617	SA1017	ACE® Iron Reagent (240 Tests/Kit)	2/Bx	•	•
116-2024	SA3018	ACE® LDH Reagent (330 Tests/Kit)	5 x 3mL/Bx	•	•
271-0012	SA1045	ACE® Lipase Reagent (150 Tests/Kit)	150/Bx	•	•
120-0927	SA2019	ACE® Magnesium Reagent (200 Tests/Kit)	1/Bx	•	•
120-4325	SA1008	ACE® Total Bilirubin Reagent (300 Tests/Kit)	3/Bx	•	•


ALFA WASSERMANN
Diagnostic Technologies, LLC
Features:

- Up to 165 photometric tests per hour; 285 per hour with ISE
- Internet connectivity allows for technical support, remote access, and laboratory integration
- Optional barcode scanner automates sample identification and streamlines workflow
- Run panels, profiles, or single tests in any combination at any time
- Create custom panels to meet the needs of your practice
- Stat interrupt allows users to load and prioritize stat samples during operation for immediate diagnosis and treatment
- Liquid, ready-to-use reagents eliminate prep time, improving operator productivity
- Automatic reagent inventory management

(120-1033)ea

ACE Axcel™ Clinical Chemistry System

This next-generation clinical chemistry system sets the standard for excellence in today's physician office laboratory (POL). With touch-screen functionality, an intuitive user interface, and Internet connectivity, the ACE Axcel offers the latest advances to make in-office laboratory testing fast, convenient, and cost-effective. Designed to meet the needs of the POL, it has been cleared by the US Food and Drug Administration specifically for use in physician office labs.

ACE Alera and ACE Axcel Reagents/Calibrators/Controls
Σ = Not available for Analyzer

Item Code	Manuf Part No.	Item Description	Size	ACE Alera	ACE Axcel
102-9352	SA1022	ACE® Total Protein Reagent (900 Tests/Kit)	6/Bx	•	•
120-4359	RX1023	ACE® Triglycerides Reagent (267 Tests/Kit)	Ea	•	•
120-5396	SA1023	ACE® Triglycerides Reagent (800 Tests/Kit)	Ea	•	•
120-0189	RX1025	ACE® Uric Acid Reagent (267 Tests/Kit)	Ea	•	•
120-9881	SA1025	ACE® Uric Acid Reagent (800 Tests/Kit)	Ea	•	•
115-4509	ACI-25	Apolipoprotein A1 Reagent	1/Bx	•	•
115-4513	S2-78	Apolipoprotein A1 Calibrator	1/Bx	Σ	•
115-4510	ACI-26	Apolipoprotein B Reagent (200 Tests/Kit)	1/Bx	•	•
108-7656	S2-76	Direct TIBC Calibrators	6/Kt	•	•
120-6806	S1-32	Gemcal Reference Serum	Pk	•	•
271-0001	S2-86	HbA1c Calibrator	1/Bx	•	•
337-2509	C2-86	HbA1c Control	1/Bx	•	•
110-5708	ACI-29	Microalbumin Reagent (200 Tests/Kit)	Ea	•	Σ
110-5709	S2-81	Microalbumin Calibrator	Ea	•	Σ
399-1398	S2-71	LDL Calibrator	Ea	•	•
120-2062	C1-1	Level 1 Chemistry Control	5/Kt	•	•
120-1319	C1-2	Level 2 Chemistry Control	5/Kt	•	•
115-4512	ACI-27	Lipoprotein(a) Reagent (40 Tests/Kit)	60/Bx	•	•
115-4515	S2-79	Lipoprotein(a) Calibrator	Ea	•	•
112-4522	C2-79	Lipoprotein Control	Ea	•	•
153-9764	SA3007	Total Thyroxine (T4) Reagent (100 Tests/Kit)	100/Bx	•	•
399-9673	S3-4	T4 Calibrator	Ea	•	•
113-3419	ACI-28	Transferrin Reagent (200 Tests/Kit)	200/Bx	•	•
113-3422	S2-80	Transferrin Calibrator Kit	6/Bx	Σ	•
112-4526	C2-80	Transferrin Control	Ea	•	Σ
271-0020	C2-85	Urine Control	Ea	•	•
118-7492	S2-85	Urine Standard	Ea	Σ	•

ACE Alera and ACE Axcel Sytem Accessories

118-7348	908962	Line Conditioner	Ea	•	•
115-4701	701303	Maintenance Log	Ea	•	Σ
108-1972	402995	ACE Alera® System Cap, Cal A Bottle Assembly	Ea	•	Σ
120-3986	402377	Probe Assembly, Fluid Arm (Green / Red Dot)	Ea	•	Σ
120-4236	E1-3	Reference Electrode Solution, 90 mL	3/Bx	•	•

ALFA WASSERMANN

Diagnostic Technologies, LLC

ACE Alera and ACE Axcel Sytem Accessories

Σ = Not available for Analyzer

Item Code	Manuf Part No.	Item Description	Size	Ace Alera	ACE Axcel
120-2470	903393	Electrode, Sodium	Ea	•	•
498-2545	906182	Pump Tubing Set	Ea	•	Σ
108-8904	907183	Waste Bottle	Ea	•	•
106-4468	403063	ACE Axcel™, ACE Alera® System Bar Code Reader	Ea	•	Σ
621-6176	27-11283	Cap, Cal B Bottle Assembly	Ea	•	•
108-0384	402994	Cap, Waste Assembly	Ea	•	Σ
743-0722	907123	Cleaning Wire	10/Pk		• •
628-6341	906078	Diskettes	10/Bx	•	Σ
426-5275	27-11347	Fan Filter	Ea	•	•
642-9771	402770	Tubing Kit, ISE Silicone	Ea	•	•
271-5808	903625	O-Ring, Electrode	Ea	•	•
271-1007	903626	O-Ring, Reference Electrode	Ea	•	•
612-0646	39-10042	Reagent Tray	Ea	•	•
271-7606	903396	Electrode Housing, Reference	Ea	•	•
290-0947	402771	ISE Reference & Waste Tubing Kit	Ea	•	•
271-3484	27-11282	Cap, Reference Bottle Assembly	Ea	•	•
488-0536	903624	Service Syringe	Ea	•	•
255-5523	903410	Wire Stylet	Ea	•	•
120-5327	906178	Electrode, Chloride	Ea	•	•
120-4730	903394	Electrode, Potassium	Ea	•	•
399-2949	906179	Electrode, Reference	Ea	•	•
108-0414	907318-1	Quick Disconnect	Ea	•	Σ
337-5867	908950	Deck Latch Key	Ea	•	•
109-2798	402996	Cap, Diluent Assembly (ACE Alera)	Ea	•	Σ
108-8179	907353	Diluent Bottle Coupling	Ea	•	Σ
825-3164	906204	ISE Electrode Holder	Ea	•	Σ
118-8095	908818	ISE Peristaltic Pump Tubing	Ea	•	•

Consumables ACE Alera and ACE Axcel

271-4651	28-10421	Segment, Sample Cup	Ea
962-0720	504675	Segment, Well (Seglet)	Ea
153-3389	ACD-4	ACE® Bio-Bags	50/Bx
120-6121	ACD-5	ACE® Evap Cap	500/Pk
120-0233	E1-2	ISE Cal B	3 x 90 mL
120-7110	OC1R1	Open Channel 1 / Reagent 1, 6 Coded 30 mL Bottles	6/Bx
120-9392	OC1R2	Open Channel 1 / Reagent 2, 6 Coded 12 mL Bottles	6/Bx
120-2957	OC2R1	Open Channel 2 / Reagent 1, 6 Coded 30 mL Bottles	6/Bx
120-6041	OC2R2	Open Channel 2 / Reagent 2, 6 Coded 12 mL Bottles	6/Bx
255-6374	OC3R1	Open Channel 3 / Reagent 1, 6 Coded 30 mL Bottles	6/Bx
255-9428	OC3R2	Open Channel 3 / Reagent 2, 6 Coded 12 mL Bottles	6/Bx
271-1862	OC4R1	Open Channel 4 / Reagent 1, 6 Coded 30 mL Bottles	6/Bx
271-7753	OC4R2	Open Channel 4 / Reagent 2, 6 Coded 12 mL Bottles	6/Bx
255-9307	OC5R2	Open Channel 5 / Reagent 2, 6 Coded 12 mL Bottles	6/Bx
984-0684	OC6R1	Open Channel 6 / Reagent 1, 6 Coded 30 mL Bottles	6/Bx
120-5832	OC6R2	Open Channel 6 / Reagent 2, 6 Coded 12 mL Bottles	6/Bx
997-7756	OC7R1	Open Channel 7 / Reagent 1, 6 Coded 30 mL Bottles	6/Bx
120-9833	OC7R2	Open Channel 7 / Reagent 2, 6 Coded 12 mL Bottles	6/Bx
271-4580	OC8R1	Open Channel 8 / Reagent 1, 6 Coded 30 mL Bottles	6/Bx
436-3956	OC8R2	Open Channel 8 / Reagent 2, 6 Coded 12 mL Bottles	6/Bx
882-3780	OC9R1	Open Channel 9 / Reagent 1, 6 Coded 30 mL Bottles	6/Bx
104-2205	OC9R2	Open Channel 9 / Reagent 2, 6 Coded 12 mL Bottles	6/Bx
106-4165	OC10R1	Open Channel 10/Reagent 1, 6 Coded 30 mL Bottles (amber)	6/Bx
655-0906	OC10R2	Open Channel 10 / Reagent 2, 6 Coded 12 mL Bottles (amber)	6/Bx



Consumables ACE Alera and ACE Axcel

Item Code	Manuf Part No.	Item Description	Size
106-4168	OC11R1	Open Channel 11 / Reagent 1, 6 Coded 30 mL Bottles (amber)	6/Bx
802-0467	OC11R2	Open Channel 11 / Reagent 2, 6 Coded 12 mL Bottles (amber)	6/Bx
120-0456	SA1020	NS Diluent	6 x 12 mL
120-5987	SA1021	System Diluent	3 x 750 mL
120-3995	SA1027	ISE Wash	6 x 30 mL
120-2576	SA1028	ISE Cleaning & Sodium Electrode Conditioner Kit	3 x 12 mL
120-1813	SA1029	ISE Calibrator A	3 x 300 mL
117-9822	SA1048	DiluQuik, 3 mL, Reconstitution of GEMCAL or Lipase Calibrator	20 x 3 mL
117-9827	SA1049	DiluQuik, 5 mL Reconstitution of Alfa Wassermann Chemistry Controls	20 x 5 mL
120-4447	ACD-2	Reaction Cuvettes	1000/Bx
120-4499	ACD-3	Seglet Wells	100/Bx
120-2618	ACD-6	Sample Cups	1000/BG

Lab Information Systems

337-0413	404169	Alfa Laboratory Information System (Alfa LIS)	Ea
337-0414	908914-1	Laboratory Information System, Management Console License	Ea
487-0000	908914-12	Laboratory Information System, Apex Integrator with Workstation	Ea
337-3083	908914-14	Laboratory Information System, Apex Connector without Workstation (2 waived devices & HL7 results)	Ea
337-0416	908914-2	Laboratory Information System, HL7 EHR Orders Interface	Ea
117-5860	908914-20	Laboratory Information System (Alera Console and Apex Workstation (Alfa LIS on a Dell Mini-tower)	Ea
118-5723	908914-22	Laboratory Information System, Onsite install and training 1st day (includes expenses)	Ea
337-3080	908914-23	Laboratory Information System, Onsite install and training 2nd day (includes expenses)	Ea
271-0018	908914-25	Laboratory Information System, Zebra Barcode Printer	Ea
337-0417	908914-3	Laboratory Information System, HL7 EHR Results Interface	Ea
337-0418	908914-4	Laboratory Information System, Uni-directional Analyzer Interface	Ea
337-0419	908914-5	Laboratory Information System, Bi-directional Analyzer Interface	Ea
337-0422	908914-6	Bar Code Software	Ea
337-0420	908914-7	Inventory Management Software	Ea
337-0412	908914-8	Allergy Management System	Ea
337-0421	908914-9	Bar Code Package Includes Software, Printer & Scanner	Ea
337-3087	908914-11	Apex Integrator without Workstation	Ea
337-3082	908914-13	Apex Connector with Workstation (2 waived devices & HL7 results)	Ea
337-3084	908914-15	Web Deployment Service	Ea
337-3085	908914-16	Apex Integrator 2nd WS	Ea
337-3086	908914-17	Bi-directional Host Query Analyzer Interface	Ea
271-0015	908914-18	HL7 EHR Demographics	Ea
271-0016	908914-19	HL7 Billing Export	Ea

Starlyte II Consumables

271-6047	E1-5	Electrode Cleaner, 90 mL	1/Ea
271-9235	E1-6	Sodium Electrode Conditioner, 100mL	1/Ea
120-6493	E1-8	ISE Reagent Fluid Pack	Ea

Starlyte System Accessories

498-5733	906186	Starlyte III Paper	5/Pk
673-9002	906203	ISC Pump Roller f/ISE Analyzer	Ea
389-8684	903408	Transport Housing, Reference Electrode	Ea

HORIBA

Medical

ABX Pentra 400 Analyzer Clinical Chemistry System

A compact, benchtop, clinical chemistry analyzer. Its autonomy and continuous loading provide enhanced productivity in a user-friendly environment. This very precise optical system meets the requirements of highly accurate photometric measurement. An extensive panel of liquid reagents simplifies workflow. Horiba presents a range of ready-to-use reagents in a single format.

- Up to 420 tests per hour with the ISE module
- 55 on-board assays with backup availability
- Save valuable time: control all of your operations from the main screen (system, reagent, and sample status), and easily upload target values for calibrators and controls

- Secure data management: Westgard multirule QC procedures with up to 3 default controls, full QC program stores up to 12 months of control history, Levey-Jennings, and statistics
- Conduct high-quality analysis: single-use reaction cuvette performs sensitive tests and all latex assays, high-performance mixer design mixes in minimal time without contaminating, and test-sequence setting program features an incompatibility function



- Control cost: optimized number of tests per cassette and limited calibrations

#P4000000NUA1
(117-1710).....ea

ABX Pentra 400 & Pentra C200 Test Menu

Σ = Not available for Analyzer

Item Code	Manuf Part No.	Item Description	Test Menu	Size	Pentra 400 Analyzer w/ISE	Pentra C200 Analyzer Chemistry
113-3181	A11A01901	ABX Pentra Reference	Reagent	Ea	•	•
109-2216	A11A01664	Albumin Cassette Reagent	Reagent	Ea	•	•
690-0015	A11A01626	ALP CP	Reagent	Ea	•	•
690-0016	A11A01627	ALT CP	Reagent	Ea	•	•
690-0017	A11A01628	Amylase CP	Reagent	Ea	•	•
690-0018	A11A01629	AST CP	Reagent	Ea	•	•
690-0023	A11A01633	Calcium CP	Reagent	Ea	•	•
118-8518	A11A01933	Cassette Creatinine CP	Reagent	Ea	•	•
690-0024	A11A01634	Cholesterol CP	Reagent	Ea	•	•
690-0094	A11A01786	CK Control	Control	4 x 3mL	•	•
690-0022	A11A01632	CK NAC CP	Reagent	Ea	•	•
690-0085	A11A01755	Clean-Chem CP	Reagent	Ea	•	•
110-4929	A11A01648	CO2 Calibrator Cassette	Calibrator	Ea	•	•
110-4170	A11A01890	CO2 Control Pentra	Control	8/Bx	•	•
110-4928	A11A01645	CO2 RTU Cassette	Reagent	Ea	•	•
690-0005	A11A01616	CRP Calibrator	Calibrator	5 x 1 mL	•	•
690-0000	A11A01611	CRP CP	Reagent	Ea	•	•
690-0084	A11A01754	Deproteinizer CP	Reagent	Ea	•	•
690-0025	A11A01635	Direct Bilirubin CP	Reagent	Ea	•	•
690-0091	A11A01769	Etching CP	Reagent	Ea	•	•
690-0019	A11A01630	GGT CP	Reagent	Ea	•	•
690-0044	A11A01667	Glucose HK CP	Reagent	Ea	•	•
690-0071	A11A01702	HbA1C WB 2x23/1x110/2x21	Reagent	Box	•	•
690-0072	A11A01703	HbA1c WB Cal 1x8mL/5x2mL	Calibrator	Box	•	•
690-0035	A11A01647	HDL Calibrator	Calibrator	2 x 1 mL	•	•
690-0026	A11A01636	HDL Direct CP	Reagent	Ea	•	•
690-0010	A11A01621	Immuno I Control L/H	Control	Ea	•	•
690-0027	A11A01637	Iron CP	Reagent	Ea	•	•
110-4171	A11A01871	LDH CP Reagent Pentra	Reagent	1/Bx	•	•
690-0049	A11A01678	LDL Calibrator	Calibrator	2 x 1 mL	•	•
690-0028	A11A01638	LDL Direct CP	Reagent	Ea	•	•
690-0021	A11A01631	Lipase CP	Reagent	Ea	•	Σ
690-0034	A11A01646	Magnesium RTU	Reagent	2 x 25mL	•	•



ABX Pentra C200 Clinical Chemistry System

Offers a unique level of automation in a compact system.

- Fixed cycles provide continuous throughput of 90 colorimetric tests/hour
- Optional high-throughput ISE module (270 tests/hour)
- Reagent management with unique Pentra cassette
- Fully automatic cuvette management with 2 racks of 96 single-use cuvettes
- Dedicated waste compartment for used cuvettes and contaminated liquids
- Peltier integrated computer and cooling unit
- Spectrophotometry: colorimetry and turbidimetry
- Potentiometry: direct (serum or plasma) and indirect (urine)
- 15 parameters per reagent tray with optional ISE module (Na⁺, K⁺, and Cl⁻)
- Reagent handling: positive reagent identification features an integrated bar-code reader; reagent needle detects level and shock; and reagent management with cassettes features an automated reagent on-board stability check plus continuous-test countdown check by capacitive effect



HORIBA
Medical

#PC200ISEOONUA1
(120-9231)ea

#PC200ISECONUA1
(690-9819)ea

ABX Pentra 400 & Pentra C200 Test Menu

Σ = Not available for Analyzer

Item Code	Manuf Part No.	Item Description	Test Menu	Size	Pentra 400 Analyzer w/ISE	Pentra C200 Analyzer Chemistry
690-0013	A11A01624	Micro-Alb Cal	Calibrator	5 x 1 mL	•	•
690-0014	A11A01625	Micro-Alb Control L/H	Control	2 x 1mL	•	•
690-0012	A11A01623	Micro-Albumin CP f/Pentra 400	Reagent	Ea	•	•
555-0036	A11A01652	Multical 3ml	Calibrator	10/Bx	•	•
690-0037	A11A01653	N Control	Control	10 x 5 mL	•	•
690-0038	A11A01654	P Control	Control	10 x 5 mL	•	•
117-8540	A11A01934	PC200 HDL Direct CP	Reagent	Ea	Σ	•
116-3956	A11A01898	Pentra TPU Calibrator	Calibrator	Ea	•	•
690-0042	A11A01665	Phosphorus CP	Reagent	Ea	•	•
690-0082	A11A01740	Potassium-E	Reagent	Ea	•	•
690-0086	A11A01758	Qualitest Solution	Reagent	Ea	•	•
117-6041	D01A00143	Reagent Grade Water 20 Liter	Reagent	20 Liter	•	•
690-0076	A11A01719	Reference	Reagent	Ea	•	•
690-0083	A11A01741	Reference-E	Reagent	Ea	•	•
690-0080	A11A01738	Sodium-E	Reagent	Ea	•	•
690-0074	A11A01717	Standard 1	Reagent	Ea	•	•
690-0075	A11A01718	Standard 2	Reagent	Ea	•	•
690-0029	A11A01639	Total Bilirubin CP	Reagent	Ea	•	•
110-4172	A11A01932	Total Protein CP Pentra	Reagent	1/Bx	•	•
690-0030	A11A01640	Triglycerides CP	Reagent	Ea	•	•
690-0031	A11A01641	Urea CP	Reagent	Ea	•	•
690-0046	A11A01670	Uric Acid CP	Reagent	Ea	•	•
690-0032	A11A01642	Urinary Proteins CP	Reagent	Ea	•	•
113-5338	A11A01674	Urine Control L/H	Control	Ea	•	•

ABX Pentra 400 & Pentra C200 Parts

108-5456	XEA935AS	Chloride Electrode Cleaning Kit	Part	Ea	•	•
690-0081	A11A01739	Chloride-Electrode	Part	Ea	•	•
119-1979	DZZ046F	Cooling Unit f/Pentra 400	Part	Ea	•	•
117-9574	A11A01891	Cuvette Segments f/Pentra 400	Part	450/Ca	•	•
118-5118	A11A01914	Cuvettes f/Pentra C200 60 Racks	Part	1/Ca	Σ	•
108-9268	A11A01851	Dummy Electrode f/Pentra 400	Part	Ea	•	Σ

ABX Pentra 400 & Pentra C200 Parts

Item Code	Manuf Part No.	Item Description	Test Menu	Size	Σ = Not available for Analyzer	
					Pentra 400 Analyzer w/ISE	Pentra C200 Analyzer Chemistry
114-6772	LAY008A	Ethylene Glycol	Part	Ea	•	•
117-5934	DZZ055A	Filter Cooling Unit	Part	Ea	•	•
119-5557	F10B40253	Filter Inline Water Tank	Part	Ea	•	•
690-0164	GBM0738	Lamp	Part	Ea	•	•
690-0146	B8089078	Liq Elem & Gasket Ring	Part	Ea	•	•
110-6458	XBA564DT	Probe f/Reagent Pentra 400	Part	Ea	•	Σ
120-2189	XBA563CT	Probe Sample Pentra 400	Part	Ea	•	Σ
690-0137	B8086583	O'Rings Replacement f/Reag Syr	Part	Ea	•	Σ
690-0115	B1034626	Reag Cup 10 mL	Part	100/Bag	•	•
690-0118	B1037307	Reag Cup 15 mL	Part	100/Bag	•	•
690-0116	B1034634	Reag Cup 4 mL	Part	100/Bag	•	•
690-0135	B8078955	Reagent Syringe (1000uL)	Part	Ea	•	•
690-0087	A11A01765	Sample Cup, Blue	Part	1000/Bg	•	•
690-0089	A11A01767	Sample Cup, White	Part	1000/Bg	•	•
111-9659	XEA778CS	Sample Rack Kit w/Adaptor	Part	Ea	•	•
690-0122	B8076812	Sample Syringe (100uL)	Part	Ea	•	•
690-0143	B8088454	Seal Replacement Tool	Part	Ea	•	•
690-0138	B8086605	Teflon Seal 1000uL	Part	Ea	•	•
690-0139	B8086648	Teflon Seal 100uL	Part	Ea	•	•
118-4803	A11A01915	Waste Bags Cuvette Pentra C200	Part	100/Bx	Σ	•

Is an antibiotic necessary?



Make a rapid determination...

VIRAL or **BACTERIAL**

MicrosCRP200

*A single finger stick of whole blood gets you
CBC and CRP results in 4.5 minutes!*

Accurate diagnosis avoids unnecessary antibiotic prescriptions.

Reimbursement codes are applicable for each CBC and CRP test. The Micros CRP200 analyzer provides STAT testing, producing quantitative results that lead to a more accurate diagnosis.

IMMEDIATE / ACCURATE / EFFECTIVE



MEDTEST

Spotchem EZ Analyzer

Automated Dry Chemistry Analyzer

- Easy, walk-away operation
- Calibration in seconds
- Whole blood or serum analysis within 3 minutes
- 13 CLIA-waived tests (lipids, diabetes, and general chemistries)
- Single strips, not just panels
- Individual sealed tests
- Reagent strips
- Whole blood sample
- CRMLN certification
- Onboard centrifuge
- Minimal waste
- Minimal maintenance
- Easy installation
- Small, compact footprint
- Mobile
- Internal computer
- Integral printer
- LIS capability

#SP-4430-FT, Analyzer
(117-2463)ea

Consumables

#10712, Centrifuge Tubes
(499-8158)100/box

#10743, Tip Set
(499-8159)500/box



#TP-5, Thermal Paper for Spotchem
(117-2462)5 rolls/box

Diabetes Product

#77424, Glucose Single Strip Test*
(499-8194)25 strips/box

General Chemistries

#77431, Lipid Panel*
(499-8157)25 strips/kit

Contains: HDL, Total Chl, Trigs & LDL-C.

#77433, Spotchem Uric Acid
(117-2461)25 strips/kit

#77434, Spotchem EZ BUN*
(117-2457)25 strips/kit

#77445, Spotchem Creatinine*
(117-2460)25 strips/kit

Panels

#77447, Spotchem Magnesium
(117-2494)10/case

#77484, Spotchem Basic Panel I*
(111-6779)ea

Contains: ALB, BUN, GLU, CRE, CAL;
ISE separate.

#77485, Spotchem Basic Panel II*
(117-2518)25/box

Contains: ALP, T-BIL, ALT, AST & TP.

Controls

#FRUQ400, Fructosamine Control Low*
(175-9710)3 x 1 mL/box

#FRUQ401, Fructosamine Control High*
(175-9711)3 x 1 mL/box

#HC-DIL, Diluent and Dropper Tops
(111-7321)20 x 5 mL/box

#HNC200, Serum Control Level I
(108-7422)10 x 5 mL

#HEC200, Serum Control Level II
(108-7424)10 x 5 mL

*CLIA-Waived Test.



Ortho Clinical Diagnostics

a **Johnson & Johnson** company

Vitros® 350 Chemistry System

Offers easy operation, maintenance, and training, whether you use it as a primary, STAT, or backup to the powerful Vitros® 5,1 FS. The Vitros® 350 comes complete with enhanced throughput, a broad accessible menu, and new ergonomic design.

Features:

- Labor-saving ease of use: easy to train and cross-train operators
- Walk-away operation: quality-control procedures are required just once a day and calibration intervals up to 6 months
- Operator interface: ergonomically designed, flat, low-glare LCD touch-screen monitor
- Onboard documentation and help function
- Versatile and efficient: up to 300 results per hour
- Onboard capacity for all 45 assays
- Primary tube sampling
- Optional positive identification via barcode scanner
- Secure results: self-monitoring alerts operators to prompt action, avoiding testing delays

#6802153
(327-3877)ea

Ortho Clinical Diagnostics

a **Johnson & Johnson** company

Vitros® ECIQ Immunodiagnostic System

Features:

- Advanced processing center: Ensures result integrity for quality result reporting
- Patient samples and reagents automatically metered into MicroWells
- Intellicheck® Technology verification at each step
- Dual-logic processing optimizes assay performance and throughput
- Reagent and Management Center: Optimizes labor by reducing operator intervention
- Up to 20 integrated reagent packs loaded any time and stored onboard for continuous, random access operation
- Automatic inventory management monitors and verifies inventory, expiration, and calibration
- Sampling center: Reduces errors for worry-free result reporting
- Unique Save-the-Sample Clot/Bubble Detection preserves the specimen and minimizes redraw
- Universal sample trays accommodate primary, secondary, and microcollection containers and sample cups without adapters or pour-off containers
- e-Connectivity® Center: Maximizes system uptime and increases staff productivity
- Highly secure, real-time 2-way communication between your system and our support teams
- Automatic data exchange and remote connectivity for immediate diagnostic and performance analysis
- Immediate resolution of technical issues without a service visit

#1922814
(822-7677)ea

Ortho Clinical Diagnostics

a **Johnson & Johnson** company

Vitros® 5600

Overview:

- Increases capacity with 150 reagent positions that allow over 100 assays to be on-board at once
- Provides greater efficiency, lower maintenance, and minimal calibration by integrating our proven MicroSlide, MicroTip, MicroWell, Intellicheck, and MicroSensor technologies
- Optimizes turnaround time and productivity through innovative Sample-Centered Processing
- Verifies system processing and result reporting while reducing errors
- **Microsensor Technology:**
- Reagent-less determination of sample quality indices without additional sample volume or impact on throughput
- Reduces the time, labor, and error potential of manual interpretation
- **Sample Centered Processing:**
- Chemistry and immunoassay sample are processed in parallel
- Eliminates the need to move trays or aliquot samples
- Improves turnaround time and increases efficiency
- Single-use tips for sample and reagent metering
- Eliminates sample and reagent carryover
- Automatically detects clots and bubbles to avoid reporting erroneous results
- **e-Connectivity®/Proactive Monitoring**
- Highly secure, real-time, 2-way communication between your system and our support teams for timely resolution of technical issues

#6802413
(801-0000)ea

SIEMENS

Advia® 1800 Clinical Chemistry System

Pushes the limit of productivity. Its extensive menu allows integration of testing for clinical chemistry, drugs-of-abuse, therapeutic drug monitoring (TDM), and specific proteins like Cystatin C and CardioPhase® hsCRP. This is an ideal solution for mid- to high-volume chemistry analysis. Featuring intuitive new software and capabilities, it is designed to deliver power and speed to meet the most demanding performance needs.

System offers:

- Throughput of up to 1,800 tests per hour provides speed to keep pace with workload peaks
- 200 basic metabolic panels per hour meet turnaround-time demands
- Large onboard reagent capacity and optional concentrated reagents reduce interruptions
- Optimal standardization of reagents and systems for efficient consumable usage
- Automation ready design: no additional hardware required
- **Automation-ready design:**
- Point-in-space aspiration offers connectivity to a track, VersaCell, or rack handler without robotics or special hardware
- Direct sampling from track prevents detaining tubes at analyzer and reduces aliquoting
- **Peak productivity:**
- Automatic Sample Retain technology provides the ability to perform reruns and dilutions without the need to retrieve primary tubes



- Consolidation power with a comprehensive menu, consistent cycle times, speed, and open-application capability
- **Total system reliability:**
- Robust hardware and simplified software design, with automatic calibration validation, and tracking
- Remote diagnostics and ILQC (Inter-Laboratory Quality Control) options

Easy efficiency:

- Reagent pause capability and onboard Sample Pause button make loading easy
- Sample integrity checking, including serum indices reporting, short sample flagging, and clot detection reduce manual intervention

#10315183

(116-1530).....ea

ADVIA 1800 Test Menu

Item Code	Manuf Part No.	Item Description	Size
117-0433	10311893	Acid Phosphatase ADVIA-200 Tests	1/Cr
118-6535	10491408	ADVIA A1c_3 Calibrator	Ea
118-6538	10485591	ADVIA Hemoglobin A1c_3 Large Kit-800 Tests	Ea
118-6537	10379673	ADVIA Hemoglobin A1c_3 Small Kit-160 Tests	Ea
117-0395	10311832	Albumin ADVIA-2220 Tests	1/Cr
119-6339	7501976	Alanine Aminotransferase-2520 Tests	Ea
116-2794	10378758	Amphetamines_2-760 Tests	1/Cr
117-0439	10311898	APO-A1/B Calibrator	1/Cr
119-6343	7499718	Aspartate Aminotransferase-2520 Tests	Ea
116-2805	10378759	Barbiturate_2-760 Tests	1/Cr
116-2808	10378760	Benzodiazepine_2-760 Tests	1/Cr
119-6358	2189494	Bilirubin, Direct-520 Tests	Ea
119-6345	2189346	Bilirubin, Total-1424 Tests	Ea
119-6334	2189915	Calcium_2 2320 Tests	Ea
116-2809	10378766	Cannabinoid_2-760 Tests	1/Cr
119-6332	7506781	Carbon Dioxide Liquid-1200 Tests	Ea
117-0412	10311854	Chem Cuvette Conditioner -1000 Tests	1/Cr
116-8188	10376501	Cholesterol_2-2450 Tests	Ea
117-0635	10445277	Chromium (VI) Validity Calibrator 50 (14 mL)	1/Cr
116-2810	10378761	Cocaine Metabolite_2-760 Tests	1/Cr
117-0399	10341111	Creatine Kinase Large Kit-1180 Tests	1/Cr
119-6364	7498541	Creatine Kinase-980 Tests	Ea



SIEMENS

ADVIA 1800 Test Menu

Item Code	Manuf Part No.	Item Description	Size
117-0627	3T149UL	Creatinine Validity Calibrator 100 (14 mL)	1/Cr
117-0633	3T249UL	Creatinine Validity Calibrator 2 (14 mL)	1/Cr
117-0626	3T139UL	Creatinine Validity Calibrator 20	1/Cr
117-0631	3T159UL	Creatinine Validity Calibrator 400 (14 mL)	1/Cr
116-2811	10490244	Creatinine Validity Test-720 Tests	1/Cr
117-0510	073-0023-02	Cuvette RRV 1 x 17	1 x 17
117-0411	B01-4178-51	Cuvette Wash Solution-2000 Tests	1/Cr
117-0509	073-0022-01	Cuvette, DTT, 1 x 20	1 x 20
119-7712	073-0571-01	Cuvette, RRV, 1 x 17	17/Bx
117-0530	094-0028-01	Cuvette, RRV, 1 x 21-Tests	1 x 21
114-3909	9A509UL	Emit Calibrator/Control Level 0-14 Tests	1/Bx
114-3910	9A529UL	Emit Calibrator/Control Level 1-14 Tests	1/Bx
114-3911	9A549UL	Emit Calibrator/Control Level 2-14 Tests	1/Bx
114-3912	9A569UL	Emit Calibrator/Control Level 3-14 Tests	1/Bx
114-3913	9A589UL	Emit Calibrator/Control Level 4-14 Tests	1/Bx
114-3914	9A609UL	Emit Calibrator/Control Level 5-14 Tests	1/Bx
117-0421	B01-4769-01	Digoxin ADVIA-800 Tests	1/Cr
119-6351	7511947	Direct HDL Cholesterol-1304 Tests	Ea
117-0420	B01-4760-01	Direct LDL Cholesterol-752 Tests	1/Cr
119-7723	10376771	Drug Calibrator II	Ea
116-2793	10378735	Ecstasy-760 Tests	1/Cr
114-3915	9X529UL	ADVIA Chemistry Emit II Plus Ecstasy Cal/Cont Level 1 Kit	1/Bx
114-3916	9X549UL	ADVIA Chemistry Emit II Plus Ecstasy Cal/Cont Level 2 Kit	1/Bx
114-3917	9X569UL	ADVIA Chemistry Emit II Plus Ecstasy Cal/Cont Level 3 Kit	1/Bx
114-3918	9X589UL	ADVIA Chemistry Emit II Plus Ecstasy Cal/Cont Level 4 Kit	1/Bx
116-2792	10379635	Ethanol-600 Tests	1/Cr
119-6360	4862277	Ferritin Small Kit-800 Tests	Ea
117-0401	B01-4138-01	Gamma-Glutamyl Transferase Flex Reagent Cartridge-2700 Tests	1/Cr
119-6365	7498649	Gamma-Glutamyl Transferase-980 Tests	Ea
116-8187	10492319	Glucose Oxidase, Revised-2450 Tests	Ea
117-0431	B01-4797-01	HbA1c Manual Pretreatment-200 Tests	1/Cr
117-0442	B03-4763-01	HDL/LDL Cholesterol Calibrator	1/Cr
117-0413	B01-4180-01	Incubation Bath Oil	Ea
119-6337	7506099	Inorganic Phosphorus-2506 Tests	Ea
119-6355	2194838	Iron_2-1015 Tests	Ea
117-4606	10377510	Iron 1800/2400-2100 Tests	Ea
119-6372	B01-4171-51	ISE Buffer-2000 Tests	Ea
117-0408	B01-4174-52	ISE Detergent Solution-100 Tests	1/Cr
117-0409	B01-4177-51	ISE Internal Filling Solution	1/Cr
117-0514	073-0052-01	ISE Reference Electrode-30000 Tests	1/Cr
117-0405	B01-4173-51	ISE Reference Solution-Tests	1/Cr
117-0440	B03-4175-01	ISE Standard Set Serum-100 Tests	1/Cr
117-0441	B03-4176-01	ISE Standard Set Urine-100 Tests	1/Cr
117-0511	073-0049-01	ISE, Electrode - CL-30000 Tests	1/Cr
117-0512	073-0050-01	ISE, Electrode - K-30000 Tests	1/Cr
117-0513	073-0051-01	ISE, Electrode - NA-30000 Tests	1/Cr
117-0520	073-0653-01	ISE, Electrode-Reference-30000 Tests	1/Cr
119-6387	073-0078-02	ISE, Perista-Tube No.	Ea
117-0436	B01-4818-01	Lactate ADVIA-300 Tests	1/Cr
117-0417	B01-4496-51	Lamp Coolant Additive	1/Cr
117-0516	073-0099-01	Lamp, Halogen	1/Cr
117-0438	B01-4840-01	Lipase ADVIA-640 Tests	1/Cr
117-0446	B03-4845-01	Liquid Specific Protein Cal	1/Cr
117-0403	B01-4148-01	Magnesium Large Kit-1120 Tests	1/Cr
116-2812	10378762	Methadone_2-760 Tests	1/Cr

SIEMENS

ADVIA 1800 Test Menu

Item Code	Manuf Part No.	Item Description	Size
117-0634	3T309UL	Nitrite Validity Calibrator 500 (14 mL)	1/Cr
116-2795	10490248	Nitrite Validity Test-720 Tests	1/Cr
116-2796	10378763	Opiate_2-760 Tests	1/Cr
116-2797	10490247	Oxidant Validity Test-720 Tests	1/Cr
117-0642	3T459UL	pH Validity Cal 11.0	1/Cr
117-0636	3T389UL	pH Validity Cal 2.0	1/Cr
117-0641	3T399UL	pH Validity Cal 3.0	1/Cr
117-0653	3T499UL	pH Validity Cal 4.5	1/Cr
116-2804	10490245	pH Validity Test-720 Tests	1/Cr
116-2800	10378764	Phencyclidine_2-760 Tests	1/Cr
116-2802	10378765	Propoxyphene_2-760 Tests	1/Cr
117-0414	B01-4181-55	Reagent Wash Probe 1 5 x 250 mL -1250 Tests	5 x 250 mL
117-0415	B01-4182-55	Reagent Wash Probe 2 5 x 250 mL-1250 Tests	5 x 250 mL
117-0416	B01-4183-51	Reagent Wash Probe 3 1 x 500 mL-500 Tests	1 x 500 mL
117-0435	B01-4806-01	Rheumatoid Factor-200 Tests	1/Cr
116-3809	073-0157-02	Sample Cups-1000 Tests	1000/Ca
119-6378	T23-1291-52	SetPoint CO2 CAL/DIL	Ea
119-6382	T03-1291-62	SetPoint Calibrator	Ea
117-0445	B03-4821-01	Special Chemistry Calibrator	1/Cr
117-0655	3T619UL	Specific Gravity Validity Cal 1.003	1/Cr
117-0661	3T629UL	Specific Gravity Validity Cal 1.020	1/Cr
116-2803	10490246	Specific Gravity Validity Test-600 Tests	1/Cr
117-0518	073-0273-01	Syringe Seal Set	1/Kt
117-0642	3T459UL	pH Validity Calibrator 11.0 (14 mL)	1/Cr
117-0636	3T389UL	pH Validity Calibrator 2.0 (14 mL)	1/Cr
117-0641	3T399UL	pH Validity Calibrator 3.0 (14 mL)	1/Cr
117-0653	3T499UL	pH Validity Calibrator 4.50 (14 mL)	1/Cr
117-0655	3T619UL	Specific Gravity 1.003 Validity Calibrator (14 mL)	1/Cr
117-0661	3T629UL	Specific Gravity 1.020 Validity Calibrator (14 mL)	1/Cr
117-0443	B03-4790-01	TDM Calibrators	1/Cr
117-0419	B01-4601-01	Total Protein_2-3400 Tests	1/Cr
117-0398	B01-4133-01	Triglycerides ADVIA-2300 Tests	1/Cr
119-7724	B03-4847-55	Urine Total Protein Calibrator	Ea
117-0620	3T129UL	Validity Negative Calibrator/Control	1/Cr
119-7721	10377509	Vancomycin_2-400 Tests	Ea
117-0444	B03-4815-01	WR-CRP Calibrator (6 Level)	1/Cr
117-0432	B01-4800-01	WR-CRP Small Kit-440 Tests	1/Cr

ADVIA 1800 Accessories

Item Code	Manuf Part No.	Item Description	Size
117-0484	040-K087-01	ADVIA Automation - 1m Spacer Kit, ADVIA WorkCell	1/Kt
117-0508	073-K047-04	ADVIA Automation - ADVIA Chemistry Riser Skirt Kit	1/Kt
117-0528	092-A001-03	ADVIA Automation - Centrifuge Module w decapper covers 60Hz	Ea
117-0481	040-A059-07	ADVIA Automation - Decapper Module @ Centrifuge Module	Ea
117-0491	040-K179-04	ADVIA Automation - Decapper/Gate Kit @ Centrifuge Module	1/Kt
117-0529	092-K016-01	ADVIA Automation - Hettich Centrifuge	1/Kt
117-0485	040-K097-02	ADVIA Automation - IMMULITE SMS Interface Kit	1/Kt
117-0496	040-0003-02	ADVIA Automation Box of 25 Sample Packs	25/Bx
117-0495	040-0002-03	ADVIA Automation Box of 5 Sample Trays	5/Bx
117-0497	040-0598-02	ADVIA Automation Box of 5 Tray Covers	5/Bx
117-0479	040-A019-05	ADVIA Automation Centaur Computer Interface Kit	1/Kt
117-0521	078-K231-02	ADVIA Automation Centaur XP Riser/Skirt Kit	1/Kt
117-0478	040-A017-05	ADVIA Automation Chemistry Interface Kit	1/Kt



ADVIA 1800 Accessories

SIEMENS

Item Code	Manuf Part No.	Item Description	Size
117-0482	040-B208-01	ADVIA Automation Decapper Cabinet	Ea
117-0490	040-K179-03	ADVIA Automation Decapper Gate Kit	1/Kt
117-0487	040-K159-01	ADVIA Automation Decapper Side Car	Ea
117-0486	040-K124-01	ADVIA Automation Dimension RxL ADVIA Interface Kit	1/Kt
117-0781	400909-01	ADVIA Automation Direct Water Feed Kit	1/Kt
117-0498	040-0722-01	ADVIA Automation LabCell Compressor	Ea
117-0883	615-4294-01	ADVIA Automation Pilot x 32 CPU	Ea
117-0884	615-4294-02	ADVIA Automation Pilot x 64 CPU	Ea
117-0477	040-A012-15	ADVIA Automation Sample Manager Module	Ea
117-0859	422458-0001	ADVIA Automation Track Connection - Left Side (VC S/N: > G0801)	Ea
117-0789	422240-03	ADVIA Automation Track Connection - Left Side (VC S/N: B0001 - G0800)	Ea
117-0861	422458-0002	ADVIA Automation Track Connection - Rear (VC S/N: > G0801)	Ea
117-0792	422240-04	ADVIA Automation Track Connection - Rear (VC S/N: B0001 - G0800)	Ea
117-0493	040-K198-06	ADVIA Automation Vista 1500 Connection Kit	1/Kt
117-0543	10469522	ADVIA Automation Vista Riser and Trim Kit	1/Kt
117-0483	040-K052-01	ADVIA Automation WorkCell Expansion 2 Stations	Ea
116-3125	664-1962-01	CentraLink - Additional 8 Port Terminal Server	Ea
116-3100	084-0021-01	CentraLink - Annual Maintenance - Instrument Connection	Ea
116-3123	664-0679-01	CentraLink - Barcode Scanner: USB Port CCD Scanner LGP6125RU1S-000	Ea
116-3095	084-0004-02	CentraLink - Client User Interface License	Ea
116-3124	664-1782-01	CentraLink - Printer: Lexmark E460DN Laser Jet, 120VAC	Ea
116-3121	615-4188-01	CentraLink - RJ45 - RJ45 Cable (10) Qty 2	2/Pk
116-3122	615-4188-02	CentraLink - RJ45 - RJ45 Cable (25')	Ea
116-3118	615-4187-01	CentraLink - RJ45 to 9 pin DB adapter (LIS)	Ea
116-3115	515-7188-03	CentraLink - UPS: PowerVar #ABCE1100-11, 120VAC, 1100VA	Ea
117-0507	073-B005-04	Rack Handler	Ea
117-5920	422761-02	VersaCell ADVIA Centaur/ ADVIA 1800/ ExL/RxL 800 down wheel- back	Ea
117-3511	422711-01	VersaCell Instrument Connection Kit A1800	Ea
117-3556	422711-02	VersaCell Instrument Connection Kit Centaur	Ea
119-7009	422711-04	VersaCell Instrument Connection Kit EXL	Ea
117-2061	422711-03	VersaCell Instrument Connection Kit RxL	Ea
117-2063	422761-05	VersaCell Interface Wheel, Back, G801 Up	Ea
117-3510	422761-06	VersaCell Interface Wheel, Left, G801 Up	Ea
117-3557	422761-04	VersaCell Interface Wheel, Right, G801 Up	Ea
117-2185	422729-02	VersaCell Kit, IMM 2000/2500, Back	Ea
117-2042	422729-03	VersaCell Kit, IMM 2000/2500, Left, G800 Down	Ea
117-2121	422729-01	VersaCell Kit, IMM 2000/2500, Right	Ea

ADVIA 1800 Instruments

Item Code	Manuf Part No.	Item Description	Size
116-1530	073-A021-01	ADVIA 1800 SYSTEM	Ea
117-0474	040-A001-08	ADVIA Automation WorkCell A - 60Hz	Ea
117-0475	040-A001-09	ADVIA Automation WorkCell B - 60Hz	Ea
117-0476	040-A001-10	ADVIA Automation WorkCell C - 60Hz	Ea
116-3109	084-K029-02	CentraLink - ADVIA 2120/120 Connection Kit	Ea
116-3113	084-K041-01	CentraLink - Lab Automation Kit	Ea
116-3114	084-K042-01	CentraLink - Non ADVIA 2120/120 Connection Kit	Ea
116-3107	084-K026-05	CentraLink ADVIA CentraLink Client Kit	Ea
118-5140	10624091	syngo Barcode Scanner	Ea
118-5146	10624078	syngo Country Kit - North America	Ea
118-5143	10624082	syngo Lab Data Manager Instrument Connection License	Ea
118-5141	10638933	syngo Server - Standard Capacity Server	Ea
117-0473	030102-03	Versacell	1/Cr
117-2061	422711-03	VersaCell Instr Connection Kit RxL	1each

SIEMENS

Dimension® EXL™ 200 Integrated Chemistry System

This system offers revolutionary LOCI® advanced chemiluminescent technology and automated, productivity-enhancing features for the smaller-sized laboratory.

Immunoassay:

- 167 tests per hour
- Homogeneous chemiluminescent technology
- Short reaction times: no washing or separation steps
- Excellent STAT testing capability
- Low sample volume required
- **Trusted, proven technologies:**
- LOCI® advanced chemiluminescence, with a maximum throughput of 624 tests per hour (437 photometric tests per hour, 187 ISE tests per hour)
- Heterogeneous module
- Photometry
- QuikLYTE® IMT
- User-defined serum index (HIL) testing
- Automation capabilities
- **Fewer manual tasks:**
- No manual reagent preparation
- No manual sample pretreatment
- No limit or dedicated STAT positions required

#778041.911

(117-8345)ea



Siemens Dimension Test Menu

Item Code	Manuf Part No.	Item Description	Size
114-3156	DC23	Acid Phosphatase Verifier, 2 Tests	1/Bx
114-3193	DF13	Albumin Flex Reagent Cartridge, 480 Tests	1/Bx
114-3165	DC37	Ethyl Alcohol Calibrator, 2 Tests	1/Bx
114-3143	DC131	ALDL Calibrator, 1 Test	1/Bx
114-3291	D830	Aliquot Wheels, 8 Tests	8/Bg
119-8562	DF150	Alk Phos rev, 360 Tests	Ea
119-8564	DC150	ALPI Calibrator, 2 Tests	Ea
114-3201	DF17A	Amylase Flex Reagent Cartridge, 240 Tests	1/Bx
114-7531	DC48B	Automated HDL Calibrator, 3 Tests	Ea
114-7548	DF48B	AHDL Cholesterol, 240 Tests	1/Bx
114-3266	DF96A	Urine Barbiturates Screen Flex Reagent Cartridge, 80 Tests	1/Bx
114-3267	DF97A	Urine Benzodiazapines Screen Flex Reagent Cartridge, 80 Tests	1/Bx
114-3207	DF23A	Calcium Flex Reagent Cartridge, 480 Tests	1/Bx
114-3297	G905	Calomel Electrodes, 1 Test	Ea
114-3255	DF87	Carbamazepine Flex Reagent Cartridge, 80 Tests	1/Bx
114-3579	RC434	DM CCRP Calibrator, 2 Test	1/Bx
114-3599	RF434	DM CCRP Flex, 120 Tests	1/Bx
114-3585	RD701	Chemistry Wash, 1 Test	Ea
114-3299	G908	Chloride Electrode, 1 Test	Ea
114-6658	DC32	CKI/MBI Calibrator, 2 Tests	Ea
114-3220	DF37	C-Reactive Protein Flex Reagent Cartridge, 120 Tests	1/Bx
114-6640	DF32	Creatine Kinase MB Revised, 120 Tests	Ea
114-6638	DF38	Creatine Kinase Revised, 480 Tests	Ea
114-3593	RF421C	Cardiac Troponin-I Flex Reagent Cartridge, 120 Tests	1/Bx
114-3572	RC421C	Cardiac Troponin-I Calibrator, 2 Tests	1/Bx



SIEMENS

Siemens Dimension Test Menu

Item Code	Manuf Part No.	Item Description	Size
114-3289	D828	Cuvette Cartridge, 12000 Tests	1/Bx
114-3184	DF108	Cyclosporine Extended Range Flex Reagent Cartridge, 80 Tests	1/Bx
114-3140	DC108A	Dade Behring CSAE Calibrators, 1 Test	1/Bx
114-6659	DC108	Dade Behring CSAE Calibrators, 1 Test	Ea
114-3909	9A509UL	Emit Calibrator/Control Level 0, 14 Tests	1/Bx
114-3910	9A529UL	Emit Calibrator/Control Level 1, 14 Tests	1/Bx
114-3911	9A549UL	Emit Calibrator/Control Level 2, 14 Tests	1/Bx
114-3912	9A569UL	Emit Calibrator/Control Level 3, 14 Tests	1/Bx
114-3913	9A589UL	Emit Calibrator/Control Level 4, 14 Tests	1/Bx
114-3914	9A609UL	Emit Calibrator/Control Level 5, 14 Tests	1/Bx
114-3616	STM2	STM Aliquot Segments, 5040 Tests	5040/Bx
114-3615	STM1	STM Pipette Tips, 5000 Tests	5000/Bx
117-9093	DF179	System Check Material, 16 Tests	Ea
118-8306	KD694	BNP/TSH Sample Diluent, 6 Tests	1/Kt
114-3606	RF523A	Low Volume NT-proBNP Flex, 72 Tests	1/Bx
114-3575	RC423A	2 Calibrator Sets Per Carton, 2 Tests	1/Bx
114-3596	RF423A	Dimension Flex, 120 Tests	1/Bx
114-3172	DC49D	Drug Calibrator II, 2 Tests	1/Bx
114-3915	9X529UL	Emit® II Plus Ecstasy Calibrator/Control Level I (10 mL), 1 Test	1/Bx
114-3916	9X549UL	Emit® II Plus Ecstasy Calibrator/Control Level 2 (10 mL), 1 Test	1/Bx
114-3917	9X569UL	Emit® II Plus Ecstasy Calibrator/Control Level 3 (10 mL), 1 Test	1/Bx
114-3918	9X589UL	Emit® II Plus Ecstasy Calibrator/Control Level 4 (10 mL), 1 Test	1/Bx
114-3196	DF137	Enzymatic Carbonate Flex Reagent Cartridge, 360 Tests	1/Bx
114-3144	DC137	Enzymatic Carbonate Calib, 2 Tests	1/Bx
119-3535	DF270B	Enzymatic Creatinine Flex Reagent Cartridge, 448 Tests	Ea
114-6660	DC35	Enzyme 1 Calibrator, 2 Tests	Ea
118-8256	KC360	Enzyme 6 Calibrator, 3 Tests	1/Kt
114-3292	ED	Enzyme Diluent, 10 Tests	1/Bx
117-8058	DC143	ENZ II Calibrator, 2 Tests	Ea
114-3583	RC621	DM EXL LOCI Troponin Calibrator, 2 Tests	1/Bx
116-8900	778021.911	EXL with LM, 1 Test	Ea
114-3600	RF440	Ferritin Flex Reagent Cartridge, 120 Tests	1/Bx
114-3580	RC440	Ferritin Calibrator, 2 Tests	1/Bx
114-3602	RF452	Free Prostate Specific Antigen Flex Reagent Cartridge, 120 Tests	1/Bx
114-3589	RF410	Free Thyroxine Flex Reagent Cartridge, 120 Tests	1/Bx
114-3225	DF45A	Gamma-Glutamyl Transferase Flex Reagent Cartridge, 288 Tests	1/Bx
114-3598	RF430	Human Chorionic Gonadotropin Flex Reagent Cartridge, 120 Tests	1/Bx
114-3578	RC430	Human Chorionic Gonadotropin Calibrator, 2 Tests	1/Bx
116-8354	DF105A	Hemoglobin A1c - rev, 120 Tests	Ea
114-3246	DF76	Immunoglobulin G Flex Reagent Cartridge, 120 Tests	1/Bx
114-3588	RD704	IMT Probe Cleaner, 12 Tests	Ea
114-3285	D130	ISE Cleaning Solution D, 2 Tests	Ea
114-3283	D129	ISE Deproteinizing Solution, 2 Tests	1/Bx
114-3280	D124	ISE External Verifier Kit, 1 Test	1/Bx
114-6641	DF54	Lactate Dehydrogenase Revised, 480 Tests	Ea
114-3607	RF530	Low Human Chorionic Gonadotropin Flex Reagent Cartridge, 60 Tests	1/Bx
114-3175	DC56	Liquid Lipase Calibrator, 2 Levels	1/Bx
114-3603	RF520	LMMB, 80 Tests	1/Bx
118-8304	KD692	LOCI CTNI Sample Diluent, 6 Tests	1/Kt
116-9000	RF616	LOCI Free T3, 120 Tests	1/Bx
114-3608	RF610	DM LOCI FT4 Flex, 120 Tests	1/Bx
119-9214	RF723	DM LOCI LV NT-proBNP, 72 Tests	Ea
114-3611	RF623	DM LOCI NT-proBNP, 120 Tests	1/Bx
114-3584	RC623	DM LOCI NTP Calibrator, 2 Tests	1/Bx
114-3610	RF621	DM LOCI TNI Flex, 144 Tests	1/Bx

SIEMENS

Siemens Dimension Test Menu

Item Code	Manuf Part No.	Item Description	Size
114-3610	RF621	DM LOCI TNI Flex, 144 Tests	1/Bx
114-3609	RF612	DM LOCI TSH Flex, 200 Tests	1/Bx
118-8303	KD691	LOCI TSH Sample Diluent, 6 Tests	1/Kt
114-3604	RF521	Troponin LTNI, 72 Tests	1/Bx
114-3141	DC114	MALB Calibrator 1 ml fill 2 sets, 2 Tests	1/Bx
114-3259	DF90A	Urine Methadone Screen Flex Reagent Cartridge, 80 Tests	1/Bx
114-3592	RF420	Mass Creatine Kinase MB Isoenzyme Flex Reagent Cartridge, 160 Tests	1/Bx
114-3571	RC420	Mass Creatine Kinase MB Isoenzyme Calibrator, 2 Tests	1/Bx
114-3576	RC425	DM MPO Calibrator, 2 Tests	1/Bx
114-3577	RC426	DM MPO Control, 5 x 2 mL	1/Bx
114-3624	S570	Multiply Ext Verifier Kit, 1 Test	1/Kt
114-3617	S500A	Multiply IMT Cartridge – Blue, 8 Tests	1/Bx
114-3622	S560	Multiply Standard C, 4 Tests	1/Bx
114-3618	S530	Multiply Urine Diluent, 1 Test	1/Bx
114-3597	RF425	DM Myeloperoxidase, 80 Tests	1/Bx
114-3573	RC422	Myoglobin Calibrator, 2 Tests	1/Bx
114-3294	G214	Na&K/Cl Electrode Fill Sol, 2 Tests	Ea
114-3160	DC28	Pseudocholinesterase Verifier, 2 Tests	1/Bx
116-5770	9J129UL	Phencyclidine, 1 L, 1 Tests Assay Emit II Plus PCP	1/Kt
114-3296	G902	Potassium Electrode, 1 Test	Ea
114-3290	D829	Printer Paper - RxL, 4 Tests	1/Bx
114-3586	RD702	Reagent Probe Cleaner, 1 Test	Ea
114-3634	WATER-4	Purified Water Diluent, 1 Test	Ea
114-3630	S640	QuikLYTE Dilution Check, 1 Test	1/Bx
114-3628	S630	QuikLYTE Flush Solution, 3 Tests	3/Bx
114-3625	S600	QuikLYTE IMT Cartridge, 4 Tests	1/Bx
114-3629	S635	QuikLYTE Sample Diluent, 6 Tests	6/Bx
114-3626	S620	QuikLYTE Standard A, 3 Tests	3/Bx
114-3627	S625	QuikLYTE Standard B, 3 Tests	3/Bx
114-3613	RXV1A	Revised Heteogeneous Immunoassay Module Reaction Vessels, 1000 Tests	Ea
114-3238	DF69A	REV TGL, 480 Tests	1/Bx
114-3200	DF167	DM DBI Flex Reagent, 480 Tests	1/Bx
114-3594	RF422A	Revised Myoglobin Flex Reagent Cartridge, 120 Tests	1/Bx
114-3167	DC38	Salicylate Calibrator, 2 Tests	1/Bx
114-3276	D105	Salt Bridge Solution, 3 Tests	1/Bx
114-3614	SDPLUS	Sample Diluent, SD PLUS, 2 Tests	1/Bx
114-3587	RD703	Sample Probe Cleaner, 1 Test	Ea
114-3270	DSC5	Small Sample Cups, 1000 Tests	1000/Bg
114-3269	DSC4	Sample Cups w/Lids (1.5mL), 1000 Tests	1000/Bg
114-3295	G901	Sodium Electrode, 1 Test	Ea
114-3174	DC51	Special Protein Calibrator, 2 Tests	1/Bx
114-3287	D202	Standard A - AR, 4 Tests	1/Bx
114-3273	D102	Standard A - ES, 12 Tests	Ea
114-3286	D200	Standard C - AR, 4 Tests	1/Bx
118-8455	TSFOIL	StreamLAB Tube Sealer Foil Cartridge, 3 Tests	1/Kt
114-3581	RC452	Total FPSA Calibrator, 2 Tests	1/Bx
114-3148	DC17	Bilirubin Calibrator, 2 Tests	1/Bx
114-3278	D107	TCO2 Electrode Fill Soln, 1 Test	1/Bx
114-3298	G907	TCO2 Electrode Sensor, 1 Test	Ea
114-3275	D104A	TCO2 Release Agent, 6 Tests	6/Bx
114-3569	RC410	Thyroid Calibrator, 2 Tests	1/Bx
114-3245	DF75A	Thyronine Uptake Flex Reagent Cartridge, 160 Tests	1/Bx
114-3197	DF14	Tobramycin Flex Reagent Cartridge, 80 Tests	1/Bx
114-3601	RF451	Total Prostate Specific Antigen Flex Reagent Cartridge, 120 Tests	1/Bx
114-3181	DF103	Transferrin Flex reagent cartridge, 120 Tests	1/Bx
114-3590	RF412	Thyroid Stimulating Hormone Flex Reagent Cartridge, 200 Tests	1/Bx



SIEMENS

Siemens Dimension Test Menu

Item Code	Manuf Part No.	Item Description	Size
114-3247	DF77	Uric Acid Flex Reagent Cartridge, 480 Tests	1/Bx
114-3277	D106	Urine Diluent, 1 Test	Ea
118-8304	KD692	CTNI Sample Diluent, 6 Tests	1/Kt
118-8303	KD691	TSH Sample Diluent, 6 Tests	1/Kt
114-3623	S565	XL Standard C, 4 Tests	1/Bx

Siemens Dimension Accessories

Item Code	Manuf Part No.	Item Description	Size
117-0484	040-K087-01	ADVIA Automation - 1m Spacer Kit, ADVIA WorkCell	1/Kt
117-0486	040-K124-01	ADVIA Automation Dimension RxL ADVIA Interface Kit	1/Kt
117-0492	040-K198-05	ADVIA Automation Vista Connection Interface Kit	1/Kt
116-3125	664-1962-01	CentraLink - Additional 8 Port Terminal Server	Ea
116-3123	664-0679-01	CentraLink - Barcode Scanner: USB Port CCD Scanner LGP6125RU1S-000	Ea
116-3095	084-0004-02	CentraLink - Client User Interface License	Ea
116-3124	664-1782-01	CentraLink - Printer: Lexmark E460DN Laser Jet, 120VAC	Ea
116-3121	615-4188-01	CentraLink - RJ45 - RJ45 Cable (10')	2/Pk
116-3122	615-4188-02	CentraLink - RJ45 - RJ45 Cable (25')	Ea
116-3118	615-4187-01	CentraLink - RJ45 to 9 pin DB adapter (LIS)	Ea
116-3115	515-7188-03	CentraLink - UPS: PowerVar #ABCE1100-11, 120VAC, 1100VA	Ea
114-4427	753000.901	RMS	Ea
116-3936	753000.501	RMS (Reconditioned)	Ea
114-4440	752025.901	RxLMax QCC PowerPak® Field Upgrade Kit	Ea
114-3817	263100.202	StreamLAB Centrifuge UPS	Ea
117-5920	422761-02	VersaCell ADVIA Centaur/ ADVIA 1800/ ExL/RxL 800 down wheel-back	Ea
117-3511	422711-01	VersaCell Instrument Connection Kit A1800 Wheel	Ea
117-3556	422711-02	VersaCell Instrument Connection Kit Centaur	Ea
119-7009	422711-04	VersaCell Instrument Connection Kit EXL	Ea
117-2061	422711-03	VersaCell Instrument Connection Kit RxL	Ea
117-2063	422761-05	VersaCell Interface Wheel, Back, G801 Up	Ea
117-3510	422761-06	VersaCell Interface Wheel, Left, G801 Up	Ea
117-3557	422761-04	VersaCell Interface Wheel, Right, G801 Up	Ea
117-2185	422729-02	VersaCell Kit, IMM 2000/2500, Back	Ea
117-2042	422729-03	VersaCell Kit, IMM 2000/2500, Left, G800 Down	Ea
117-2121	422729-01	VersaCell Kit, IMM 2000/2500, Right	Ea
114-4411	765025.901	XpandPlus QCC PowerPak® Field Upgrade Kit	Ea
116-3109	084-K029-02	CentraLink - ADVIA 2120/120 Connection Kit	Ea
116-3106	084-B003-09	CentraLink - ADVIA CentraLink High End Server Assembly	Ea
116-3111	084-K040-09	CentraLink - Configured High End Server Kit	Ea
116-3117	515-7245-05	CentraLink - Dell 22 Flat Panel LCD Wide Screen Monitor	Ea
116-3113	084-K041-01	CentraLink - Lab Automation Kit	Ea
116-3114	084-K042-01	CentraLink - Non ADVIA 2120/120 Connection Kit	Ea
116-3107	084-K026-05	CentraLink ADVIA CentraLink Client Kit	Ea
117-2061	422711-03	VersaCell Instr Connection Kit RxL	Ea
114-4405	765000.901	Xpand Plus Basic	Ea
114-4449	765000.511	Xpand Plus w/HM Reconditioned	Ea

HORIBA

Medical



ABX Micros 60 Hematology Analyzer

- Complete hematology profile:
16 parameters, 60 samples/hour (open tube) or 55 samples/hour (closed tube), and microsampling: 10 µL of whole blood per cycle
- Concept and technology: impedance and selective lysing (LMG), no compressor (stepper motor), and low noise
- Automatic maintenance: automatic cleaning of the sampling probe and 16 parameters: HGB, HCT, MCV, MCH, MCHC, PLT, MPV, RDW, LYM (% AND #), MON (% AND #), GRAN (% AND #), RBC, WBC, and PLT histograms
- List of options: closed-tube sampling mode (4-position tube holder), bar-code reader (optional), and standalone smart-card reader (stores 78 results with graphics on rewritable memory smart card)

#MIC60002NUA1, Analyzer, Flexible
(109-8690)ea

#MIC60023NUA1B, Analyzer Only
(119-0988)ea

#MIC60021NUA1B, Analyzer with Lite DM Software, Laptop,
and Bar-Code Wand
(117-9276)ea



ABX Micros CRP 200 Analyzer

- Rapid patient diagnosis: 16 CBC parameters with simultaneous CRP option, 3-part WBC differential, whole-blood microsampling, information management, and STAT results
- Complete system designed to fit your needs: minimal daily maintenance, ideal for capillary collection, 4-position tube holder (QC, 500 µL, 3 mL, and 5 mL), and automatic sampling probe cleaning
- Information management: patient history log, cumulative reporting, and CLIA/HIPAA compliant (with easy-to-use tutorial guide)
- Therapeutic decision making and monitoring:
 - Pediatrics: determination of viral or bacterial origin
 - Oncology: patient-monitoring chemotherapy
 - Rheumatology: monitoring anti-inflammation therapy
 - Intensive care: postoperative monitoring
 - Emergency services: rapid patient diagnosis

#MCRP2001NUA1
(110-8430)ea

Micros 60 Test Menu

Item Code	Manuf Part No.	Item Description	Size	Test Menu	Σ = Not available for Analyzer	
					Micros 60 Analyzer Open/Closed	Micros CRP 200 Analyzer
110-5018	501015	ABX CRP Reagent 100/Bx	Ea	Reagent	Σ	•
435-1352	403010	ABX Mini Clean 1ltr	Ea	Reagent	•	•
106-5492	0906014	Alphalyse 360 F/ABX Micros	1/Bt	Reagent	•	Σ
266-5033	0602050	Cobas Minipack Lmg ABX Minidil LMG	150/Pk	Reagent	•	Σ
111-3999	0501016	CRP Calibrator	Ea	Calibrator	Σ	•
111-3997	D01A00146	Micros 45 Minotrol 12 x Tri Level	Ea	Control	•	Σ
110-1813	D01A00145	Micros 60 Minotrol 12 x Tri Level	12/Kt	Control	•	Σ
806-8290	802010	Mini Dil Lmg 10ltr Cube	1 Ea/Ca	Reagent	•	•
102-1792	702010	Mini Lyse 1ltr	Ea	Reagent	•	•
690-8400	2043745	Minocal Calibrator 2 x 2 mL	Kt	Calibrator	•	•
403-5757	2043756	Minotrol Tri Level 12 x 2 mL	Bx	Control	•	Σ
101-1198	2043752	Minotrol Tri Level 6 x 2 mL	6/Bx	Control	•	Σ
110-5908	2042207	Minotrol CRP Level 3, 2 x 2 mL	2/Pk	Control	Σ	•
110-8785	2042205	Minotrol CRP (Low), 2 x 1 mL	2/Pk	Control	Σ	•
110-8786	2042206	Minotrol CRP (Normal), 2 x 2 mL	2/Pk	Control	Σ	•
453-3069	2043748	Minotrol Tri-level Control 3 x 2 mL	Ea	Control	•	Σ
118-2961	D01A00161	Micros 60 Lite Minotrol w/CD 12 x 2.5	Ea	Part	•	Σ
403-8940	2052111	Micros/45 Micros/60 Memory Card	Ea	Part	•	Σ
690-0904	2042017	Smart Card	Ea	Part	•	Σ


ABX Pentra 60 C+ Hematology Analyzer

- Horiba Medical holds 2 patents, Double Hydrodynamic Sequential System (DHSS™) and Multi-Distribution Sampling System (MDSS™)
- Concept and technology: cytochemistry, impedance (real cell volume measurement), optical DHSS technology (analysis of the internal cellular structure by measuring light absorbance), and MDSS technology (precise homogenization of blood samples with reagents)
- No compressor, no shear valve, and less maintenance
- Throughput to 60 samples/hour
- Perfect differentiation of the 5 WBC subpopulations with DHSS technology
- Three histograms for RBC, BAS/WBC, and PLT, together with the 5 DIFF matrix
- Basophils counted through specific channel
- Data management on stand-alone PC
- Microsampling from whole blood (CBC: 30 µL and DIFF: 53 µL)

#P60CP001NUA1

(111-3703)ea

ABX Pentra XL 80 Hematology Analyzer

- HORIBA Medical holds 2 patents, Double Hydrodynamic Sequential System (DHSS™) and Multi-Distribution Sampling System (MDSS™)
- Cytology platform performance: 80 tests/hour, large-capacity autoloader(100 tubes), stat sampling on open or closed tubes, 26 parameters (CBC: 10 and DIFF: 10), microsampling on whole blood (30 µL in CBC mode and 53 µL in CBC + DIFF mode), automatic sample rerun, integrated validation station, and compatible with ABX Pentra ML (Multilink System) to centralize hematology operations
- System memory stores up to 10,000 patient results: precise patient reports show tests results, demographics, graphs, flags, specific dilution ratios in CDR mode, and remarks, with the option to export data
- Quality assurance: 3 active control levels identified by bar-code scanner and results displayed on Levey-Jennings graphs
- Precise, reliable results from DHSS and MDSS technologies

#P80XL001NUA1

(114-7359)ea


ABX Pentra 60 & ABX Pentra XL 80 Hematology Test Menu

Σ = Not available for Analyzer

Item Code	Manuf Part No.	Item Description	Size	Test Menu	Pentra 60	Pentra 80
338-8145	906003	ABX Lysebio 400 mL	Ea	Reagent	•	•
108-9175	906013	Lysebio F/ABX Pentra Hema 400 mL	Ea	Reagent	•	•
348-4269	D01A00112	Diffrol Level Tri-Level, 3 x 3 mL	Ea	Control	•	•
110-1239	D01A00142	Diffrol Tri-Level, 12 x 3 mL	Ea	Control	•	•
338-4034	0206010	Eosinofix Plus 1 L	Ea	Reagent	•	•
690-8400	2043745	Minocal Calibrator	Ea	Calibrator	•	•
338-9153	903010	ABX Cleaner	Ea	Reagent	•	•
338-5479	901020BX	ABX 20 L Diluent	Ea	Reagent	•	•
324-2019	2062015	Diffrol QC Disk	Ea	Part	•	•



QBC STAR Dry Hematology Analyzer

The QBC STAR dry hematology analyzer can deliver a 9-part CBC with the touch of a single button. The STAR's unique microprocessor performs self-calibration, ensuring the accuracy of results.

#429001, STAR Analyzer
(108-4403)ea

Contains: the QBC STAR instrument, 2 rolls of thermal printer paper, screwdriver, forceps & power cord for 120V & 60Hz operation. Patient ID Keyboard, 10-1.8MM black lancets & 10-2.3MM black lancets.



#424304, Hematology Controls
(117-4049)2 x 6/pkg
Contains: 6 vials ea of 2 different blood samples.



#429580, Internal Printer Paper
(932-7453)3/pkg



QBC AUTOREAD™ Plus Centrifugal Hematology System

- Results in 7 minutes; 9 parameters, and a 2-part differential are displayed and printed
- Improved readability and accuracy
- Totally dry procedure needs no user maintenance, no pumps or tubes to clean, no bulky wet reagents, no circulating blood, and no user calibration

#428605, AutoRead Plus System
(168-3721)ea

Contains: Analyzer, software, QBC centrifuge, printer, operator's manual, AccuTube user guide, AccuTube starter pack, universal power pack, pipetter & spacer, workstation, & space-saver stand.

#423406, AccuTube
(168-6070)100/box

#429576, QBC Autoread Plus Analyzer Without Printer
(800-7690)ea

Contains: Analyzer only. Does not include accessories, centrifuge, or printer for 120/230 V, 50/60 Hz operation.

#424613, QBC Autoread Plus Calibration Check Rod
(101-3663)ea

Contains: specially coded metal carrier rod, used to calibrate the QBC Autoread Plus.



Ortho Clinical Diagnostics
a **JOHNSON & JOHNSON** company

Ortho ProVue®

Overview:

- Immunohematology: The first fully automated blood bank instrument for the North American market to be based on the ID-MTS™ Gel Test™
- Operational improvement: Eliminate labor-intensive, time-consuming processes while increasing productivity via random access, discrete testing, and continuous throughput
- Risk management: Enhances process control with dispense verification, bar-code tracking from start to finish, and multilevel password protection
- Quality control features monitor reagent expiration, gel card integrity, and can reduce transcription errors due to a bidirectional interface
- ID-MTS™ Gel Test™ technology: With over 2,300 users in North America, this patented technology, based on column agglutination, can standardize your blood bank's complete menu of testing
- Flexible, powerful automation to meet your changing laboratory needs
- Flexible, powerful automation to meet your changing laboratory needs

Regardless of the type or size of your facility, Ortho offers a variety of affordable and flexible instrument platforms. Since the ID-MTS™ Gel Test™ requires no wash steps or multiple reagent additions, it is easily automated with a choice of affordable and flexible instrument platforms. These solutions secure positive patient identification, increase efficiencies, and maximize the use of resources.

#MTS213784
(117-9664)ea



QBC™ Capillary Tubes

Designed for use with the QBC™ AutoRead Plus System.

#424241, Capillary Tubes
(168-5256)100/box





Ac•T diff 2™ Hematology Analyzer

Each analyzer brings the added convenience and safety of a closed-vial sampling as standard. This fully automated 16-parameter or 18-parameter system offers a CBC and COULTER® histogram differential, providing all of the functions that make producing results fast, accurate, and easy.

Ac•T™ diff2-Tainer with diff Ac•T™ Tainer Tubing (#6605541), 1 Yr. Bus. Hrs. 5-Days/Week Warranty
(331-4254)ea

Ac•T diff™ with diff Ac•T Pak™ Tubing #6605540, 1 Yr. Bus. Hrs. 5-Days/Week Warranty
(331-2293)ea

- 16-parameter CBC with three-part differential
- Target volume: ideal for performing 1–100 CBCs a day
- Closed-vial sampling accommodates a variety of collection tube sizes
- Automated probe wipe safely cleans probe
- Small 18-μL sample size in either closed-vial or open-vial modes
- Compact: only 2 square feet of counter space required
- Efficient data management meets laboratory needs
- Automated quality assurance functions improve productivity
- Controls and calibrators designed for closed-vial sampling
- Zero routine maintenance requirements and high system reliability help contain cost

Ac•T diff 2™ and Ac•T diff™ Reagents, Controls, Calibrator

diff Ac•T™ Tainer, 1 liter
(331-4642)ea

diff Ac•T™ Tainer, 4 x 4 liters
(331-6645)ea

diff Ac•T™ Pak, 15 liters
(823-4044)ea

Ac•T™ Pak, 15 liters
(371-3817)ea

Ac•T™ Rinse, Shut-Down Diluent, 500 mL
(357-3531)ea

4C Plus Tri-Pack Cell Control, 3 x 3.3 mL
(331-7048)3/pkg



Ac•T 5diff™ CP Hematology System

Excellent Accuracy, Increased Safety
Fully automated system that delivers five-part differential testing for complete results in an affordable, easy-to-use package. With the convenience of cap piercing and an extensive data management system, the CP delivers reliability, accuracy, and more consistent patient care every time.

- 20-parameter CBC with five-part differential
- Target volume: ideal for performing 10–100 CBCs/day
- Closed-vial sampling accommodates a variety of collection tube sizes
- Automatic cap piercing with multiple-tube capacity
- Extensive data management capabilities for optimal record keeping
- 53-μL sample size for CBC
- 60-samples-per-hour throughput ensures productivity
- Zero routine maintenance requirement and high system reliability help contain costs

#6605705, Analyzer with Data Management Printer
(331-3976)ea

Reagents, Controls, Calibrator, and Diluent, 20 liters
(686-0122)ea

Reagents, Controls, Calibrator, and Hgb Lyse, 400 mL
(102-4036)ea

Reagents, Controls, Calibrator, and Fix, 1 liter
(102-4038)ea

Reagents, Controls, Calibrator, WBC Lyse, 1 L
(102-4037)ea

Reagents, Controls, Calibrator, and Rinse, 1 liter
(102-4035)ea

Reagents, Controls, and Calibrator, 2 x 2 mL
(331-9685)ea

Ac•T Control Plus 5diff
(861-1263)ea



Access® 2 Immunoassay System

With an intuitive operator interface, enhanced sample handling and innovative menu, the state-of-the-art Access 2 system delivers new levels of laboratory efficiency, standardization and cost effectiveness. A broad menu of high-quality assays from Beckman Coulter and productivity-enhancing features on the Access 2 offer laboratories a flexible option for immunoassay testing. Maximum throughput of 100 tests per hour for one-step assays and 50 tests per hour for two-step assays.

Product features:

- Touch-panel user interface improves data and sample management
- View rack and run completion time on screen, in real time
- Quickly view information on requested test, pending work and test results
- Perform user-defined reflex testing to eliminate repetitive, manual processes
- Link up to 4 systems through a single LIS interface for greater efficiency
- Instruments share information such as test requests and test results



- Streamline work flow with the system's flexible sample handling, using bar-coded primary and/or aliquot tubes
- Magnetic particle separation with chemiluminescent detection technology achieves wide analytical range and fast assay processing

- 24 resident assays with continuous random access and STAT capability
- Onboard refrigeration enhances reagent stability

#973111, Access® 2
(104-7912)ea

Σ = Not available for Analyzer

Access 2 and UniCel Dxl Consumables

Item Code	Manuf Part No.	Item Description	Size	UniCel® Dxl 600 Access® Immunoassay System	UniCel® Dxl 800 Access® Immunoassay System	Access 2 Immunoassay System
331-0152	80769	Aspirate Probe Cleaning Kit, Access® 2	Kit	•	•	•
331-0153	81606	Sample Racks, 13 mm, 1 package, 6/package	6/Pk	•	•	•
331-0154	81608	Sample Racks, 16 mm x 100 mm, 1 package, 6/package	6/Pk	•	•	•
331-0155	81609	Elevated Sample Racks, 16 mm x 75 mm, 1 package, 6/package	6/Pk	•	•	•
331-0156	81901	Reaction/Aliquot Vessel Cartridges, Access®, 1 box, 16 x 98/box	16x98/Bx	•	•	•
331-0161	81908	Sample Diluent A, Access®, 1 vial, 4.0 mL	Bx	•	•	•
331-0162	81910	System Check Solution, Access®, 6 bottles, 4.0 mL/bottle	6/Bx	•	•	•
331-0163	81911	Contrad 70, 1.0 L	Ea	•	•	•
331-0164	81912	Citranox, 1 gallon	1gal	•	•	•
331-0165	81914	Sample Container, 3.0 mL, Access®, 1 bag, 500/bag	500/Bg	•	•	•
331-0166	81915	Insert Cups, 1.0 mL/13 mm, Access®, 1 bag, 1000/bag	1000/Bg	•	•	•
331-0167	81916	Pediatric Insert Cups, 1.0 mL, 1 bag, 1000/bag	1000/Bg	•	•	•
331-0168	81917	Insert Cups, 2.0 mL/16 mm, Access®, 1 bag, 1000/bag	1000/Bg	•	•	•
331-0169	81920	Insert Cup Caps, 1.0 mL/13 mm, Access®, 1 bag, 1000/bag	1000/Bg	•	•	•
331-0171	81922A	Sample Container Caps, 3.0 mL, Access®, 1 bag, 1000/bag	1000/Bx	•	•	•
331-0172	8409B	Aspirate probe, Access® 2	Ea	•	•	•
331-0173	973001	Disposable aspirate probe brushes	10/Pk	•	•	•
331-0174	973052	Sample Tray Barcode Label Kit	57/Kt	•	•	•
117-8884	A16793	Wash Buffer II, UniCel® Dxl, 1 container, 10 liters	1/Bx	•	•	Σ
203-0683	104838	Swab Applicators, Fiber-Free Polyester, 1 package, 100/package	100/Pk	•	•	•
331-0038	386167	Reaction Vessels, UniCel® Dxl 800, 1000/bag, 10 bags/box	10bg/Bx	•	•	Σ
331-0186	81902	Sample Cups, 2.0 mL/13 mm, Access®, 1 bag, 1000/bag	1000/Bg	•	•	•
331-0187	81906	Substrate, Access®, 4 bottles, 130 mL/bottle	4/Bx	•	•	•
331-0188	A16792	Wash Buffer II, Access®, 4 bottles, 1.9 L/bottle	4/Bx	•	•	•
331-0202	81904	Waste Bags, Access®, 1 box, 20 bags/box	20/Bx	•	•	•

**UniCel® Dxl 600 Access®
Immunoassay System**

With up to 50 reagents onboard, the UniCel Dxl 600 Access Immunoassay System enables you to have the flexibility to run both routine and specialized tests. The Dxl 600 allows laboratories to decrease process steps and improve turnaround time, with astounding ease of use. Uses the same reagents as other UniCel Access systems, so your laboratory can deliver consistent patient results, no matter which Beckman Coulter immunoassay system performs the test. A broad menu across a wide variety of disease states enables your laboratory to consolidate workstations and recapture testing currently being sent out. Load reagents and all consumables on the fly, without using the console computer or pausing the instrument, minimizing operator intervention. Large onboard capacity of up to 60 samples at once with continuous loading sample racks allow samples to be loaded directly from the centrifuge, reducing manual processes.

Product features:

- Onboard capacity of 50 refrigerated reagent packs enhance flexibility
- Large onboard reagent capacity simplifies reagent management
- Rapid assay kinetics improve turnaround times
- System supports multiple assay protocols
- Magnetic particle separation with enhanced chemiluminescent detection delivers exceptional clinical performance
- Routine sample loading can be interrupted to load STAT samples quickly

- "Point-in-space" pipetting technology aligned with CLSI guidelines allows integration with laboratory automation systems
- Automatically creates aliquots based on the volume of requested tests
- System can be configured to aspirate additional volume required for user-defined reflex or rerun test scenarios



**Why have 50
Reagent spaces
on-board?**

**To provide the
flexibility you
want.**

- System stores aliquots for up to 3 hours in a refrigerated, low-evaporation compartment
- Auto-matched pipettors work in parallel to deliver sample and reagent into the reaction vessels
- Achieves a peak throughput of 200 tests per hour

#A45867, UniCel® Dxl 600
(864-0017)ea

Σ = Not available for Analyzer

Access 2 and UniCel Dxl Test Menu

Item Code	Manuf Part No.	Item Description	Size	UniCel® Dxl 600 Access® Immunoassay System	UniCel® Dxl 800 Access® Immunoassay System	Access 2 Immunoassay System
109-2668	A12985	Access TPO Antibody Reagent, 100 Determinations, 2 x 50 tests	2x50/Bx	•	•	Σ
109-2669	A18227	Access TPO Antibody Calibrators (S0-S5), 1 set, 6 x 2.0 mL	1/Bx	•	•	•
109-3044	A16365	EPO Calibrators (S0-S5), Access®, 1 set, S0 - 10.0 mL; S1-S5 - 2.5 mL	6/Kt	•	•	•
109-3392	A16364	Access EPO Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	•
109-7646	A10827	Access DHEA-S Calibrators (S0-S5), 1 set, 6 x 2.0 mL	Bx	•	•	•
110-7415	A32898	Access Thyroglobulin Antibody II Reagent, 100 Determinations, 2 x 50 tests	100/Bx	•	•	•
110-7489	A36920	Access Thyroglobulin Antibody II Calibrators (S0-S5), 1 set, S0 - 4.0 mL; S1-S5 - 2.5 mL	Ea	•	•	•
114-2975	A48617	Access SHBG Reagent, 100 Determinations, 2 x 50 tests	2x50/Bx	•	•	•
114-4521	A48618	Access SHBG Calibrators (S0-S5), 1.0 mL per vial per level	6/Bx	•	•	•
115-4627	A48619	Access SHBG Quality Control (QC) (QC1-QC2), 3 vials/level, 2 x 2.0 mL	Ea	•	•	•
115-4628	A36100	Access Inhibin A QC (QC1-QC3), 2 vials/level, 6 x 2.5 mL	Ea	•	•	•
115-4630	A36097	Access Inhibin A Reagent, 100 Determinations, 2 x 50 tests	Ea	•	•	•
117-3886	A98032	Access Folate Reagent, 100 Determinations, 2 x 50 tests	2x50/Ea	•	•	•
117-3889	A98033	Access Folate Calibrator (S0-S5), 4.0 mL per vial	5/Bx	•	•	•
117-9902	A32493	Access sTfR Reagent, 100 Determinations, 2 x 50 tests	Ea	•	•	•
331-0054	33006	Access Vitamin B12 Calibrator S0, 1 x 4.0 mL	Kit	•	•	•
331-0060	33025	Access Ferritin Calibrators (S0-S5), 1 set, 6 x 4.0 mL	6/Kt	•	•	•
331-0061	33200	Access CEA Reagent, 100 determinations, 2 x 50 tests	100/Kt	•	•	•



Access 2 and UniCel DxI Test Menu

Σ = Not available for Analyzer

Item Code	Manuf Part No.	Item Description	Size	UniCel® DxI 600 Access® Immunoassay System	UniCel® DxI 800 Access® Immunoassay System	Access 2 Immunoassay System
331-0062	33205	Access CEA Calibrators (S0-S5) 1 set, 6 x 2.5 mL	6/Kt	•	•	•
331-0063	33206	Access CEA Diluent, 1 vial, 4.0 mL	Kit	•	•	•
331-0064	33209	Access CEA QC (QC1-QC2), 3 vials/level, 6 x 2.5 mL	6/Kt	•	•	•
331-0065	33211	Access AFP Reagent (300 Test Kit), 6 x 50 Tests, and Calibrators (S0-S6), 1 calibrator set, 7 x 2.5 mL	300/Kt	•	•	•
331-0066	33216	Access AFP Sample Diluent, 1 vial, 14.0 mL	Kit	•	•	•
331-0075	33500	Access Total Beta-hCG Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	•
331-0076	33505	Access Total Beta-hCG Calibrators (S0-S5), 1 set, 6 x 4.0 mL	6/Kt	•	•	•
331-0079	33520	Access hFSH Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	•
331-0080	33525	Access hFSH Calibrators (S0-S5), 1 set, 6 x 4.0 mL	6/Kt	•	•	•
331-0081	33530	Access Prolactin Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	•
331-0082	33535	Access Prolactin Calibrators (S0-S5), S0 - 4.0 mL; S1-S5 - 2.5 mL	6/Kt	•	•	•
331-0083	33540	Access Estradiol Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	•
331-0084	33545	Access Estradiol Calibrators (S0-S5), S0 - 4.0 mL; S1-S5 - 2.5 mL	6/Kt	•	•	•
331-0085	33546	Access Estradiol Calibrator S0, 1 x 4.0 mL	Kit	•	•	•
331-0086	33550	Access Progesterone Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	•
331-0087	33555	Access Progesterone Calibrators (S0-S5), S0 - 4.0 mL; S1-S5 - 2.5 mL	6/Kt	•	•	•
331-0088	33556	Access Progesterone Calibrator S0, 1 x 4.0 mL	Kit	•	•	•
331-0091	33570	Access Unconjugated Estriol Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	•
331-0092	33575	Access Unconjugated Estriol Calibrators (S0-S6), S0 4.0 mL; S1-S6 2.5 mL	7/Kt	•	•	•
331-0093	33580	Access Ultrasensitive hGH Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	Σ
331-0094	33585	Access Ultrasensitive hGH Calibrators (S0-S5), 1 set, 6 x 2.0 mL	6/Kt	•	•	•
331-0095	33600	Access Cortisol Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	Σ
331-0096	33605	Access Cortisol Calibrators (S0-S5), 1 set, 6 x 4.0 mL	6/Kt	•	•	•
331-0097	33606	Access Cortisol Calibrator S0, 1 x 4.0 mL	Kit	•	•	•
331-0100	33710	Access Digoxin Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	•
331-0101	33715	Access Digoxin Calibrators (S0-S5), 1 set, 6 x 4.0 mL	6/Kt	•	•	•
331-0102	33716	Access Digoxin Calibrator S0, 1 x 4.0 mL	Kit	•	•	•
331-0103	33800	Access Total T4 Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	Σ
331-0104	33805	Access Total T4 Calibrators (S0-S5), 1 set, 6 x 4.0 mL	6/Kt	•	•	•
331-0109	33830	Access Total T3, 100 Determinations, 2 x 50 tests	100/Kt	•	•	•
331-0110	33835	Access Total T3 Calibrators (S0-S5), 1 set 6 x 4.0 mL	6/Kt	•	•	•
331-0113	33860	Access Thyroglobulin Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	Σ
331-0114	33865	Access Thyroglobulin Calibrators (S0-S5), 1 set, 6 x 2.0 mL	6/Kt	•	•	•
331-0115	33866	Access Thyroglobulin Sample Diluent, 1 x 14.0 mL	Kit	•	•	Σ
331-0116	33880	Access Free T4 Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	•
331-0117	33885	Access Free T4 Calibrators (S0-S5), 6 x 2.5 mL	6/Kt	•	•	•
331-0122	34430	Access Rubella IgG Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	•
331-0123	34435	Access Rubella IgG Calibrators (S0-S5), 1 set, 6 x 1.0 mL	6/Kt	•	•	•
331-0124	34439	Access Rubella IgG QC (QC1-QC2), 3 vials per level; 6 x 2.5 mL	6/Kt	•	•	•
331-0128	34470	Access Toxo IgM II Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	•
331-0129	34475	Access Toxo IgM II Calibrators (C0-C1), 1 set, 2 x 1.5 mL	2/Kt	•	•	•
331-0130	34479	Access Toxo IgM II QC (QC1-QC2), 3 vials/level, 6 x 3.5 mL	6/Kt	•	•	•
331-0131	35000	Access Total IgE Reagent, 100 determinations, 2 x 50 tests	100/Kt	•	•	Σ

Access 2 and UniCel DxI Test Menu

Σ = Not available for Analyzer

Item Code	Manuf Part No.	Item Description	Size	UniCel® DxI 600 Access® Immunoassay System	UniCel® DxI 800 Access® Immunoassay System	Access 2 Immunoassay System
331-0132	35005	Access Total IgE Calibrators (S0-S6), 1 set, S0 - 6.0 mL; S1-S6-4.0 mL	6/Kt	•	•	•
331-0133	35006	Access Total IgE Calibrator S0, 1 vial, 6.0 mL	Kit	•	•	•
331-0134	37200	Access Hybritech PSA Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	•
331-0135	37205	Access Hybritech PSA Calibrators (S0-S5), 1 set, 6 x 2.5 mL	6/Kt	•	•	•
331-0136	37206	Access Hybritech PSA Sample Diluent, 1 x 14.0 mL	Kit	•	•	•
331-0137	37209	Access Hybritech PSA QC (QC1-QC3), 3 levels, 5.0 mL each	3/Kt	•	•	•
331-0140	37219	Access Hybritech free PSA QC (QC1-QC2), 1 vial/level, 2 x 5.0 mL	2/Kt	•	•	•
331-0141	37300	Access Ostase Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	•
331-0142	37305	Access Ostase Calibrators (S0-S5), 1 set, 6 x 2.5 mL	6/Kt	•	•	•
331-0143	37309	Access Ostase QC (QC1-QC2), 1 set, 2 x 4.0 mL	2/Kt	•	•	•
331-0144	386357	Access OV Monitor Reagent, 100 determinations, 2 x 50 tests	100/Kt	•	•	•
331-0145	386358	Access OV Monitor Calibrators (S0-S5), 1 set, 6 x 2.5 mL	6/Kt	•	•	•
331-0146	386371	Access CK-MB Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	•
331-0148	387620	Access BR Monitor Reagent, 100 determinations, 2 x 50 tests	100/Kt	•	•	•
331-0149	387647	Access BR Monitor Calibrators (S0-S5), 1 set, 6 x 1.5 mL	6/Kt	•	•	•
331-0150	387687	Access GI Monitor Reagent, 100 determinations, 2 x 50 tests	100/Kt	•	•	•
331-0151	387688	Access GI Monitor Calibrators (S0-S5), 1 set, 6 x 2.5 mL	6/Kt	•	•	•
331-0176	973244	Access Myoglobin Calibrators (S0-S5), 1 set, 6 x 1.0 mL	6/Kt	•	•	•
331-0177	33005	Access Vitamin B12 Calibrators (S0-S5), 6 x 4.0 mL	6/Kt	•	•	•
331-0190	33000	Access Vitamin B12 Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	•
331-0192	A16972	Access Intact PTH Reagent, 100 determinations, 2 x 50 tests	1/Bx	•	•	•
331-0193	A16953	Access Intact PTH Calibrators (S0-S5, RB), 1 set, S0-S5 - 1.0 mL, RB - 2 x 4.0 mL	1/Bx	•	•	•
331-0196	A99250	Access Folate Calibrator S0, 4 mL per vial	1/Kt	•	•	•
331-0201	386372	Access CK-MB Calibrators (S0-S5), 1 set, 6 x 2.0 mL	6/Kt	•	•	•
331-1081	33815	Access Thyroid Uptake Calibrator (S0), 1 set, 6 x 1.0 mL	6/Kt	•	•	•
331-1929	33825	Access HYPERsensitive hTSH Calibrators (S0-S5), 1 set, 6 x 4.0 mL	6/Kt	•	•	•
331-4528	33810	Access Thyroid Uptake Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	•
331-6439	33415	Access Ultrasensitive Insulin Calibrators (S0-S5), 6 x 2.0 mL	6/Kt	•	•	•
331-7037	33410	Access Ultrasensitive Insulin Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	•
331-7403	A10826	Access DHEA-S Reagent, 100 Determinations, 2 x 50 tests	2x50/Bx	•	•	•
331-7405	33020	Access Ferritin Reagent, 100 determinations, 2 x 50 tests	100/Kt	•	•	•
331-7408	A13430	Access Free T3 Calibrators (S0-S5), 1 set, 6 x 2.5 mL	6/Kt	•	•	•
331-7409	A13422	Access Free T3 Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	•
331-7412	37210	Access Hybritech Free PSA Reagent, 100 Determinations, 2 x 50 tests	2x50/Kt	•	•	•
331-7414	37215	Access Hybritech Free PSA Calibrators(S0-S5), 1 set, S0 5.0mL; S1-S5 2.5mL	6x2.5/Kt	•	•	•
331-9133	33515	Access hLH Calibrators (S0-S5), 1 set, 6 x 4.0 mL	6/Kt	•	•	•
331-9501	33565	Access Testosterone Calibrators (S0-S5), 1 set, 6 x 2.5 mL	6/Kt	•	•	•
331-9890	33560	Access Testosterone Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	Σ
331-9959	33510	Access hLH Reagent, 100 Determinations, 2 x 50 tests	100/Kt	•	•	•
372-6559	33820	Access HYPERsensitive hTSH/Fast TSH Reagent, 100 Determinations, 2 x 50 tests	1/Kt	•	•	Σ
861-0546	A32494	Access sTfR Calibrators (S0-S5), 1 set, S0 - 4.0 mL; S1-S5 - 2.5 mL	Ea	•	•	•
861-0563	973243	Access Myoglobin Reagent, 100 determinations, 2 x 50 tests	100/Kt	•	•	•

CLINICAL LABORATORY



Provide convenient and quality patient care in your practice while you also increase profitability.

Henry Schein Medical's LabTeam™ is a single source for advice within the laboratory setting. With specific expertise in lab consultation and analysis, the LabTeam™ improves your patients' experiences in the office, while navigating profitability and regulatory requirements.

“ Our Clinical Lab Specialist has been very helpful to our bottom line by assisting in the research and selection of proper billing codes for us to obtain reimbursement. The analyzers he has introduced to our practice are state-of-the-art. I would not hesitate to recommend Henry Schein and their Clinical Lab Specialists to any physician looking to set up a lab or upgrade a current one. ”

Larry W. Mayfield
CEO
Ocala Family Medical Center



FEATURES:

- Consultation
- Design
- Financing
- Selection
- Implementation
- Support/Training

HENRY SCHEIN®
LabTeam™
SOLUTIONS FOR CLINICAL LABS

BENEFITS:

- Improve patient satisfaction with convenient onsite testing
- Expedite ability to diagnose with rapid turnaround of results
- New or expanded sources of revenue
- Stay current on technology/reimbursement/regulatory changes
- Custom laboratory management solutions for cost-effective operations

Contact Your Lab Consultant

1.800.428.9200

E-mail: labquestions@henryschein.com

Visit: www.henryschein.com/labteam

HENRY SCHEIN®
MEDICAL



FINANCIAL SOLUTIONS



Take ownership of business decisions that impact the financial health of your practice.

Henry Schein Financial Services (HSFS) provides avenues for optimizing financial health and establishing a platform for ongoing growth.



FEATURES:

- Financing, Leasing and Tools
- Risk Management
- Tax Strategy
- Practice Start Up/Expansion



FINANCIAL SERVICES

Business Solutions for Healthcare Practitioners

BENEFITS:

- Experts in medical practice financials
- Industry-leading approval rating
- Low interest rates
- Accelerated timeline – fast to close
- Flexible terms
- Streamlined process due to strong relationships with medical-focused lenders

1.800.443.2756

Visit: www.henryschein.com/mhsfs



SIEMENS

Immolute® 1000 Immunoassay System

Small benchtop immunoassay analyzer. The extensive menu, low cost of operation, reliability, and ease of use makes this system a great choice for laboratories with low-volume immunoassay tests.

Dependable operation:

- Throughput up to 120 tests per hour
- Reliability allows laboratories confidence in assay turnaround times
- Cost of ownership is low due to minimal service, operator interventions, and operator time

Improved laboratory efficiency:

- Few consumables
- Daily maintenance in less than 10 minutes
- Sample and reagent level sensors
- Easy entry of assay parameters, expiration dates, and calibration adjustment intervals
- Integration of the computer, keyboard, and monitor provide a consolidated footprint

Flexible solutions:

- Turbo assays (Troponin I, CK-MB, Myoglobin, Intact PTH, and HCG) allow results in 15 minutes for critical tests
- 12-position reagent carousel

#10380058

(114-3635).....ea



Immolute Test Menu

Item Code	Manuf Part No.	Item Description	Size	Product Type
114-3406	L2UNZ	3G Allergy Sample Diluent	Ea	Control
114-3405	L2UNCVM	3G Allergy Specific IGE CVM	Ea	Control
393-0577	LTSZ	3rd Generation TSH Sample Diluent	Ea	Control
400-7098	LACCM	ACTH Control Module	2/Bx	Control
114-3300	LACCVM	ACTH CVM	Ea	Control
114-3303	LAPCVM	AFP CVM	Ea	Control
358-3099	LAPZ	AFP Sample Diluent	Bt	Control
114-3317	LHACVM	Albumin (Urinary) CVM	Ea	Control
719-4838	LHACM	Albumin Control Module	2/Bx	Control
670-0735	LHAZ	Albumin Sample Diluent	Ea	Control
114-3342	L2ATC	Allergen Cap	1000/Bx	Control
114-3567	PL2	Allergen Panel 20T	Ea	Reagent
114-3568	PL4	Allergen Panel 40T	Ea	Reagent
111-6649	L2ATS2	Allergen Septa	250/Bg	Control
114-3341	L2	Allergens 20T	Ea	Reagent
114-3409	L4	Allergens 40T	Ea	Reagent
114-3302	LAOCVM	Androstenedione CVM	Ea	Control
115-7036	L2MCCVM	Anti-HBC IGM CVM	Ea	Control
115-7034	LAHCVM	Anti-HBS CVM	Ea	Control
114-0754	LTGCVM	Anti-TG AB CVM	Ea	Control
114-0753	LTOCVM	Anti-TPO AB (Im1000) CVM	Ea	Control
114-7293	10485435	Ascorbic Acid 2pk	Ea	Control
357-9907	BMCM	B2-Microglobulin Control Module	Ea	Control
114-6652	LBMCM	Beta-2 Microglobulin CVM Cal Verifier BMG 4 Levels	Ea	Control
414-0937	LBMZ	Beta-2 Microglobulin Sample Diluent	Bx	Control
114-3305	LBRCVM	BR-MA (Im1000) CVM	Ea	Control
358-3100	LBRZ	BR-MA Sample Diluent	Ea	Control

Immolute® 2000 XPi Immunoassay System

This highly reliable system is easy to use and has one of the largest automated immunoassay menus available. Medium- to high-volume laboratories can incorporate specialty and allergy testing into routine workloads to improve productivity and efficiency.

Enhanced laboratory workflow:

- Maximum throughput of 200 tests per hour
- Rerun, dilution, and reflex capabilities minimize technologist's time
- Auto-start functionality—including automated daily maintenance and quality control scheduling—reduces hands-on operations

Flexible solutions:

- Cost-effective automation available by connecting to VersaCell® solutions
- Scalable connectivity to Siemens track-based automation solutions

Improved efficiency:

- Intuitive, easy-to-use software
- 19" touch screen monitor for easy navigation
- Color-coded consumable levels always available at a glance
- User configurable dilution and reflex testing to meet laboratory needs

(116-8979).....ea



Immolute Test Menu

Item Code	Manuf Part No.	Item Description	Size	Product Type
114-3307	LCLCM	Calcitonin Control Module	2/VI	Control
114-3308	LCLCVM	Calcitonin CVM	Ea	Control
833-6025	K9CON	Canine Control Module	Ea	Control
116-2873	LCTCVM	Canine T4 CVM	Ea	Control
833-6026	K9TCM	Canine Thyroid Control Module	2x3 mL/Ea	Control
611-7524	LLICM	Canine TLI Control Module	Ea	Control
116-5391	LLICVM	Canine TLI CVM	Ea	Control
116-5393	LKTCVM	Canine TSH CVM	Ea	Control
116-5912	LCBCVM	Carbamazepine CVM	Ea	Control
108-4541	LCECVM	CEA CVM	4/Kt	Control
114-3306	LCEZ	CEA Sample Diluent	Ea	Control
116-9013	LMBCVM	CK-MB (IM1000 / 2000) CVM 4 x 2 mL	Ea	Control
114-3310	LCOCVM	Cortisol CVM	Ea	Control
366-5989	LCOZ	Cortisol Sample Diluent 1 x 25 mL	24/Ca	Control
611-1048	PECM	C-Peptide Control Module	Ea	Control
114-3332	LPEPCVM	C-Peptide CVM	Ea	Control
114-3397	L2PEZ	C-Peptide Sample Diluent	Ea	Control
114-5951	L2CRZ	CRP Sample Diluent	Ea	Control
114-3312	LDSCVM	Dhea-So4 CVM	Ea	Control
115-7033	LDSZ	Dhea-So4 Sample Diluent	Ea	Control
114-3311	LDICVM	Digoxin CVM	Ea	Control
114-3407	L2ZC	Diluent Caps	250/Pk	Control
114-3408	L2ZT	Diluent Tubes	250/Pk	Control
645-8455	LEPCM	EPO Control Module	3/Bx	Control
114-3314	LEPCVM	EPO CVM	Ea	Control
350-0208	E2TC1	Estradiol Control Module	Bx	Control
114-3345	L2E2Z	Estradiol Sample Diluent	Ea	Control

SIEMENS

Immolute Test Menu

Item Code	Manuf Part No.	Item Description	Size	Product Type
112-4460	LE2Z	Estradiol Sample_ Diluent	Ea	Control
110-7206	LFECVM	Ferritin CVM	Ea	Control
780-0516	LFEZ	Ferritin Sample Diluent	1/Pk	Control
114-1009	L2FVZ	Folic Acid / Vitamin B12 Sample Diluent	Ea	Control
110-8340	LFOCVM	Folic Acid CVM	Ea	Control
781-9617	LFOZ	Folic Acid Sample Diluent	Ea	Control
114-3333	LPFCVM	Free PSA CVM	Ea	Control
108-4550	LF3CVM	Free T3 CVM	4/Kt	Control
114-1105	LFT4CVM	Free T4 CVM	Ea	Control
806-5122	LFSZ	FSH Sample Diluent (IM 1000)	Ea	Control
114-3346	L2FSZ	FSH Sample Diluent (IM 2000)	Ea	Control
117-2689	LGACM	Gastrin Control Module	Ea	Control
110-8565	LRXT	Gear Driven Reaction Tubes	1000/Bg	Control
114-3316	LGHCV	Growth Hormone CVM	Ea	Control
115-7040	LHPGCM	H Pylori Control Module	Ea	Control
110-8341	LHPGCV	H Pylori IGG CVM	Ea	Control
114-3343	L2CGZ	HCG (IM 2000/2500) Sample Diluent (1x50 mL)	Ea	Control
888-7720	LCGZ4	HCG (IM1000) Sample Diluent (1x100 mL)	Ea	Control
851-4579	LCGZ	HCG (IM1000) Sample Diluent (1x25 mL)	Ea	Control
114-7218	LCGCM	HCG Control Module	Ea	Control
108-4545	LCGCV	HCG CVM	4/Kt	Control
114-5934	L2HOZ	Homocystein Sample Diluent	Ea	Control
110-7208	LHOCVM	Homocysteine CVM	Ea	Control
613-0262	LCRPCVM	HS-CRP CVM	4/Kt	Control
988-5311	LGCOCM	IGF Control Module	2/St	Control
114-3315	LGFCVM	IGF-1 CVM	Ea	Control
114-5670	LGBCVM	IGF-BP3 CVM	Ea	Control
117-6827	L2GFZ	IGF-I Sample Diluent	20mL/Bt	Control
119-3731	L1KIGW1	IGG/IGM Sample Diluent	Ea	Control
116-2025	LKUE31	IML 1000 Unconjugated Estriol 100 Tests	Ea	Reagent
114-3387	L2KTS6M	IML 2000 3rd Gen. Tsh Kit 6000 Tests	Ea	Reagent
114-0827	L2KTS6	IML 2000 3rd Gen. Tsh Kit 600 Tests	Ea	Reagent
110-7444	L2KTS2	IML 2000 3rd Gen.Tsh Kit 200 Tests	Ea	Reagent
113-7863	L2KAC2	IML 2000 ACTH Kit 200 Tests	Ea	Reagent
113-8532	L2KAP2	IML 2000 AFP Kit 200 Tests	Ea	Reagent
114-3349	L2KAP6	IML 2000 AFP Kit 600 Tests	1/Kt	Reagent
110-7438	L2KHA2	IML 2000 Albumin Kit 200 Tests	Ea	Reagent
113-8533	L2KA02	IML 2000 Androstenedione Kit 200 Tests	Ea	Reagent
113-8524	L2KMC2(D)	IML 2000 Anti-HBC LGM Kit 200 Tests Domestic	Ea	Reagent
113-8523	L2KHC2(D)	IML 2000 Anti-HBC Kit 200 Tests Domestic	Ea	Reagent
113-8525	L2KAB2(D)	IML 2000 Anti-HBS KIT 200 Tests - Rev	Ea	Reagent
114-3348	L2KAH6(D)	IML 2000 Anti-HBS Kit 600 Tests Domestic	1/Kt	Reagent
110-7443	L2KTG2	IML 2000 Anti-TG Kit 200 Tests	Ea	Reagent
114-3382	L2KTG6	IML 2000 Anti-TG Kit 600 Tests	Ea	Reagent
110-7445	L2KT02	IML 2000 Anti-TPO Kit 200 Tests	Ea	Reagent
114-3385	L2KT06	IML 2000 Anti-TPO Kit 600 Tests	1/Kt	Reagent
113-8555	L2KBM2	IML 2000 Beta-2 Microglobulin 200 Tests	Ea	Reagent
113-8548	L2KBR2	IML 2000 BR-MA Kit 200 Tests	Ea	Reagent
114-3350	L2KBR6	IML 2000 BR-MA Kit 600 Tests	1/Kt	Reagent
113-8442	L2KCL2	IML 2000 Calcitonin Kit 200 Tests	Ea	Reagent
114-2306	L2KCT2	IML 2000 Canine T4 Kit 200 Tests	1/Bx	Reagent
114-3354	L2KCT6	IML 2000 Canine T4 Kit 600 Tests	1/Kt	Reagent
114-2307	L2KL12	IML 2000 Canine TLI Kit 200 Tests	1/Bx	Reagent
114-3370	L2KKT6	IML 2000 Canine TSH Kit 600 Tests	1/Kt	Reagent
113-3064	L2KCB2	IML 2000 Carbamazepine Kit 200 Tests	200/Bx	Reagent



Immulite Test Menu

Item Code	Manuf Part No.	Item Description	Size	Product Type
113-8549	L2KCE2	IML 2000 CEA Kit 200 Tests	Ea	Reagent
113-8520	L2KCVG2	IML 2000 CMV Igg Kit 200 Tests	Ea	Reagent
114-3355	L2KCVG6	IML 2000 CMV Igg Kit 600 Tests	1/Kt	Reagent
110-7433	L2KCO2	IML 2000 Cortisol Kit 200 Tests	Ea	Reagent
114-3352	L2KCO6	IML 2000 Cortisol Kit 600 Tests	1/Kt	Reagent
114-3374	L2KPEP2	IML 2000 C-Peptide Kit 200 Tests	1/Kt	Reagent
114-3375	L2KPEP6	IML 2000 C-Peptide Kit 600 Tests	1/Kt	Reagent
113-8535	L2KDS2	IML 2000 Dhea-SO4 Kit 200 Tests	Ea	Reagent
114-3356	L2KDS6	IML 2000 Dhea-SO4 Kit 600 Tests	1/Kt	Reagent
113-3074	L2KDI2	IML 2000 Digoxin Kit 200 Tests	200/Bx	Reagent
114-3359	L2KEP6	IML 2000 EPO Kit 600 Tests	1/Kt	Reagent
117-9021	L2KE22	IML 2000 Estradiol Kit 200 Tests	1/Bx	Reagent
113-3066	L2KFE2	IML 2000 Ferritin Kit 200 Tests	200/Bx	Reagent
113-0735	L2KFE6	IML 2000 Ferritin Kit 600 Tests	Ea	Reagent
112-4548	L2KF02	IML 2000 Folic Acid Kit 200 Tests	Ea	Reagent
210-7153	L2KF06	IML 2000 Folic Acid Kit 600 Tests	1/Bx	Reagent
111-9414	L2KPF2(D)	IML 2000 Free PSA Kit 200 Tests Domestic	Ea	Reagent
110-7436	L2KF32	IML 2000 Free T3 Kit 200 Tests	Ea	Reagent
114-3364	L2KF33M	IML 2000 Free T3 Kit 3000 Tests	1/Kt	Reagent
114-0828	L2KF36	IML 2000 Free T3 Kit 600 Tests	Ea	Reagent
110-7435	L2KFS2	IML 2000 FSH Kit 200 Tests	Ea	Reagent
116-6981	L2KFS6	IML 2000 FSH Kit 600 Tests	Ea	Reagent
113-8556	L2KGA2	IML 2000 Gastrin Kit 200 Tests	Ea	Reagent
116-8016	L2KGRH2	IML 2000 Growth Hormone Kit 200 Tests Rev	Ea	Reagent
113-8527	L2KHGP2	IML 2000 H. Pylori Semi-Quant. 200 Tests	Ea	Reagent
114-3367	L2KHGP6	IML 2000 H. Pylori Semi-Quant. 600 Tests	Ea	Reagent
113-0837	L2KHB2(D)	IML 2000 HBSAG Kit 200 Tests Domestic	Ea	Reagent
113-8536	L2KCG2	IML 2000 HCG Kit 200 Tests	Ea	Reagent
114-3351	L2KCG6	IML 2000 HCG Kit 600 Tests	1/Kt	Reagent
115-2861	L2KHVG6	IML 2000 Herpes L&I IGG Kit 600 Tests	1/Bx	Reagent
110-7432	L2KCRP2	IML 2000 High Sensitivity CRP 200 Tests	Ea	Reagent
114-3353	L2KCRP6	IML 2000 High Sensitivity CRP 600 Tests	1/Kt	Reagent
113-8518	L2KGB2	IML 2000 IGFBP-3 Kit 200 Tests	Ea	Reagent
113-8519	L2KGF2	IML 2000 IGF-I Kit 200 Tests	Ea	Reagent
113-8515	L2KIN2	IML 2000 Insulin Kit 200 Tests	Ea	Reagent
110-7439	L2KLH2	IML 2000 LH Kit 200 Tests	Ea	Reagent
114-3371	L2KLH6	IML 2000 LH Kit 600 Tests	1/Kt	Reagent
113-8513	L2KMY2	IML 2000 Myoglobin Kit 200 Tests	Ea	Reagent
113-8550	L2KOP2	IML 2000 OM-Ma Kit 200 Tests	Ea	Reagent
113-8551	L2KPA2	IML 2000 Pap Kit 200 Tests	Ea	Reagent
114-3373	L2KPA6	IML 2000 Pap Kit 600 Tests	1/Kt	Reagent
113-3061	L2KPB2	IML 2000 Phenobarbital Kit 200 Tests	200/Bx	Reagent
113-3077	L2KPN2	IML 2000 Phenytoin Kit 200 Tests	200/Bx	Reagent
113-8540	L2KPW2	IML 2000 Progesterone Kit 200 Tests	Ea	Reagent
114-3378	L2KPW6	IML 2000 Progesterone Kit 600 Tests	1/Kt	Reagent
110-7442	L2KPR2	IML 2000 Prolactin Kit 200 Tests	Ea	Reagent
114-3377	L2KPR6	IML 2000 Prolactin Kit 600 Tests	1/Kt	Reagent
117-7376	L2KPTS6(D)	IML 2000 PSA Kit 600 Tests Domestic	1/Bx	Reagent
110-7441	L2KPP2	IML 2000 PTH (M/P) Kit 200 Tests	Ea	Reagent
114-3376	L2KPP6	IML 2000 PTH (M/P) Kit 600 Tests	1/Kt	Reagent
113-8508	L2KPD2	IML 2000 Pylilinks-D Kit 200 Tests	Ea	Reagent
113-3076	L2KRT2	IML 2000 Rapid TSH Kit 200 Tests	200/Bx	Reagent
114-3380	L2KRT3M	IML 2000 Rapid TSH Kit 3000 Tests	1/Kt	Reagent
113-0833	L2KRT6	IML 2000 Rapid TSH Kit 600 Tests	Ea	Reagent
113-8528	L2KRM2(D)	IML 2000 Rubella IGM Kit 200 Tests Domestic	Ea	Reagent

SIEMENS

Immulite Test Menu

Item Code	Manuf Part No.	Item Description	Size	Product Type
113-8529	L2KRUB2	IML 2000 Rubella Quant. IGG Kit 200 Tests	Ea	Reagent
113-8541	L2KSH2	IML 2000 SHBG Kit 200 Tests	Ea	Reagent
114-8357	L2KSY2(D)	IML 2000 Syphilis Screen 200 Tests Domestic	Ea	Reagent
114-8360	L2KSY6	IML 2000 Syphilis Screen 600 Tests	Ea	Reagent
113-8543	L2KTB2	IML 2000 TBG Kit 200 Tests	Ea	Reagent
113-3075	L2KTN2	IML 2000 Theophylline Kit 200 Tests	200/Bx	Reagent
114-3384	L2KTN6	IML 2000 Theophylline Kit 600 Tests	1/Kt	Reagent
113-8547	L2KTY2(D)	IML 2000 Thyroglobulin Kit 200 Tests Domestic	Ea	Reagent
114-3381	L2KTC6	IML 2000 Tobramycin Kit 600 Tests	Ea	Reagent
111-6224	L2KIE2	IML 2000 Total IGE Kit 200 Tests	200/Bx	Reagent
114-3368	L2KIE6	IML 2000 Total IGE Kit 600 Tests	1/Kt	Reagent
110-7449	L2KT32	IML 2000 Total T3 Kit 200 Tests	Ea	Reagent
114-3392	L2KT36M	IML 2000 Total T3 Kit 6000 Tests	1/Kt	Reagent
114-3391	L2KT36	IML 2000 Total T3 Kit 600 Tests	Ea	Reagent
110-7450	L2KT42	IML 2000 Total T4 Kit 200 Tests	Ea	Reagent
114-3393	L2KT43M	IML 2000 Total T4 Kit 3000 Tests	1/Kt	Reagent
114-3394	L2KT46	IML 2000 Total T4 Kit 600 Tests	1/Kt	Reagent
110-7446	L2KTW2	IML 2000 Total Testosterone Kit 200 Tests	Ea	Reagent
114-3389	L2KTW6	IML 2000 Total Testosterone Kit 600 Tests	Ea	Reagent
113-8531	L2KTXP2	IML 2000 Toxo. Quant. IGG Kit 200 Tests	Ea	Reagent
114-3390	L2KTXP6	IML 2000 Toxo. Quant. IGG Kit 600 Tests	Ea	Reagent
113-8530	L2KTM2(D)	IML 2000 Toxo. IGM (U-Cap)Kit 200 Tests Domestic	Ea	Reagent
114-3383	L2KTM6(D)	IML 2000 Toxo. IGM (U-Cap)Kit 600 Tests Domestic	Ea	Reagent
113-0729	L2KTU2	IML 2000 T3-Uptake Kit 200 Tests	Ea	Reagent
113-8542	L2KUE32	IML 2000 Unconjugated Estriol 200 Tests	Ea	Reagent
114-3358	L2KEF6	IML 2000 Unconjugated Estriol 600 Tests	Ea	Reagent
113-3134	L2KVA2	IML 2000 Valproic Acid Kit 200 Tests	200/Bx	Reagent
110-7451	L2KVB2	IML 2000 Vitamin B12 Kit 200 Tests	Ea	Reagent
926-6861	L2KVB6	IML 2000 Vitamin B12 Kit 600 Tests	1/Bx	Reagent
115-2406	FP24L4	Imm. Food Panel # 24 ASR 40 Tests	Kt	Reagent
115-2417	FP6L4	Imm. Food Panel # 6 40 Tests	Kt	Reagent
115-2428	GP3L4	Imm. Grass Panel # 3 40 Tests	Kt	Reagent
114-3007	F207L2	Imm. Clam	20/Bx	Reagent
358-0126	E88L2	Imm. Mouse	20/Bx	Reagent
358-0021	F13L4	Imm. Peanut	40/Bx	Reagent
358-0188	F210L2	Imm. Pineapple	20/Bx	Reagent
358-0077	T19L2	Imm. Acacia	20/Bx	Reagent
358-0172	D70L2	Imm. Acarus Siro	20/Bx	Reagent
114-4727	M12L2	Imm. Aerobasidium Pullulans	20/Pk	Reagent
114-3038	T2L4	Imm. Alder	40/Bx	Reagent
358-0025	F20L2	Imm. Almond	20/Bx	Reagent
358-0269	F76L4	Imm. Alpha Lactalbumin	40/Bx	Reagent
358-0152	M6L4	Imm. Alternaria Tenuis	40/Bx	Reagent
358-0151	I206L2	Imm. American Cockroach	20/Bx	Reagent
115-2396	EP1L4	Imm. Animal Panel # 1 40 Tests	Kt	Reagent
115-2397	EP2L4	Imm. Animal Panel # 2 ASR 40 Tests	Kt	Reagent
115-2398	EP70L4	Imm. Animal Panel # 70 ASR 40 Tests	Kt	Reagent
115-2399	EP71L4	Imm. Animal Panel # 71 ASR 40 Tests	Kt	Reagent
115-2400	EP72L4	Imm. Animal Panel # 72 ASR 40 Tests	Kt	Reagent
117-8741	F49L4	Imm. Apple	40/Bx	Reagent
358-0206	F237L2	Imm. Apricot	20/Bx	Reagent
358-0066	M3L4	Imm. Asp. Fumigatus	40/Bx	Reagent
358-0215	F261L2	Imm. Asparagus	20/Bx	Reagent
114-4726	M207L2	Imm. Aspergillus Niger	20/Bx	Reagent
358-0358	T73L2	Imm. Australian Pine	20/Bx	Reagent
358-0285	F96L2	Imm. Avocado	20/Bx	Reagent



SIEMENS

Immulite Test Menu

Item Code	Manuf Part No.	Item Description	Size	Product Type
115-2463	W67L2	Imm. Baccharis	20/Kt	Reagent
358-0046	G17L2	Imm. Bahia Grass	20/Bx	Reagent
358-0034	F45L2	Imm. Baker'S Yeast	20/Bx	Reagent
114-3020	F92L2	Imm. Banana	20/Bx	Reagent
358-0038	F6L2	Imm. Barley	20/Bx	Reagent
358-0218	F269L2	Imm. Basil ASR	20/Bx	Reagent
358-0349	T218L2	Imm. Bayberry/Sweet Gale	20/Bx	Reagent
114-3042	T5L2	Imm. Beech	20/Bx	Reagent
358-0030	F27L2	Imm. Beef	20/Bx	Reagent
358-0047	G2L4	Imm. Bermuda Grass	40/Bx	Reagent
358-0270	F77L4	Imm. Beta Lactoglobulin	40/Bx	Reagent
114-3039	T3L4	Imm. Birch	40/Bx	Reagent
114-3012	F280L2	Imm. Black Pepper	20/Bx	Reagent
358-0170	D201L2	Imm. Blomia Tropicalis	20/Bx	Reagent
358-0254	F37L2	Imm. Blue Mussel	20/Bx	Reagent
358-0231	F288L2	Imm. Blueberry	20/Bx	Reagent
358-0135	F18L2	Imm. Brazil Nut	20/Bx	Reagent
358-0353	T401L2	Imm. Brazilian Peppertree	20/Bx	Reagent
114-3015	F403L2	Imm. Brewer's Yeast	20/Bx	Reagent
358-0214	F260L2	Imm. Broccoli	20/Bx	Reagent
358-0132	F11L2	Imm. Buckwheat	20/Bx	Reagent
358-0117	E78L2	Imm. Budgerigar Feather ASR	20/Bx	Reagent
358-0192	F216L2	Imm. Cabbage	20/Bx	Reagent
358-0105	E201L2	Imm. Canary Feathers	20/Bx	Reagent
358-0296	G71L2	Imm. Canary Grass	20/Bx	Reagent
358-0337	M5L4	Imm. Candida Albicans	40/Bx	Reagent
358-0177	F102L2	Imm. Cantaloupe	20/Bx	Reagent
358-0378	W82L2	Imm. Careless Weed	20/Bx	Reagent
358-0243	F31L4	Imm. Carrot	40/Bx	Reagent
358-0271	F78L4	Imm. Casein	40/Bx	Reagent
358-0182	F202L2	Imm. Cashew	20/Bx	Reagent
358-0016	E1L4	Imm. Cat Dander Epithelium	40/Bx	Reagent
358-0234	F291L2	Imm. Cauliflower	20/Bx	Reagent
358-0276	F85L2	Imm. Celery	20/Bx	Reagent
114-5270	M202L2	Imm. Cephalo. Acrem.	20/Bx	Reagent
358-0273	F81L2	Imm. Cheddar Cheese	20/Bx	Reagent
358-0274	F82L2	Imm. Cheese-Mold Type	20/Bx	Reagent
358-0208	F242L2	Imm. Cherry	20/Bx	Reagent
358-0237	F299L2	Imm. Chestnut	20/Bx	Reagent
358-0123	E85L2	Imm. Chicken Feathers	20/Bx	Reagent
358-0042	F83L2	Imm. Chicken Meat	20/Bx	Reagent
358-0226	F279L2	Imm. Chili-Pepper	20/Bx	Reagent
358-0178	F105L2	Imm. Chocolate	20/Bx	Reagent
358-0261	F50L2	Imm. Chub Mackerel	20/Bx	Reagent
115-2422	F220L2	Imm. Cinnamon	20/Kt	Reagent
358-0065	M2L4	Imm. Cladosporium Herba	40/Bx	Reagent
358-0155	W13L2	Imm. Cocklebur	20/Bx	Reagent
358-0058	I6L2	Imm. Cockroach	20/Bx	Reagent
114-3014	F36L2	Imm. Coconut	20/Bx	Reagent
358-0145	F3L4	Imm. Codfish	40/Bx	Reagent
358-0139	F221L2	Imm. Coffee	20/Bx	Reagent
358-0087	W1L4	Imm. Common Ragweed	40/Bx	Reagent
358-0050	G7L2	Imm. Common Reed Grass	20/Bx	Reagent
358-0373	W43L2	Imm. Common Sagebrush	20/Bx	Reagent
358-0150	F8L4	Imm. Corn	40/Bx	Reagent
358-0074	T14L2	Imm. Cottonwood	20/Bx	Reagent
358-0107	E4L2	Imm. Cow Dander	20/Bx	Reagent

SIEMENS

Immulite Test Menu

Item Code	Manuf Part No.	Item Description	Size	Product Type
114-5189	F23L2	Imm. Crab	20 Tests	Reagent
358-0209	F244L2	Imm. Cucumber	20/Bx	Reagent
358-0290	G14L2	Imm. Cultivated Oat Grass	20/Bx	Reagent
358-0288	G12L4	Imm. Cultivated Rye Grass	40/Bx	Reagent
358-0291	G15L2	Imm. Cultivated Wheat Grass	20/Bx	Reagent
114-5272	M16L2	Imm. Curvularia Lunata	20/Bx	Reagent
358-0095	W8L2	Imm. Dandelion	20/Bx	Reagent
358-0171	D3L2	Imm. Derm. Microceras	20/Bx	Reagent
358-0015	D2L4	Imm. Derma. Farinae	40/Bx	Reagent
358-0014	D1L4	Imm. Derma. Pteronyssinus	40/Bx	Reagent
358-0019	E5L4	Imm. Dog Dander	40/Bx	Reagent
358-0017	E2L4	Imm. Dog Epithelium	40/Bx	Reagent
358-0374	W46L2	Imm. Dog Fennel	20/Bx	Reagent
358-0124	E86L2	Imm. Duck Feathers	20/Bx	Reagent
114-3009	F245L2	Imm. Egg	20/Bx	Reagent
358-0020	F1L4	Imm. Egg White	40/Bx	Reagent
358-0040	F75L4	Imm. Egg Yolk	40/Bx	Reagent
358-0085	T8L2	Imm. Elm	20/Bx	Reagent
358-0096	W9L4	Imm. English Plantain	40/Bx	Reagent
114-4607	M14L2	Imm. Epicoccum Purpura	20/Bx	Reagent
358-0076	T18L2	Imm. Eucalyptus	20/Bx	Reagent
358-0176	D74L2	Imm. Euroglyphus Maynei	20/Bx	Reagent
114-4606	W4L2	Imm. False Ragweed	20/Bx	Reagent
358-0365	W17L2	Imm. Firebush	20/Bx	Reagent
358-0179	F147L2	Imm. Flounder	20/Bx	Reagent
115-2402	FP1L4	Imm. Food Panel # 1	40/Kt	Reagent
115-2403	FP13L4	Imm. Food Panel # 13 ASR	40/Kt	Reagent
115-2404	FP15L4	Imm. Food Panel # 15 ASR	40/Kt	Reagent
115-2405	FP2L4	Imm. Food Panel # 2	40/Kt	Reagent
115-2407	FP25L4	Imm. Food Panel # 25 ASR	40/Kt	Reagent
115-2409	FP26L4	Imm. Food Panel # 26 ASR	40/Kt	Reagent
115-2410	FP27L4	Imm. Food Panel # 27 ASR	40/Kt	Reagent
115-2412	FP28L4	Imm. Food Panel # 28 ASR	40/Kt	Reagent
115-2413	FP3L4	Imm. Food Panel # 3	40/Kt	Reagent
115-2414	FP5L4	Imm. Food Panel # 5	40/Kt	Reagent
115-2415	FP50L4	Imm. Food Panel # 50 ASR	40/Kt	Reagent
115-2416	FP51L4	Imm. Food Panel # 51 ASR	40/Kt	Reagent
115-2418	FP7L4	Imm. Food Panel # 7	40/Kt	Reagent
115-2419	FP73L4	Imm. Food Panel # 73 ASR	40/Kt	Reagent
115-2420	FP8L4	Imm. Food Panel # 8 ASR	40/Kt	Reagent
114-4728	M9L2	Imm. Fusarium Moniliforme	20/Bx	Reagent
358-0035	F47L2	Imm. Garlic	20/Bx	Reagent
114-3047	W3L2	Imm. Giant Ragweed	20/Bx	Reagent
358-0143	F270L2	Imm. Ginger	20/Bx	Reagent
114-3017	F79L4	Imm. Gluten	40/Bx	Reagent
358-0175	D73L2	Imm. Glycyphagus Domesticus	20/Bx	Reagent
358-0118	E80L2	Imm. Goat Epithelium	20/Bx	Reagent
358-0363	W12L2	Imm. Goldenrod	20/Bx	Reagent
358-0110	E70L2	Imm. Goose Feathers	20/Bx	Reagent
358-0213	F259L2	Imm. Grape	20/Bx	Reagent
358-0187	F209L2	Imm. Grapefruit	20/Bx	Reagent
115-2426	GP1L4	Imm. Grass Panel # 1	40/Kt	Reagent
115-2427	GP2L4	Imm. Grass Panel # 2	40/Kt	Reagent
115-2429	GP4L4	Imm. Grass Panel # 4	40/Kt	Reagent
358-0246	F315L2	Imm. Green Bean	20/Bx	Reagent
358-0133	F12L2	Imm. Green Pea	20/Bx	Reagent



SIEMENS

Immulite Test Menu

Item Code	Manuf Part No.	Item Description	Size	Product Type
358-0217	F263L2	Imm. Green Pepper	20/Bx	Reagent
358-0108	E6L2	Imm. Guinea Pig Epithelium	20/Bx	Reagent
358-0260	F42L2	Imm. Haddock	20/Bx	Reagent
358-0241	F303L2	Imm. Halibut	20/Bx	Reagent
358-0134	F17L4	Imm. Hazelnut	40/Bx	Reagent
358-0352	T4L4	Imm. Hazelnut Tree	40/Bx	Reagent
114-4602	M8L2	Imm. Helminthosporium Halodes	20/Pk	Reagent
358-0299	I1L4	Imm. Honey Bee Venom	40/Bx	Reagent
114-4604	M45L2	Imm. Hormodendrum Horde	20/Pk	Reagent
358-0018	E3L4	Imm. Horse Dander	40/Bx	Reagent
114-3028	H2L4	Imm. House Dust Hollister-Stier	40/Bx	Reagent
115-2430	HP1L4	Imm. House Dust Panel # 1	Kt	Reagent
358-0059	I70L2	Imm. Imported Fire Ant	20/Bx	Reagent
115-2431	IP1L4	Imm. Inhalant Panel # 1 ASR	40/Bx	Reagent
115-2432	IP10L4	Imm. Inhalant Panel # 10 ASR	40/Bx	Reagent
115-2433	IP2L4	Imm. Inhalant Panel # 2 ASR	40/Bx	Reagent
115-2434	IP3L4	Imm. Inhalant Panel # 3 ASR	40/Bx	Reagent
115-2435	IP4L4	Imm. Inhalant Panel # 4 ASR	40/Bx	Reagent
115-2436	IP5L4	Imm. Inhalant Panel # 5 ASR	40/Bx	Reagent
115-2437	IP6L4	Imm. Inhalant Panel # 6 ASR	40/Bx	Reagent
115-2438	IP7L4	Imm. Inhalant Panel # 7 ASR	40/Bx	Reagent
115-2439	IP8L4	Imm. Inhalant Panel # 8 ASR	40/Bx	Reagent
115-2440	IP9L4	Imm. Inhalant Panel # 9 ASR	40/Bx	Reagent
114-5996	T23L4	Imm. Italian Cypress	40/Bx	Reagent
358-0346	T17L4	Imm. Japanese Cedar	40/Bx	Reagent
358-0045	G10L2	Imm. Johnson Grass	20/Bx	Reagent
358-0051	G8L4	Imm. June Grass - Kentucky Blue	40/Bx	Reagent
358-0275	F84L2	Imm. Kiwi Fruit	20/Bx	Reagent
358-0278	F88L2	Imm. Lamb	20/Bx	Reagent
358-0088	W10L2	Imm. Lamb's Quarters	20/Bx	Reagent
358-0060	K82L4	Imm. Latex	40/Bx	Reagent
358-0186	F208L2	Imm. Lemon	20/Bx	Reagent
358-0173	D71L2	Imm. Lepidoglyphus Destructor	20/Bx	Reagent
358-0191	F215L2	Imm. Lettuce	20/Bx	Reagent
358-0181	F182L2	Imm. Lima Bean	20/Bx	Reagent
115-2423	F306L2	Imm. Lime	20/Bx	Reagent
114-3031	T103L2	Imm. Live Oak Tree	20/Bx	Reagent
114-3041	T43L2	Imm. Loblolly Pine	20/Bx	Reagent
358-0272	F80L2	Imm. Lobster	20/Bx	Reagent
114-3019	F90L2	Imm. Malt	20/Bx	Reagent
358-0280	F91L2	Imm. Mango	20/Bx	Reagent
358-0069	T1L2	Imm. Maple	20/Bx	Reagent
114-8958	G4L2	Imm. Meadow Fescue	20/Bx	Reagent
358-0292	G16L2	Imm. Meadow Foxtail Grass	20/Bx	Reagent
114-3035	T21L2	Imm. Melaleuca	20/Bx	Reagent
358-0043	F87L2	Imm. Melon	20/Bx	Reagent
358-0078	T20L2	Imm. Mesquite	20/Bx	Reagent
358-0136	F2L4	Imm. Milk	40/Bx	Reagent
115-2395	DP1L4	Imm. Mite Panel # 1 ASR	40/Bx	Reagent
115-2443	MP1L4	Imm. Mold Panel # 1	40/Bx	Reagent
358-0111	E71L2	Imm. Mouse Epithelium	20/Bx	Reagent
358-0116	E76L2	Imm. Mouse Serum Protein	20/Bx	Reagent
358-0112	E72L2	Imm. Mouse Urine	20/Bx	Reagent
358-0067	M4L2	Imm. Mucor Racemosus	20/Bx	Reagent
358-0094	W6L4	Imm. Mugwort	40/Bx	Reagent
358-0189	F212L2	Imm. Mushroom	20/Bx	Reagent

SIEMENS

Immulite Test Menu

Item Code	Manuf Part No.	Item Description	Size	Product Type
358-0279	F89L2	Imm. Mustard	20/Bx	Reagent
358-0093	W20L2	Imm. Nettle	20/Bx	Reagent
114-3046	T7L4	Imm. Oak	40/Bx	Reagent
358-0084	T77L2	Imm. Oak Mix - Red, White, Black	20/Bx	Reagent
358-0039	F7L2	Imm. Oat	20/Bx	Reagent
115-2442	KP1L4	Imm. Occupational Panel # 1 ASR	40/Bx	Reagent
358-0086	T9L4	Imm. Olive	40/Bx	Reagent
358-0036	F48L2	Imm. Onion	20/Bx	Reagent
114-3013	F33L4	Imm. Orange	40/Bx	Reagent
114-3022	G3L4	Imm. Orchard Grass	40/Bx	Reagent
358-0228	F283L2	Imm. Oregano ASR	20/Bx	Reagent
358-0377	W7L2	Imm. Ox-Eye Daisy	20/Bx	Reagent
358-0233	F290L2	Imm. Oyster	20/Bx	Reagent
358-0057	I4L2	Imm. Paper Wasp Venom	20/Bx	Reagent
358-0369	W21L4	Imm. Parietaria Judaica	40/Bx	Reagent
358-0366	W19L4	Imm. Parietaria Officinalis	40/Bx	Reagent
358-0128	E91L2	Imm. Parrot Feather ASR	20/Bx	Reagent
358-0277	F86L2	Imm. Parsley ASR	20/Bx	Reagent
358-0284	F95L2	Imm. Peach	20/Bx	Reagent
358-0282	F94L2	Imm. Pear	20/Bx	Reagent
358-0079	T22L2	Imm. Pecan	20/Bx	Reagent
358-0026	F201L2	Imm. Pecan Nut	20/Bx	Reagent
358-0332	M305L2	Imm. Penicillium Brevi	20/Bx	Reagent
358-0063	M1L4	Imm. Penicillium Notatum	40/Bx	Reagent
358-0267	F65L2	Imm. Perch	20/Bx	Reagent
358-0048	G5L4	Imm. Perennial Rye Grass	40/Bx	Reagent
114-4730	M13L2	Imm. Phoma Betae	20/Bx	Reagent
358-0109	E7L2	Imm. Pigeon Droppings	20/Bx	Reagent
358-0210	F253L2	Imm. Pine Nut ASR	20/Bx	Reagent
358-0146	F300L2	Imm. Pinto Bean	20/Bx	Reagent
358-0183	F203L2	Imm. Pistachio	20/Bx	Reagent
358-0339	M70L2	Imm. Pityrosporum Orbiculare	20/Bx	Reagent
358-0212	F255L2	Imm. Plum	20/Bx	Reagent
358-0362	T96L2	Imm. Poplar	20/Bx	Reagent
358-0029	F26L2	Imm. Pork	20/Bx	Reagent
358-0031	F35L2	Imm. Potato	20/Bx	Reagent
114-3032	T210L2	Imm. Privet	20/Bx	Reagent
358-0199	F225L2	Imm. Pumpkin	20/Bx	Reagent
358-0083	T72L2	Imm. Queen Palm	20/Bx	Reagent
358-0371	W36L2	Imm. Rabbit Bush	20/Bx	Reagent
358-0120	E82L4	Imm. Rabbit Epithelium	40/Bx	Reagent
358-0113	E73L2	Imm. Rat Epithelium	20/Bx	Reagent
358-0115	E75L2	Imm. Rat Serum Proteins	20/Bx	Reagent
358-0125	E87L2	Imm. Rat Urine	20/Bx	Reagent
358-0114	E74L2	Imm. Rat Urine Proteins	20/Bx	Reagent
114-3033	T219L2	Imm. Red Cedar	20/Bx	Reagent
358-0230	F287L2	Imm. Red Kidney Bean	20/Bx	Reagent
114-3037	T27L2	Imm. Red Maple	20/Bx	Reagent
114-3044	T71L2	Imm. Red Mulberry	20/Bx	Reagent
115-2425	F381L2	Imm. Red Snapper	20/Bx	Reagent
114-4731	M11L2	Imm. Rhizopus Nigricans	20/Bx	Reagent
358-0044	F9L4	Imm. Rice	40/Bx	Reagent
358-0090	W16L2	Imm. Rough Marsh Elder	20/Bx	Reagent
358-0156	W14L2	Imm. Rough Pigweed	20/Bx	Reagent
358-0089	W11L2	Imm. Russian Thistle	20/Bx	Reagent
358-0149	F5L2	Imm. Rye	20/Bx	Reagent



SIEMENS

Immulite Test Menu

Item Code	Manuf Part No.	Item Description	Size	Product Type
114-5190	F41L2	Imm. Salmon	20/Bx	Reagent
358-0372	W37L2	Imm. Saltbush	20/Bx	Reagent
358-0364	W15L2	Imm. Scale	20/Bx	Reagent
358-0250	F338L2	Imm. Scallop	20/Bx	Reagent
358-0129	F10L2	Imm. Sesame Seed	20/Bx	Reagent
358-0119	E81L2	Imm. Sheep Epithelium	20/Bx	Reagent
358-0091	W18L2	Imm. Sheep Sorrel	20/Bx	Reagent
358-0207	F24L4	Imm. Shrimp	40/Bx	Reagent
358-0249	F337L2	Imm. Sole	20/Bx	Reagent
358-0022	F14L4	Imm. Soybean	40/Bx	Reagent
358-0190	F214L2	Imm. Spinach	20/Bx	Reagent
358-0064	M10L2	Imm. Stemphylium Botryosum	20/Bx	Reagent
115-2444	M88L2	Imm. Stemphylium Solani	20/Bx	Reagent
358-0033	F44L2	Imm. Strawberry	20/Bx	Reagent
114-3008	F21L2	Imm. Sugar Cane	20/Bx	Reagent
358-0324	K84L2	Imm. Sunflower	20/Bx	Reagent
358-0262	F54L2	Imm. Sweet Potato	20/Bx	Reagent
358-0286	G1L2	Imm. Sweet Vernal Grass	20/Bx	Reagent
358-0121	E83L2	Imm. Swine Epithelium	20/Bx	Reagent
358-0072	T11L2	Imm. Sycamore	20/Bx	Reagent
358-0049	G6L4	Imm. Timothy Grass	40/Bx	Reagent
358-0341	O201L2	Imm. Tobacco	20/Bx	Reagent
358-0028	F25L4	Imm. Tomato	40/Bx	Reagent
358-0052	G9L2	Imm. Top Grass	20/Bx	Reagent
115-2445	TP1L4	Imm. Tree Panel # 1	40/Bx	Reagent
115-2446	TP2L4	Imm. Tree Panel # 2	40/Bx	Reagent
115-2447	TP3L4	Imm. Tree Panel # 3 ASR	40/Bx	Reagent
115-2448	TP4L4	Imm. Tree Panel # 4	40/Bx	Reagent
115-2449	TP5L4	Imm. Tree Panel # 5 ASR	40/Bx	Reagent
115-2450	TP6L4	Imm. Tree Panel # 6	40/Bx	Reagent
115-2451	TP7L4	Imm. Tree Panel # 7	40/Bx	Reagent
115-2452	TP8L4	Imm. Tree Panel # 8 ASR	40/Bx	Reagent
115-2453	TP9L4	Imm. Tree Panel # 9	40/Bx	Reagent
358-0325	M15L2	Imm. Trichoderma Viride	20/Bx	Reagent
358-0184	F204L2	Imm. Trout	20/Bx	Reagent
114-3016	F40L2	Imm. Tuna	20/Bx	Reagent
358-0127	E89L2	Imm. Turkey Feathers	20/Bx	Reagent
358-0229	F284L2	Imm. Turkey Meat	20/Bx	Reagent
358-0174	D72L2	Imm. Tyrophagus Putrecen	20/Bx	Reagent
358-0203	F234L2	Imm. Vanilla ASR	20/Bx	Reagent
358-0289	G13L2	Imm. Velvet Grass	20/Bx	Reagent
114-3010	F256L2	Imm. Walnut	20/Bx	Reagent
358-0070	T10L2	Imm. Walnut Tree	20/Bx	Reagent
358-0147	F329L2	Imm. Watermelon	20/Bx	Reagent
115-2455	WP1L4	Imm. Weed Panel # 1	40/Bx	Reagent
115-2456	WP2L4	Imm. Weed Panel # 2	40/Bx	Reagent
115-2457	WP3L4	Imm. Weed Panel # 3	40/Bx	Reagent
115-2459	WP5L4	Imm. Weed Panel # 5 ASR	40/Bx	Reagent
115-2460	WP6L4	Imm. Weed Panel # 6 ASR	40/Bx	Reagent
115-2462	WP7L4	Imm. Weed Panel # 7 ASR	40/Bx	Reagent
358-0092	W2L2	Imm. Western Ragweed	20/Bx	Reagent
358-0032	F4L4	Imm. Wheat	40/Bx	Reagent
358-0205	F236L2	Imm. Whey	20/Bx	Reagent
358-0075	T15L2	Imm. White Ash	20/Bx	Reagent
115-2454	T37L2	Imm. White Bald Cypress	20/Bx	Reagent
358-0055	I2L2	Imm. White Faced Hornet	20/Bx	Reagent
114-3040	T41L2	Imm. White Hickory	20/Bx	Reagent

SIEMENS

Immulite Test Menu

Item Code	Manuf Part No.	Item Description	Size	Product Type
358-0082	T70L2	Imm. White Mulberry	20/Bx	Reagent
358-0154	T16L2	Imm. White Pine	20/Bx	Reagent
358-0295	G70L2	Imm. Wild Rye Grass	20/Bx	Reagent
358-0073	T12L2	Imm. Willow	20/Bx	Reagent
358-0157	W75L2	Imm. Wing Scale	20/Bx	Reagent
358-0375	W5L2	Imm. Wormwood	20/Bx	Reagent
358-0302	I5L2	Imm. Yellow Hornet	20/Bx	Reagent
358-0056	I3L2	Imm. Yellow Jacket Venom	20/Bx	Reagent
358-0148	F360L2	Imm. Yogurt	20/Bx	Reagent
358-0351	T280L2	Imm. Locust Tree	20/Bx	Reagent
358-0348	T211L2	Imm. Sweet Gum	20/Bx	Reagent
358-0023	F15L2	Imm. White Bean	20/Bx	Reagent
325-4745	LKAC1	Immulite® 1000 ACTH-100 Tests	1/Bx	Reagent
102-6776	LKAP1	Immulite® 1000 AFP-100 Tests	Ea	Reagent
114-3319	LKA01	Immulite® 1000 Androstenedione-100 Tests	1/Kt	Reagent
943-5814	LKHC5	Immulite® 1000 Anti-HBC 500 Tests	1/Bx	Reagent
606-6494	LKHC1(D)	Immulite® 1000 Anti-HBC Domestic-100 Tests	1/Bx	Reagent
110-6205	LKMC1(D)	Immulite® 1000 Anti-HBC IGM Domestic-100 Tests	1/Bx	Reagent
984-9469	LKMC5	Immulite® 1000 Anti-HBC IGM 500 Tests	1/Kt	Reagent
984-1141	LKAH1(D)	Immulite® 1000 Anti-HBS Domestic-100 Tests	1/Kt	Reagent
984-9289	LKAH5(D)	Immulite® 1000 Anti-HBS Domestic-500 Tests	1/Kt	Reagent
885-2474	LKTG1	Immulite® 1000 Anti-TG AB 100 Tests	Ea	Reagent
885-2475	LKT01	Immulite® 1000 Anti-TPO AB 100 Tests	Ea	Reagent
358-3101	LKBM1	Immulite® 1000 Beta-2 Microglobulin 100 Tests	Ea	Reagent
120-0760	LKBR1	Immulite® 1000 BR-MA (Ca15-3) 100 Tests	Ea	Reagent
693-6628	LKCL1	Immulite® 1000 Calcitonin 100 Tests	1/Bx	Reagent
231-4394	LKLI1	Immulite® 1000 Canine TLI 100 Tests	1/Bx	Reagent
106-1742	LKCT5	Immulite® 1000 Canine Total T4 500 Tests	500/Bx	Reagent
833-9213	LKCT1	Immulite® 1000 Canine Total T4 100 Tests	100/Ea	Reagent
833-9216	LKKT1	Immulite® 1000 Canine TSH 100 Tests	Ea	Reagent
967-8077	LKCB1	Immulite® 1000 Carbamazepine-100 Tests	Ea	Reagent
358-3102	LKCE1	Immulite® 1000 CEA-100 Tests	100/Bx	Reagent
868-0311	LKMB1	Immulite® 1000 CK-MB-100 Tests	1/Bx	Reagent
111-3509	LSKCP1	Immulite® 1000 CK-MB Turbo 100 Tests	100/Box	Reagent
114-6280	LKCV1	Immulite® 1000 CMV IGG 100 Tests	Ea	Reagent
833-6028	LKCO1	Immulite® 1000 Cortisol-100 Tests	Ea	Reagent
104-2751	LKPE1	Immulite® 1000 C-Peptide-100 Tests	1/Bx	Reagent
114-3325	LKPEP1	Immulite® 1000 C-Peptide-100 Tests	1/Kt	Reagent
611-1046	LKDS1	Immulite® 1000 DHEA-SO4-100 Tests	100/Bx	Reagent
116-8543	LKDG1	Immulite® 1000 Digoxin-100 Tests	1/Bx	Reagent
968-1265	LKDI1	Immulite® 1000 Digoxin 100 Tests	Ea	Reagent
110-5133	LKE25	Immulite® 1000 Estradiol 500 Tests	500/Bx	Reagent
968-4452	LKFE1	Immulite® 1000 Ferritin-100 Tests	Ea	Reagent
114-3326	LKPF1(D)	Immulite® 1000 FPSA Domestic-100 Tests	1/Kt	Reagent
831-6902	LKF31	Immulite® 1000 Free T3 100 Tests	100/Bx	Reagent
113-2841	LKF35	Immulite® 1000 Free T3 500 Tests	Ea	Reagent
113-3741	LKFT41	Immulite® 1000 Free T4 100 Tests	1/Bx	Reagent
611-1047	LKFS1	Immulite® 1000 FSH-100 Tests	Ea	Reagent
104-8487	LKHPG1	Immulite® 1000 H. Pylori IGG 100 Tests	1/Bx	Reagent
758-3743	LKHB1(D)	Immulite® 1000 HBSAG Domestic-100 Tests	1/Bx	Reagent
969-7213	LKCH(D)	Immulite® 1000 HBSAG Confirm Kit Domestic-20 Tests	1/Bx	Reagent
609-8308	LKCG1	Immulite® 1000 HCG (Total HCG)-100 Tests	Ea	Reagent
113-9787	LKCG5	Immulite® 1000 HCG (Total HCG) 500 Tests	1/Bx	Reagent
114-6281	LKHS1	Immulite® 1000 Herpes I & II IGG-100 Tests	Ea	Reagent
786-7458	LKGRH1	Immulite® 1000 High Rev-100 Tests	100/Kt	Reagent



Immulite Test Menu

Item Code	Manuf Part No.	Item Description	Size	Product Type
114-3320	LKCRP1	Immulite® 1000 High Sensitivity Crp-100 Tests	1/Kt	Reagent
767-9373	LKH01	Immulite® 1000 Homocysteine-100 Tests	Ea	Reagent
114-6283	LKGB1	Immulite® 1000 IGFBP-3-100 Tests	Ea	Reagent
609-8226	LKGF1	Immulite® 1000 IGF-I-100 Tests	1/Bx	Reagent
114-3322	LKGF5	Immulite® 1000 IGF-I 500 Tests	Ea	Reagent
468-6054	LKIN1	Immulite® 1000 Insulin-100 Tests	Ea	Reagent
470-2005	LKIN5	Immulite® 1000 Insulin 500 Tests	Ea	Reagent
234-3045	LKPP1	Immulite® 1000 Intact PTH-100 Tests	1/Bx	Reagent
114-6274	LSKPTZ	Immulite® 1000 Intact PTH Turbo-50 Tests	Ea	Reagent
120-4584	LKLH1	Immulite® 1000 LH-100 Tests	1/Bx	Reagent
228-8873	LKMY1	Immulite® 1000 Myoglobin-100 Tests	1/Bx	Reagent
120-4412	LKOP1	Immulite® 1000 OM-MA (Ca125) 100 Tests	Ea	Reagent
833-6031	LKPA1	Immulite® 1000 Pap 100 Tests	1/Bx	Reagent
833-9217	LKPB1	Immulite® 1000 Phenobarbital-100 Tests	Ea	Reagent
969-0828	LKPN1	Immulite® 1000 Phenytoin-100 Tests	Ea	Reagent
114-3327	LKPG5	Immulite® 1000 Progesterone 500 Tests	Ea	Reagent
949-3217	LKPR1	Immulite® 1000 Prolactin-100 Tests	Ea	Reagent
117-5509	LKPTS1(D)	Immulite® 1000 PSA Domestic-100 Tests	Ea	Reagent
110-8352	LKPD1	Immulite® 1000 Pyrilinks®-D 100 Tests	1/Kt	Reagent
103-0870	LKRT1	Immulite® 1000 Rapid TSH-100 Tests	1/Bx	Reagent
114-3328	LKRT5	Immulite® 1000 Rapid TSG 500t-500 Tests	1/Kt	Reagent
114-2344	LKRUB1	Immulite® 1000 Rubella Quantitative IGG-100 Tests	100/Bx	Reagent
114-3329	LKSH1	Immulite® 1000 SHBG-100 Tests	1/Kt	Reagent
969-7203	LKTN1	Immulite® 1000 Theophylline-100 Tests	Ea	Reagent
115-6477	LKUP5	Immulite® 1000 Third Generation PSA 500 Tests	Ea	Reagent
600-3088	LKTS5	Immulite® 1000 Third Generation TSH 500 Tests	500/Bx	Reagent
104-6370	LKTY1(D)	Immulite® 1000 Thyroglobulin Domestic-100 Tests	100/Bx	Reagent
970-3583	LKTU1	Immulite® 1000 Thyroid Uptake 100 Tests	Ea	Reagent
427-4859	LKIE1	Immulite® 1000 Total IGE-100 Tests	1/Bx	Reagent
112-8510	LKT31	Immulite® 1000 Total T3-100 Tests	Ea	Reagent
970-3585	LKT41	Immulite® 1000 Total T4-100 Tests	Ea	Reagent
984-7023	LKT45	Immulite® 1000 Total T4 500 Tests	1/Bx	Reagent
320-6957	LKTW1	Immulite® 1000 Total Testosterone-100 Tests	100/Kt	Reagent
832-3277	LKHA1	Immulite® 1000 Urinary Albumin-100 Tests	100/Bx	Reagent
355-1212	LKHA5	Immulite® 1000 Urinary Albumin 500 Tests	Ea	Reagent
609-8309	LKVA1	Immulite® 1000 Valproic Acid-100 Tests	100/Bx	Reagent
117-1820	LKVF41	Immulite® 1000 Veterinary Free T4-100 Tests	Ea	Reagent
970-3589	LKVB1	Immulite® 1000 Vitamin B12-100 Tests	Ea	Reagent
114-1451	L2KFT42	Immulite® 2000 Free T4 Kit 200 Tests	Ea	Reagent
114-3361	L2KFT43M	Immulite® 2000 Free T4 Kit 3000 Tests	1/Kt	Reagent
114-3362	L2KFT46	Immulite® 2000 Free T4 Kit 600 Tests	Ea	Reagent
115-7043	L2INZ	Insulin (IM 2000/2500) Sample Diluent	Ea	Control
106-8580	LINZ	Insulin (IM 1000) Sample Diluent	Ea	Control
115-7042	LINCM	Insulin Control Module	Ea	Control
601-0254	LINCVM	Insulin CVF	4/Kt	Control
110-7209	LPPCVM	Intact PTH CVM	Ea	Control
110-1145	LLHCVM	LH (IM1000) CVM	4/Bx	Control
114-3395	L2LHCVM	LH (IM2000 / 2500) CVM	Ea	Control
114-3330	LMH15	Microsampling Insert Set Of 1	1/Kt	Control
114-3337	LSMC	Microsampling Tube Caps	1000/Pk	Control
114-3331	LMST	Microsampling Tubes	1000/Pk	Control
114-0172	L2M1Z	Multi Diluent 1	Ea	Control
117-1455	LMYCVM	Myoglobin CVM	4/Bx	Control
114-4946	LOMZ	OM-MA (IM1000) Sample Diluent	1/Bt	Control
114-5669	L2OPZ	OM-MA (IM2000/2500) Sample Diluent	Ea	Control
108-4553	LOPCVM	OM-MA CVM	4/Kt	Control

SIEMENS

Immulite Test Menu

Item Code	Manuf Part No.	Item Description	Size	Product Type
358-6281	PACM	Pap Control Module	Ea	Control
116-5915	LPACVM	Pap CVM	Ea	Control
116-5285	LPBCVM	Phenobarbital CVM	Ea	Control
114-3334	LPNCVM	Phenytoin CVM	Ea	Control
113-9755	L2KPM	Probe Cleaning Kit	Ea	Control
833-9220	LKPM	Probe Cleaning Kit (Immulite 1000)	Ea	Control
361-4974	LPWS2	Probe Wash (Immulite 1000)	Ea	Control
111-6650	L2PWSM	Probe Wash Module	2 Bt/Bx	Control
110-1146	LPRCVM	Prolactin CVM	4/Bx	Control
114-3399	L2PRZ	Prolactin Sample Diluent	Ea	Control
345-8798	LPRZ	Prolactin Sample Diluent	Ea	Control
472-4320	LPSZ	PSA Sample Diluent	1/Bx	Control
358-3105	LPHCM	PTH Control Module	Bx	Control
114-3398	L2PHZ	PTH Sample Diluent	Ea	Control
489-0065	LPHZ	PTH Sample Diluent 1 x 2 mL	Ea	Control
844-1236	LPDCM	Pyrilinks-D Control Module	Ea	Control
110-8355	LPDCVM	Pyrilinks-D CVM	Ea	Control
110-8353	LPDZ	Pyrilinks-D Sample Diluent	Ea	Control
115-7325	LRUBCM	Rubella IGG Control Module	Ea	Control
114-3335	LRUBCVM	Rubella Quant. IGG CVM	Ea	Control
970-3590	LSCP	Sample Cups (Immulite 1000)	1000/Pk	Control
834-2401	SDCM	Serum Drug Control Module	Ea	Control
114-6803	LSHCM	SHBG Control Module	Ea	Control
114-3336	LSHCVM	SHBG CVM	Ea	Control
116-6721	L2SNCCM	Specific IGE Negative Con	Ea	Control
833-9224	LSUBX	Substrate (Immulite 1000)	1000/Bx	Control
111-6646	L2SUBM	Substrate Module 2000	2 Bt/Bx	Control
799-4964	LT3Z	T3 Sample Diluent	1/Bx	Control
358-3107	LT4Z	T4 Sample Diluent	Ea	Control
114-7337	LTBCM	TBG Control Module	Ea	Control
114-3338	LTBCVM	TBG CVM	Ea	Control
115-7044	L2TBZ	TBG Sample Diluent	Ea	Control
116-5918	LTNCVM	Theophylline CVM	Ea	Control
106-6154	LTGCM	Third Gen. TSH Control Module	2/Bx	Control
108-4563	LUPCVM	Third Generation PSA CVM	4/Kt	Control
114-3339	LYCVM	Thyroglobulin CVM	Ea	Control
114-0751	L2TYZ	Thyroglobulin Sample Diluent 1 x 2.5 mL	Ea	Control
114-0750	LYTZ	Thyroglobulin Sample Diluent 1 x 2.5 mL	Ea	Control
611-1044	LAACM	Thyroid Autoab Control Module	2/Bx	Control
967-4889	IECM	Total IGE Control Module	Ea	Control
114-3318	LIECVM	Total IGE CVM	Ea	Control
113-9754	LT3CVM	Total T3 CVM	Ea	Control
110-7210	LT4CVM	Total T4 CVM	Ea	Control
111-7433	LTWCVM	Total Testosterone CVM	Ea	Control
115-7323	LTXPCM	Toxoplasma IGG Control Module	Ea	Control
115-7039	LTXPCVM	Toxoplasma Quant. IGG CVM	Ea	Control
116-9923	LTICVM	Troponin I (IM1000) CVM	Ea	Control
111-3513	LTIZ	Troponin I Sample Diluent	1/Bx	Control
114-0575	L2TSZ	TSH Sample Diluent	Ea	Control
114-3340	LVACVM	Valproic Acid CVM	Ea	Control
116-1927	L2KVF42	Veterinary Free T4	Ea	Reagent
110-8345	LVBCVM	Vitamin B12 (IM1000/2000) CVM	Ea	Control
627-3555	LVBZ	Vitamin B12 Sample Diluent	Ea	Control
833-9222	LKWT	Water Test	Ea	Control



SIEMENS

Immolute Accessories

Item Code	Manuf Part No.	Item Description	Size	Product Type
117-0781	400909-01	ADVIA Automation Direct Water Feed Kit	1/Kt	Accs
116-3125	664-1962-01	Centralink - Additional 8 Port Terminal Server	Ea	Accs
116-3104	084-0028-01	Centralink - Backup License - ADVIA 120 Instrument Connection	Ea	Accs
116-3102	084-0026-01	Centralink - Backup License - Automation Module	Ea	Accs
116-3101	084-0024-01	Centralink - Backup License - Client User Interface	Ea	Accs
116-3103	084-0027-01	Centralink - Backup License - Non-ADVIA 120 Instrument Connection	Ea	Accs
116-3100	084-0021-01	Centralink - Backup License - Server With 1 Client User Interface	Ea	Accs
116-3123	664-0679-01	Centralink - Barcode Scanner: USB Port CCD Scanner LGP6125RU1S-000	Ea	Accs
116-3095	084-0004-02	Centralink - Client User Interface License	Ea	Accs
116-3124	664-1782-01	Centralink - Printer: Lexmark E460DN Laser Jet, 120VAC	Ea	Accs
116-3121	615-4188-01	Centralink - RJ45 - RJ45 Cable (10)	2/Pk	Accs
116-3122	615-4188-02	Centralink - RJ45 - RJ45 Cable (25')	Ea	Accs
116-3118	615-4187-01	Centralink - RJ45 To 9 Pin DB Adapter (LIS)	Ea	Accs
116-3115	515-7188-03	Centralink - UPS: PowerVar #ABCE1100-11, 120VAC, 1100VA	Ea	Accs
115-2466	400930-02	Imm. Allergen Wedges Bar Coded 1 - 33	33/Bx	Accs
115-2467	400930-03	Imm. Allergen Wedges Bar Coded 34 - 66	33/Bx	Accs
115-2468	400930-04	Imm. Allergen Wedges Bar Coded 67 - 99	33/Bx	Accs
114-3844	422710-06	Kit, INT Wheel, XP, Left, G801 Up	Ea	Accs
114-3840	422710-02	Kit, INT Wheel, XP, Rear, G800 Down	Ea	Accs
114-3843	422710-05	Kit, INT Wheel, XP, Rear, G801 Up	Ea	Accs
114-3839	422710-01	Kit, INT Wheel, XP, Right, G800 Down	Ea	Accs
114-3842	422710-04	Kit, INT Wheel, XP, Right, G801 Up	Ea	Accs
114-3841	422710-03	Kit, INT Wheel, XP, Left, G800 Down	Ea	Accs
117-5920	422761-02	Versacell ADVIA Centaur/ ADVIA 1800/ Exl/Rxl 800 Down Wheel- Back	Ea	Accs
117-3511	422711-01	Versacell Instrument Connection Kit A1800	Ea	Accs
117-3556	422711-02	Versacell Instrument Connection Kit Centaur	Ea	Accs
119-7009	422711-04	Versacell Instrument Connection Kit EXL	Ea	Accs
117-2061	422711-03	Versacell Instrument Connection Kit RXL	Ea	Accs
117-2063	422761-05	Versacell Interface Wheel, Back, G801 Up	Ea	Accs
117-3510	422761-06	Versacell Interface Wheel, Left, G801 Up	Ea	Accs
117-3557	422761-04	Versacell Interface Wheel, Right, G801 Up	Ea	Accs
117-2185	422729-02	Versacell Kit, IMM 2000/2500, Back	Ea	Accs
117-2042	422729-03	Versacell Kit, IMM 2000/2500, Left, G800 Down	Ea	Accs
117-2121	422729-01	Versacell Kit, IMM 2000/2500, Right	Ea	Accs
114-3835	422709-03	Versacell Kit, INT Wheel, A1800, Left, G800 Down	Ea	Accs
114-3838	422709-06	Versacell Kit, INT Wheel, A1800, Left, G801 UP	Ea	Accs
114-3834	422709-02	Versacell Kit, INT Wheel, A1800, Rear, G800 Down	Ea	Accs
114-3837	422709-05	Versacell Kit, INT Wheel, A1800, Rear, G801 UP	Ea	Accs
114-3833	422709-01	Versacell Kit, INT Wheel, A1800, Right, G800 Down	Ea	Accs
114-3836	422709-04	Versacell Kit, INT Wheel, A1800, Right, G801 UP	Ea	Accs

Immolute Instruments

Item Code	Manuf Part No.	Item Description	Size	Product Type
116-3109	084-K029-02	Centralink - ADVIA 2120/120 Connection Kit	Ea	Ins
116-3106	084-B003-09	Centralink - ADVIA Centralink High-end Server Assembly	Ea	Ins
116-3111	084-K040-09	Centralink - Configured High-end Server Kit	Ea	Ins
116-3117	515-7245-05	Centralink - Dell 22 Flat Panel LCD Wide Screen Monitor	Ea	Ins
116-3113	084-K041-01	Centralink - Lab Automation Kit	Ea	Ins
116-3114	084-K042-01	Centralink - Non ADVIA 2120/120 Connection Kit	Ea	Ins
116-3107	084-K026-05	Centralink ADVIA Centralink Client Kit	Ea	Ins
114-3635	030006-03	Immolute 1000	Ea	Ins
114-3636	030006-03R	Immolute 1000 (Requalified)	Ea	Ins
840-9390	030002-13	Immolute 2000	Ea	Ins
840-9399	030002-13R	Immolute 2000 (Requalified)	Ea	Ins
116-8979	030002-33	Immolute 2000 XPI	Ea	Ins
114-3637	030007-13	Immolute 2500	Ea	Ins

SIEMENS

Advia Centaur® XP Immunoassay System

This high-performance system has more onboard reagents and dedicated STAT capabilities to maximize productivity, regardless of volume or types of tests. Offers "always ready" continuous operation without interruption.

System offers:

- High throughput of up to 240 tests per hour
- Intuitive software for ease of operation
- STAT port guarantees priority sampling at any time
- Data archive reduces administrative tasks
- Extensive menu with 30 onboard reagents

Powerful productivity:

- Data archive feature reduces time for administrative tasks
- Universal sample rack design eliminates manual tasks

Excellent patient care:

- Automatic dilution, repeat, and reflex testing
- Smart Algorithm Software automatically repeats and confirms reactive testing

Flexible solutions:

- Connects to Siemens Aptio™ Automation, Advia® Automation, StreamLAB® Automation, and VersaCell® Solutions
- Siemens Remote Service (SRS) proactively maintains Siemens systems

#078-A011-03

(116-9224).....ea



Advia Centaur Test Menu

Item Code	Manuf Part No.	Item Description	Size	Product Type
116-8008	156-B849-01	13MM EZEE Nest Sample Cups	1000/Bx	Control
114-3704	112219	Acid/Base RGT 1 & 2	1/Cr	Control
114-3689	110763	AFP 100 Tests	1/Bx	Reagent
114-3690	110764	AFP 500 Tests	1/Bx	Reagent
114-3763	124975	Ascorbic Acid/DLL	1/Cr	Control
117-8258	10630917	ATG - Optimized Control Set	Ea	Control
117-4902	125747	ATG 100T Kit	100/Bx	Reagent
117-8249	10492398	ATG 100T Kit Optimized	Ea	Reagent
114-3775	125748	ATG 500T Kit	1/Bx	Reagent
117-8250	10492399	ATG 500T Kit Optimized	Ea	Reagent
114-3777	127636	ATG Diluent Ancillary	2/Cr	Control
114-3780	128280	ATG Diluent, LBLD	2/Cr	Control
117-6383	10630936	ATPO - optimized 1, 2 QC Material	Ea	Control
117-4891	125745	ATPO 100 Tests	100/Bx	Reagent
117-8244	10630886	ATPO 100 Tests Optimized	Ea	Reagent
114-3773	125746	ATPO 500 Tests	1/Bx	Reagent
117-6381	10630887	ATPO 500 Tests Optimized	1/Bx	Reagent
114-3770	125733	ATPO Diluent	2/Cr	Control
117-4898	126076	ATPO QC Kit	Ea	Control
114-3712	116735	BR 250 Tests	1/Bx	Reagent
114-3711	116734	BR 50 Tests	1/Bx	Reagent
114-3713	116746	BR PRETRT	2/Cr	Control
114-3787	128532	CA 125 II 100 Tests	1/Bx	Reagent
114-3788	128533	CA 125 II 500 Tests	1/Bx	Reagent
114-3789	128534	CA 15-3 100 Tests	1/Bx	Reagent



SIEMENS

Advia Centaur Test Menu

Item Code	Manuf Part No.	Item Description	Size	Product Type
114-3790	128535	CA 15-3 500 Tests	1/Bx	Reagent
114-3786	128531	CA 15-3 MCM	1/Cr	Control
114-3751	123521	CA 19-9 250 Tests	1/Bx	Reagent
117-0545	10491244	CA 19-9 250 Tests Revised	1/Cr	Reagent
114-3750	123519	CA 19-9 50 Tests	1/Bx	Reagent
117-0546	10491379	CA 19-9 50 Tests Revised	1/Cr	Reagent
114-3752	124061	CA 19-9 Diluent	1/Cr	Control
117-0550	10491974	CA 19-9 Diluent 25 mL Vial	1/VI	Control
114-3759	124615	CA 19-9 Diluent	2/Cr	Control
117-0549	10491972	CA 19-9 Diluent	2/Pk	Control
114-3791	128536	CAL CA 125	2/Cr	Control
117-8255	10630915	CAL 1 Optimized	2/Pk	Control
117-4904	127657	CAL 1	2/Pk	Control
114-3781	128304	CAL 27	2/Cr	Control
114-3654	109562	CAL 6 (E2-6 II)	2/Cr	Control
114-3655	109563	CAL 6 (E2-6 II)	6/Cr	Control
114-3784	128433	CAL 66 (INS)	2/Cr	Control
114-3793	129021	CAL 67 (CPEP-S)	2/Cr	Control
118-4468	10492491	Cal 80	2/Pk	Control
114-3757	124579	CAL 88 2-Pack	2/Cr	Control
117-1120	672172000	CAL C 2-Pack	2/Pk	Control
114-3899	672182000	CAL C 6-Pack	6/Cr	Control
119-3052	128538	CAL CA 15-3 2-Pack	2/Pk	Control
116-4502	672173000	CAL D 2-Pack	Ea	Control
114-3900	672183000	CAL D 6-Pack	6/Cr	Control
114-3708	112776	CAL G (BR) 2-Pack	2/Cr	Control
114-3709	112777	CAL G (BR) 6-Pack	6/Cr	Control
114-3896	672174000	CAL K (CKMB) 2-Pack	2/Pk	Control
114-3897	672179000	CAL K (CKMB) 6-Pack	6/Cr	Control
114-3884	570403	CAL N (PHTN) 2-Pack	2/Cr	Control
114-3885	570404	CAL N (PHTN) 6-Pack	6/Cr	Control
117-6382	10630932	CAL O - Optimized 2-Pack	2/Pk	Control
117-4895	125737	CAL O 2-Pack	2/Pk	Control
114-3728	118220	CAL Q (PSA) 2-Pack	2/Cr	Control
114-3729	118221	CAL Q (PSA) 6-Pack	6/Cr	Control
114-3703	111083	CAL U 2-Pack	2/Cr	Control
114-3881	570061	CAL V (VANC) 2-Pack	2/Pk	Control
114-3882	570062	CAL V (VANC) 6-Pack	6/Cr	Control
114-3764	125014	CAL Y (CPSA) 2-Pack	2/Cr	Control
117-0541	10379810	Calibrator 30 2-Pack	2/Pk	Control
117-0542	10379811	Calibrator 30 6-Pack	6/Pk	Control
114-3649	104474	Cardiac Markers 1, 2, 3 F/G	1/Cr	Control
114-3687	110761	CEA 100 Tests	1/Bx	Reagent
114-3688	110762	CEA 500 Tests	1/Bx	Reagent
117-1123	672255000	CEA Diluent 10 mL vial	1/VI	Control
114-3665	110322	CEA Diluent 2-Pack	2/Cr	Control
117-1126	672427000	CEA MCM Set	Ea	Control
114-3693	110773	CKMB 100 Tests	1/Kt	Reagent
114-3694	110774	CKMB 500 Tests	1/Kt	Reagent
114-3903	672230000	CKMB Diluent - 10 mL Vial	2/Cr	Control
114-3706	112748	CLN SOLN 12-Pack	12/Cr	Control
114-3794	129026	CPEP-S 100 Tests	1/Kt	Reagent
114-3760	124830	CPSA 100 Tests	1/Bx	Reagent
114-3768	125425	CPSA Control Set	1/Cr	Control
114-3765	125308	CPSA PRTRT RGT 1-Pack	1/Cr	Control

SIEMENS

Advia Centaur Test Menu

Item Code	Manuf Part No.	Item Description	Size	Product Type
114-3714	116993	CTNI 100 Tests	1/Bx	Reagent
114-3734	119714	DGTN 250 Tests	1/Bx	Reagent
114-3731	118453	DGTN 50 Tests	1/Bx	Reagent
114-3692	110772	DIG 250 Tests	1/Bx	Reagent
114-3691	110771	DIG 50 Tests	1/Bx	Reagent
114-3702	110785	E2-6 500 Tests	1/Bx	Reagent
117-0548	10491910	eE2 Diluent 1 x 25 mL	Ea	Control
117-3133	10491878	eE2 Diluent 2 x 5 mL Pack	Ea	Control
117-0544	10490889	Enhanced Estradiol 100 Tests	1/Kt	Reagent
117-0547	10491445	Enhanced Estradiol 500 Tests	1/Kt	Reagent
114-3678	110746	FER 250 Tests	1/Bx	Reagent
114-3677	110745	FER 50 Tests	1/Bx	Reagent
114-3733	118551	FOL BA 100 Tests	1/Bx	Reagent
114-3761	124838	FOL BA 500 Tests	1/Bx	Reagent
114-3769	125500	FOL BA Diluent 2-Pack	2/Cr	Control
114-3762	124974	FOL BA Diluent LBLD	1/Cr	Control
114-3652	104533	FRT4 MCM Set	1/Cr	Control
114-3683	110755	FSH 100 Tests	1/Bx	Reagent
114-3684	110756	FSH 500 Tests	1/Bx	Reagent
114-3737	119781	FT3 V2 250 Tests	1/Bx	Reagent
114-3735	119780	FT3 V2 50 Tests	1/Bx	Reagent
117-8290	6981117	HBeAg QC Pack	Ea	Control
114-3754	124553	HCY DIL LBLD	1/Cr	Control
114-3755	124577	Homocysteine 100 Tests	1/Bx	Reagent
114-3756	124578	Homocysteine 500 Tests	1/Bx	Reagent
114-3758	124581	Homocysteine Diluent	1/Cr	Control
114-3696	110778	IGE 250 Tests	1/Bx	Reagent
114-3695	110777	IGE 50 Tests	1/Bx	Reagent
114-3666	110324	IGE Diluent 2-Pack	2/Cr	Control
114-3785	128434	Insulin 100 Tests	1/Bx	Reagent
114-3798	129883	Insulin Diluent 2-Pack	2/Cr	Control
117-0551	10492368	Intact PTH2 100 Tests	1/Bx	Reagent
117-0552	10492369	Intact PTH2 500 Tests	1/Bx	Reagent
117-0553	10492378	IPTH 1,2,3 QC material	1/Bx	Control
114-3795	129501	IPTH 100 Tests	1/Bx	Reagent
114-3796	129502	IPTH 500 Tests	1/Bx	Reagent
114-3797	129863	IPTH Control Set	1/Cr	Control
117-0555	10492394	IPTH2 Calibrator 2-Pack	2/Pk	Control
117-0554	10492387	IPTH2 Calibrator 6-Pack	6/Pk	Control
114-3643	078-K139-01	Kit Sample Tips (6480/Pkg)	6480/Pk	Control
114-3642	078-K138-01	Kit, Cuvette 3000 Pack	3000/Pk	Control
114-3641	078-K137-01	Kit, Sample Cup 1500 Pack	1500/Pk	Control
114-3653	106382	Label, BC, Cardiac Marker	1/Cr	Control
114-3930	986400	Label, BC, Ligand Plus	1/Cr	Control
114-3710	114925	LBL, BC, TMR Marker Plus 1, 2 Intl	1/Cr	Control
114-3682	110754	LH 300 Tests	1/Bx	Reagent
114-3681	110753	LH 60 Tests	1/Bx	Reagent
114-3929	986000	LIGAND PLUS 1, 2, 3 F/G	1/Cr	Control
117-1121	672177000	Multi-Diluent 1 50 ml vial	Ea	Control
114-3656	110312	Multi-Diluent 1 2-Pack	2/Cr	Control
114-3657	110313	Multi-Diluent 1 6-Pack	6/Cr	Control
117-0570	112370	Multi-Diluent 10 10 mL Vial	Ea	Control
114-3719	117227	Multi-Diluent 10 2-Pack	2/Cr	Control
114-3658	110314	Multi-Diluent 2 2-Pack	2/Cr	Control



SIEMENS

Advia Centaur Test Menu

Item Code	Manuf Part No.	Item Description	Size	Product Type
114-3659	110315	Multi-Diluent 3	2/Cr	Control
114-3660	110316	Multi-Diluent	2/Cr	Control
114-3661	110317	Multi-Diluent 5	1/Cr	Control
114-3883	570099	Multi-Diluent 7	1/Cr	Control
114-3662	110318	Multi-Diluent 7	2/Cr	Control
114-3648	100438	Multi-Diluent 9	1/Cr	Control
114-3720	117228	Multi-Diluent 11	2/Cr	Control
114-3716	117031	MYO 250 Tests	1/Bx	Reagent
114-3715	117030	MYO 50 Tests	1/Bx	Reagent
114-3743	120539	PHTN 250 Tests	1/Bx	Reagent
114-3742	120538	PHTN 50 Tests	1/Bx	Reagent
114-3766	125362	PRGE 250 Tests	1/Bx	Reagent
114-3732	118530	PRGE 50 Tests	1/Bx	Reagent
114-3686	110758	PRL 250 Tests	1/Bx	Reagent
114-3685	110757	PRL 50 Tests	1/Bx	Reagent
114-3725	118156	PSA 100 Tests	1/Bx	Reagent
114-3726	118157	PSA 500 Tests	1/Bx	Reagent
114-3721	117710	Rubella IGG 100 Tests	1/Bx	Reagent
114-3722	117711	Rubella IGG QC Kit	1/Cr	Control
114-3723	117733	Rubella IGM 50 Tests	1/Bx	Reagent
114-3724	117734	Rubella IGM QC Kit	1/Cr	Control
114-3894	672010000	Sample Cup Caps	1000/Pk	Control
114-3638	076-B007-01	Septum Caps	1/Cr	Control
114-3667	110328	T3/T4/VB12 ANC RGT	2/Cr	Control
114-3668	110329	T3/T4/VB12 ANC RGT	6/Cr	Control
114-3671	110735	T4 100 Tests	1/Bx	Reagent
114-3672	110736	T4 500 Tests	1/Bx	Reagent
114-3901	672202000	T4 Diluent	1/Cr	Control
114-3717	117143	THCG Diluent	2/Cr	Control
114-3718	117157	THCG Diluent	6/Cr	Control
117-0574	115249	THCG Diluent 50 mL Vial	1/VI	Control
114-3783	128380	THEO 2 250 Tests Kit	1/Bx	Reagent
114-3782	128379	THEO 2 50 Tests Kit	1/Bx	Reagent
114-3748	121431	TOBR 250 Tests	1/Bx	Reagent
114-3747	121430	TOBR 50 Tests	1/Bx	Reagent
114-3739	120154	TOXO IGG 100 Tests	1/Bx	Reagent
114-3740	120382	TOXO IGG QC Kit	1/Cr	Control
114-3738	120153	TOXO IGM 50 Tests	1/Bx	Reagent
114-3741	120383	TOXO IGM QC Kit	1/Cr	Control
114-3669	110732	TSH 100 Tests	1/Bx	Reagent
114-3670	110734	TSH 500 Tests	1/Bx	Reagent
114-3904	672422000	TSH MCM Set	1/Cr	Control
114-3697	110780	TSH3 100 Tests	1/Bx	Reagent
114-3698	110781	TSH3 500 Tests	1/Bx	Reagent
114-3767	125363	TSTO 250 Tests	1/Bx	Reagent
114-3699	110782	TSTO 50 Tests	1/Bx	Reagent
114-3928	976700	Tumor Marker Plus 1, 2	1/Cr	Control
114-3676	110744	TUP 250 Tests	1/Bx	Reagent
114-3675	110743	TUP 50 Tests	1/Bx	Reagent
114-3746	121197	Vancomycin 250 Tests	1/Bx	Reagent
114-3730	118451	Vancomycin 50 Tests	1/Bx	Reagent
114-3679	110747	VB12 100 Tests	1/Bx	Reagent
114-3680	110748	VB12 500 Tests	1/Bx	Reagent
117-0556	105441	VB12 Diluent 15 mL Vial	1/VI	Control
114-3664	110321	VB12 Diluent	2/Cr	Control

SIEMENS

Advia Centaur Test Menu

Item Code	Manuf Part No.	Item Description	Size	Product Type
114-3650	104523	VB12 MCM Set	1/Cr	Control
114-3744	121093	VB12 RA Kit (2 Btls)	2bt/Cr	Control
117-3328	10493589	Vitamin D Calibrator	2pk	Control
117-3327	10630911	Vitamin D Calibrator	6pk	Control
117-3325	10632114	Vitamin D Diluent 25 mL Bottle	25ml	Control
117-3324	10494100	Vitamin D Diluent	2pk	Control
117-3326	10632229	Vitamin D QC Pack	1pack	Control
117-3322	10491994	Vitamin D Total 100 Tests	100 Test	Reagent
117-3323	10631021	Vitamin D Total 500 Tests	500 Test	Reagent
114-3705	112351	Wash 1 Reagent	1/Cr	Control

Advia Centaur Accessories

Item Code	Manuf Part No.	Item Description	Size	Product Type
117-0525	086-0094-01	ADVIA Centaur CP Ancillary Tray	Ea	Accs
117-1492	086-K027-01	Centaur CP Direct Plumbing Kit	Ea	Accs
118-2605	078-K190-02	Centaur XP Direct Plumbing Kit	Ea	Accs
116-3125	664-1962-01	CentraLink - Additional 8 Port Terminal Server	Ea	Accs
116-3123	664-0679-01	CentraLink - Barcode Scanner: USB Port CCD Scanner LGP6125RU1S-000	Ea	Accs
116-3095	084-0004-02	CentraLink - Client User Interface License	Ea	Accs
116-3124	664-1782-01	CentraLink - Printer: Lexmark E460DN Laser Jet, 120VAC	Ea	Accs
116-3121	615-4188-01	CentraLink - RJ45 - RJ45 Cable (10) Qty 2	2/Pk	Accs
116-3122	615-4188-02	CentraLink - RJ45 - RJ45 Cable (25')	Ea	Accs
116-3118	615-4187-01	CentraLink - RJ45 to 9 pin DB adapter (LIS)	Ea	Accs
116-3115	515-7188-03	CentraLink - UPS: PowerVar #ABCE1100-11, 120VAC, 1100VA	Ea	Accs
114-3844	422710-06	Kit, Int Wheel, XP, Left, G801 Up	Ea	Accs
114-3840	422710-02	Kit, Int Wheel, XP, Rear, G800 Down	Ea	Accs
114-3843	422710-05	Kit, Int Wheel, XP, Rear, G801 Up	Ea	Accs
114-3839	422710-01	Kit, Int Wheel, XP, Right G800 Down	Ea	Accs
114-3842	422710-04	Kit, Int Wheel, XP, Right, G801 Up	Ea	Accs
114-3841	422710-03	Kit, Int Wheel, XP, Left, G800 Down	Ea	Accs
114-3817	263100.202	StreamLAB Centrifuge UPS	Ea	Accs
117-5920	422761-02	VersaCell ADVIA Centaur/ ADVIA 1800/ ExL/RxL 800 Down Wheel-Back	Ea	Accs
117-3511	422711-01	VersaCell Instrument Connection Kit A1800	Ea	Accs
117-3556	422711-02	VersaCell Instrument Connection Kit Centaur	Ea	Accs
119-7009	422711-04	VersaCell Instrument Connection Kit EXL	Ea	Accs
117-2063	422761-05	VersaCell Interface Wheel, Back, G801 Up	Ea	Accs
117-3510	422761-06	VersaCell Interface Wheel, Left, G801 Up	Ea	Accs
117-3557	422761-04	VersaCell Interface Wheel, Right, G801 Up	Ea	Accs
117-2185	422729-02	VersaCell Kit, IMM 2000/2500, Back	Ea	Accs
117-2042	422729-03	VersaCell Kit, IMM 2000/2500, Left, G800 Down	Ea	Accs
117-2121	422729-01	VersaCell Kit, IMM 2000/2500, Right	Ea	Accs
114-3835	422709-03	VersaCell Kit, INT Wheel, A1800, Left, G800 Down	Ea	Accs
114-3838	422709-06	VersaCell Kit, INT Wheel, A1800, Left, G801 Up	Ea	Accs
114-3834	422709-02	VersaCell Kit, INT Wheel, A1800, Rear, G800 Down	Ea	Accs
114-3837	422709-05	VersaCell Kit, INT Wheel, A1800, Rear, G801 Up	Ea	Accs
114-3833	422709-01	VersaCell Kit, INT Wheel, A1800, Right, G800 Down	Ea	Accs
114-3836	422709-04	VersaCell Kit, INT Wheel, A1800, Right, G801 Up	Ea	Accs



SIEMENS

Advia Centaur Instruments

Item Code	Manuf Part No.	Item Description	Size	Product Type
114-3646	086-0287-01	7' Mobile Table / Workbench	Ea	Ins
117-0474	040-A001-08	ADVIA Automation WorkCell A - 60Hz	Ea	Ins
117-0475	040-A001-09	ADVIA Automation WorkCell B - 60Hz	Ea	Ins
117-0476	040-A001-10	ADVIA Automation WorkCell C - 60Hz	Ea	Ins
114-3645	086-A002-01	ADVIA Centaur CP	Ea	Ins
117-5460	078-A011R03	ADVIA Centaur XP (refurb)	Ea	Ins
116-3109	084-K029-02	CentraLink - ADVIA 2120/120 Connection Kit	Ea	Ins
116-3106	084-B003-09	CentraLink - ADVIA CentraLink High-end Server Assembly	Ea	Ins
116-3111	084-K040-09	CentraLink - Configured High-end Server Kit	Ea	Ins
116-3117	515-7245-05	CentraLink - Dell 22 Flat Panel LCD Wide Screen Monitor	Ea	Ins
116-3113	084-K041-01	CentraLink - Lab Automation Kit	Ea	Ins
116-3114	084-K042-01	CentraLink - Non ADVIA 2120/120 Connection Kit	Ea	Ins
116-3107	084-K026-05	CentraLink ADVIA CentraLink Client Kit	Ea	Ins
117-2061	422711-03	VersaCell Instrument Connection Kit R&L	Ea	Ins

HENRY SCHEIN®
MEDICAL

MEDICAL WEB SITE

**Quick &
Easy Ordering!**

1

Login

2

Order

3

Checkout

Which dynamic features do the Henry Schein Medical Web sites offer your practice?

- 24-hour/365-day availability to supply ordering
- Track orders and shipments and process returns
- Generate comprehensive utilization reports via our dynamic report wizard
- Set budgets and electronically approve orders for greater security and accuracy
- Access purchase history and real-time pricing
- Custom shopping list capability for easy ordering
- Paperless billing statement option
- Shop by procedure, chronic conditions, or simple, everyday essentials
- Access free MSDS documents for quick OSHA compliance
- Retrieve electronic statements and invoices for a paperless environment
- Access to promotions created uniquely for your specialty



www.henryschein.com

Call:
1-800-P-SCHEIN
(1-800-772-4346)
8am-9pm (et)

Fax:
1-800-329-9109
24 Hours

www.henryschein.com/medical

Due to manufacturer price increases, prices are subject to change without notification. Offer valid on our full array of products. Upon any request by Medicare, Medicaid, Tricare or other payer, you must report the value of any benefit received from a discount program (e.g., points redemptions, gifts or other special awards). ©2013 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors. No adjustments from prior sales. Offer not to be combined with other promotions and/or special contract pricing. Void where prohibited. Local regulatory requirements may apply to use or installation of certain products. It is your responsibility to understand and comply with any such requirements prior to purchase, use or installation of products. †For informational purposes only. Source: ReimbursementCodes.com and/or www.cms.gov. Customer is responsible for verification of billing/coding in accordance with applicable specific circumstances.

13MS8011

 **HENRY SCHEIN®**
Henry Schein, Inc. • 135 Duryea Road • Melville, NY 11747

PRSRT STD
U.S. POSTAGE
PAID
HENRY SCHEIN, INC.

CLINICAL LABORATORY



Provide convenient and quality patient care in your practice while you also increase profitability.

Henry Schein Medical's LabTeam™ is a single source for advice within the laboratory setting. With specific expertise in lab consultation and analysis, the LabTeam™ improves your patients' experiences in the office, while navigating profitability and regulatory requirements.



FEATURES:

- Consultation • Design
- Financing • Selection
- Implementation • Support/Training

 **HENRY SCHEIN®**
LabTeam™
SOLUTIONS FOR CLINICAL LABS

BENEFITS:

- Improve patient satisfaction with convenient onsite testing
- Expedite ability to diagnose with rapid turnaround of results
- New or expanded sources of revenue
- Stay current on technology/reimbursement /regulatory changes
- Custom laboratory management solutions for cost-effective operations

Contact Your Lab Consultant
1.800.428.9200
E-mail: labquestions@henryschein.com
Visit: www.henryschein.com/labteam

 **HENRY SCHEIN®**
MEDICAL

 **Rely on Us**