

SINGLE-USE AEROSOL EVACUATOR

A simple solution for aerosol control to protect your patients and staff.



✓ Disposable, cost-effective device for aerosol reduction in the operator.

✓ Use during hygiene, restorative, oral surgery, and other procedures.

✓ Hands-free, ease-of-use for both the hygienist and dentist.

✓ Use in conjunction with saliva ejector.



SINGLE-USE AEROSOL EVACUATOR

NEW

The Single-Use Aerosol Evacuator is a disposable device that attaches to any standard HVE connection to provide an effective level of aerosol reduction. The extraoral product placement allows for hands-free suction of aerosols at the exit of the oral cavity, delivering comfort for the patient and ease-of-use for the hygienist or dentist. This low-profile, single-piece device provides minimal interference at the procedure location, especially when compared with intraoral options.

Advantages

- Proprietary patented design.
- Universal size with malleable "wings" that comfortably adjust to patient's mouth.
- Affordable aerosol control.
- Low-profile provides unobstructed view of the procedure area.

Specs

- Material: High-density polyethylene (HDPE)
- Color: White
- Quantity: 50/pack
- Dimensions: 4.25" x 3" (L x W)
- Use with saliva ejector



1. Internal Wings
2. External Wings
3. Lower Teeth Bumper
4. Suction Inlet
5. Valve Stem

Ordering Information

Single-Use Aerosol Evacuator
Item No.: (570-2441)

Complementary Products



Maxima PRO 2LART™ Optic (570-2094)
Maxima PRO 2SAR™ Non-Optic (570-2095)

Available through

HENRY SCHEIN®
SPECIAL MARKETS

www.henryschein.com/sm

For more information, please see product Instructions for Use.

